Kit of reagents for the execution of Voges-Proskauer test

DESCRIPTION

VP (NaOH) Reagent is a kit containing the reagents for the execution of Voges-Proskauer test used for identification of microorganisms performing 2-3 butilneglicole fermentation.

CONTENT OF THE KIT

Each package contains 2 DROPPERS and the follow reagents:

- ALPHA NAPHTHOL 5 Vials x 10 ml
- NaOH 40% 5 Vials x 10 ml

METHOD PRINCIPLE

Voges-Proskauer reaction distinguishes organisms that ferment glucose to acids from those that ferment to acetoin (2-3 butilenglicole fermentation). After 48 hours incubation ALPHA NAPHTHOL and NaOH 40% reagents are added to sample. Positive reaction to Voges-Proskauer test is highlighted by the formation of a red coloration obtained from acetyl-methyl-carbinol oxidation in the presence of alkali and oxygen.

REAGENTS TABLE

Reagents	Composition
ALPHA NAPHTHOL	Alpha Naphthol 6 g, Ethyl Alcohol 100 ml
NaOH 40%	Sodium Hydroxide 40 g, Distilled Water 100 ml

TEST PROCEDURE AND RESULTS INTERPRETATION

The VP test is performed on cultures of 48 hours of incubation prepared dissolving a pure culture in a tube of MR-VP BROTH (ref. 610032 or 620032) properly prepared, or in a tube of MR-VP MEDIUM (ref. 20149). Add 3 drops of ALPHA NAPHTHOL and 2 drops of NaOH 40% solution to 3 mL of broth culture. The positive test is given by the development of a bright red colour 15 minutes after addition of reagents.

QUALITY CONTROL FOR THE USER

Visual inspection

NaOH 40% is a colourless clear solution;
ALPHA NAPHTHOL is a beige colour clear solution.

Microbiology control

Each batch of V.P. TEST Reagent is subjected to quality control using a bacterial culture of *Escherichia coli* ATCC 25922 as negative control and one of *Enterobacter cloacae* ATCC 13047 as positive control.

PRECAUTIONS

VP (NaOH) Reagent is classified as hazardous under current legislation; for its use it is recommend to consult the MSDS. The product is a disposable device to be used only for in vitro diagnostic use, is intended for professional use and should be used in the laboratory by properly trained operator, with approved methods of aseptic and safety against pathogens.

STORAGE

Store at 2-8°C away from light in its original package, until the expiry date shown on the label. Eliminate if signs of deterioration or contamination are evident.

ELIMINATION OF THE USED MATERIAL

After use, VP (NaOH) Reagent and the material that comes in contact with the sample must be decontaminated and disposed in accordance with the techniques used in the laboratory for decontamination and disposal of potentially infectious materials.

BIBLIOGRAPHY

- Edwin H. Lenette: *Manual of Clinical Microbiology* (1995).
- Blazevic, D.J., and Ederer, G.M.: *Principles of biochemical tests in diagnostic microbiology*. 63-67. New York, John Wiley & Sons, 1975

PRESENTATION

Product	REF	
VP (NaOH) Reagent	80280	100 tests

TABLE OF SYMBOLS

LOT	Batch code	IVD	In vitro Diagnostic Medical Device	Manufacturer	Use by	Fragile, handle with care
REF	Catalogue number	Temperature limitation	Contains Sufficient for <n> tests	Consult instructions for use	Do not reuse	

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LIOFILCHEM® S.r.l.