iQ Precision Kit

REF 800-3014

TABLE OF CONTENTS

English	2
Français (FR)	4
Deutsch (DE)	7
Español (ES)	10
Português Portugal (PT-PT)	12
中文 (ZH-CN)	15
Magyar (HU)	17
Polski (PL)	20
Slovenčina (SK)	23
Türkçe (TR)	26
Hrvatski (HR)	28
Български (BG)	31
Србија (SR)	34
Português Brasil (PT-BR)	36
Tiếng Việt (VI)	39

iQ Precision Kit

REF 800-3014

Rx Only in the U.S.A.

INTENDED USE

For the verification of the iQ200 analyzer at the time of installation or during routine service visits.

For *in vitro* Diagnostic Use with iQ200 Series Analyzers.

PRINCIPLE

iQ Mid Control and iQ Positive Control are suspensions of fixed human red blood cells in a buffered, isotonic balanced solution. The iQ Mid Control is used in conjunction with the iQ Negative Control and iQ Positive Control in separate tubes arranged in the QC tube rack prior to iQ200 Automated Urinalysis System aspiration. Within run and total imprecision using the Positive, Negative and Mid controls are based upon CLSI EP-10 Guidelines which uses the mean, variance and coefficient of variation of the three level control runs to evaluate imprecision and carryover of each individual instrument.

WARNING AND PRECAUTIONS

- The iQ Precision Kit has been determined to be nonhazardous under guidelines issued by OSHA in 29 CFR 1910.1200(d).
- Appropriate precautions in the collection, handling, storage and disposal of specimens should be exercised accordingly.
- Wear appropriate personal protective equipment when handling the materials.
- Avoid skin and eye contact. Refer to the Safety Data Sheet for first aid measures.
- Consult local or state authorities with regard to recommended method of disposal.

GHS HAZARD CLASSIFICATION

iQ Positive Control
iQ Mid Control

WARNING

H317
H361
P201
P280

P308+P313
P333+P313

P362+P364

May cause an allergic skin reaction.
Suspected of damaging fertility or the unborn child.
Obtain special instructions before use.
Wear protective gloves, protective clothing and eye/face protection.
IF exposed or concerned: Get medical advice/attention.
If skin irritation or rash occurs: Get medical advice/attention.
Take off contaminated clothing and wash it before use.
N-(3-Chloroallyl)hexaminium Chloride 0.1 - 0.5%

SDS Safety Data Sheet is available at techdocs.beckmancoulter.com

NON-REACTIVE INGREDIENTS

iQ Negative Control: Water, Salts and Preservatives
iQ Positive Control, iQ Mid Control: Human erythrocytes, water, salts, dyes and preservatives. The stabilized human erythrocytes are <1% in aqueous suspension and are non-reactive for HbsAg, HCV 3.0 Ab, HIV 1/2, HIV-1 p24 Ag, RPR and ALT <82 IU/mL. Stabilized human erythrocytes are chemically fixed, which minimizes any hazardous characteristics.

STORAGE AND STABILITY

- Store iQ Positive Control and iQ Mid Control at 2°C to 8°C (35.6°F to 46.4°F).
- Store iQ Negative Control at 2°C to 28°C (35°F to 82.4°F).
- iQ Precision Kit is stable for eight months from the date of manufacture.

- Open expiration: 30 days from bottle open date.
- Do not freeze.
- Do not store in direct sunlight.
- Bring to room temperature before use.
- Do not save leftover material or mix control material from different bottles.
- Do not use iQ Precision Kit after the expiration date printed on the bottle labels.

MATERIALS PROVIDED

- 2 x 125 mL iQ Positive Control
- 1 x 125 mL iQ Mid Control
- 1 x 125 mL iQ Negative Control

MATERIALS REQUIRED, BUT NOT PROVIDED

- Approved tubes listed in the iQ200 Instructions for Use. Do not use preservative tubes.
- Validation Analyzer software.
- iQ200 Series Automated Urine Microscopy Analyzer.
- Personal Protective Equipment.

PROCEDURE

Follow the directions listed in the iQ200 Series Service Manual.

PRECISION KIT NOMINAL ASSAY VALUES AND RANGES

CONTROL	NOMINAL ASSAY VALUE	RANGE
POS	1,000	± 100
MID	500	± 100
NEG	0	+ 20

Important Note: The actual assay values for the precision kit control material may change from LOT number to LOT number.

TECHNICAL SUPPORT

For additional information, or if damaged product is received, call Beckman Coulter Customer Service at 800-854-3633 (USA or Canada) or contact your local Beckman Coulter Representative.

A Glossary of Symbols (PN C27648) is available at techdocs.beckmancoulter.com.

TRADEMARKS

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

May be covered by one or more pat. - see www.beckmancoulter.com/patents.

Revision History

Revision AA, 10/2018

Initial Release

Revision AB, 03/2019

Add Bulgarian, Hungarian, Polish, Portuguese and Vietnamese languages

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Phiên bản AB, 03/2019
Beckman Coulter Ireland Inc., Lismeehan, O'Callaghan's Mills Co. Clare, Ireland (001) 703-527-3887