

MP Biomedicals, LLC

29525 Fountain Parkway Solon, Ohio 44139 Telephone: 440/337-1200 Toll Free: 800/854-0530 Fax: 440/337-1180 mailto: biotech@mpbio.com

mailto: <u>biotech@mpbio.com</u> web: <u>http://www.mpbio.com</u>

## TECHNICAL INFORMATION

Catalog Number: 195043, 1672346, 1672348

Amphotericin B

## Structure:

Molecular Formula: C<sub>47</sub>H<sub>73</sub>NO<sub>17</sub>

Molecular Weight: 924.11

CAS #: 1397-89-3

Synonyms: Fungizone; Amphozone; Fungilin; Ampho-Moronal

Source: Streptomyces species.

**Physical Description:** Yellow to orange powder. Depending on pH, amphotericin in solution is yellow to dark yellow.

**Description:** An antifungal agent. This product contains no animal-derived components.

Pka values: 5.5 and 10.0.3

**Solubility:** Soluble in acidic water (pH 2) or basic water (pH 11) (about 0.1 mg/ml). Water solubility can be increased with the addition of deoxycholic acid, sodium salt (increasing solubility with increasing concentration of deoxycholic acid). Soluble in Dimethylformamide (DMF) (2 to 4 mg/ml), DMF + 1 M HCl (3:1) (60 to 80 mg/ml), DMSO (30 to 40 mg/ml), propylene glycol; slightly soluble in methyl alcohol; insoluble in water, anhydrous alcohol, ether, benzene, and toluene. Solutions are stable for long periods between pH 4 and 10.1

Cat# 16723 is a concentrated solution of amphotericin at neutral to slightly basic pH, which is slightly cloudy due to the low solubility of amphotericin in aqueous buffer. The concentrated solution is stable at 2-8°C for up to 3 weeks. When diluted to the working concentration (see below for recommended concentration), the amphotericin will completely go into solution. Diluted amphotericin is stable for up to 7 days in culture at 37°C. Store diluted solutions aliquoted at -20°C or below for up to approximately 6 months.

Sterilization of solutions should be by filtration through teflon membrane filters; Amphotericin B should not be gamma irradiated or autoclaved.

Formulation (for 16723):

Tormulation for 10725).					
Component	mg/liter	Mol. Wt.	Mol. (mM)		
Amphotericin B	0.25	924.1	0.00027		
Sodium Deoxycholate	250	414.6	0.60		

**Description:** An antimicrobial/antifungal agent effect against yeast and other fungi at a recommended concentration of 2.5 ug/L (prepared solution use 10 ml/L). It is a mixture of polyenes. It interrupts cell membrane permeability by binding sterols resulting in the loss of low molecular weight compounds from the cell. Amphotericin B may be toxic to some insect cell types. It is inactive against bacteria, rickettsia and viruses. It is typically stable in media at 37°C in the dark for approximately 3 days.<sup>21</sup> **Availability:** 

Catalog Number	Description	Size

195043		100 mg 250 mg 500 mg 1 g 5 g
1672348	Amphotericin B Solution, 250 ug/ml; formulated with deoxycholate and prepared in water for increased solubility, pH = 7-8	20 ml 50 ml

## References:

- Merck Index, 12th Ed., No. 627.
- Pandey, R.C. and Rinehart, K.L., J. Antibiot., v. 29, 1035 (1976).
- Mechlinski, et al., Tetrahedron Letters, 3873 (1970).
- Holz, R.W., Antibiotics (N.Y.), v. 5, 313 (1979).
- Medoff, G., et al., Ann. Rev. Pharmacol. Toxicol., v. 23, 303 (1983).
- Clarkes's Isolation and Identification of Drugs, 2nd Ed., Moffatt, M.C., et al. (eds.); Pharmaceutical Press, p. 351 (1988).
- Martindale, The Extra Pharmacopoeia, 30th Ed.; Pharmaceutical Press, pp. 315-318 (1993).
- The Pharmacological Basis of Therapeutics, 8th Ed., Goodman and Gilman (eds.); Pergamon Press, p. 1165 (1990).
- Data for Biochemical Research, 3rd Ed.; Oxford Press, pp. 296-297 (1989).
- Analytical Profiles of Drug Substances, 6th Ed., Florey, K. (ed.); Academic Press, pp. 1-42 (1977).
  Advances in Lipid Research, v. 14, 127 (1976).
- Wildfeuer, A., et al., "In vitro activity of voriconazole against yeasts, moulds and dermatophytes in comparison with fluconazole, amphotericin B and griseofulvin." Arzneimittelforschung, v. 47:11, 1257-1263 (1997).
- Kretschmar, M., et al., "Effects of amphotericin B incorporated into liposomes and in lipid suspensions in the treatment of murine candidiasis." Arzneimittelforschung, v. 46:7, 711-715 (1996).
- Villani, P., et al., "Clinical and pharmacokinetic evaluation of a new lipid-based delivery system of amphotericin B in AIDS patients." Arzneimittelforschung, v. 46:4, 445-449 (1996).
- Espinel-Ingroff, A., "Clinical utility of in vitro antifungal susceptibility testing." Rev. Esp. Quimioter., v. 13:2, 161-166 (2000).
  Borrelli, P., et al., "Liposomal amphotericin B as first line and secondary prophylactic treatment for visceral leishmaniasis in a patient infected with HIV." Ann. Ital. Med. Int., v. 15:2, 169-171 (2000).
- Eldem, T. and Arican-Cellat, N., "High-performance liquid chromatographic determination of amphotericin B in a liposomal pharmaceutical product and validation of the assay." J. Chromatogr. Sci., v. 38:8, 338-344 (2000).
- Hahn, R.C. and Hamdan, J.S., "Effects of amphotericin B and three azole derivatives on the lipids of yeast cells of Paracoccidioides brasiliensis." Antimicrob. Agents Chemother., v. 44:7, 1997-2000 (2000).
- Younsi, M., et al., "Amphotericin B resistance and membrane fluidity in Kluyveromyces lactis strains." Antimicrob. Agents Chemother., v. 44:7, 1911-1916 (2000).
- Zager, R.A., "Polyene antibiotics: relative degrees of in vitro cytotoxicity and potential effects on tubule phospholipid and ceramide content." Am. J. Kidney Dis., v. 36:2, 238-249 (2000).
- Perlman, D., "Use of Antibiotics in cell culture media," Methods in Enzymology; Cell Culture, Jakoby, W.B. and Pastan, I.H. (eds.), Academic Press: New York, p. 112 (1979).