Toutes les
photos (1)

Documents

 FDS

 COA/CoQ

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1.09204 ► Sigma-Aldrich.

Giemsa's azur eosin methylene blue solution

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

[Ajouter au panier](#)

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Sigma-Aldrich
1.01424
May-Grünwald's eosine-methylene blue
solution modified
for microscopy

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DESCRIPTION

Description générale

Giemsa's Azure Eosin Methylene Blue solution - for microscopy, can be used for the staining of blood and bone marrow smears, paraffin sections and clinical-cytological specimens. Giemsa staining is one of the standard procedures in histology, demonstrating various cells with their morphological features better than with H&E stain. When used in hematological applications, Giemsa's stain is frequently used in combination with other staining solutions, e. g. with May-Grünwald's solution for Pappenheim (MGG) overview staining. This ready-to-use staining solution generally stains the nuclei red, based on the molecular interaction between the eosin Y dye and an Azure B-DNA complex. Both dyes assemble to an Eosin Y - Azure B-DNA complex and the intensity of the resulting stain depends on the content of Azure B and the ratio of Azure B : Eosin Y. Furthermore, the resulting stain can vary depending on the influence of fixation, staining times, pH-value of the solutions or buffer substances. In histology and clinico-cytological applications, Giemsa's staining without additional stainings is used as an extended overview staining method. In this

method, the color of the various cell components is influenced by pretreatment of the specimen material. Here, chromatin-containing structures (e. g. cell nuclei) appear in various blue shades, while the acidophilic components show up in a variety of red shades. The 500 ml bottle is sufficient for 3500-5000 stainings. This product is registered as IVD and CE marked. For more details, please see instructions for use (IFU). The IFU can be downloaded from this webpage.

Application

For microscopy. For the staining of blood and bone-marrow smears, paraffin sections and clinical-cytological specimens

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For staining results, please refer to the Instructions for Use link under the Documents Section below.

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


Plasmodium vivax Liver and Blood Stages Recruit the Druggable Host Membrane Channel Aquaporin-3.

Dora Posfai et al.

Cell chemical biology, 27(6), 719-727 (2020-04-25)

Plasmodium vivax infects hepatocytes to form schizonts that cause blood infection, or dormant hypnozoites that can persist for months in the liver before leading to relapsing blood infections. The molecular processes that drive *P. vivax* schizont and hypnozoite survival remain

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
Sigma-Aldrich

1.09093

AFB-Fluor

Staining kit for fluorescence-microscopic detection of acid fast bacteria

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Contactez notre Service technique

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© 2022 Merck KGaA, Darmstadt, Allemagne et/ou ses sociétés affiliées. Tous droits réservés.