

# ECMatrix™-511 Silk E8 Laminin Substrate

## Stem Cell Reagent

Cat. # CC161-100ML

FOR RESEARCH USE ONLY.  
NOT FOR USE IN DIAGNOSTIC PROCEDURES.  
NOT FOR HUMAN OR ANIMAL CONSUMPTION.

Pack size: 100 mL

Store at 2-8°C



## Data Sheet

page 1 of 2

### Description

Human pluripotent stem cells (ES and iPS cells) express  $\alpha 6 \beta 1$  as the major integrin species and therefore can be maintained stably and expanded efficiently in feeder-free conditions on culture vessels coated with its binding partner laminin-511. However, laminin-511 is not suitable for large-scale production because of its large molecular weight and heterotrimeric nature. Professor Kiyotoshi Sekiguchi's group (Matrixome, Inc.) have solved this problem by producing a recombinant E8 fragment of laminin-511 at large-scale while retaining the full integrin binding activity.<sup>1</sup> The ECMatrix™-511 E8 Laminin Substrate can be used to culture pluripotent stem cells in feeder-free conditions with numerous added benefits over traditional methods including:

#### Features and Benefits

- **Animal-free, xeno-free format:** Consistent from lot-to-lot with no prescreening required
- **No plate precoat required:** Save time by simply adding to media while passaging cells
- **Supports single cell passaging w/out ROCKi:** Great for CRISPR editing or clonal isolation
- **Higher adhesion and growth rates:** Get to your experiments faster
- **Easy to handle:** No chilling of cell culture consumables required

### Storage and Handling

ECMatrix™-511 E8 Laminin Substrates should be stored at 2-8°C. Avoid multiple freeze-thaw cycles and protect from light.

### Presentation

1) 1 X 100 mL ECMatrix™-511 E8 Laminin Substrate (1.6  $\mu$ g/mL in recombinant HSA/PBS). Expressed in transgenic silkworm cocoon.

### Quality Control Testing

- Purity (SDS-Page): > 95%
- Endotoxin Test:  $\leq$  750 EU/mg
- Mycoplasma Test: Negative
- Sterility Test: Negative
- Integrin Binding Assay (kDa):  $\leq$  10 nM

### Protocol

Depending on application, either a precoat or non-precoat method can be used to culture pluripotent stem cells.

#### Non-Precoat Method

1. Detach cells into small clumps or single cells using Accutase.
2. Add ECMatrix™-511 to fresh media at a final concentration of 0.25  $\mu$ g/cm<sup>2</sup> (for example: for one well of a 6-well plate add 5  $\mu$ L of the 0.5 mg/mL stock solution).
3. Add cells to the ECMatrix-511™/Media and plate the cells at desired density.

#### Precoat Method

1. Dilute the 0.5 mg/mL stock solution with sterile PBS to achieve a 2.5  $\mu$ g/mL working solution.
2. Coat dishes with ECMatrix™-511 at 0.25  $\mu$ g/cm<sup>2</sup> (for example, for one well of a 6-well plate add 1 mL of the 2.5  $\mu$ g/mL working solution).
3. Incubate for 1 hour at 37°C, 3 hours at room temperature or overnight at 4°C.
4. Before use, remove remaining fluid from the coated surface (do not rinse).
5. Detach cells into small clumps using Accutase.
6. Plate the cells at desired density.

*Note: Do not allow the plates to dry, briefly spin down all liquids in the tube before use, avoid repeated freeze-thaw cycles.*


### References

1. Sekiguchi K et al. Laminin E8 fragments support efficient adhesion and expansion of dissociated human pluripotent stem cells. *Nature Commun.* 2012;3:1236.
2. Yamanaka S, et al. A novel efficient feeder-free culture system for the derivation of human induced pluripotent stem cells. *Sci Rep.* 2014 Jan 8;4:3594.
3. Takashima Y, et al. Resetting transcription factor control circuitry toward ground-state pluripotency in human. *Cell.* 2014 Sep 11;158(6):1254-1269.
4. Miyazaki T, et al. Efficient Adhesion Culture of Human Pluripotent Stem Cells Using Laminin Fragments in an Uncoated Manner. *Sci Rep.* 2017 Jan 30;7:41165.


**SPECIES LEGEND:** H Human Ca Canine M Mouse R Rat Rb Rabbit B Bovine P Porcine WR Most Common Vertebrates

Please visit [www.millipore.com](http://www.millipore.com) for additional product information and references.

Submit your published journal article, and earn credit toward future purchases. Visit [www.millipore.com/publicationrewards](http://www.millipore.com/publicationrewards) to learn more!



**Figure 1. Growth of human induced pluripotent stem cells on Silk vs. CHO derived ECMatrix™-511 E8 Laminin Substrates.** Epithelial or PBMC derived human iPSCs grown on silkworm or CHO derived ECMatrix™-511 E8 Laminin Substrates have similar proliferation rates and pluripotent morphologies over a 3-day culture period.



**Figure 1. Crystal violet staining of human induced pluripotent stem cells.** Epithelial derived human iPSCs passaged using a split ratio of 1:5 to 1:20 grown on silkworm or CHO derived ECMatrix™-511 E8 Laminin Substrates have similar proliferation rates over a 3-day culture period. Cells stained blue with crystal violet.

ECMatrix™ is a trademark of Merck KGaA  
Matrigel is a registered trademark of Corning Incorporated

■ antibodies ■ Multiplex products ■ biotools ■ cell culture ■ enzymes ■ kits ■ proteins/peptides ■ siRNA/cDNA products

Please visit [www.millipore.com](http://www.millipore.com) for additional product information, test data and references

EMD Millipore Corporation, 28820 Single Oak Drive, Temecula, CA 92590, USA 1-800-437-7500

Technical Support: T: 1-800-MILLIPORE (1-800-645-5476) • F: 1-800-437-7502

**FOR RESEARCH USE ONLY.** Not for use in diagnostic procedures. Not for human or animal consumption. Purchase of this Product does not include any right to resell or transfer, either as a stand-alone product or as a component of another product. Any use of this Product for purposes other than research is strictly prohibited.

EMD Millipore®, the M mark, Upstate®, Chemicon®, Linco® and all other registered trademarks, unless specifically identified above in the text as belonging to a third party, are owned by Merck KGaA, Darmstadt, Germany. Copyright ©2008-2019 Merck KGaA, Darmstadt, Germany. All rights reserved.



We Buy 100% Certified  
Renewable Energy