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
COO/COA

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DN25 ► Sigma-Aldrich.

Deoxyribonuclease I from bovine pancreas

★★★★★ 5.0 (2)

lyophilized powder, Protein ≥85 %, ≥400 Kunitz units/mg protein

Synonyme(s):

Deoxyribonuclease, DNase I, Deoxyribonuclease 5'-oligonucleotide-hydrolase

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
Sigma-Aldrich

D5025

Deoxyribonuclease I from bovine pancreas

Type IV, lyophilized powder, ≥2,000 Kunitz units/mg protein

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DESCRIPTION

Description générale

DNase I, or Deoxyribonuclease I, is an endonuclease isolated from bovine pancreas.

- Our DNase I Digests double str and single stranded DNA into oligo and mononucleotides.
- Bovine pancreatic DNase exists as four isozymes, having isoelectric points for A, B, C and D: 5.22, 4.96, 5.06 and 4.78.3. The predominant form is A, with smaller amounts of B and C, and only minor amount of D.
- DNase I structure resembles the structure of exonuclease III. It includes two central β sheets. Each β sheet is composed of six β -strands. This complex of β sheets is surrounded by extensive loop and α -helical regions. This enzyme shares structural similarity to exonuclease III. [\[1\]](#)

Application

- Decreases viscosity providing better yields by removing DNA in primary cell isolation:
- Incorporating labelled bases into DNA: DNA nick
- Radioactive labelling
- Bioprocessing applications: DNA removal
- Eliminating genomic DNA from RNA preparations before RT-PCR
- *In vitro* transcription
- Nick translation
- DNase footprinting
- Actin regulation of actin polymerization in cells, and cell apoptosis
- UV crosslinking of proteins to nucleic acids
- DNase play a role in the regulation of actin polymerization in cells and is involved in apoptosis process [\[1\]](#)

Used for the removal of DNA from protein samples.

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DNase I is an endonuclease that acts on phosphodiester bonds adjacent to pyrimidines to produce polynucleotides with terminal 5'-phosphates. In the presence of Mg^{2+} , DNase I cleaves each strand of DNA independently and the cleavage sites are random. Both DNA strands are cleaved at approximately the same site in the presence of Mn^{2+} . Divalent cations such as Mn^{2+} , Ca^{2+} , Co^{2+} , and Zn^{2+} are activators of the enzyme. A concentration of 5 mM Ca^{2+} stabilizes the enzyme against proteolytic digestion. The pH optimum is found to be between 7 and 8. DNase I from bovine pancreas consists of four chromatographically distinguishable components, A, B, C, and D, with their molar ratios being 4:1:1 respectively. Only minor amounts of D are found. 2-Mercaptoethanol, chelators, sodium dodecyl sulfate (SDS) and actin are known to inhibit the enzyme activity.

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Our Deoxyribonuclease DNase I, is essentially RNase free product, which support the product application.

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One Kunitz unit will produce a change in A_{260} of 0.001 per min per mL at pH 5.0 at 25 °C, using DNA, Type I or III, as substrate.

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Crude preparation, contains calcium chloride

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10 mg/mL solution of DNase I in 0.15 M NaCl may lose <10% of its activity when stored for a week in aliquots at -20 °C. The same solutions stored in aliquots at 2-8 °C can lose approximately 20% activity. It remains active for up to five hours at 60 °C between pH 5 and 7, and loses activity in <10 minutes at 68 °C. It loses activity at the rate of 6%/hour in acetate buffer (pH 5.0) and tris buffer ((pH 7.2) at 1 mg/mL concentration.

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Protein determined by biuret.

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Inhibiteur

N9653

Netropsin dihydrochloride, from *Streptomyces netropsis*, ≥98% (HPLC and TLC), powder

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

Sigma-Aldrich

C5138**Collagenase from *Clostridium histolyticum***

suitable for release of physiologically active rat hepatocytes, Type IV, 0.5-5.0 FALGPA units/mg solid,...

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Enzymes of Molecular Biology

Weir, A.F.*Methods in Molecular Biology*, 16(2), 2-16 null

Purification of rat spermatogenic cells and preliminary biochemical analysis of these cells.

ML Meistrich et al.

Biology of reproduction, 25(5), 1065-1077 (1981-12-01)

Comparison of the multiple forms of bovine pancreatic deoxyribonuclease.

J Salnikow et al.

The Journal of biological chemistry, 245(21), 5685-5690 (1970-11-10)

The neurospora plasma membrane ATPase is an electrogenic pump.

G A Scarborough

Proceedings of the National Academy of Sciences of the United States of America, 73(5), 1485-1488 (1976-05-01)

Biochemical evidence is presented which demonstrates that the *Neurospora crassa* plasma membrane ATPase (ATP phosphohydrolase, EC 3.6.1.3) is an electrogenic pump. The electrical potential across the *Neurospora* plasma membrane, as monitored by [¹⁴C]SCN⁻ uptake by isolated *Neurospora* plasma membrane vesicles

Kaposi Sarcoma-Associated Herpesvirus Glycoprotein H Is Indispensable for Infection of Epithelial, Endothelial, and Fibroblast Cell Types.

Murali Muniraju et al.

Journal of virology, 93(16) (2019-05-31)

Kaposi sarcoma-associated herpesvirus (KSHV) is an emerging pathogen and is the causative infectious agent of Kaposi sarcoma and two malignancies of B cell origin. To date, there is no licensed KSHV vaccine. Development of an effective vaccine against KSHV continues

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
Protocoles

Enzymatic Assay of Deoxyribonuclease I (EC 3.1.21.1)

To standardize a procedure for the enzymatic assay of Deoxyribonuclease I.

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
★★★★★ · il y a 3 mois

DNASE 1

This reagent works well for our purpose in our cell culture prep for cardiomyocytes

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© 2022 Merck KGaA, Darmstadt, Allemagne et/ou ses sociétés affiliées. Tous droits réservés.