Bilirubin Auto Direct FS* (Bilirubine Auto Directe FS*)

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| Référence | Composition du kit | | | |
|------------------|--------------------|------------|---|----------------|
| 1 0821 99 10 021 | R1 | 5 x 20 mL | + | R2 1 x 25 mL |
| 1 0821 99 10 026 | R1 | 5 x 80 mL | + | R2 1 x 100 mL |
| 1 0821 99 10 023 | R1 | 1 x 800 mL | + | R2 1 x 200 mL |
| 1 0821 99 10 704 | R1 | 8 x 50 mL | + | R2 8 x 12,5 mL |
| 1 0821 99 10 930 | R1 | 4 x 20 mL | + | R2 2 x 10 mL |

Emploi Prévu

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4. Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. ClinChemLabMed 2007;45(9):1240-1243.
5. Guder WG, Zawta B et al. The Quality of Diagnostic Samples. 1st ed. Darmstadt: GIT Verlag; 2001; p. 18-9.
6. Young DS. Effects of Drugs on Clinical Laboratory Tests. 5th ed. Volume 1 and 2. Washington, DC: The American Association for Clinical Chemistry Press 2000.
7. Young DS. Effects on Clinical Laboratory Tests - Drugs Disease, Herbs & Natural Products, <https://clinf.wiley.com/aaccweb/aacc/>, accessed in December 2020. Published by AACC Press and John Wiley and Sons, Inc.

DiaSys Diagnostic Systems GmbH
Alte Strasse 9 65558 Holzheim
Allemagne
www.diasys-diagnostics.com

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1. Thomas L ed. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft, 1998: p. 192-202.
2. Tolman KG, Rej R. Liver function. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 1125-77.
3. Rand RN, di Pasqua A. A new diazo method for the determination of bilirubin. Clin Chem 1962;6:570-8.
4. Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. ClinChemLabMed 2007;45(9):1240-1243.
5. Guder WG, Zawta B et al. The Quality of Diagnostic Samples. 1st ed. Darmstadt: GIT Verlag; 2001; p. 18-9.
6. Young DS. Effects of Drugs on Clinical Laboratory Tests. 5th ed. Volume 1 and 2. Washington, DC: The American Association for Clinical Chemistry Press 2000.
7. Young DS. Effects on Clinical Laboratory Tests - Drugs Disease, Herbs & Natural Products, <https://clinfex.wiley.com/aaccweb/aacc/>, accessed in December 2020. Published by AACC Press and John Wiley and Sons, Inc.

DiaSys Diagnostic Systems GmbH
Alte Strasse 9 65558 Holzheim
Allemagne
www.diasys-diagnostics.com

* Fluid Stable = Liquide & Stable

Bilirubin Auto Direct FS

Application for serum and plasma samples

This application was set up and evaluated by DiaSys. It is based on the standard equipment at that time and does not apply to any equipment modifications undertaken by unqualified personnel.

| Identification | |
|-------------------------------------|------|
| This method is usable for analysis: | Yes |
| Twin reaction: | No |
| Name: | DBIL |
| Shortcut: | |
| Reagent barcode reference: | 018 |
| Host reference: | |

| Technic | |
|---------------------------------------|------------|
| Type: | End point |
| First reagent:[μL] | 180 |
| Blank reagent | Yes |
| Sensitive to light | |
| Second reagent:[μL] | 45 |
| Blank reagent | Yes |
| Sensitive to light | |
| Main wavelength:[nm] | 546 |
| Secondary wavelength:[nm] | 660 |
| Polychromatic factor: | 1.000 |
| 1 st reading time [min:sec] | (04:24) |
| Last reading time [min:sec] | 10:00 |
| Reaction way: | Increasing |
| Linear Kinetics | |
| Substrate depletion: Absorbance li | |
| Linearity: Maximum deviation [%] | |
| Fixed Time Kinetics | |
| Substrate depletion: Absorbance limit | |
| Endpoint | |
| Stability: Largest remaining slope | |
| Prozone Limit [%] | |

| Reagents | |
|----------|--|
| Decimals | |
| Units | |

| Sample | |
|--------------------------------------|------------------|
| Diluent | DIL A (NaCl) |
| Hemolysis: | |
| Agent [μL] | 0 (no hemolysis) |
| Cleaner | |
| Sample [μL] | 0 |
| Technical limits | |
| Concentration technical limits-Lower | 0.1 |
| Concentration technical limits-Upper | 7 |
| SERUM | |
| Normal volume [μL] | 8 |
| Normal dilution (factor) | 1 |
| Below normal volume [μL] | 15 |
| Below normal dilution (factor) | 1 |
| Above normal volume [μL] | 2 |
| Above normal dilution (factor) | 1 |
| URIN | |
| Normal volume [μL] | 8 |
| Normal dilution (factor) | 1 |
| Below normal volume [μL] | 15 |
| Below normal dilution (factor) | 1 |
| Above normal volume [μL] | 2 |
| Above normal dilution (factor) | 1 |
| PLASMA | |
| Normal volume [μL] | 8 |
| Normal dilution (factor) | 1 |
| Below normal volume [μL] | 15 |
| Below normal dilution (factor) | 1 |
| Above normal volume [μL] | 2 |
| Above normal dilution (factor) | 1 |
| CSF | |
| Normal volume [μL] | 8 |
| Normal dilution (factor) | 1 |
| Below normal volume [μL] | 15 |
| Below normal dilution (factor) | 1 |
| Above normal volume [μL] | 2 |
| Above normal dilution (factor) | 1 |
| Whole blood | |
| Normal volume [μL] | 8 |
| Normal dilution (factor) | 1 |
| Below normal volume [μL] | 15 |
| Below normal dilution (factor) | 1 |
| Above normal volume [μL] | 2 |
| Above normal dilution (factor) | 1 |

| Results | |
|---------------------------|-------|
| Decimals | 2 |
| Units | mg/dL |
| Correlation factor-Offset | 0.000 |
| Correlation factor-Slope | 1.000 |

| Range | |
|-------------|----------|
| Gender | All |
| Age | |
| SERUM | >= <=0.2 |
| URINE | |
| PLASMA | >= <=0.2 |
| CSF | |
| Whole blood | |
| Gender | |
| Age | |
| SERUM | |
| URINE | |
| PLASMA | |
| CSF | |
| Whole blood | |

| Contaminants | |
|---|--|
| Please refer to r910 Carryover Pair Table | |

| Calibrators details | |
|---------------------|---------------|
| Calibrator list | Concentration |
| Cal. 1/Blank | 0 |
| Cal. 2 | * |
| Cal. 3 | |
| Cal. 4 | |
| Cal. 5 | |
| Cal. 6 | |
| Max delta abs. | |
| Cal. 1 | 0.002 |
| Cal. 2 | 0.005 |
| Cal. 3 | |
| Cal. 4 | |
| Cal. 5 | |
| Cal. 6 | |
| Drift limit [%] | 0.8 |

| Calculations | |
|--------------|---|
| Model | X |
| Degree | 1 |

* Enter calibrator value