Glucose GOD FS *

CODE CQN : HD

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1. Thomas L. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998. p. 131-7.
2. Sacks DB. Carbohydrates. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 750-808.
3. Barham D, Trinder P. An improved color reagent for the determination of blood glucose by the oxidase system. Analyst 1972; 97: 142-5.
4. Guder WG, Zawta B et al. The Quality of Diagnostic Samples. 1st ed. Darmstadt: GIT Verlag; 2001; p. 30-1.
5. Sacks DB, Bruns DE, Goldstein DE, Mac Laren NK, Mc Donald JM, Parrott M. Guidelines and recommendations for laboratory analysis in the diagnosis and management of diabetes mellitus. Clin Chem 2002; 48: 436-72.
6. Young DS. Effects of Drugs on Clinical Laboratory Tests. 5th ed. Volume 1 and 2. Washington, DC: The American Association for Clinical Chemistry Press 2000.
7. Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. ClinChemLabMed 2007;45(9):1240-1243.

Fabricant

DiaSys Diagnostic Systems GmbH
Alte Strasse 9 65558 Holzheim (Allemagne)

Glucose Standard FS* (Standard Glucose FS*)

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1. Barham D, Trinder P. An improved color reagent for the determination of blood glucose by the oxidase system. *Analyst* 1972; 97: 142-5.
2. Guder WG, Zawta B et al. The Quality of Diagnostic Samples. 1st ed. Darmstadt: GIT Verlag; 2001; p. 30-1.
3. Sacks DB. Carbohydrates. In: Burtis CA, Ashwood ER, editors. *Tietz Textbook of Clinical Chemistry*. 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 750-808.
4. Sacks DB, Bruns DE, Goldstein DE, Mac Laren NK, Mc Donald JM, Parrott M. Guidelines and recommendations for laboratory analysis in the diagnosis and management of diabetes mellitus. *Clin Chem* 2002; 48: 436-72.
5. Young DS. *Effects of Drugs on Clinical Laboratory Tests*. 5th. ed. Volume 1 and 2. Washington, DC: The American Association for Clinical Chemistry Press, 2000.
6. Thomas L. *Clinical Laboratory Diagnostics*. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998. p. 131-7.
7. Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. *ClinChemLabMed* 2007;45(9):1240-1243.

Fabricant

DiaSys Diagnostic Systems GmbH
Alte Strasse 9 65558 Holzheim Allemagne

Glucose GOD FS 10'

Application for serum and plasma samples

This application was set up and evaluated by DiaSys. It is based on the standard equipment at that time and does not apply to any equipment modifications undertaken by unqualified personnel.

Identification	
This method is usable for analysis:	Yes
Twin reaction:	No
Name:	GLUC GOD
Shortcut:	
Reagent barcode reference:	036
Host reference:	036

Technic	
Type:	End point
First reagent:[μL]	180
Blank reagent	Yes
Sensitive to light	
Second reagent:[μL]	
Blank reagent	
Sensitive to light	
Main wavelength:[nm]	508
Secondary wavelength:[nm]	700
Polychromatic factor:	1.0000
1 st reading time [min:sec]	(-00:12)
Last reading time [min:sec]	09:48
Reaction way:	Increasing
Linear Kinetics	
Substrate depletion: Absorbance limit	
Linearity: Maximum deviation [%]	
Fixed Time Kinetics	
Substrate depletion: Absorbance limit	
Endpoint	
Stability: Largest remaining slope	
Prozone Limit [%]	

Reagents	
Decimals	
Units	

Sample	
Diluent	DIL A (NaCl)
Hemolysis:	
Agent [μL]	0 (no hemolysis)
Cleaner	
Sample [μL]	0
Technical limits	
Concentration technical limits-Lower	3.0000
Concentration technical limits-Upper	400.0000
SERUM	
Normal volume [μL]	2.0
Normal dilution (factor)	1
Below normal volume [μL]	
Below normal dilution (factor)	
Above normal volume [μL]	2.0
Above normal dilution (factor)	6
URINE	
Normal volume [μL]	2.0
Normal dilution (factor)	1
Below normal volume [μL]	
Below normal dilution (factor)	
Above normal volume [μL]	2.0
Above normal dilution (factor)	6
PLASMA	
Normal volume [μL]	2.0
Normal dilution (factor)	1
Below normal volume [μL]	
Below normal dilution (factor)	
Above normal volume [μL]	2.0
Above normal dilution (factor)	6
CSF	
Normal volume [μL]	2.0
Normal dilution (factor)	1
Below normal volume [μL]	
Below normal dilution (factor)	
Above normal volume [μL]	2.0
Above normal dilution (factor)	6
Whole blood	
Normal volume [μL]	2.0
Normal dilution (factor)	1
Below normal volume [μL]	
Below normal dilution (factor)	
Above normal volume [μL]	2.0
Above normal dilution (factor)	6

Results	
Decimals	1
Units	mg/dL
Correlation factor-Offset	0.0000
Correlation factor-Slope	1.0000

Range	
Gender	All
Age	
SERUM	>=70.0 <=115.0
URINE	
PLASMA	>=70.0 <=115.0
CSF	
Whole blood	
Gender	
Age	
SERUM	
URINE	
PLASMA	
CSF	
Whole blood	

Contaminants	
Please refer to r910 Carryover Pair Table	

Calibrators details	
Calibrator list	Concentration
Cal. 1/Blank	0
Cal. 2	*
Cal. 3	
Cal. 4	
Cal. 5	
Cal. 6	
Max delta abs.	
Cal. 1	0.006
Cal. 2	0.040
Cal. 3	
Cal. 4	
Cal. 5	
Cal. 6	
Drift limit [%]	0.80

Calculations	
Model	X
Degree	1

* Enter calibrator value