ASAT (GOT) FS* (IFCC mod.)

Avec/Sans Pyridoxal-5-Phosphate FS (P-5-P)

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Référence	Composition du kit			
1 2601 99 10 021	R1	5 x 20 mL	+	R2 1 x 25 mL
1 2601 99 10 026	R1	5 x 80 mL	+	R2 1 x 100 mL
1 2601 99 10 023	R1	1 x 800 mL	+	R2 1 x 200 mL
1 2601 99 10 704	R1	8 x 50 mL	+	R2 8 x 12,5 mL
1 2601 99 10 917	R1	8 x 60 mL	+	R2 8 x 15 mL
1 2601 99 90 314	R1	10 x 20 mL	+	R2 2 x 30 mL

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- Thomas L. Alanine aminotransferase (ALT), Aspartate aminotransferase (AST). In: Thomas L, editor. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998. p. 55-65.
- Moss DW, Henderson AR. Clinical enzymology. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 617-721.
- Bergmeyer HU, Horder M, Rej R. Approved Recommendation (1985) on IFCC Methods for the Measurement of Catalytic Concentration of Enzymes. L.Clin. Chem. Clin. Biochem 1986; 24: 497-510.
- Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. ClinChemLabMed 2007;45(9):1240-1243.
- Guder WG, Zawta B et al. The Quality of Diagnostic Samples. 1st ed. Darmstadt: GIT Verlag; 2001; p. 18-9.
- Young DS. Effects of Drugs on Clinical Laboratory Tests. 5th ed. Volume 1 and 2. Washington, DC: The American Association for Clinical Chemistry Press 2000.
- Young DS. Effects on Clinical Laboratory Tests - Drugs Disease, Herbs & Natural Products, <https://clinf.wiley.com/aaccweb/aacc/>, accessed in September 2021. Published by AACC Press and John Wiley and Sons, Inc.
- Schumann G, Bonora R, Ceriotti F, Féraud G et al. IFCC primary reference procedure for the measurement of catalytic activity concentrations of enzymes at 37°C. Part 5: Reference procedure for the measurement of catalytic concentration of aspartate aminotransferase. Clin Chem Lab Med 2002;40:725-33.
- Lorentz K, Röhle G, Siekmann L. Einführung der neuen Standardmethoden 1994 zur Bestimmung der katalytischen

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1. Thomas L. Alanine aminotransferase (ALT), Aspartate aminotransferase (AST). In: Thomas L, editor. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998. p. 55-65.
2. Moss DW, Henderson AR. Clinical enzymology. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 617-721.
3. Bergmeyer HU, Horder M, Rej R. Approved Recommendation (1985) on IFCC Methods for the Measurement of Catalytic Concentration of Enzymes. L.Clin. Chem. Clin. Biochem 1986; 24: 497-510.
4. Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. ClinChemLabMed 2007;45(9):1240-1243.
5. Guder WG, Zawta B et al. The Quality of Diagnostic Samples. 1st ed. Darmstadt: GIT Verlag; 2001; p. 18-9.
6. Young DS. Effects of Drugs on Clinical Laboratory Tests. 5th ed. Volume 1 and 2. Washington, DC: The American Association for Clinical Chemistry Press 2000.
7. Young DS. Effects on Clinical Laboratory Tests - Drugs Disease, Herbs & Natural Products, <https://clinfo.wiley.com/aaccweb/aacc/>, accessed in September 2021. Published by AACC Press and John Wiley and Sons, Inc.
8. Schumann G, Bonora R, Ceriotti F, Féraud G et al. IFCC primary reference procedure for the measurement of catalytic activity concentrations of enzymes at 37°C. Part 5: Reference procedure for the measurement of catalytic concentration of aspartate aminotransferase. Clin Chem Lab Med 2002;40:725-33.
9. Lorentz K, Röhle G, Siekmann L. Einführung der neuen Standardmethoden 1994 zur Bestimmung der katalytischen Enzymkonzentrationen bei 37 °C. DG Klinische Chemie Mitteilungen 26; 1995; Heft 4.
10. Zawta B, Klein G, Bablok W. Temperature Conversion in Clinical Enzymology? Klin. Lab. 1994; 40: 33-42.

DiaSys Diagnostic Systems GmbH
Alte Strasse 9 65558 Holzheim
Germany
www.diasys-diagnostics.com

* Fluid Stable = Liquide & Stable

ASAT (GOT) FS (IFCC mod.)

Application for serum and plasma samples

This application was set up and evaluated by DiaSys. It is based on the standard equipment at that time and does not apply to any equipment modifications undertaken by unqualified personnel.

Identification	
This method is usable for analysis:	Yes
Twin reaction:	No
Name:	AST
Shortcut:	
Reagent barcode reference:	011
Host reference:	

Technic	
Type:	Linear Kinetic
First reagent:[μL]	160
Blank reagent	Yes
Sensitive to light	
Second reagent:[μL]	40
Blank reagent	Yes
Sensitive to light	
Main wavelength:[nm]	340
Secondary wavelength:[nm]	405
Polychromatic factor:	1.000
1 st reading time [min:sec]	5:48
Last reading time [min:sec]	8:48
Reaction way:	Decreasing
Linear Kinetics	
Substrate depletion: Absorbance li	0.2700
Linearity: Maximum deviation [%]	100
Fixed Time Kinetics	
Substrate depletion: Absorbance limit	
Endpoint	
Stability: Largest remaining slope	
Prozone Limit [%]	

Reagents	
Decimals	
Units	

Sample	
Diluent	DIL A (NaCl)
Hemolysis:	
Agent [μL]	0 (no hemolysis)
Cleaner	
Sample [μL]	0
Technical limits	
Concentration technical limits-Lower	2.00
Concentration technical limits-Upper	700
SERUM	
Normal volume [μL]	12
Normal dilution (factor)	1
Below normal volume [μL]	20
Below normal dilution (factor)	1
Above normal volume [μL]	2
Above normal dilution (factor)	1
URIN	
Normal volume [μL]	12
Normal dilution (factor)	1
Below normal volume [μL]	20
Below normal dilution (factor)	1
Above normal volume [μL]	2
Above normal dilution (factor)	1
PLASMA	
Normal volume [μL]	12
Normal dilution (factor)	1
Below normal volume [μL]	20
Below normal dilution (factor)	1
Above normal volume [μL]	2
Above normal dilution (factor)	1
CSF	
Normal volume [μL]	12
Normal dilution (factor)	1
Below normal volume [μL]	20
Below normal dilution (factor)	1
Above normal volume [μL]	2
Above normal dilution (factor)	1
Whole blood	
Normal volume [μL]	12
Normal dilution (factor)	1
Below normal volume [μL]	20
Below normal dilution (factor)	1
Above normal volume [μL]	2
Above normal dilution (factor)	1

Results	
Decimals	1
Units	U/L
Correlation factor-Offset	0.000
Correlation factor-Slope	1.000

Range	
Gender	Male
Age	
SERUM	>= <=35.0
URINE	
PLASMA	>= <=35.0
CSF	
Whole blood	
Gender	Female
Age	
SERUM	>= <=31.0
URINE	
PLASMA	>= <=31.0
CSF	
Whole blood	

Contaminants	
Please refer to r910 Carryover Pair Table	

Calibrators details	
Calibrator list	Concentration
Cal. 1/Blank	0
Cal. 2	*
Cal. 3	
Cal. 4	
Cal. 5	
Cal. 6	
Max delta abs.	
Cal. 1	0.002
Cal. 2	0.005
Cal. 3	
Cal. 4	
Cal. 5	
Cal. 6	
Drift limit [%]	0.8

Calculations	
Model	X
Degree	1

* Enter calibrator value

ASAT (GOT) FS (IFCC mod.) with P-5-P activation

Application for serum and plasma samples

This application was set up and evaluated by DiaSys. It is based on the standard equipment at that time and does not apply to any equipment modifications undertaken by unqualified personnel.

Identification	
This method is usable for analysis:	Yes
Twin reaction:	No
Name:	AST
Shortcut:	
Reagent barcode reference:	64
Host reference:	

Technic	
Type:	Linear Kinetic
First reagent:[μL]	160
Blank reagent	Yes
Sensitive to light	
Second reagent:[μL]	40
Blank reagent	Yes
Sensitive to light	
Main wavelength:[nm]	340
Secondary wavelength:[nm]	405
Polychromatic factor:	1.000
1 st reading time [min:sec]	5:48
Last reading time [min:sec]	8:48
Reaction way:	Decreasing
Linear Kinetics	
Substrate depletion: Absorbance li	0.3500
Linearity: Maximum deviation [%]	100
Fixed Time Kinetics	
Substrate depletion: Absorbance limit	
Endpoint	
Stability: Largest remaining slope	
Prozone Limit [%]	

Reagents	
Decimals	
Units	

Sample	
Diluent	DIL A (NaCl)
Hemolysis:	
Agent [μL]	0 (no hemolysis)
Cleaner	
Sample [μL]	0
Technical limits	
Concentration technical limits-Lower	2.00
Concentration technical limits-Upper	675
SERUM	
Normal volume [μL]	12
Normal dilution (factor)	1
Below normal volume [μL]	20
Below normal dilution (factor)	1
Above normal volume [μL]	2
Above normal dilution (factor)	1
URIN	
Normal volume [μL]	12
Normal dilution (factor)	1
Below normal volume [μL]	20
Below normal dilution (factor)	1
Above normal volume [μL]	2
Above normal dilution (factor)	1
PLASMA	
Normal volume [μL]	12
Normal dilution (factor)	1
Below normal volume [μL]	20
Below normal dilution (factor)	1
Above normal volume [μL]	2
Above normal dilution (factor)	1
CSF	
Normal volume [μL]	12
Normal dilution (factor)	1
Below normal volume [μL]	20
Below normal dilution (factor)	1
Above normal volume [μL]	2
Above normal dilution (factor)	1
Whole blood	
Normal volume [μL]	12
Normal dilution (factor)	1
Below normal volume [μL]	20
Below normal dilution (factor)	1
Above normal volume [μL]	2
Above normal dilution (factor)	1

Results	
Decimals	1
Units	U/L
Correlation factor-Offset	0.000
Correlation factor-Slope	1.000

Range	
Gender	Male
Age	
SERUM	>= <=35.0
URINE	
PLASMA	>= <=35.0
CSF	
Whole blood	
Gender	Female
Age	
SERUM	>= <=31.0
URINE	
PLASMA	>= <=31.0
CSF	
Whole blood	

Contaminants	
Please refer to r910 Carryover Pair Table	

Calibrators details	
Calibrator list	Concentration
Cal. 1/Blank	0
Cal. 2	*
Cal. 3	
Cal. 4	
Cal. 5	
Cal. 6	
Max delta abs.	
Cal. 1	0.002
Cal. 2	0.005
Cal. 3	
Cal. 4	
Cal. 5	
Cal. 6	
Drift limit [%]	0.8

Calculations	
Model	X
Degree	1

* Enter calibrator value