# HbA1c FS\*

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1. Thomas L. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998. p. 142-48.
2. Sacks DB. Carbohydrates. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3<sup>rd</sup> ed. Philadelphia: W.B. Saunders Company; 1999. p. 790-6.
3. Sacks DB. Carbohydrates. In: Burtis CA, Ashwood ER, Bruns DE, editors. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. 4th edition St. Louis Missouri: Elsevier Saunders; 2006; p. 878-884.
4. Jeppsson JO, Kobold U, Barr J, Finke A et al. Approved IFCC reference method for the measurement of HbA1c in human blood. Clin Chem Lab Med 2002; 40: 78-89.
5. Hoelzel W, Weykamp C et al. IFCC Reference System for Measurement of Hemoglobin A1c in Human Blood and the National Standardization Schemes in the United States, Japan, and Sweden: A Method-Comparison Study. Clin Chem 2004; 50 (1): 166-74.
6. The Diabetes Control and Complications Trial Research Group. The effect of intensive treatment of diabetes in the development and progression of long-term complications in insulin-dependent diabetes mellitus. N Engl J Med. 1993; 329: 977-86.
7. Little RR, Rohlfing CL, Wiedmeyer HM, Myers GL et al. The National Glycohemoglobin Standardization Program: A Five-Years Progress Report. Clin Chem 2001; 47: 1985-92.
8. Data on file at DiaSys Diagnostic Systems GmbH.
9. Pantheghini M, John WG on behalf of the IFCC Scientific Division. Implementation of haemoglobin A1c results traceable to the IFCC reference system: the way forward. Clin Chem Lab Med 2007; 45(8): 942-4.
10. Nordin G., Dybkær R. Recommendation for term and measurement unit for "HbA1c". Clin Chem Lab Med 2007; 45(8): 1081-2.
11. Sacks DB. Translating Hemoglobin A1c into Average Blood Glucose: Implications for Clinical Chemistry. Clinical Chemistry 2008; 54: 1756-8.
12. Weykamp C. Carbamylated Hemoglobin Interference in Glycohemoglobin Assays. Clin Chem 1999; 45: 438-9.
13. Young DS. Effects of Drugs on Clinical Laboratory Tests. 5th ed. Volume 1 and 2. Washington, DC: The American Association for Clinical Chemistry Press 2000.
14. Sacks DB, Arnold M, Bakris GL, Bruns DE, AR Horvath et al. Guidelines and recommendations for laboratory analysis in the diagnosis and management of diabetes mellitus. Clin Chem 2011; 57(6): e1-e47.
15. Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. ClinChemLabMed 2007; 45(9): 1240-1243.
16. Ferri S, Kim S, Tsugawa W, Sode K. Review of Fructosyl Amino Acid Oxidase Engineering Research: A Glimpse into the Future of Hemoglobin A1c Biosensing. Journal of Diabetes Science and Technology 2009; 3(3): 585-592.

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1. Jeppsson JO, Kobold U, Barr J, Finke A et al. Approved IFCC reference method for the measurement of HbA1c in human blood. Clin Chem Lab Med 2002; 40: 78-89.
2. Hoelzel W, Weykamp C et al. IFCC Reference System for Measurement of Hemoglobin A1c in Human Blood and the National Standardization Schemes in the United States, Japan, and Sweden: A Method-Comparison Study. Clin Chem 2004; 50 (1): 166-74.
3. The Diabetes Control and Complications Trial Research Group. The effect of intensive treatment of diabetes in the development and progression of long-term complications in insulin-dependent diabetes mellitus. N Engl J Med.1993; 329: 977-86.
4. Little RR, Rohlfing CL, Wiedmeyer HM, Myers GL et al. The National Glycohemoglobin Standardization Program: A Five-Years Progress Report. Clin Chem 2001; 47: 1985-92.
5. Data on file at DiaSys Diagnostic Systems GmbH.
6. Panteghini M, John WG on behalf of the IFCC Scientific Division. Implementation of haemoglobin A1c results traceable to IFCCreference system: the way forward. Clin Chem Lab Med 2007; 45(8): 942-4
7. Nordin G, Dybkær R. Recommendation for term and measurement unit for "HbA1c". Clin Chem Lab Med 2007; 45(8): 1081-2.
8. Sacks DB. Translating Hemoglobin A1c into Average Blood Glucose: Implications for Clinical Chemistry. Clinical Chemistry 2008; 54: 1756-8.
9. Weykamp C. Carbamylated Hemoglobin Interference in Glyco-hemoglobin Assays. Clin Chem 1999; 45: 438-9.
10. Young DS. Effects of Drugs on Clinical Laboratory Tests. 5th. ed. Volume 1 and 2. Washington, DC: The American Association for Clinical Chemistry Press, 2000.
11. Thomas L. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998. p. 142-48.
12. Sacks DB. Carbohydrates. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3<sup>rd</sup> ed. Philadelphia: W.B. Saunders Company; 1999. p. 790-6.
13. Sacks DB. Carbohydrates. In: Burtis CA, Ashwood ER, Bruns DE, editors. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. 4th edition St. Louis Missouri: Elsevier Saunders; 2006; p. 878-884.
14. Sacks DB, Arnold M, Bakris GL, Bruns DE, AR Horvath et al. Guidelines and recommendations for laboratory analysis in the diagnosis and management of diabetes mellitus. Clin Chem 2011; 57(6): e1-e47.
15. Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. ClinChemLabMed 2007;45(9):1240–1243.
16. Ferri S, Kim S, Tsugawa W, Sode K. Review of Fructosyl Amino Acid Oxidase Engineering Research: A Glimpse into the Future of Hemoglobin A1c Biosensing. Journal of Diabetes Science and Technology 2009; 3(3): 585-592.

## Fabricant

DiaSys Diagnostic Systems GmbH  
Alte Strasse 9 65558 Holzheim (Allemagne)

# HbA1c net FS (reaction 1 hemoglobin)

## Application for whole blood samples

This application was set up and evaluated by DiaSys. It is based on the standard equipment at that time and does not apply to any equipment modifications undertaken by unqualified personnel

Identification	
This method is usable for analysis:	Yes
Twin reaction:	Yes
Name:	HbA1c net
Shortcut:	
Reagent barcode reference:	723
Host reference:	

Technic	
Type:	End point
First reagent:[μL]	150
Blank reagent	Yes
Sensitive to light	Yes
Second reagent:[μL]	
Blank reagent	
Sensitive to light	
Main wavelength:[nm]	570
Secondary wavelength:[nm]	800
Polychromatic factor:	1.000
1 st reading time [min:sec]	(04:24)
Last reading time [min:sec]	04:24
Reaction way:	Increasing
Linear Kinetics	
Linearity: Maximum deviation [%]	
Fixed Time Kinetics	
Endpoint	
Prozone Limit [%]	

Reagents	
Decimals	2
Units	g/dL

Sample	
Diluent	System water
Hemolysis:	
Agent [μL]	NET(R951) 200
Cleaner	CLN A (R900)
Sample [μL]	10
Technical limits	HbA1c [mmol/mol]: 20 -150
Concentration technical limits-Lower	6
Concentration technical limits-Upper	30
SERUM	
Normal volume [μL]	25
Normal dilution (factor)	1
Below normal volume [μL]	25
Below normal dilution (factor)	1
Above normal volume [μL]	25
Above normal dilution (factor)	1
URIN	
Normal volume [μL]	25
Normal dilution (factor)	1
Below normal volume [μL]	25
Below normal dilution (factor)	1
Above normal volume [μL]	25
Above normal dilution (factor)	1
PLASMA	
Normal volume [μL]	25
Normal dilution (factor)	1
Below normal volume [μL]	25
Below normal dilution (factor)	1
Above normal volume [μL]	25
Above normal dilution (factor)	1
CSF	
Normal volume [μL]	25
Normal dilution (factor)	1
Below normal volume [μL]	25
Below normal dilution (factor)	1
Above normal volume [μL]	25
Above normal dilution (factor)	1
Whole blood	
Normal volume [μL]	25
Normal dilution (factor)	1
Below normal volume [μL]	25
Below normal dilution (factor)	1
Above normal volume [μL]	25
Above normal dilution (factor)	1

Results	
Decimals	1
Units	mmol/mol
Correlation factor-Offset	0.000
Correlation factor-Slope	1.000

Range	
Gender	All
Age	
SERUM	
URINE	
PLASMA	
CSF	
Whole blood	>=20 <=42
Gender	
Age	
SERUM	
URINE	
PLASMA	
CSF	
Whole blood	

Contaminants	
Contaminant 1	
Wash with	
Cycle	
Volume [μL]	
Contaminant 2	
Wash with	
Cycle	
Volume [μL]	
Contaminant 3	
Wash with	
Cycle	
Volume [μL]	

Calibrators details		
Calibrator list		Concentration
Cal. 1/Blank		*
Cal. 2		*
Cal. 3		
Cal. 4		
Cal. 5		
Cal. 6		
	Max delta abs.	
Cal. 1	0.005	
Cal. 2	0.005	
Cal. 3		
Cal. 4		
Cal. 5		
Cal. 6		
Drift limit [%]	0.8	
Calculations		
Model		X
Degree		1

\* Enter calibrator value

Calculation of HbA1c /Hb ratio is done automatically.

For values in percent according to DCCT/NGSP please enter 2.15 offset and a slope of 0.0915.

# HbA1c net FS (reaction 2 HbA1c)

## Application for whole blood samples

This application was set up and evaluated by DiaSys. It is based on the standard equipment at that time and does not apply to any equipment modifications undertaken by unqualified personnel

Identification	
This method is usable for analysis:	Yes
Twin reaction:	Yes
Name:	HbA1c net
Shortcut:	
Reagent barcode reference:	723
Host reference:	

Technic	
Type:	Fixed time kinetic
First reagent:[ $\mu$ L]	
Blank reagent	
Sensitive to light	
Second reagent:[ $\mu$ L]	50
Blank reagent	Yes
Sensitive to light	Yes
Main wavelength:[nm]	660
Secondary wavelength:[nm]	800
Polychromatic factor:	1.000
1 st reading time [min:sec]	5:00
Last reading time [min:sec]	10:00
Reaction way:	Increasing
Linear Kinetics	
Linearity: Maximum deviation [%]	
Fixed Time Kinetics	
Endpoint	
Prozone Limit [%]	

Reagents	
Decimals	3
Units	g/dL

Sample	
Diluent	System water
Hemolysis:	
Agent [ $\mu$ L]	NET(R951) 200
Cleaner	CLN A (R900)
Sample [ $\mu$ L]	10
Technical limits	HbA1c [mmol/mol]: 20 -150
Concentration technical limits-Lower	0.3
Concentration technical limits-Upper	2.0
SERUM	
Normal volume [ $\mu$ L]	25
Normal dilution (factor)	1
Below normal volume [ $\mu$ L]	25
Below normal dilution (factor)	1
Above normal volume [ $\mu$ L]	25
Above normal dilution (factor)	1
URIN	
Normal volume [ $\mu$ L]	25
Normal dilution (factor)	1
Below normal volume [ $\mu$ L]	25
Below normal dilution (factor)	1
Above normal volume [ $\mu$ L]	25
Above normal dilution (factor)	1
PLASMA	
Normal volume [ $\mu$ L]	25
Normal dilution (factor)	1
Below normal volume [ $\mu$ L]	25
Below normal dilution (factor)	1
Above normal volume [ $\mu$ L]	25
Above normal dilution (factor)	1
CSF	
Normal volume [ $\mu$ L]	25
Normal dilution (factor)	1
Below normal volume [ $\mu$ L]	25
Below normal dilution (factor)	1
Above normal volume [ $\mu$ L]	25
Above normal dilution (factor)	1
Whole blood	
Normal volume [ $\mu$ L]	25
Normal dilution (factor)	1
Below normal volume [ $\mu$ L]	25
Below normal dilution (factor)	1
Above normal volume [ $\mu$ L]	25
Above normal dilution (factor)	1

Results	
Decimals	1
Units	mmol/mol
Correlation factor-Offset	0.000
Correlation factor-Slope	1.000

Range	
Gender	All
Age	
SERUM	
URINE	
PLASMA	
CSF	
Whole blood	>=20 <=42
Gender	
Age	
SERUM	
URINE	
PLASMA	
CSF	
Whole blood	

Contaminants	
Contaminant 1	
Wash with	
Cycle	
Volume [ $\mu$ L]	
Contaminant 2	
Wash with	
Cycle	
Volume [ $\mu$ L]	
Contaminant 3	
Wash with	
Cycle	
Volume [ $\mu$ L]	

Calibrators details		
Calibrator list		Concentration
Cal. 1/Blank		*
Cal. 2		*
Cal. 3		
Cal. 4		
Cal. 5		
Cal. 6		
	Max delta abs.	
Cal. 1	0.005	
Cal. 2	0.015	
Cal. 3		
Cal. 4		
Cal. 5		
Cal. 6		
Drift limit [%]	0.8	
Calculations		
Model		X
Degree		1

\* Enter calibrator value

Calculation of HbA1c /Hb ratio is done automatically.

For values in percent according to DCCT/NGSP please enter 2.15 offset and a slope of 0.0915.

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IFCC: International Federation of Clinical Chemistry [1,2,7]  
DCCT: Diabetes Control and Complications Trial [3]  
NGSP: National Glycohemoglobin Standardization Program [4]

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1. Jeppsson JO, Kobold U, Barr J, Finke A et al. Approved IFCC reference method for the measurement of HbA1c in human blood. Clin Chem Lab Med 2002; 40: 78-89.
2. Hoelzel W, Weykamp C et al. IFCC Reference System for Measurement of Hemoglobin A1c in Human Blood and the National Standardization Schemes in the United States, Japan, and Sweden: A Method-Comparison Study. Clin Chem 2004; 50 (1): 166-74.
3. The Diabetes Control and Complications Trial Research Group. The effect of intensive treatment of diabetes in the development and progression of long-term complications in insulin-dependent diabetes mellitus. N Engl J Med.1993; 329: 977-86.
4. Little RR, Rohlfing CL, Wiedmeyer HM, Myers GL et al. The National Glycohemoglobin Standardization Program: A Five-Years Progress Report. Clin Chem 2001; 47: 1985-92.
5. Data on file at DiaSys Diagnostic Systems GmbH.
6. Panteghini M, John WG on behalf of the IFCC Scientific Division. Implementation of haemoglobin A1c results traceable to IFCCreference system: the way forward. Clin Chem Lab Med 2007; 45(8): 942-4.Nordin G., Dybkær R. Recommendation for term and measurement unit for "HbA1c". Clin Chem Lab Med 2007; 45(8): 1081-2.
7. Sacks DB. Translating Hemoglobin A1c into Average Blood Glucose: Implications for Clinical Chemistry. Clinical Chemistry 2008; 54: 1756-8.
8. Weykamp C. Carbamylated Hemoglobin Interference in Glyco-hemoglobin Assays. Clin Chem 1999; 45: 438-9.
9. Young DS. Effects of Drugs on Clinical Laboratory Tests. 5th. ed. Volume 1 and 2. Washington, DC: The American Association for Clinical Chemistry Press, 2000.
10. Thomas L. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998. p. 142-48.
11. Sacks DB. Carbohydrates. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3<sup>rd</sup> ed. Philadelphia: W.B. Saunders Company; 1999. p. 790-6.
12. Sacks DB. Carbohydrates. In: Burtis CA, Ashwood ER, Bruns DE, editors. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. 4th edition St. Louis Missouri: Elsevier Saunders; 2006; p. 878-884.
13. Sacks DB, Arnold M, Bakris GL, Bruns DE, AR Horvath et al. Guidelines and recommendations for laboratory analysis in the diagnosis and management of diabetes mellitus. Clin Chem 2011; 57(6): e1-e47.
14. Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. ClinChemLabMed 2007;45(9):1240-1243.
15. Ferri S, Kim S, Tsugawa W, Sode K. Review of Fructosyl Amino Acid Oxidase Engineering Research: A Glimpse into the Future of Hemoglobin A1c Biosensing. Journal of Diabetes Science and Technology 2009; 3(3): 585-592.

## Fabricant

DiaSys Diagnostic Systems GmbH  
Alte Strasse 9 65558 Holzheim (Allemagne)

## Hemoglobin TWIN test Application for HbA1c net FS

Reagent code 10 334

### Application for whole blood samples

This application was set up and evaluated by DiaSys. It is based on the standard equipment at that time and does not apply to any equipment modifications undertaken by unqualified personnel.

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Analysis item condition setting (M) [Subwindow button]		
Sample Type *	Serum	Urine
Reac. sample vol.	15	15
Diluent method	With dilution	With dilution
Undil. sample vol.	5	5
Diluent volume	100	100
Diluent position	**	**

\*Has to be defined for controls and Blank as well

\*\*Position number must be between 1 and 45 respectively 1 and 50 (depends on rotor)

Standards setting					
FV	#	BLK H	9.999	STD H	9.999
		BLK L	-9.999	STD L	-9.999

Multi-STD setting (subwindow button)								
Formula: Linear correction		Blank: Any value		Axis conv. : No convert.			Points: 2	
FV		Reac. smp. vol.	Dil. method	Undiluted smp. vol.	Diluent volume	Diluent position	STD-H	STD-L
BLK/0	#	15	With dilution	1.6	100	**	9.999	-9.999
1	#	15	With dilution	5.0	100	**	9.999	-9.999
2								
3								
4								
5								

#entered by user

For adjustment and maintenance of the sample probe, please contact your local distributor/service technician.

Individual results for hemoglobin and HbA1c in g/dL must not be used for diagnostic purposes.

## HbA1c TWIN test Application for HbA1c net FS

Reagent code 10 334

### Application for whole blood samples

This application was set up and evaluated by DiaSys. It is based on the standard equipment at that time and does not apply to any equipment modifications undertaken by unqualified personnel.

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Analysis test condition setting (M) [Subwindow button]		
Sample Type *	Serum	Urine
Reac. sample vol.	15	15
Diluent method	With dilution	With dilution
Undil. sample vol.	5	5
Diluent volume	100	100
Diluent position	**	**

\* Has to be defined for controls and Blank as well

\*\*Position number must be between 1 and 45 respectively 1 and 50 (depends on rotor)

Standards setting					
FV	#	BLK H	9.999	STD H	9.999
		BLK L	-9.999	STD L	-9.999

Multi-STD setting (subwindow button)								
Formula: Linear correction		Blank: Any value		Axis conv. : No convert.		Points: 2		
FV		Reac. smp. vol.	Dil. method	Undiluted. smp. vol.	Diluent volume	Diluent position	STD-H	STD-L
BLK/0	#	15	With dilution	1.6	100	**	9.999	-9.999
1	#	15	With dilution	5.0	100	**	9.999	-9.999
2								
3								
4								
5								

#entered by user

For adjustment and maintenance of the sample probe, please contact your local distributor/service technician.

Individual results for hemoglobin and HbA1c in g/dL must not be used for diagnostic purposes.