# CRP FS\*

## CODE CQN : HT

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1. Thompson D, Milford-Ward A, Whicher JT. The value of acute phase protein measurements in clinical practice. Ann Clin Biochem 1992; 29: 123-31.
2. Gabay C, Kushner I. Acute-phase proteins and other systemic responses to inflammation. N Engl J Med 1999; 340: 448-54.
3. Hansson LO, Lindquist L. C-reactive protein: its role in the diagnosis and follow-up of infectious diseases. Curr Opin Infect Diseases 1997; 10: 196-201.
4. Sipe JD. Acute-phase proteins in osteoarthritis. Semin Arthritis Rheum 1995; 25: 75-86.
5. Guder WG, Zawta B et al. The Quality of Diagnostic Samples. 1st ed. Darmstadt: GIT Verlag; 2001. p. 24-5.
6. Dati F, Schumann G, Thomas L, Aguzzi F, Baudner S, Bienvenu J et al. Consensus of a group of professional societies and diagnostic companies on guidelines for interim reference ranges for 14 proteins in serum based on the standardization against the IFCC/BCR/CAP reference material (CRM 470). Eur J Clin Chem Clin Biochem 1996; 34: 517-20.
7. Schlebusch H, Liappis N, Klein G. High sensitive CRP and creatinine: reference intervals from infancy to childhood. Poster präsentiert am AACC/CSCC; Juli/August 2001, Chicago, Illinois.
8. Young DS. Effects of Drugs on Clinical Laboratory Tests. 5th ed. Volume 1 and 2. Washington, DC: The American Association for Clinical Chemistry Press 2000.
9. Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. ClinChemLabMed 2007;45(9):1240-1243.

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1. Guder WG, Zawta B et al. The Quality of Diagnostic Samples. 1st ed. Darmstadt: GIT Verlag; 2001. p. 24 -5.
2. Young DS. Effects of Drugs on Clinical Laboratory Tests. 5th. ed. Volume 1 and 2. Washington, DC: The American Association for Clinical Chemistry Press, 2000.
3. Dati F, Schumann G, Thomas L, Aguzzi F, Baudner S, Bienvenu J et al. Consensus of a group of professional societies and diagnostic companies on guidelines for interim reference ranges for 14 proteins in serum based on the standardization against the IFCC/BCR/CAP reference material (CRM 470). Eur J Clin Chem Clin Biochem 1996; 34: 517-20.
4. Schlebusch H, Liappis N, Klein G. High sensitive CRP and creatinine: reference intervals from infancy to childhood. Poster presented at AACC/CSCC; July/August 2001, Chicago, Illinois.
5. Thompson D, Milford-Ward A, Whicher JT. The value of acute phase protein measurements in clinical practice. Ann Clin Biochem 1992; 29: 123-31.
6. Gabay C, Kushner I. Acute-phase proteins and other systemic responses to inflammation. N Engl J Med 1999; 340: 448-54.
7. Hansson LO, Lindquist L. C-reactive protein: its role in the diagnosis and follow-up of infectious diseases. Curr Opin Infect Diseases 1997; 10: 196-201.
8. Sipe JD. Acute-phase proteins in osteoarthritis. Semin Arthritis Rheum 1995; 25: 75-86.
9. Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. ClinChemLabMed 2007;45(9):1240-1243.

## Fabricant

DiaSys Diagnostic Systems GmbH  
Alte Strasse 9 65558 Holzheim Allemagne

## CRP FS

### Application for serum and plasma samples

This application was set up and evaluated by DiaSys. It is based on the standard equipment at that time and does not apply to any equipment modifications undertaken by unqualified personnel.

Identification	
This method is usable for analysis:	Yes
Twin reaction:	No
Name:	CRP
Shortcut:	
Reagent barcode reference:	706
Host reference:	706

Technic	
Type:	End point
First reagent:[μL]	180
Blank reagent	Yes
Sensitive to light	
Second reagent:[μL]	36
Blank reagent	No
Sensitive to light	
Main wavelength:[nm]	340
Secondary wavelength:[nm]	700
Polychromatic factor:	1.0000
1 st reading time [min:sec]	(04:24)
Last reading time [min:sec]	10:00
Reaction way:	Increasing
Linear Kinetics	
Substrate depletion: Absorbance limit	
Linearity: Maximum deviation [%]	
Fixed Time Kinetics	
Substrate depletion: Absorbance limit	
Endpoint	
Stability: Largest remaining slope	
Prozone Limit [%]	

Reagents	
Decimals	
Units	

Sample	
Diluent	DIL A (NaCl)
Hemolysis:	
Agent [μL]	0 (no hemolysis)
Cleaner	
Sample [μL]	0
Technical limits	
Concentration technical limits-Lower	2.0000
Concentration technical limits-Upper	250.0000
SERUM	
Normal volume [μL]	11.0
Normal dilution (factor)	1
Below normal volume [μL]	
Below normal dilution (factor)	
Above normal volume [μL]	2.0
Above normal dilution (factor)	1
URINE	
Normal volume [μL]	11.0
Normal dilution (factor)	1
Below normal volume [μL]	
Below normal dilution (factor)	
Above normal volume [μL]	2.0
Above normal dilution (factor)	1
PLASMA	
Normal volume [μL]	11.0
Normal dilution (factor)	1
Below normal volume [μL]	
Below normal dilution (factor)	
Above normal volume [μL]	2.0
Above normal dilution (factor)	1
CSF	
Normal volume [μL]	11.0
Normal dilution (factor)	1
Below normal volume [μL]	
Below normal dilution (factor)	
Above normal volume [μL]	2.0
Above normal dilution (factor)	1
Whole blood	
Normal volume [μL]	11.0
Normal dilution (factor)	1
Below normal volume [μL]	
Below normal dilution (factor)	
Above normal volume [μL]	2.0
Above normal dilution (factor)	1

Results	
Decimals	2
Units	mg/L
Correlation factor-Offset	0.0000
Correlation factor-Slope	1.0000

Range	
Gender	All
Age	
SERUM	>= <=5.00
URINE	
PLASMA	>= <=5.00
CSF	
Whole blood	
Gender	
Age	
SERUM	
URINE	
PLASMA	
CSF	
Whole blood	

Contaminants	
Please refer to r910 Carryover Pair Table	

Calibrators details	
Calibrator list	Concentration
Cal. 1/Blank	0
Cal. 2	*
Cal. 3	*
Cal. 4	*
Cal. 5	*
Cal. 6	*
Max delta abs.	
Cal. 1	0.0100
Cal. 2	0.0100
Cal. 3	0.0100
Cal. 4	0.0100
Cal. 5	0.0150
Cal. 6	0.0300
Drift limit [%]	2.00

Calculations	
Model	Cubic Spline
Degree	

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1. Guder WG, Zawta B et al. The Quality of Diagnostic Samples. 1st ed. Darmstadt: GIT Verlag; 2001. p. 24 -5.
2. Dati F, Schumann G, Thomas L, Aguzzi F, Baudner S, Bienvenu J et al. Consensus of a group of professional societies and diagnostic companies on guidelines for interim reference ranges for 14 proteins in serum based on the standardization against the IFCC/BCR/CAP reference material (CRM 470). Eur J Clin Chem Clin Biochem 1996; 34: 517-20.
3. Thompson D, Milford-Ward A, Whicher JT. The value of acute phase protein measurements in clinical practice. Ann Clin Biochem 1992; 29: 123-31.
4. Gabay C, Kushner I. Acute-phase proteins and other systemic responses to inflammation. N Engl J Med 1999; 340: 448-54.
5. Hansson LO, Lindquist L. C-reactive protein: its role in the diagnosis and follow-up of infectious diseases. Curr Opin Infect Diseases 1997; 10: 196-201.
6. Sipe JD. Acute-phase proteins in osteoarthritis. Semin Arthritis Rheum 1995; 25: 75-86.
7. Young DS. Effects of Drugs on Clinical Laboratory Tests. 5th ed. Volume 1 and 2. Washington, DC: The American Association for Clinical Chemistry Press 2000.
8. Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. ClinChemLabMed 2007;45(9):1240-1243.

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Analysis Test Condition Setting (M)		
Sample Type	Serum	Urine
Reac. sample vol.	4.8	4.8
Diluent method	No dil	No dil
Undil. sample vol.	0	0
Diluent volume	0	0
Diluent position	0	0

Prozone	
Prozone form	No
Prozone limit	9.999
Prozone judge	Upper limit
Judge limit	9.999
M-DET.P.m	0
M-DET.P.n	0
S-DET.P.p	0
S-DET.P.r	0

MULTI-STD Setting								
Formula	Logit Log 3	Axis Conv	No conv					
Blank	Blank is 0	Points	6					
	FV	Reac. smp. vol.	Dil. method	Dil. smp. vol.	Diluent vol.	Diluent pos.	STD H	STD L
BLK	#	4.8	No dil	0	0	0	9.999	-9.999
1	#	4.8	No dil	0	0	0	9.999	-9.999
2	#	4.8	No dil	0	0	0	9.999	-9.999
3	#	4.8	No dil	0	0	0	9.999	-9.999
4	#	4.8	No dil	0	0	0	9.999	-9.999
5	#	4.8	No dil	0	0	0	9.999	-9.999

# entered by user