# Myoglobine FS \*

CODE CQN : HT

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1. Stone MJ, Willerson JT, Gomez-Sanchez CE, Waterman MR. Radioimmunoassay of myoglobin in human serum. Results in patients with acute myocardial infarction. J Clin Invest 1975; 56: 1334-9.
2. Bhayana V, Henderson AR. Biochemical markers of myocardial damage. Clin Biochem 1995; 28: 1-29.
3. Mair J, Artner-Dworzak E, Lechleitner P, Morass B, Smidt J, Wagner I et al. Early diagnosis of acute myocardial infarction by a newly developed rapid immunoturbidimetric assay for myoglobin. Br Heart J 1992; 68: 462-8.
4. Zaninotto M, Altinier S, Lachin M, Celegon L, Plebani M. Strategies for the early diagnosis of acute myocardial infarction using biochemical markers. Am J Pathol 1999; 111: 399-405.
5. De Winter RJ, Koster RW, Sturk A, Sanders GT. Value of myoglobin, troponin T and CK-MB mass in ruling out myocardial infarction in the emergency room. Circulation 1995; 92: 3401-7.
6. Laperche T, Steg PG, Dehoux M, Benessiano I, Grollier G, Aliot E et al. A study of biochemical markers of reperfusion early after thrombolysis for acute myocardial infarction. Circulation 1995; 92: 2079-86.
7. Guder WG, Zawta B et al. The Quality of Diagnostic Samples. 1st ed. Darmstadt: GIT Verlag; 2001. p. 38-9.
8. Young DS. Effects of Drugs on Clinical Laboratory Tests. 5th ed. Volume 1 and 2. Washington, DC: The American Association for Clinical Chemistry Press 2000.
9. Baum H, Booksteegers P, Steinbeck G, Neumeier D. A rapid assay for the quantification of myoglobin: evaluation and diagnostic relevance in the diagnosis of acute myocardial infarction. Eur J Clin Chem Biochem 1994; 32: 853-8.
10. Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: Mechanisms, detection and prevention. Clin Chem Lab Med 2007; 45(9): 1240-1243.

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- Guder WG, Zawta B et al. The Quality of Diagnostic Samples. 1st ed. Darmstadt: GIT Verlag; 2001. p. 38-9.
- Young DS. Effects of Drugs on Clinical Laboratory Tests. 5th. ed. Volume 1 and 2. Washington, DC: The American Association for Clinical Chemistry Press, 2000

3. Mair J, Artner-Dworzak E, Lechleitner P, Morass B, Smidt J, Wagner I et al. Early diagnosis of acute myocardial infarction by a newly developed rapid immunoturbidimetric assay for myoglobin. *Br Heart J* 1992; 68: 462-8.
4. Stone MJ, Willerson JT, Gomez-Sanchez CE, Waterman MR. Radioimmunoassay of myoglobin in human serum. Results in patients with acute myocardial infarction. *J Clin Invest* 1975; 56: 1334-9.
5. Bhayana V, Henderson AR. Biochemical markers of myocardial damage. *Clin Biochem* 1995; 28: 1-29.
6. Zaninotto M, Altinier S, Lachin M, Celegon L, Plebani M. Strategies for the early diagnosis of acute myocardial infarction using biochemical markers. *Am J Pathol* 1999; 111: 399-405.
7. De Winter RJ, Koster RW, Sturk A, Sanders GT. Value of myoglobin, troponin T and CK-MB mass in ruling out myocardial infarction in the emergency room. *Circulation* 1995; 92: 3401-7.
8. Laperche T, Steg PG, Dehoux M, Benessiano I, Grollier G, Aliot E et al. A study of biochemical markers of reperfusion early after thrombolysis for acute myocardial infarction. *Circulation* 1995; 92: 2079-86.
9. Baum H, Booksteegers P, Steinbeck G, Neumeier D. A rapid assay for the quantification of myoglobin: evaluation and diagnostic relevance in the diagnosis of acute myocardial infarction. *Eur J Clin Chem Biochem* 1994; 32: 853-8.
10. Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: Mechanisms, detection and prevention. *Clin Chem Lab Med* 2007; 45(9): 1240–1243.

#### Fabricant

DiaSys Diagnostic Systems GmbH  
 Alte Strasse 9 65558 Holzheim Allemagne

## Myoglobin FS

### Application for serum and plasma samples

This application was set up and evaluated by DiaSys. It is based on the standard equipment at that time and does not apply to any equipment modifications undertaken by unqualified personnel.

Identification	
This method is usable for analysis:	Yes
Twin reaction:	No
Name:	MYO
Shortcut:	
Reagent barcode reference:	711
Host reference:	711

Technic	
Type:	Fixed time kinetic
First reagent:[μL]	150
Blank reagent	Yes
Sensitive to light	
Second reagent:[μL]	50
Blank reagent	No
Sensitive to light	
Main wavelength:[nm]	570
Secondary wavelength:[nm]	
Polychromatic factor:	
1 st reading time [min:sec]	05:00
Last reading time [min:sec]	08:00
Reaction way:	Increasing
Linear Kinetics	
Substrate depletion: Absorbance limit	
Linearity: Maximum deviation [%]	
Fixed Time Kinetics	
Substrate depletion: Absorbance limit	
Endpoint	
Stability: Largest remaining slope	
Prozone Limit [%]	

Reagents	
Decimals	
Units	

Sample	
Diluent	DIL A (NaCl)
Hemolysis:	
Agent [μL]	0 (no hemolysis)
Cleaner	
Sample [μL]	0
Technical limits	
Concentration technical limits-Lower	6.0000
Concentration technical limits-Upper	600.0000
SERUM	
Normal volume [μL]	5.0
Normal dilution (factor)	1
Below normal volume [μL]	
Below normal dilution (factor)	
Above normal volume [μL]	5.0
Above normal dilution (factor)	6
URINE	
Normal volume [μL]	5.0
Normal dilution (factor)	1
Below normal volume [μL]	
Below normal dilution (factor)	
Above normal volume [μL]	5.0
Above normal dilution (factor)	6
PLASMA	
Normal volume [μL]	5.0
Normal dilution (factor)	1
Below normal volume [μL]	
Below normal dilution (factor)	
Above normal volume [μL]	5.0
Above normal dilution (factor)	6
CSF	
Normal volume [μL]	5.0
Normal dilution (factor)	1
Below normal volume [μL]	
Below normal dilution (factor)	
Above normal volume [μL]	5.0
Above normal dilution (factor)	6
Whole blood	
Normal volume [μL]	5.0
Normal dilution (factor)	1
Below normal volume [μL]	
Below normal dilution (factor)	
Above normal volume [μL]	5.0
Above normal dilution (factor)	6

Results	
Decimals	1
Units	μg/L
Correlation factor-Offset	0.0000
Correlation factor-Slope	1.0000

Range	
Gender	All
Age	
SERUM	>= <=70
URINE	
PLASMA	>= <=70
CSF	
Whole blood	
Gender	
Age	
SERUM	
URINE	
PLASMA	
CSF	
Whole blood	

Contaminants	
Please refer to r910 Carryover Pair Table	

Calibrators details	
Calibrator list	Concentration
Cal. 1/Blank	0
Cal. 2	*
Cal. 3	*
Cal. 4	*
Cal. 5	*
Cal. 6	
Max delta abs.	
Cal. 1	0.0050
Cal. 2	0.0100
Cal. 3	0.0100
Cal. 4	0.0100
Cal. 5	0.0100
Cal. 6	
Drift limit [%]	2.00

Calculations	
Model	Akima Spline
Degree	

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1. Bhatnagar D, Durrington PN. Measurement and clinical significance of apolipoproteins A-I and B. In: Rifai N, Warnick GR, Dominiczak MH, eds. Handbook of lipoprotein testing. Washington: AACC Press, 1997: p. 177-98.
2. Rifai N, Bachorik PS, Albers JJ. Lipids, lipoproteins and apolipoproteins. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 809-61.
3. Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. ClinChemLabMed 2007;45(9):1240-1243.
4. Guder WG, Zawta B et al. The Quality of Diagnostic Samples. 1st ed. Darmstadt: GIT Verlag; 2001; p. 18-9.
5. Young DS. Effects of Drugs on Clinical Laboratory Tests. 5th ed. Volume 1 and 2. Washington, DC: The American Association for Clinical Chemistry Press 2000.
6. Young DS. Effects on Clinical Laboratory Tests - Drugs Disease, Herbs & Natural Products, <https://clinfx.wiley.com/aaccweb/aacc/>, accessed on December 2020. Published by AACC Press and John Wiley and Sons, Inc.
7. Jungner I, Marcovina SM, Walldius G, Holme I, Kolar W, Steiner E. Apolipoprotein B and A-I values in 147576 Swedish males and females, standardized according to the World Health Organization-International Federation of Clinical Chemistry First International Reference Materials. Clin Chem 1998; 44: 1641-9.

DiaSys Diagnostic Systems GmbH  
Alte Strasse 9 65558 Holzheim Allemagne  
[www.diasys-diagnostics.com](http://www.diasys-diagnostics.com)

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Analysis Test Condition Setting (M)		
Sample Type	Serum	Urine
Reac. sample vol.	1.0	1.0
Diluent method	No dil	No dil
Undil. sample vol.	0	0
Diluent volume	0	0
Diluent position	0	0

Prozone	
Prozone form	No
Prozone limit	9.999
Prozone judge	Upper limit
Judge limit	9.999
M-DET.P.m	0
M-DET.P.n	0
S-DET.P.p	0
S-DET.P.r	0

MULTI-STD Setting								
Formula	Logit Log 3	Axis Conv	No conv					
Blank	Blank is 0	Points	6					
	FV	Reac. smp. vol.	Dil. method	Dil. smp. vol.	Diluent vol.	Diluent pos.	STD H	STD L
BLK	#	1.0	No dil	0	0	0	9.999	-9.999
1	#	1.0	With dil	10	40	0	9.999	-9.999
2	#	1.0	With dil	20	30	0	9.999	-9.999
3	#	1.5	With dil	20	30	0	9.999	-9.999
4	#	2.0	With dil	20	30	0	9.999	-9.999
5	#	1.0	No dil	0	0	0	9.999	-9.999

# entered by user