# Lp(a) 21 FS \*

CODE CQN : HT

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\*Northwest Lipid Research Laboratories

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1. Rifai N, Bachorik PS, Albers JJ. Lipids, lipoproteins and apolipoproteins. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 809-61.
2. Marcovina SM, Koschinsky ML. Lipoprotein (a): Structure, measurement and clinical significance. In: Rifai N, Warnick GR, Dominiczak MH, eds. Handbook of lipoprotein testing. Washington: AACC Press; 1997. p.283-313.
3. Guder WG, Zawta B et al. The Quality of Diagnostic Samples. 1st ed. Darmstadt: GIT Verlag; 2001. p. 36-37.
4. Young DS. Effects of Drugs on Clinical Laboratory Tests. 5th ed. Volume 1 and 2. Washington, DC: The American Association for Clinical Chemistry Press 2000.
5. Riesen WF. Lipid metabolism. In: Thomas L, editor. Clinical laboratory diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998. p. 174-5.
6. Marcovina SM, Koschinsky ML et al. Report of the national heart, lung, and blood institute workshop of Lipoprotein(a) and cardiovascular disease: recent advances and future directions. Clin Chem 2003; 49(11): 1785-96.
7. Nordestgaard BG, Chapman MJ, Ginsberg HN. Lipoprotein (a): EAS Recommendations for Screening, Desirable Levels and Management. The European Atherosclerosis Society (EAS) Consensus Panel 2012.
8. Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: Mechanisms, detection and prevention. Clin Chem Lab Med 2007; 45(9): 1240-1243.

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1. Guder WG, Zawta B et al. The Quality of Diagnostic Samples. 1<sup>st</sup> ed. Darmstadt: GIT Verlag; 2001; p. 36-7.
2. Young DS. Effects of Drugs on Clinical Laboratory Tests. 5th. ed. Volume 1 and 2. Washington, DC: The American Association for Clinical Chemistry Press, 2000.
3. Nordestgaard BG, Chapman MJ, Ginsberg HN. Lipoprotein (a): EAS Recommendations for Screening, Desirable Levels and Management. The European Atherosclerosis Society (EAS) Consensus Panel 2012.
4. Riesen WF. Lipid metabolism. In: Thomas L, editor. Clinical laboratory diagnostics. 1<sup>st</sup> ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998. p. 174-5.
5. Rifai N, Bachorik PS, Albers JJ. Lipids, lipoproteins and apolipoproteins. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3<sup>rd</sup> ed. Philadelphia: W.B Saunders Company; 1999. p. 809-61.
6. Marcovina SM, Koschinsky ML. Lipoprotein (a): Structure, measurement and clinical significance. In: Rifai N, Warnick GR, Dominiczak MH, eds. Handbook of lipoprotein testing. Washington: AACC Press; 1997. p. 283-313.
7. Marcovina SM, Koschinsky ML et al. Report of the national heart, lung, and blood institute workshop of Lipoprotein(a) and cardiovascular disease: recent advances and future directions. Clin Chem 2003; 49(11): 1785-96.
8. Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: Mechanisms, detection and prevention. Clin Chem Lab Med 2007; 45(9): 1240-1243.

#### Fabricant

**IVD** **CE** DiaSys Diagnostic Systems GmbH  
Alte Strasse 9 65558 Holzheim Allemagne

## Lp(a) 21 FS

### Application for serum and plasma samples

This application was set up and evaluated by DiaSys. It is based on the standard equipment at that time and does not apply to any equipment modifications undertaken by unqualified personnel.

| Identification                      |     |
|-------------------------------------|-----|
| This method is usable for analysis: | Yes |
| Twin reaction:                      | No  |
| Name:                               | LPA |
| Shortcut:                           |     |
| Reagent barcode reference:          | 710 |
| Host reference:                     | 710 |

| Technic                               |                    |
|---------------------------------------|--------------------|
| Type:                                 | Fixed time kinetic |
| First reagent:[μL]                    | 120                |
| Blank reagent                         | Yes                |
| Sensitive to light                    |                    |
| Second reagent:[μL]                   | 60                 |
| Blank reagent                         | No                 |
| Sensitive to light                    |                    |
| Main wavelength:[nm]                  | 700                |
| Secondary wavelength:[nm]             |                    |
| Polychromatic factor:                 |                    |
| 1 st reading time [min:sec]           | 04:48              |
| Last reading time [min:sec]           | 10:00              |
| Reaction way:                         | Increasing         |
| Linear Kinetics                       |                    |
| Substrate depletion: Absorbance limit |                    |
| Linearity: Maximum deviation [%]      |                    |
| Fixed Time Kinetics                   |                    |
| Substrate depletion: Absorbance limit |                    |
| Endpoint                              |                    |
| Stability: Largest remaining slope    |                    |
| Prozone Limit [%]                     |                    |

| Reagents |  |
|----------|--|
| Decimals |  |
| Units    |  |

| Sample                               |                  |
|--------------------------------------|------------------|
| Diluent                              | DIL A (NaCl)     |
| Hemolysis:                           |                  |
| Agent [μL]                           | 0 (no hemolysis) |
| Cleaner                              |                  |
| Sample [μL]                          | 0                |
| Technical limits                     |                  |
| Concentration technical limits-Lower | 2.0000           |
| Concentration technical limits-Upper | 110.0000         |
| SERUM                                |                  |
| Normal volume [μL]                   | 4.0              |
| Normal dilution (factor)             | 1                |
| Below normal volume [μL]             |                  |
| Below normal dilution (factor)       |                  |
| Above normal volume [μL]             | 2.0              |
| Above normal dilution (factor)       | 1                |
| URINE                                |                  |
| Normal volume [μL]                   | 4.0              |
| Normal dilution (factor)             | 1                |
| Below normal volume [μL]             |                  |
| Below normal dilution (factor)       |                  |
| Above normal volume [μL]             | 2.0              |
| Above normal dilution (factor)       | 1                |
| PLASMA                               |                  |
| Normal volume [μL]                   | 4.0              |
| Normal dilution (factor)             | 1                |
| Below normal volume [μL]             |                  |
| Below normal dilution (factor)       |                  |
| Above normal volume [μL]             | 2.0              |
| Above normal dilution (factor)       | 1                |
| CSF                                  |                  |
| Normal volume [μL]                   | 4.0              |
| Normal dilution (factor)             | 1                |
| Below normal volume [μL]             |                  |
| Below normal dilution (factor)       |                  |
| Above normal volume [μL]             | 2.0              |
| Above normal dilution (factor)       | 1                |
| Whole blood                          |                  |
| Normal volume [μL]                   | 4.0              |
| Normal dilution (factor)             | 1                |
| Below normal volume [μL]             |                  |
| Below normal dilution (factor)       |                  |
| Above normal volume [μL]             | 2.0              |
| Above normal dilution (factor)       | 1                |

| Results                   |        |
|---------------------------|--------|
| Decimals                  | 2      |
| Units                     | mg/dL  |
| Correlation factor-Offset | 0.0000 |
| Correlation factor-Slope  | 1.0000 |

| Range       |            |
|-------------|------------|
| Gender      | All        |
| Age         |            |
| SERUM       | >= <=30.00 |
| URINE       |            |
| PLASMA      | >= <=30.00 |
| CSF         |            |
| Whole blood |            |
| Gender      |            |
| Age         |            |
| SERUM       |            |
| URINE       |            |
| PLASMA      |            |
| CSF         |            |
| Whole blood |            |

| Contaminants                              |  |
|---|--|
| Please refer to r910 Carryover Pair Table |  |

| Calibrators details |               |
|---------------------|---------------|
| Calibrator list     | Concentration |
| Cal. 1/Blank        | 0             |
| Cal. 2              | *             |
| Cal. 3              | *             |
| Cal. 4              | *             |
| Cal. 5              | *             |
| Cal. 6              | *             |
| Max delta abs.      |               |
| Cal. 1              | 0.0100        |
| Cal. 2              | 0.0100        |
| Cal. 3              | 0.0100        |
| Cal. 4              | 0.0100        |
| Cal. 5              | 0.0250        |
| Cal. 6              | 0.0400        |
| Drift limit [%]     | 2.00          |

| Calculations |              |
|--------------|--------------|
| Model        | Cubic Spline |
| Degree       |              |

\* Enter calibrator value