Cystatin C FS* (Cystatine C FS*)

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4. Le Bricon T., Leblanc I et al. Evaluation of renal function in
intensive care: plasma cystatin C vs. creatinine and derived
glomerular filtration rates Clin Chem Lab Med 2005;
43(9):953-957.
5. Le Bricon T., Thervet E., Benlakehal M. et al. Changes in
Plasma cystatin C after renal transplantation and acute
rejection in adults. Clin Chem 1999; 45(12):2243-9
6. Ustundag Y., Samsar U. et al. Analysis of glomerular filtration
rate, serum cystatin C levels, and resistive index values in
cirrhosis patients. Clin Chem Lab Med 2007; 45(7):890-94.
7. Stabuc B., Vrhovec L. et al. Improved prediction of decreased
creatinine clearance by serum cystatin C : Use in cancer
patients before and during chemotherapy. Clin Chem 2000;
46(2):193-7.
8. Pucci L., Triscornia S., Lucchesi D. et al. Cystatin C and
estimates of renal function: searching for a better measure of
kidney function in diabetic patients. Clin Chem 2007;
53(3):480-8.
9. Strevens H, Wide-Svensson D, Grubb A. Serum cystatin C is
a better marker for preeclampsia than serum creatinine or
serum urate. Scand J Clin Lab Invest 2001; 61:575-80.
10. Bakker AJ, Mücke M. Gammopathy interference in clinical
chemistry assays: mechanisms, detection and prevention.
ClinChemLabMed 2007;45(9):1240-1243.
11. Guder WG, Zawta B. et al. The Quality of Diagnostic
Samples. 1st ed. Darmstadt: GIT Verlag; 2001; p. 24-5.
12. Wiesli P., Schwegler B. et al. Serum cystatin C is sensitive to
small changes in thyroid function. Clinica Chimica Acta 2003;
338: 87-90.
13. Young DS. Effects of Drugs on Clinical Laboratory Tests. 5th
ed. Volume 1 and 2. Washington, DC: The American
Association for Clinical Chemistry Press 2000.
14. Young DS. Effects on Clinical Laboratory Tests - Drugs
Disease, Herbs & Natural Products, <https://clinf.wiley.com/aaccweb/aacc/>, accessed in August 2021. Published by
AACC Press and John Wiley and Sons, Inc.
15. Bökenkamp A, Grabensee A, Stoffel-Wagner B, Hasan C,
Henne T, Offner G, Lentze MJ. The β 2-microglobulin/cystatin
C ratio: a potential marker of post-transplant
lymphoproliferative disease. Clin Nephrol 2001; 58: 417-22.
16. Randers E, Krue S, Erlandsen EJ, Danielsen H, Hansen H,
Hansen LG. Reference interval for serum cystatin C in
children. ClinChem 1999; 45:1856-60.
17. Finney H, Newman DJ, Price CP. Adult reference ranges for
serum cystatin C, creatinine and predicted creatinine
clearance. Ann Clin Biochem 2000; 37: 49-59.

DiaSys Diagnostic Systems GmbH
Alte Strasse 9 65558 Holzheim
Germany
www.diasys-diagnostics.com

* Fluid Stable = Liquide & Stable