# Immunoglobuline G FS\*

CODE CQN : HT

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1. Thomas L. Clinical Laboratory Diagnostics. 1<sup>st</sup> ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998. p. 667-78.
2. Johnson AM, Rohlfes EM, Silverman LM. Proteins. In: Burtis CA, Ashwood ER. editors. Tietz textbook of clinical chemistry. 3<sup>rd</sup> ed. Philadelphia: W. B. Saunders Company; 1999. p. 507-12.
3. Bartl R, Hoechtlen-Vollmar W, Thomas L. Monoclonal immunoglobulins. In: Thomas L. Clinical Laboratory Diagnostics. 1<sup>st</sup> ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998. p. 742-58.
4. Guder WG, Narayanan S et al. List of Analytes; Preanalytical Variables. 1<sup>st</sup> ed. Darmstadt: Git Verlag, 1996: 16-7.
5. Young DS. Effects of Drugs on Clinical Laboratory Tests. 5th ed. Volume 1 and 2. Washington, DC: The American Association for Clinical Chemistry Press 2000.
6. Dati F, Schumann G, Thomas L, Aguzzi F, Baudner S, Bienvenu J et al. Consensus of a group of professional societies and diagnostic companies on guidelines for interim reference ranges for 14 proteins in serum based on the standardization against the IFCC/BCR/CAP reference material (CRM 470). Eur J Clin Chem Clin Biochem 1996; 34: 517-20.
7. Heil R, Koberstein R, Zawta B. Referenzbereiche für Kinder und Erwachsene. Roche Diagnostics 2004. p. 46–47.
8. Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: Mechanisms, detection and prevention. Clin Chem Lab Med 2007; 45(9): 1240–1243.

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1. Guder WG, Narayanan S et al. List of Analytes; Preanalytical Variables. 1st ed. Darmstadt: Git Verlag, 1996: 16-7.
2. Young DS. Effects of Drugs on Clinical Laboratory Tests. 5th. ed. Volume 1 and 2. Washington, DC: The American Association for Clinical Chemistry Press, 2000.
3. Dati F, Schumann G, Thomas L, Aguzzi F, Baudner S, Bienvenu J et al. Consensus of a group of professional societies and diagnostic companies on guidelines for interim reference ranges for 14 proteins in serum based on the standardization against the IFCC/BCR/CAP reference material (CRM 470). Eur J Clin Chem Clin Biochem 1996;34:517-20.
4. Heil R, Koberstein R, Zawta B. Referenzbereiche für Kinder und Erwachsene. Roche Diagnostics 2004. p. 46 - 47.
5. Thomas L. Clinical Laboratory Diagnostics. 1<sup>st</sup> ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998. p. 667-78.
6. Johnson AM, Rohlfes EM, Silverman LM. Proteins. In: Burtis CA, Ashwood ER. editors. Tietz textbook of clinical chemistry. 3<sup>rd</sup> ed. Philadelphia: W. B. Saunders Company; 1999. p. 507-12.
7. Bartl R, Hoechtlen-Vollmar W, Thomas L. Monoclonal immunoglobulins. In: Thomas L. Clinical Laboratory Diagnostics. 1<sup>st</sup> ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998. p. 742-58.
8. Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: Mechanism, detection and prevention. Clin Chem Lab Med 2007; 45(9): 1240–1243.

## Fabricant

DiaSys Diagnostic Systems GmbH  
Alte Strasse 9 65558 Holzheim  
Allemagne

# Immunoglobulin G FS

## Application for serum and plasma samples

This application was set up and evaluated by DiaSys. It is based on the standard equipment at that time and does not apply to any equipment modifications undertaken by unqualified personnel.

| Identification                      |     |
|-------------------------------------|-----|
| This method is usable for analysis: | Yes |
| Twin reaction:                      | No  |
| Name:                               | IGG |
| Shortcut:                           |     |
| Reagent barcode reference:          | 715 |
| Host reference:                     | 715 |

| Technic                               |            |
|---------------------------------------|------------|
| Type:                                 | End point  |
| First reagent:[μL]                    | 200        |
| Blank reagent                         | Yes        |
| Sensitive to light                    |            |
| Second reagent:[μL]                   | 40         |
| Blank reagent                         | No         |
| Sensitive to light                    |            |
| Main wavelength:[nm]                  | 570        |
| Secondary wavelength:[nm]             |            |
| Polychromatic factor:                 |            |
| 1 st reading time [min:sec]           | (04:24)    |
| Last reading time [min:sec]           | 06:00      |
| Reaction way:                         | Increasing |
| Linear Kinetics                       |            |
| Substrate depletion: Absorbance limit |            |
| Linearity: Maximum deviation [%]      |            |
| Fixed Time Kinetics                   |            |
| Substrate depletion: Absorbance limit |            |
| Endpoint                              |            |
| Stability: Largest remaining slope    |            |
| Prozone Limit [%]                     |            |

| Reagents |  |
|----------|--|
| Decimals |  |
| Units    |  |

| Sample                               |                  |
|--------------------------------------|------------------|
| Diluent                              | DIL A (NaCl)     |
| Hemolysis:                           |                  |
| Agent [μL]                           | 0 (no hemolysis) |
| Cleaner                              |                  |
| Sample [μL]                          | 0                |
| Technical limits                     |                  |
| Concentration technical limits-Lower | 6.0000           |
| Concentration technical limits-Upper | 3282.0000        |
| SERUM                                |                  |
| Normal volume [μL]                   | 2.0              |
| Normal dilution (factor)             | 1                |
| Below normal volume [μL]             |                  |
| Below normal dilution (factor)       |                  |
| Above normal volume [μL]             | 2.0              |
| Above normal dilution (factor)       | 6                |
| URINE                                |                  |
| Normal volume [μL]                   | 2.0              |
| Normal dilution (factor)             | 1                |
| Below normal volume [μL]             |                  |
| Below normal dilution (factor)       |                  |
| Above normal volume [μL]             | 2.0              |
| Above normal dilution (factor)       | 6                |
| PLASMA                               |                  |
| Normal volume [μL]                   | 2.0              |
| Normal dilution (factor)             | 1                |
| Below normal volume [μL]             |                  |
| Below normal dilution (factor)       |                  |
| Above normal volume [μL]             | 2.0              |
| Above normal dilution (factor)       | 6                |
| CSF                                  |                  |
| Normal volume [μL]                   | 2.0              |
| Normal dilution (factor)             | 1                |
| Below normal volume [μL]             |                  |
| Below normal dilution (factor)       |                  |
| Above normal volume [μL]             | 2.0              |
| Above normal dilution (factor)       | 6                |
| Whole blood                          |                  |
| Normal volume [μL]                   | 2.0              |
| Normal dilution (factor)             | 1                |
| Below normal volume [μL]             |                  |
| Below normal dilution (factor)       |                  |
| Above normal volume [μL]             | 2.0              |
| Above normal dilution (factor)       | 6                |

| Results                   |        |
|---------------------------|--------|
| Decimals                  | 1      |
| Units                     | mg/dL  |
| Correlation factor-Offset | 0.0000 |
| Correlation factor-Slope  | 1.0000 |

| Range       |                  |
|-------------|------------------|
| Gender      | All              |
| Age         |                  |
| SERUM       | >=700.0 <=1600.0 |
| URINE       |                  |
| PLASMA      | >=700.0 <=1600.0 |
| CSF         |                  |
| Whole blood |                  |
| Gender      |                  |
| Age         |                  |
| SERUM       |                  |
| URINE       |                  |
| PLASMA      |                  |
| CSF         |                  |
| Whole blood |                  |

| Contaminants                              |  |
|---|--|
| Please refer to r910 Carryover Pair Table |  |

| Calibrators details |               |
|---------------------|---------------|
| Calibrator list     | Concentration |
| Cal. 1/Blank        | 0             |
| Cal. 2              | *             |
| Cal. 3              | *             |
| Cal. 4              | *             |
| Cal. 5              | *             |
| Cal. 6              | *             |
| Max delta abs.      |               |
| Cal. 1              | 0.0100        |
| Cal. 2              | 0.0100        |
| Cal. 3              | 0.0150        |
| Cal. 4              | 0.0200        |
| Cal. 5              | 0.0300        |
| Cal. 6              | 0.0500        |
| Drift limit [%]     | 2.00          |

| Calculations |              |
|--------------|--------------|
| Model        | Akima Spline |
| Degree       |              |

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- Guder WG, Narayanan S et al. List of Analytes; Preanalytical Variables. 1st ed. Darmstadt: Git Verlag, 1996: 16-7.
- Dati F, Schumann G, Thomas L, Aguzzi F, Baudner S, Bienvenu J et al. Consensus of a group of professional societies and diagnostic companies on guidelines for interim reference ranges for 14 proteins in serum based on the standardization against the IFCC/BCR/CAP reference material (CRM 470). Eur J Clin Chem Clin Biochem 1996;34:517-20.
- Heil R, Koberstein R, Zawta B. Referenzbereiche für Kinder und Erwachsene. Roche Diagnostics 2004. p. 46 - 47.
- Thomas L. Clinical Laboratory Diagnostics. 1<sup>st</sup> ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998. p. 667-78.
- Johnson AM, Rohlfis EM, Silverman LM. Proteins. In: Burtis CA, Ashwood ER. editors. Tietz textbook of clinical chemistry. 3<sup>rd</sup> ed. Philadelphia: W. B. Saunders Company; 1999. p. 507-12.
- Bartl R, Hoechtlen-Vollmar W, Thomas L. Monoclonal immunoglobulins. In: Thomas L. Clinical Laboratory Diagnostics. 1<sup>st</sup> ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998. p. 742-58.
- Young DS. Effects of Drugs on Clinical Laboratory Tests. 5th ed. Volume 1 and 2. Washington, DC: The American Association for Clinical Chemistry Press 2000.
- Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: Mechanism, detection and prevention. Clin Chem Lab Med 2007; 45(9): 1240–1243.

#### Fabricant

DiaSys Diagnostic Systems GmbH  
Alte Strasse 9 65558 Holzheim Allemagne

## Immunoglobulin G FS

Chemistry code 10 721

### Application for serum and plasma samples

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| Analysis Test Condition Setting (M) |        |        |
|-------------------------------------|--------|--------|
| Sample Type                         | Serum  | Urine  |
| Reac. sample vol.                   | 1.0    | 1.0    |
| Diluent method                      | No dil | No dil |
| Undil. sample vol.                  | 0      | 0      |
| Diluent volume                      | 0      | 0      |
| Diluent position                    | 0      | 0      |

| Prozone       |             |
|---------------|-------------|
| Prozone form  | No          |
| Prozone limit | 9.999       |
| Prozone judge | Upper limit |
| Judge limit   | 9.999       |
| M-DET.P.m     | 0           |
| M-DET.P.n     | 0           |
| S-DET.P.p     | 0           |
| S-DET.P.r     | 0           |

| MULTI-STD Setting |             |                 |             |                |              |              |       |        |
|-------------------|-------------|-----------------|-------------|----------------|--------------|--------------|-------|--------|
| Formula           | Logit Log 3 | Axis Conv       | No conv     |                |              |              |       |        |
| Blank             | Blank is 0  | Points          | 6           |                |              |              |       |        |
|                   | FV          | Reac. smp. vol. | Dil. method | Dil. smp. vol. | Diluent vol. | Diluent pos. | STD H | STD L  |
| BLK               | #           | 1.0             | No dil      | 0              | 0            | 0            | 9.999 | -9.999 |
| 1                 | #           | 1.0             | No dil      | 0              | 0            | 0            | 9.999 | -9.999 |
| 2                 | #           | 1.0             | No dil      | 0              | 0            | 0            | 9.999 | -9.999 |
| 3                 | #           | 1.0             | No dil      | 0              | 0            | 0            | 9.999 | -9.999 |
| 4                 | #           | 1.0             | No dil      | 0              | 0            | 0            | 9.999 | -9.999 |
| 5                 | #           | 1.0             | No dil      | 0              | 0            | 0            | 9.999 | -9.999 |

# entered by user