

Sudan Black

IVD In-vitro diagnostic medical device **CE**

CND Code: W01030708

Catalog number	Unit size
05-B27002	150 ml

Packaging Primary container: white bottle in High Density Polyethylene (HDPE). Useful capacity 150 ml. HDPE cap. Tamper evident cap.

Wear, water, alcohol and solvents resistant PVC label. Scratchproof ink resistant to water and alcohol.

Expected aim Product for the preparation of cyto-histological samples for optical microscopy.

Application Cryostatic sections of 10 micron. Staining solution for fat.

Method

- 1) Sections in ethanol 95°
- 2) Sudan Black
- 3) Differentiate rapidly in ethanol 95°
- 4) Wash in distilled water
- 5) Mount with aqueous mount

7-10 minutes

Results Fatsgreenish - black

Components

Components	CAS	CE	Index
Sudan Black B C.I. 26150	4197-25-5	2240871	-
Ethylene Glycerol	107-21-1	2034733	603-027-00-1
Ethanol 95°	64-17-5	2005786	603-002-00-5

Warning and precaution

The product must be used exclusively by specialized technical operators.
Carefully read the information on the classification of dangerous substances on the label. Always refer to the safety data sheet where are available the information on the risks presented by the mixture, the precautionary measures during use, the measures first aid and the intervention in the event of accidental release.
Do not use if the primary container is damaged.

Storage Store the preparation at room temperature. Keep the containers tightly closed.

Stability After the first opening, the product is usable until the expiry date, if correctly stored. Product validity: 2 years.

Disposal Hazardous preparation: observe all state and local environmental regulations regarding waste disposal.

Date of issue: April 2021