

## Scott's water

**IVD** In-vitro diagnostic medical device **CE**

**CND Code: W01030799**

Catalog number	Unit size
05-M05023	500 ml

**Packaging** Primary container: white bottle in High Density Polyethylene (HDPE). Useful capacity 500 ml. HDPE cap. Tamper evident cap.

Wear, water, alcohol and solvents resistant PVC label. Scratchproof ink resistant to water and alcohol.

**Expected aim** Product for the preparation of cyto-histological samples for optical microscopy.

**Application** Colour change solution for hematoxylin. Concentrated buffer 10x.

**Method** Hematoxylin-eosin

- |                                    |           |
|------------------------------------|-----------|
| 1) Sections to water               |           |
| 2) Hematoxylin according to Harris | 3 minutes |
| 3) Scott's water                   | 1 minute  |
| 4) Aqueous eosin solution 1%       | 5 minutes |
| 5) Wash in tap water               | 5 minutes |
| 6) Dehydrate rapidly               |           |
| 7) Clearing agent and mount        |           |

### Components

Components	CAS	CE	Index
Sodium bicarbonate	144-55-8	-	-
Magnesium sulfate	7487-88-9	-	-
Deionized water			

### Warning and precaution

The product must be used exclusively by specialized technical operators.  
Carefully read the information on the classification of dangerous substances on the label. Always refer to the safety data sheet where are available the information on the risks presented by the mixture, the precautionary measures during use, the measures first aid and the intervention in the event of accidental release.  
Do not use if the primary container is damaged.

**Storage** Store the preparation at 2-8°C. Keep the containers tightly closed.

**Stability** After the first opening, the product is usable until the expiry date, if correctly stored. Product validity: 2 years.

**Disposal** Hazardous preparation: observe all state and local environmental regulations regarding waste disposal.

Date of issue: January 2017