

# RBC Lysis Buffer (10X)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                     |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog# / Size</b>   | 420301 / 100 mL<br>420302 / 500 mL                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Regulatory Status</b> | RUO                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Other Names</b>       | Red blood cell lysis buffer                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Description</b>       | Red Blood Cell (RBC) Lysis Buffer has been designed, formulated, and tested to ensure optimal lysis of RBCs in single cell suspensions with minimal effects on leukocytes. RBC Lysis Buffer is supplied as a 10X solution containing ammonium chloride, potassium carbonate, and EDTA, and should be diluted in deionized water prior to use. Nucleated RBCs are not effectively lysed with ammonium chloride. Test |

## Product Details

|                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
|-------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Formulation</b>            | This red blood cell (RBC) lysis buffer is supplied as a 10X solution and should be diluted to 1X in deionized water. 100 ml of 10X concentrate will yield a quantity of 1X solution that is sufficient to lyse 500 samples.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Preparation</b>            | The pH of the 1X solution should fall within the range of pH 7.1-7.4. Adjust the pH if necessary. Warm the 1X solution to room temperature prior to use.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Storage &amp; Handling</b> | Store RBC lysis buffer between 2°C and 8°C.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Application Notes</b>      | <p><b>A. Lysis of Mouse Spleen RBCs:</b></p> <ol style="list-style-type: none"> <li>1. Harvest mouse spleen and prepare a single cell suspension.</li> <li>2. Pellet the cells by centrifugation (350 x g); aspirate the supernatant.</li> <li>3. Dilute the 10X Red Blood Cell Lysis Buffer to 1X working concentration with deionized water and resuspend the pellet in 5 ml of 1X Lysis Buffer.</li> <li>4. Incubate on ice for 4-5 minutes with occasional shaking.</li> <li>5. Stop the reaction by diluting the Lysis Buffer with 20-30 ml of 1X PBS.</li> <li>6. Spin the cells (350 x g) and discard the supernatant.</li> <li>7. Resuspend the pellet in the appropriate buffer (e.g., BioLegend Cell Staining Buffer Cat. No. 420201), wash 1X.</li> <li>8. Count cells, adjust density, and proceed with cell staining procedures.</li> </ol> <p><b>B. Lysis of Human Peripheral Blood RBCs:</b></p> <ol style="list-style-type: none"> <li>1. Dilute the 10X RBC Lysis Buffer to 1X working concentration with deionized water. Warm the 1X solution to room temperature prior to use.</li> <li>2. Add 2.0 ml of 1X RBC Lysis Buffer to each tube containing up to 100 µl of whole blood.</li> <li>3. Gently vortex each tube immediately after adding the lysing solution. Incubate at room temperature, protected from light, for 10-15 minutes.</li> <li>4. Centrifuge 350 x g for 5 minutes. Aspirate supernatant without disturbing pellet.</li> <li>5. Resuspend the pellet in the appropriate buffer (e.g., BioLegend Cell Staining Buffer, Cat. No. 420201), wash 1X.</li> <li>6. Resuspend and proceed with further procedures.</li> </ol> |

## Application References

(PubMed link indicates BioLegend citation)

1. Adams RA, *et al.* 2007. *J. Exp. Med.* 204:571.
2. Feild-Corbett C, *et al.* 2009. *Stem Cells*. [PubMed](#)
3. de Resende MM, *et al.* 2010. *Physiol Genomics*. 42:437. [PubMed](#)
4. Roybal JD, *et al.* 2011. *Mol Cancer Res*. 9:25. [PubMed](#)
5. Boldin MP, *et al.* 2011. *J Exp Med*. [PubMed](#)
6. Kahlenberg JM, *et al.* 2011. *J. Immunol.* 187:6143. [PubMed](#)
7. Deiluiis JA, *et al.* 2012. *Am J Physiol Lung Cell Mol Physiol*. 302:399. [PubMed](#)
8. Wu CY, *et al.* 2012. *PLoS One*. 7:e42365. [PubMed](#)
9. Hsia BJ, *et al.* 2012. *J Allergy Clin Immunol*. 130:205. [PubMed](#)
10. Rico C, *et al.* 2012. *Carcinogenesis*. 33:2283. [PubMed](#)
11. Jin J, *et al.* 2012. *Arterioscler Thromb Vasc Biol*. 32:2901. [PubMed](#)
12. Saraiva AM, 2013. *Hum Immunol*. 74:207. [PubMed](#)
13. Koon HW, *et al.* 2013. *J Vis exp*. 68:4208. [PubMed](#)
14. Ma Dy, *et al.* 2013. *J Virol*. 87:3361. [PubMed](#)
15. Han W, 2013. *J. Immunol*. 190:4786. [PubMed](#)
16. Garcia-Acros I, *et al.* 2013. *J Biol Chem*. 288:14046. [PubMed](#)
17. Zhao JL, *et al.* 2013. *Elife*. 21:537. [PubMed](#)
18. Guo W, *et al.* 2013. *J Gerontol A Biol Sci Med Sci*. [PubMed](#)
19. Yamasaki S, *et al.* 2013. *Clin Vaccine immunol*. 20:1508. [PubMed](#)
20. Kraft P, *et al.* 2013. *Stroke*. 44:3202. [PubMed](#)
21. Liu X, *et al.* 2013. *Invest Ophthalmol Vis Sci*. 54:7386. [PubMed](#)
22. Fullerton AM, *et al.* 2013. *Toxicol Sci*. 136:72. [PubMed](#)
23. Cekic C, *et al.* 2013. *J Exp Med*. 210:2693. [PubMed](#)
24. Varga T, *et al.* 2013. *J. Immunol*. 191:5695. [PubMed](#)
25. Ma Y, *et al.* 2014. *Cancer Res*. 74:436. [PubMed](#)
26. Chang YJ, *et al.* 2014. *Cytotherapy*. 16:560. [PubMed](#)
27. Kim MJ, *et al.* 2014. *Mol Pharmacol*. 86:284. [PubMed](#)
28. Cekic C, *et al.* 2014. *Cancer Res*. 74:7250. [PubMed](#)
29. Szulzewsky F, *et al.* 2015. *PLoS One*. 10:116644. [PubMed](#)

## Product Citations

1. Dobie R, *et al.* 2020. Cell Reports. 29(7):1832-1847.e8.. [PubMed](#)
2. Strickley JD, *et al.* 2019. Nature. 575:519. [PubMed](#)
3. Shah D, *et al.* 2021. Oncoimmunology. 10:1939601. [PubMed](#)
4. Hutton C, *et al.* 2021. Cancer Cell. 39:1227. [PubMed](#)
5. Sandovici I, *et al.* 2022. Dev Cell. 57:63. [PubMed](#)
6. Deiluiis J, *et al.* 2012. Am J Physiol Lung Cell Mol Physiol. 302:L399. [PubMed](#)
7. Jin J, *et al.* 2012. Arterioscler Thromb Vasc Biol. 32:2901. . [PubMed](#)
8. Barros-Martins J, *et al.* 2016. J Immunol. 196: 3642 - 3652. [PubMed](#)
9. Zhou B, *et al.* 2020. Nat Immunol. 1.868055556. [PubMed](#)
10. Tan J, *et al.* 2021. iScience. 24(8):102835. [PubMed](#)
11. Asem M, *et al.* 2020. Cancer Res. 1156:80. [PubMed](#)
12. Lin C, *et al.* 2020. Cancer Immunol Res. 632:8. [PubMed](#)

## Antigen Details

---

|           |       |
|-----------|-------|
| Cell Type | Tregs |
| Gene ID   | NA    |

## Related Protocols

---

[Cell Surface Flow Cytometry Staining Protocol](#)

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

\*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, [www.biolegend.com/ordering#license](http://www.biolegend.com/ordering#license)). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

8999 BioLegend Way, San Diego, CA 92121 [www.biolegend.com](http://www.biolegend.com)  
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587