
Microstream™ Advance

Neonatal-Infant Intubated CO₂ Filter Line

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




EN Instructions for Use

Intended Use

The Microstream™ Advance Neonatal-Infant Intubated CO₂ Filter Line is intended to conduct a sample of the patient's breathing from a ventilator or anesthesia machine to a gas measurement device for measuring the percentage of CO₂ in the patient's exhaled breath. The set is intended for single patient use only. Intended population: Intubated Neonatal-Infant patients.

Intended Users

The Microstream™ Advance Neonatal-Infant Intubated CO₂ Filter Line is intended to be used by licensed healthcare providers only.

Indications for Use

The Microstream™ Advance Neonatal-Infant Intubated CO₂ Filter Line is indicated for measurement of exhaled carbon dioxide (CO₂) for monitoring purposes. These devices are limited by the indications for use of the connected Philips monitoring equipment in healthcare facilities. These devices are intended to interact with patient intact skin.

Contraindications

There are no known contraindications.

WARNINGS

- Loose or damaged connections may compromise ventilation or cause an inaccurate measurement of respiratory gases. Securely connect all components, turning the sampling line connector clockwise into the monitor CO₂ port until it can no longer be turned. Check connections for leaks according to standard clinical procedures, and ensure CO₂ values appear.
- Do not attempt to clean, disinfect, sterilize or flush any part of the sampling line. Re-use of single-use accessories could pose a cross-contamination risk to the patient or damage the functioning of the monitor.
- When used with a suction catheter, do not place the airway adapter between the suction catheter and endotracheal tube. This is to ensure that the airway adapter does not interfere with the functioning of the suction catheter, and that the suction process does not damage the airway adapter. Damage of the airway adapter may harm the patient.
- Ensure no kinks are present in CO₂ tubing, as kinked tubing may cause inaccurate CO₂ sampling.
- Do not use 989803204531, 989803204311, 989803204341, or 989803204661 with products with protruding inner connector pieces (i.e., the Hamilton flow sensor), as this may cause breakage of the airway adapter.

Caution

Federal (US) law restricts this device to sale by or on the order of a medical practitioner.

Required Compatible Components

These devices are to be used with Philips monitors using Microstream™ technology.

Airway adapters for neonatal and infant patients:

- Added dead space <0.5 cc.
- To be used with an endotracheal tube ≤ 4.5mm bore.

Applying Neonatal/Infant Intubated CO₂ Filter Line

When performing the following procedure, reference the step-by-step patient application graphics located on Page i.

- ① Firmly connect the small-end (male) of the Microstream™ airway adapter to the female-end of the Wye piece (ventilation source).
- ② Firmly connect the patient's endotracheal tube connector into the large-end (female) of the Microstream™ airway adapter.
- ③ Attach CO₂ connector to monitor as shown in this graphic.
- ④ **For the 989803204311**, fit one dryer inside the incubator and the other dryer outside the incubator.
- ⑤ To remove the airway adapter from the breathing circuit, twist the airway adapter outwards by a quarter turn while firmly gripping the closest section of the ventilation source or patient's endotracheal tubing, and disconnect.

Operational Notes

- During setup, ensure that airway adapter can be easily attached and detached from the tubing before proceeding.
- During nebulization or suction, in order to avoid moisture buildup and sampling line occlusion, remove the sampling line connector from the CO₂ port on the monitor.
- The response time of the 2m sampling line is ≤3.5 seconds, 4m sampling lines ≤5.9 seconds.
- Sampling lines with Extended Duration or High Humidity in their names include a moisture reduction component (Nafion®* or its equivalent) for use in higher humidity environments where long duration use of CO₂ sampling is required.
- Ensure CO₂ values appear on the screen following connection of the sampling line.

TM* Trademark of its respective owner.

Reordering Information

Part Number	Description	Quantity
989803204531	Neonatal-Infant Intubated CO2 Filter Line, Extended Duration Note: Length: 2m	25
989803204311	Neonatal-Infant Intubated CO2 Filter Line, High Humidity Note: Length: 2m	25
989803204341	Neonatal-Infant Intubated CO2 Filter Line, Extended Duration Note: Length: 4m	25
989803204661	Neonatal-Infant Intubated CO2 Filter Line, Extended Duration Note: Length: 2m	25

Disposal

Discard any product that becomes damaged or deteriorated. Follow approved medical waste disposal methods specified by your facility or local regulations.

Incident Reporting

Any serious incident that has occurred in relation to this device should be reported to Philips and the competent authority of the European Economic Area (EEA) countries, including Switzerland and Turkey, in which the user and/or patient is established.

Environmental Specifications

The following table provides the temperature and humidity ranges that must be maintained for these products during storage and transport:

Specification	Storage	Transport
Temperature	-20°C to 70°C (-4°F to 158°F)	-20°C to 70°C (-4°F to 158°F)
Humidity	10% RH to 90% RH	10% RH to 90% RH

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Copyright © 2020
Microstream™ is a trademark of Oridion
Nederland B.V.
Medical 1987 Ltd. Oridion Medical 1987 Ltd.