

TrueMark[™] MSI Assay

USER GUIDE

For Microsatellite Instability Testing

for use with:

3500/3500xL Genetic Analyzer

SeqStudio[™] Genetic Analyzer

TrueMark[™] MSI Analysis Software

Catalog Numbers A45295

Publication Number MAN0018868

Revision A.0



Life Technologies Ltd | 7 Kingsland Grange | Woolston, Warrington WA1 4SR | United Kingdom

For descriptions of symbols on product labels or product documents, go to thermofisher.com/symbols-definition.

The information in this guide is subject to change without notice.

DISCLAIMER: TO THE EXTENT ALLOWED BY LAW, THERMO FISHER SCIENTIFIC INC. AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

Revision history: Pub. No. MAN0018868

Revision	Date	Description
A.0	13 December 2019	New document for the TrueMark™ MSI Assay.

Important Licensing Information: These products may be covered by one or more Limited Use Label Licenses. By use of these products, you accept the terms and conditions of all applicable Limited Use Label Licenses.

TRADEMARKS: All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

©2019 Thermo Fisher Scientific Inc. All rights reserved.

Contents

■	CHAPTER 1	Product information	5
	Product description	5	
	Kit overview	5	
	About the primers	5	
	Dyes used in the kit	5	
	Markers amplified by the kit	6	
	Standards and controls that are required	6	
	TrueMark™ MSI Assay Amplification Control profile	7	
	Contents and storage	7	
	Required materials not supplied	8	
■	CHAPTER 2	Prepare for capillary electrophoresis	10
	3500/3500xL Genetic Analyzer	10	
	Perform spectral calibration	10	
	Electrophoresis software setup	10	
	Create an instrument protocol	11	
	Create the size standard	11	
	Create the sizecalling protocol	12	
	Create the Assay	15	
	SeqStudio™ Genetic Analyzer	16	
	Perform a system dye calibration	16	
	Create the TrueMark™ MSI Assay run module	18	
	Create the TrueMark™ MSI Assay Size Standards	19	
	Create a new plate setup	19	
■	CHAPTER 3	Prepare and run the samples	22
	Sample naming requirements for TrueMark™ MSI Analysis Software	22	
	DNA sample preparation	24	
	Set up the PCR reactions	25	
	Run the PCR	26	
	Prepare samples for electrophoresis	26	

■	CHAPTER 4 Analyze the data with the TrueMark™ MSI Analysis Software	28
	Software access restrictions	28
	Sign in to the TrueMark™ MSI Analysis Software (all users)	29
	Import sample files	29
	Check for import errors	31
	Perform an initial QC of the imported data	31
	View and interpret the results	32
	View the electropherogram traces	33
	Representative data	34
	(Optional) Accept and approve a specimen	40
	Generate a report	40
	Export results	41
■	APPENDIX A Troubleshooting	43
	Extra peaks in the electropherogram	44
	Causes of extra peaks	44
	Stutter definition	44
	Artifact definition	45
	Artifact observation	45
■	APPENDIX B Analyze the data with GeneMapper™ Software	47
	Overview of the GeneMapper™ Software	47
	Set up the GeneMapper™ Software (v5.0 or later) for analysis (before first use of the kit)	47
	Define custom display and plot settings	47
	Define the Size Standard	48
	Create an analysis method	49
	Create an analysis method	49
	Analysis method settings	50
	Analyze and edit sample files with GeneMapper™ Software	55
■	APPENDIX C Safety	56
	Chemical safety	57
	Biological hazard safety	58
■	APPENDIX D Documentation and support	59
	Related documentation	59
	Customer and technical support	59
	Limited product warranty	60



Product information

■ Product description	5
■ Contents and storage	7
■ Required materials not supplied	8

IMPORTANT! Before using this product, read and understand the information in the “Safety” appendix in this document.

Product description

Kit overview The Applied Biosystems™ TrueMark™ MSI Assay detects the presence of microsatellite instability in DNA samples through multiplex PCR and fragment analysis. Fragment analysis is performed on the Applied Biosystems™ 3500/3500xL Genetic Analyzer or the Applied Biosystems™ SeqStudio™ Genetic Analyzer. Data is analyzed with the Applied Biosystems™ TrueMark™ MSI Analysis Software for easy reporting of results.

About the primers The TrueMark™ MSI Assay primers are manufactured using the same synthesis and purification improvements as the primers in our forensic DNA profiling kits. These improvements enhance the assay signal-to-noise ratio and simplify the interpretation of results.

Dyes used in the kit

Dye	Color	Label
6-FAM™	Blue	Samples and controls
VIC™	Green	
NED™	Yellow	
TAZ™	Red	
SID™	Purple	
LIZ™	Orange	GeneScan™ 600 LIZ™ Size Standard v2.0

Markers amplified by the kit

Table 1 TrueMark™ MSI Assay markers

Marker designation	Chromosomal location	Type	Dye label
BAT-25	4q12	MSI	6-FAM™
NR-24	2q11.1		
NR-21	14q11.2		
BAT-40	1p12		VIC™
CAT-25	7q34		
NR-22	11q24.2		NED™
NR-27	11q22.2		
ABI-19	1q42.3		
ABI-20B	1q21.3		SID™
ABI-17	17p12		
ABI-16	17p13.2		
BAT-26	2p21		
ABI-20A	12q24.13		
TH01	11p15.5	Human Identification (HID)	VIC™
PentaD	21q22.3		TAZ™

Standards and controls that are required

For the TrueMark™ MSI Assay, the panel of standards needed for PCR amplification and PCR product sizing are:

- **TrueMark™ MSI Assay Amplification Control**—A negative MSI and PCR amplification control for evaluating the efficiency of the amplification step. TrueMark™ MSI Assay Amplification Control is included in the kit. See “TrueMark™ MSI Assay Amplification Control profile” on page 7.
- **GeneScan™ 600 LIZ™ Size Standard v2.0**—Used for obtaining sizing results. This standard, which has been evaluated as an internal size standard, yields precise sizing results for PCR products. Order the GeneScan™ 600 LIZ™ Size Standard v2.0 (Cat. No. 4408399) separately.

TrueMark™ MSI Assay Amplification Control profile




Figure 1 TrueMark™ MSI Assay Amplification Control (2 ng) amplified at 29 PCR cycles with the TrueMark™ MSI Assay and analyzed on an Applied Biosystems™ SeqStudio™ Genetic Analyzer (Y-axis scale 0 to 25,000)

Contents and storage

Table 2 TrueMark™ MSI Assay [Cat. No. A45295]

Contents	Amount	Storage
TrueMark™ MSI Assay Master Mix	480 µL	-25°C to -15°C 2-8°C after first use. Store for ≤6 months or the expiration date of the kit (whichever comes first)
TrueMark™ MSI Assay Primer Mix	130 µL	
TrueMark™ MSI Assay Amplification Control	120 µL	
TrueMark™ MSI Assay No-Template Control	520 µL	

Required materials not supplied

Unless otherwise indicated, all materials are available through **thermofisher.com**.
MLS: Fisher Scientific (**fisherscientific.com**) or other major laboratory supplier.

Table 3 Materials for generating PCR products

Item	Source
Instruments and equipment	
One of the following: ^[1]	
• ProFlex™ 96-well PCR System	4484075
• ProFlex™ 2 × 96-well PCR System	4484076
• ProFlex™ 3 × 32-Well PCR System	4484073
Benchtop microcentrifuge	MLS
Vortex mixer	MLS
Adjustable micropipettors	MLS
Plates and other consumables	
RecoverAll™ Total Nucleic Acid Isolation Kit for FFPE or equivalent	AM1975
TE Buffer	12090015
MicroAmp™ 96-Well Tray	N8010541
MicroAmp™ 96-Well Tray/Retainer Set	403081
MicroAmp™ Optical 96-Well Reaction Plate with Barcode	4306737
MicroAmp™ Clear Adhesive Film, or equivalent	4306311
Aerosol-resistant pipette tips	MLS
Other plastic consumables	thermofisher.com/plastics

^[1] You can use an equivalent thermal cycler. Optimize the protocols for other thermal cyclers.

Table 4 Materials for capillary electrophoresis

Item	Source
Instruments and equipment	
One of the following genetic analyzers:	
• 3500/3500xL Genetic Analyzer with:	
– 3500 Data Collection v3 Software or later (Windows™ 7 or 10 operating system)	
– 50 cm POP-7™ Polymer array	
• SeqStudio™ Genetic Analyzer with:	
– (Recommended) SeqStudio™ Plate Manager	
	Contact your local sales office.

Item	Source
Benchtop microcentrifuge	MLS
Vortex mixer	MLS
Adjustable micropipettors	MLS
General reagents	
GeneScan™ 600 LIZ™ Size Standard v2.0	4408399
Hi-Di™ Formamide	4311320
DS-36 (Dye Set J6) Matrix Standard Kit	4425042
Reagents for the 3500/3500xL Genetic Analyzer	
Anode Buffer Container	4393927
Cathode Buffer Container	4408256
Polymer, POP-7™ (96)	A26073
Conditioning reagent	4393718
Capillary array, 8-Capillary, 50-cm	4404685
Capillary array, 24-Capillary, 50-cm	4404689
Retainer & Base Set (Standard) for 3500/3500xL Genetic Analyzers, 96 well	4410228
Septa for 3500/3500xL Genetic Analyzers, 96 well	4412614
Reagents for the SeqStudio™ Genetic Analyzer	
Septa for SeqStudio™ Genetic Analyzer, 96 well	A36541
SeqStudio™ Cartridge or SeqStudio™ Cartridge v2	A33671 or A41331
Cathode Buffer Container	A33401
Reservoir Septa	A35640
Plates and other consumables	
MicroAmp™ Optical 96-Well Reaction Plate with Barcode	4306737
MicroAmp™ Clear Adhesive Film, or equivalent	4306311

Table 5 Materials to analyze data

Item	Source
TrueMark™ MSI Analysis Software	Download the software from downloads.thermofisher.com/TrueMark MSI Analysis Software.zip
(Optional) GeneMapper™ Software (v5.0 or later)	Contact your local sales office.

Prepare for capillary electrophoresis


■ 3500/3500xL Genetic Analyzer	10
■ SeqStudio™ Genetic Analyzer	16

3500/3500xL Genetic Analyzer

Perform spectral calibration

Perform a spectral calibration using the DS-36 (Dye Set J6) Matrix Standard Kit (J6 Dye Set) (Cat. No. 4425042).

The following figure is an example of a passing 6-dye spectral calibration.



Electrophoresis software setup

(Recommended) Ensure that your genetic analyzer is running on the Windows™ operating system with 3500 Data Collection Software v3 or later.

The following instructions cover setting up an assay (ABI_MSI_Assay), instrument protocol (ABI_MSI), and sizecalling protocol (GS600[80-400]). For more detailed information on the instrument procedures, see the documents listed in Appendix D, “Documentation and support”.

Create an instrument protocol

1. Navigate to the **Instrument Protocols** library.
2. Click **New**.
3. Specify the settings listed below.

- **Capillary length**—50cm
- **Polymer**—POP7
- **Dye Set**—J6
- **Run Module**—
FragmentAnalysis50_POP7xl
- **Protocol Name**—ABI_MSI
- **Oven Temperature (°C)**—60
- **Run Time (sec)**—1,000
- **Run Voltage (kVolts)**—19.5
- **PreRun Time (sec)**—180
- **PreRun Voltage (kVolts)**—15
- **Injection Time (sec)**—15
- **Injection Voltage (kVolts)**—1.6
- **Data Delay (sec)**—1

4. Click **Save**.

Create the size standard

A size standard defines the sizes of known fragments. It is used to generate a standard curve. The standard curve is used to determine the sizing of unknown samples. When you create a sizcalling (fragment analysis) protocol, you add a size standard to the protocol.

1. Navigate to the **Size Standards** library.
2. Click **New**.

3. Specify the settings listed below.

Create New Size Standard

Setup a Size Standard

* Size Standard: GS600LIZ(80-400) ☒ Locked

Description: J6

* Dye Color: Orange

Enter sizes in the field below separated by a comma, space, or return then click the 'Add Size(s) >>' button to add them to the current size standard definition.

Enter new Size Standard definition: (e.g. 11.0, 34.2, 55)

* Current Size Standard definition: Delete Selected Sizes

80.0
100.0
114.0
120.0
140.0
160.0
180.0
200.0
214.0
220.0
240.0
250.0
260.0
280.0
300.0
314.0
320.0
340.0
360.0
380.0
400.0

Close Save

- **Size Standard name**—GS600LIZ (80-400); Select **Locked**.
- **Description**—J6
- **Dye Color**—Orange
- **Size Standards** (list in left text box)—80.0, 100.0, 114.0, 120.0, 140.0, 160.0, 180.0, 200.0, 214.0, 220.0, 240.0, 250.0, 260.0, 280.0, 300.0, 314.0, 320.0, 340.0, 360.0, 380.0, and 400.0

4. Click **Save**.

Create the
sizecalling
protocol

1. Navigate to the **Sizecalling Protocols** library.
2. Click **New**.

3. Specify the settings listed below.

Create New Sizecalling Protocol

Setup a Sizecalling Protocol

* Protocol Name: GS600(80-400) Locked

Description:

Size Standard: GS600_LIZ_(80-400)

Sizecaller: SizeCaller v1.1.0

Analysis Settings QC Settings

Analysis Range: Full Sizing Range: Full Size Calling Method: Local Southern

Analysis Start Point: 0 Sizing Start Size: 0 Primer Peak: Present

Analysis Stop Point: 1000000 Sizing Stop Size: 100000

	Blue	Green	Yellow	Red	Purple	Orange
Minimum Peak Height	200	200	200	200	200	200

Common Settings

Use Smoothing: None

Use Baselining (Baseline Window (Pts)): ☒ 51

Minimum Peak Half Width: 2

Peak Window Size: 15

Polynomial Degree: 3

Slope Threshold Peak Start: 0.0

Slope Threshold Peak End: 0.0

Close Save

- **Sizecaller**—SizeCaller v1.1.0
- **Analysis Range and Sizing Range**—Full
- **Size Calling Method**—Local Southern
- **Primer Peak**—Present
- **Minimum Peak Height**—200
- **Use Smoothing**—None
- **Use Baselining (Baseline Window) (Pts)**—51
- **Minimum Peak Half Width**—2
- **Peak Window Size**—15
- **Polynomial Degree**—3
- **Slope Threshold Peak Start/End**—0.0 (both)

4. Click the **QC Settings** tab, and specify the settings listed below.

Create New Sizecalling Protocol

Setup a Sizecalling Protocol

* Protocol Name: GS600(80-400) ☐ Locked

Description:

Size Standard: GS600_LIZ_(80-400)

Sizecaller: SizeCaller v1.1.0

Analysis Settings | **QC Settings**

Size Quality

Fail if Value is	Suspect Range	Pass if Value is
< 0.5	0.5 - 0.75	≥ 0.75

Assume Linearity from (bp): 0 To (bp): 800

Pull Up

Actuate Pull-Up flag if Pull-Up Ratio ≤ 0.1 and Pull-Up Scan ≤ 1


Close Save


- **Fail if Value is** <0.5
- **Suspect Range**—0.5 – 0.75
- **Pass if Value is** ≥0.75
- **Assume Linearity from**—0 to 800bp
- **Actuate Pull-Up flag if Pull-Up Ratio** is ≤0.1 and **Pull-Up Scan** is ≤1

5. Click **Save**.

IMPORTANT! Normalization is not applied to samples with Size Quality flags. Specify analysis settings that accurately detect and size the size standard, and QC settings with appropriate pass fail ranges. The 3500 Series Data Collection Software does not support re-analyzing data with new settings.

Create the Assay

1. Navigate to the **Assays** library.
The list of items in the library may be filtered based on the library filtering user preference. Click **Disable Filters** to show all items in the list.
2. Click  **New**.
3. Specify the settings.



- ① **Assay Name**—ABI_MSI_Assay
- ② **Application Type**—Fragment
 - **Assay Name**—ABI_MSI_Assay
 - **Color**—Black
 - **Application Type**—Fragment
 - **Do you wish to assign multiple protocols to this assay?**—No
 - **Instrument Protocol**—ABI_MSI (select from the dropdown list)
 - **Sizecalling Protocol**—GS600(80-400) (select from the dropdown list)

4. Click **Save**.


SeqStudio™ Genetic Analyzer

Perform a system dye calibration


A system dye calibration requires ~30 minutes to complete.

Prepare the dye set calibration standards and plate as described in the *DS-36 (Dye Set J6) Matrix Standard Kit Product Information Sheet* (Pub. No. 4426052).

1. In the home screen, tap **Settings** ▶ **Maintenance and Service** ▶ **Calibration** ▶ **Dye Calibration**.
2. Touch the injection group for the dye set in the plate, then tap **Dye set**.



3. Tap **Matrix Standard**, then select a system dye calibration standard provided with the instrument.



4. Tap **Calibrate**.
The calibration run starts.

IMPORTANT! If the dye calibration fails:

- The results of the calibration are not saved, and the calibration plate is not moved to **Run History**.
 - The instrument does not allow you to rerun the plate setup for a failed calibration. Close the calibration screen, then start a new calibration.
-

Spectral Quality Value

A spectral Quality Value reflects the confidence that the individual dye emission signals can be separated from the overall measured fluorescence signal. It is a measure of the consistency between the final matrix and the data from which it was computed. A Quality Value of 1.0 indicates high consistency, providing an ideal matrix with no detected pull-up/pull-down peaks.

In rare cases, a high Quality Value can be computed for a poor matrix. This can happen if the matrix standard contains artifacts, leading to the creation of one or more extra peaks. The extra peaks cause the true dye peak to be missed by the algorithm, and can lead to a higher Quality Value than would be computed with the correct peak. Therefore, it is important to visually inspect the spectral calibration profile for each capillary.

Condition number

A Condition Number indicates the amount of overlap between the dye peaks in the fluorescence emission spectra of the dyes in the dye set.

If there is no overlap in a dye set, the Condition Number is 1.0 (ideal conditions), the lowest possible value. The condition number increases with increasing peak overlap.

The ranges that the software uses to determine if a capillary passes or fails are:

Dye Set	Quality Value Minimum	Condition Number Maximum
J6	0.95	8.0

Create the TrueMark™ MSI Assay run module

1. In the SeqStudio™ Plate Manager, open a new or existing plate, then navigate to the **Plate** tab.
2. Click **Actions** ► **Manage run modules** from the top right corner of the screen.
3. Select the default **FragAnalysis** run module, then click **Copy**.
4. Enter the **Run module name**—**ABI_MSI_Assay**
5. Change the **Separation (Run) Time** to **1,000** seconds.

Copy, delete or edit run modules

Run module name ①

Injection Time Between 1 to 600 seconds	<input type="text" value="7"/>	Separation (Run) Time Between 300 to 14000 seconds	<input type="text" value="1000"/> ②
Injection Voltage Between 0 to 13000 volts	<input type="text" value="1200"/>	Separation (Run) Voltage Between 0 to 13000 volts	<input type="text" value="9000"/>

[Advanced options](#) ▾

Capillary heater temperature setting (degrees Celsius)	<input type="text" value="60"/>	Between 40 to 60 degrees Celsius
Pre-Run Voltage (volts)	<input type="text" value="13000"/>	Between 0 to 13000 volts
Length Of Pre-Run (seconds)	<input type="text" value="180"/>	Between 0 to 1000 seconds
Separation (Run) Voltage Ramp Duration (seconds)	<input type="text" value="300"/>	Between 0 to 1800 seconds

Cancel Save

- ① Run module name
② Separation (Run)Time

6. Click **Save**.

Create the TrueMark™ MSI Assay Size Standards

1. In the SeqStudio™ Plate Manager, open a new or existing plate, then navigate to the **Plate** tab.
2. Click **Actions** ► **Manage size standards** from the top right corner of the screen.
3. Select the 80-400 sizes, then save with the name **GS600_LIZ_(80-400)**.

Create a new plate setup

1. Open the SeqStudio™ Plate Manager and select **New**.
2. Enter or accept the default **Plate name**, then click **Next**.
3. Enter the plate properties.
 - a. Select **Fragment analysis** as the **Application type**.
 - b. Click **Setting Details**, then click **Copy** to make a copy of the default settings to edit.
 - c. In the Analysis Settings window, enter a unique **Name** for the MSI analysis settings.
 - d. In the Analysis Settings window, change all minimum peak heights to **200**.

- e. In the Analysis Settings window, change the Peak window size to 5.

The screenshot shows the 'Analysis settings' window. At the top is a title bar labeled 'Analysis settings'. Below it is a tab labeled 'Fragment analysis'. The main area contains several settings sections:

- Name:** A text field containing 'ABI_MSI_SeqStudio_Analysis_Settings'. A callout line labeled '1' points to this field.
- Analysis range:** A dropdown menu set to 'Full', followed by 'from point' 1 and 'to point' 10000.
- Sizing range:** A dropdown menu set to 'Full', followed by 'from size' 1 and 'to size' 12000.
- Size calling method:** A dropdown menu set to 'Local Southern'.
- Primer peak:** A dropdown menu set to 'Present'.
- Minimum peak height:** A section with six color-coded boxes (Blue, Green, Yellow, Red, Purple, Orange). Each box has a checked checkbox and a value of 200. A callout line labeled '2' points to the Orange box.
- Common settings:** A section with several controls:
 - Use smoothing:** A dropdown menu set to 'None'.
 - Use baselining:** A checked checkbox followed by a value of 51.
 - Minimum peak half width:** A text field with the value 2.
 - Polynomial degree:** A text field with the value 3.
 - Slope threshold peak start:** A text field with the value 0.
 - Peak window size:** A text field with the value 5. A callout line labeled '3' points to this field.
 - Slope threshold peak end:** A text field with the value 0.

At the bottom right of the window are two buttons: 'Cancel' and 'Save'. A callout line labeled '4' points to the 'Save' button.

- ① **Name** for the Analysis settings.
- ② **Minimum peak height**
- ③ **Peak window size**
- ④ **Save** button

4. Click **Save** to save the **Analysis Settings**, then **Close**.

5. Ensure that the new **Fragment analysis settings** are selected, then click **Next** to proceed to the Plate tab.

① MSI fragment analysis settings

6. In the Plate tab, select:
 - **Size standard—GS600_LIZ_(80-400)**
 - **Dye set—J6 (DS-36)™**
 - **Run module 1—ABI_MSI_Assay**

Note: For information on creating a run module, see “Create the TrueMark™ MSI Assay run module” on page 18.
7. Click **Next**.
8. In the **Save the plate setup** window:
 - a. Enter or accept the default **Plate name**.
 - b. Ensure that the **ABI_MSI_SeqStudio_Analysis_Settings** is selected as the **Fragment analysis settings**.
 - c. Click **Save**.

Prepare and run the samples

■ Sample naming requirements for TrueMark™ MSI Analysis Software	22
■ DNA sample preparation	24
■ Set up the PCR reactions	25
■ Run the PCR	26
■ Prepare samples for electrophoresis	26

Sample naming requirements for TrueMark™ MSI Analysis Software

To be successfully imported into the TrueMark™ MSI Analysis Software, the sample file (FSA) names must follow the correct naming conventions.

In the following examples, "SpecimenID" becomes the main name for identifying the specimen within the software and exports. SpecimenID text cannot contain an underscore (_), because only the text before the first underscore is imported as the specimen ID.

Paired samples convention:

- SpecimenID_T_*.fsa (tumor tissue sample)
- SpecimenID_N_*.fsa (normal adjacent tissue sample from same individual)

Sample file (FSA) names must meet the following conventions.

Sample type	File naming conventions	Guidelines
Specimen	<p>For paired samples, the naming convention is:</p> <ul style="list-style-type: none"> • <SpecimenID>_T_<*>.fsa—Tumor tissue sample • <SpecimenID>_N_<*>.fsa—Normal tissue sample from the same individual, adjacent to the Tumor tissue sample <p>where:</p> <ul style="list-style-type: none"> • <SpecimenID> is user-defined, but is identical in the Normal (N) and Tumor (T) tissue samples • <*> is user-defined 	<ul style="list-style-type: none"> • The <SpecimenID> prefix identifies the specimen within the TrueMark™ MSI Analysis Software and exported file names. Ensure that the <SpecimenID> text does not include an underscore (_). • Within a batch, each <SpecimenID> must be unique, unless two files will be analyzed as a Tumor-Normal sample pair. If <SpecimenID> duplicates are detected, the software imports only the last file, in alphanumeric sort order. For example, if the files are named "SpecimenIDBlue_aaa.fsa" and "SpecimenIDBlue_zzz.fsa", the software imports only "SpecimenIDBlue_zzz.fsa". • IMPORTANT! The TrueMark™ MSI Analysis Software will not import a <SpecimenID>_N_<*>.fsa file if there is no matching <SpecimenID>_T_<*>.fsa file to import. However, a <SpecimenID>_T_



Sample type	File naming conventions	Guidelines
Specimen	For unpaired samples, the naming convention is: <SpecimenID>_<*>.fsa where: <SpecimenID> and <*> are user-defined	<*>.fsa file will be imported even if there is no <SpecimenID>_N_<*>.fsa file to import. <ul style="list-style-type: none">IMPORTANT! If you are running replicate reactions on the same plate, assign the replicates a unique <SpecimenID> before the _T or _N to ensure that the replicates are processed correctly within the software. For example:<ul style="list-style-type: none">"Spec1.rep1_T" and "Spec1.rep1_N""Spec1.rep2_T" and "Spec1.rep2_N"
Negative control	NEG<*>.fsa where: <*> is user-defined Note: The TrueMark™ MSI Assay Amplification Control sample file name must begin with "NEG" to be properly analyzed. For example, "NEGAmpl.1_A12_daytimestamp.fsa".	File names that begin with "NEG" are analyzed as negative control samples. The TrueMark™ MSI Analysis Software displays a  (Warning) Review Flag if it assigns an Unstable call to any reportable marker for a negative control sample.
No template control	NTC<*>.fsa where: <*> is user-defined	File names that begin with "NTC" are no template control samples. The TrueMark™ MSI Analysis Software displays a  (Warning) Review Flag if it assigns a call other than No Call to any reportable marker for a no template control sample.

Table 6 Examples of files that will or will not import

Files selected for import	Import result
20190917.plate1.tst123_T_A01_datetime.fsa	Tumor/Normal pair of files imported as <SpecimenID> "20190917.plate1.tst123"
20190917.plate1.tst123_N_A02_datetime.fsa	
20190924.plate1.sample1.tumor_T_A02.fsa	Unpaired Tumor file imported as <SpecimenID> "20190924.plate1.sample1.tumor"
20190924.plate1.sample1.normal_N_A01.fsa	Unpaired Normal file is not imported
20190924.plate1.sample1.normal_A01.fsa	Unpaired file imported as <SpecimenID> "20190924.plate1.sample1.normal"
specimen1_injection1.fsa	Both files have the same <SpecimenID>, "specimen1". The software imports only the last file, in alphanumeric sort order. In this example, the software imports only "specimen1_injection2.fsa".
specimen1_injection2.fsa	
specimen2_something.fsa, from the 3500/3500xL Genetic Analyzer	Both files imported as <SpecimenIDs> "specimen2" and "specimen3", because each <SpecimenID> is unique. A single batch can contain specimens from different plates and instrument types.

Files selected for import	Import result
specimen3_something.fsa, from the SeqStudio™ Genetic Analyzer	Both files imported as <SpecimenIDs> "specimen2" and "specimen3", because each <SpecimenID> is unique. A single batch can contain specimens from different plates and instrument types.
specimen4_T_something.fsa, from the 3500/3500xL Genetic Analyzer	Neither file imported, because Tumor/Normal pair of files from different instrument types is not supported.
specimen4_N_something.fsa, from the SeqStudio™ Genetic Analyzer	
HiDi_something.fsa	Files with HiDi prefix in filename are rejected

Note: Files are resized upon import. The failure of one file to meet the sizing quality threshold may prevent any file in the batch from being imported. Open the **Import Manager** for more information on the files reporting errors.

DNA sample preparation

The TrueMark™ MSI Assay accepts DNA that has been extracted by various methods. This workflow has been tested on fixed formalin paraffin-embedded (FFPE) DNA extracted with the RecoverAll™ Total Nucleic Acid Isolation Kit for FFPE (Cat. No. AM1975). DNA inputs of 2–5 ng are recommended in order to optimize results, and to minimize offscale signals.

All protocols should be optimized with your standard laboratory procedures. See the *RecoverAll™ Total Nucleic Acid Isolation Kit Protocol* (Pub. No. 1975M) for more information.

Set up the PCR reactions

IMPORTANT! Perform all steps on ice.

Thaw the TrueMark™ MSI Assay Primer Mix, TrueMark™ MSI Assay Amplification Control, and TrueMark™ MSI Assay No-Template Control on ice. Thoroughly mix the components by vortexing 3–5 seconds, then centrifuge 3–5 seconds before use.

Recommended contents per 10 µL reaction

- TrueMark™ MSI Assay Master Mix — 4 µL
- TrueMark™ MSI Assay Primer Mix — 1 µL
- DNA — 2 to 5 ng

IMPORTANT! Tumor and normal sample pairs should have similar amounts of starting DNA in the PCR reactions to minimize PCR artifacts.

- TrueMark™ MSI Assay No-Template Control — to a final reaction volume of 10 µL.
1. On ice, prepare sufficient PCR reaction mix for the required number of reactions plus 1 additional reaction for overage.

Table 7 Example reaction mix (1 ng/ µL DNA)

Component	Volume per reaction	Volume (12 reactions)
TrueMark™ MSI Assay Master Mix	4 µL	48 µL
TrueMark™ MSI Assay Primer Mix	1 µL	12 µL
TrueMark™ MSI Assay No-Template Control	3 µL	36 µL
Total volume	8 µL	96 µL

2. To the labeled reaction plate, add the following components.
 - a. Add 8 µL of PCR reaction mix to each sample, amplification control, or no template control (NTC) well.
 - b. Add 2 µL of sample DNA (1 ng/ µL) to the sample wells.
 - c. Add 2 µL of TrueMark™ MSI Assay No-Template Control to NTC wells.
 - d. Add 1 µL of TrueMark™ MSI Assay Amplification Control and 1 µL of TrueMark™ MSI Assay No-Template Control to the amplification control wells.

Note: Input DNA quantity and quality affect fragment analysis results.

3. Cover the plate with adhesive film, then centrifuge 3–5 seconds to bring the mixture to the bottom of the tube and eliminate air bubbles.
4. Immediately proceed to “Run the PCR” on page 26.

Run the PCR

1. Program the thermal cycling conditions.

IMPORTANT! If you are using the ProFlex™ 96-well PCR System, select the GeneAmp™ PCR System 9700 simulation mode (**Edit ▶ Manage Steps ▶ Advanced Options ▶ Simulation Mode**).

Step	Temperature	Time	Cycles
Hot start	95°C	11 minutes	1
Denature	94°C	20 seconds	29
Anneal/ Extend	59°C	2 minutes	
Final extension	60°C	25 minutes	1
Hold	4°C	∞	

2. Set the reaction volume to 10 µL, then load the plate into the thermal cycler.
3. Close the heated cover, then start the run.
4. When the run is complete, remove the plate from the thermal cycler.

IMPORTANT! Protect the amplified DNA from light.

Amplified DNA can be stored at 2°C to 8°C for up to 2 weeks, or at –25°C to –15°C for long-term storage.

Prepare samples for electrophoresis

Prepare the samples for electrophoresis immediately before loading.

1. Prepare the mix of Hi-Di™ Formamide and GeneScan™ 600 LIZ™ Size Standard v2.0 for the required number of reactions plus 1 additional reaction for overage.

Component	Volume per reaction	Volume (12 reactions)
GeneScan™ 600 LIZ™ Size Standard v2.0	1 µL	12 µL
Hi-Di™ Formamide	17 µL	204 µL
Total volume	18 µL	216 µL

IMPORTANT! The volume of size standard is a suggested amount. Determine the appropriate amount based on your experiments and results.

2. Thoroughly mix the components by vortexing 3–5 seconds, then centrifuge 3–5 seconds before use.

3. Prepare the fragment analysis reactions. To a MicroAmp™ Optical 96-Well Reaction Plate, add the following components.

Component	Volume per reaction	Volume (12 reactions)
Hi-Di™ Formamide and GeneScan™ 600 LIZ™ Size Standard v2.0 mix	18 µL	216 µL
PCR product (see “Run the PCR” on page 26)	2 µL	24 µL
Total volume	20 µL	240 µL

Note: For blank wells, add 10 µL of Hi-Di™ Formamide.

	1	2	3	4	5	6	7	8	9	10	11	12
A	NEG	Specimen7_T	Hi-Di									
B	NTC	Specimen7_N	Hi-Di									
C	Specimen1_T	Specimen8_T	Hi-Di									
D	Specimen2_T	Specimen8_N	Hi-Di									
E	Specimen3_T	Specimen9_T	Hi-Di									
F	Specimen4	Specimen9_N	Hi-Di									
G	Specimen5	Specimen10_T	Hi-Di									
H	Specimen6	Specimen10_N	Hi-Di									
	Injection 1			Injection 2			Injection 3			Injection 4		

Figure 2 Example 3500xL Genetic Analyzer plate layout

4. Seal the reaction plate with MicroAmp™ Clear Adhesive Film.
5. Thoroughly mix the components by vortexing 3–5 seconds, then centrifuge 10–20 seconds before use.
6. Denature the DNA fragments:
 - a. Incubate the mixture at 95°C for 3 minutes.
 - b. Incubate the mixture at 4°C, or on ice, for 2 minutes.
7. Centrifuge the plate for 10–20 seconds to ensure that all sample mixtures are at the bottom of the wells.
8. Remove the MicroAmp™ Clear Adhesive Film, then seal the plate with a septa.
9. Assemble the plate with the retainer and base, then load on the instrument. Reactions can be run on the 3500/3500xL Genetic Analyzer or the SeqStudio™ Genetic Analyzer.

See the instrument user guide for specifics on setting up the run.



Analyze the data with the TrueMark™ MSI Analysis Software


■ Software access restrictions	28
■ Sign in to the TrueMark™ MSI Analysis Software (all users)	29
■ Import sample files	29
■ Check for import errors	31
■ Perform an initial QC of the imported data	31
■ View and interpret the results	32
■ View the electropherogram traces	33
■ Representative data	34
■ (Optional) Accept and approve a specimen	40
■ Generate a report	40
■ Export results	41

For information on data analysis or troubleshooting with GeneMapper™ Software, see the *TrueMark™ MSI Assay User Guide* (Pub. No. MAN0018868).

Software access restrictions



- The TrueMark™ MSI Analysis Software allows up to five users to be signed in simultaneously. When that limit is reached, no other users are allowed to sign in.
- Your access to functions in the software is based on the permissions associated with your user account. For more information, see the *TrueMark™ MSI Analysis Software User Guide* (Pub. No. MAN0018874).
- If your system is configured for password expiration, you will be periodically prompted to change your password.
- If your system is configured to monitor failed log in attempts, you will be locked out of the software if you incorrectly enter your user name or password more than the specified number of times.

Sign in to the TrueMark™ MSI Analysis Software (all users)

1. (First sign-in only) Obtain your user name and password from your TrueMark™ MSI Analysis Software Administrator.
2. On the computer desktop, double-click  (**MSI Client**) to start the client.

IMPORTANT! The MSI Server must already be running on the computer so that the client can connect to the server. The server is configured to automatically start whenever the computer is started.

If you get a connection error when you start the client, you may need to manually start the server. Try each of the following actions to resolve the connection error; perform the actions in the order listed.

1. On the computer desktop, double-click  (**MSI Server**).
 2. On the computer desktop, right-click  (**MSI Server**), then select **Run as administrator**.
 3. See the *TrueMark™ MSI Analysis Software User Guide* (Pub. No. MAN0018874).
-
3. Enter your **User Name** and **Password**, then click **Log in**.
 4. (The first time that you sign in to a new MSI Client) Accept the End User License Agreement to continue.

The first time that you sign in, the **Home** screen is empty. To begin using the software, you must import sample files (FSA). See “Import sample files” on page 29.


Import sample files

1. In the **Home** screen toolbar, click **Import Samples**.
2. Navigate to and select the sample files (FSA) to import.
You can select a folder to import all sample files in the folder, or you can select individual sample files. If you select individual sample files, ensure that you select the normal (**N**) and tumor (**T**) sample files for paired samples.

IMPORTANT! We recommend that you group sample files from different injections or plates in separate folders. If you import tumor and normal sample files that have the same <SampleID> but are from different injections or plates into the same batch, the software may pair the samples, which increases the risk of a miscall (false positive).



WARNING! Do not import more than 96 sample files in a single import operation. Importing more than 96 sample files may cause some or all specimens in the batch to be miscalled.



The default location is:

<installation drive>\Applied Biosystems\MSI Client\User Files\Import

Note: To change the default save location, see the *TrueMark™ MSI Analysis Software User Guide* (Pub. No. MAN0018874).

3. Click **Import**.
4. In the **Batch Information** dialog box, enter a **Batch Name**, (optional) enter the **Instrument ID** and **Operator**, then click **Save**.



Note: It may take several minutes to complete the import process.

The TrueMark™ MSI Analysis Software automatically analyzes the samples at import, using the current analysis settings. The results are displayed in the **Home** screen.

To change the analysis settings, see the *TrueMark™ MSI Analysis Software User Guide* (Pub. No. MAN0018874).

IMPORTANT! Changes to the analysis settings apply only to batches that are imported after you save the changes. To reanalyze samples with new analysis settings, import the samples again into a new batch. Changes do not affect existing batches.


Check for import errors

1. In the **Batches** pane of the **Home** screen, check that the number of imported specimens reflects the number of sample files (FSA) that you selected for import. Each specimen that has paired samples will have two sample files.
2. In the **Home** screen toolbar, check for a warning symbol  on the **Import Manager** button.
3. If the number of specimens and sample files do not agree, or if  appears on the **Import Manager** button, see the *TrueMark™ MSI Analysis Software User Guide* (Pub. No. MAN0018874).








IMPORTANT! We recommend that you also check the electropherogram plot title to confirm that the expected sample files were selected for each specimen (see *TrueMark™ MSI Analysis Software User Guide* (Pub. No. MAN0018874)). If you import tumor and normal sample files that have the same <SampleID> but are from different injections or plates into the same batch, the software may pair the samples, which increases the risk of a miscall (false positive).

Perform an initial QC of the imported data

1. In the **Home** tab, look for any specimens that display **Review Flag**.
2. In the **Home** tab, check the **Well Pair Status**. This column displays **Match** if the HID markers in the specimen both have similar fragment sizes between the tumor and normal files, indicating the samples likely belong to the same specimen.



The screenshot shows the TrueMark MSI Analysis Software 1.0 Home tab. The interface includes a top navigation bar with 'Assay', 'Tools', and 'Help'. The main content area is divided into a left sidebar and a central table. The sidebar contains a 'Batches' section with a dropdown menu showing 'All (1)'. Below this, there is a list of batches, with one batch selected: 'Example7-20191029091254'. The central table displays a list of specimens with the following columns: Specimen ID, Overall Call, Review Flag, Well Pair Status, Well, and Review Status. The table contains several rows of data, including specimens with 'MSS', 'No Call', and 'MSI-High' overall calls. Annotations 1, 2, and 3 point to the 'Home' tab, the 'Review Flag' column, and the 'Well Pair Status' column respectively.

Specimen ID	Overall Call	Review Flag	Well Pair Status	Well	Review Status
MSI1	MSS			B01	New
MSI49	No Call		No Call	C02,D02	New
MSI58	MSI-High		Match	E07,F07	New
NEG	MSS			H10	New
NTC	No Call			A10	New
S06-009596.2.1	MSI-High		Match	H04,G04	New
S07-004838.2.1	MSI-High		Match	B02,A02	New

- ① Home tab
- ② Review Flag
- ③ Well Pair Status

For more detailed information see the *TrueMark™ MSI Analysis Software User Guide* (Pub. No. MAN0018874).

View and interpret the results

Open the **Batch Summary** tab. The **Batch Summary** is populated when data import is complete.

The **Overall Call** for the specimen is based on the percentage of markers reporting instability.

- View the **Overall Call** for each specimen.
 - MSS**—Microsatellite stable
 - MSI-Low**—Low levels of microsatellite instability
 - MSI-High**—High levels of microsatellite instability
 - No Call**—At least one marker had no call.

Note: If every marker **No Call** is manually changed to a call of either **Unstable** or **Stable**, then the overall call will be adjusted from **No Call** to **MSS**, **MSI-Low**, or **MSI-High**.

- View the number of **Unstable**, **Stable**, or **No Call** markers.

Interpret the results

<div> <div>Batch Summary</div> <div>Specimen Data</div> <div>Approvals</div> <div>Audit Records</div> </div>						
<input type="checkbox"/>	Specimen ID	Overall Call	Unstable	Stable	No Call	Total
<input checked="" type="checkbox"/>	MSI1	MSS	0	13	0	13
<input type="checkbox"/>	MSI49	No Call	3	7	3	13
<input type="checkbox"/>	MSI58	MSI-High	13	0	0	13
<input type="checkbox"/>	Neg	MSS	0	13	0	13
<input type="checkbox"/>	NTC	No Call	0	0	13	13
<input type="checkbox"/>	S06-009596.2.1	MSI-High	13	0	0	13
<input type="checkbox"/>	S07-004838.2.1	MSI-High	13	0	0	13

- Overall Call for the specimen
- Number of **Unstable** markers
- Number of **Stable** markers
- Number of **No Call** markers
- Total number of markers used for the Overall Call

See the *TrueMark™ MSI Analysis Software User Guide* (Pub. No. MAN0018874) for more detailed information on viewing individual markers and changing calls.

View the electropherogram traces

1. Click the **Specimen ID** in the **Home** or **Batch Summary** tabs to open the relevant specimen in the **Specimen Data** tab.
2. Review the electropherogram traces in the **Specimen Data** tab.
 - a. Click the marker row in the **Marker Details** table to display the relevant dye and size range for that marker.

Note: The **Marker Details** table can be used to identify which markers have a warning in the **Review Flag** column, indicating that a review of the call is recommended. The **Marker Details** table can also be used to override the **Auto Call** with a **Manual Call**, and/or exclude the marker from the **Overall Call** assessment and from reports and exports.

Representative data

Examples of microsatellite instability in colon tumor tissue samples




Figure 3 Example electropherogram of a FFPE colon tumor tissue sample.

- ① Tumor trace—Displays extra peaks that are not present in the normal trace.
- ② Normal trace—Displays standard peaks for the marker.

The colon tumor tissue sample was amplified with the TrueMark™ MSI Assay, run on a SeqStudio™ Genetic Analyzer, and analyzed with TrueMark™ MSI Analysis Software (Y-axis scale 0–12,000 RFU). The extra peaks in the tumor trace have a smaller fragment size than the main peak for each marker, indicating microsatellite instability in the tumor sample.




Figure 4 Alternative view: tumor and normal traces are overlaid to more easily evaluate smaller differences in fragment size distribution.

Automatic calling of low frequency small deletions in synthetic constructs

Instability in colon tumor tissue microsatellites are readily identified by large base pair deletions. However, extra-colonic tumors can feature small deletions that are more difficult to distinguish. To accurately report microsatellite instability in colon and other tumor tissue types, it is important to identify these small deletions.

We designed synthetic constructs with varying lengths of small deletions for microsatellite markers in the TrueMark™ MSI Assay, then evaluated them at several allele frequencies. Examples of these small deletions are shown in Figure 5 and Figure 6.

IMPORTANT! Run the tumor and normal samples on the same plate. The software can give inaccurate results if you use tumor and normal sample files from different runs. Additionally, the signal strength of tumor and normal samples should match as closely as possible in order to detect small deletions or low allele frequencies.




Figure 5 Examples of microsatellite stability and instability in synthetic constructs.

The synthetic construct containing the wild-type allele is displayed in brown. The synthetic construct containing the mutant allele is displayed in blue.

- ① Marker BAT-25—Stable microsatellite with no deletions
- ② Marker BAT-25—Large unstable microsatellite with a 5 bp deletion at 50% allele frequency
- ③ Marker BAT-25—Unstable microsatellite with a 1 bp deletion at 50% allele frequency
- ④ Marker ABI-20A—Unstable microsatellite with a 1 bp deletion at 50% allele frequency

It is possible to detect small deletions (2–5 bp) at lower allele frequencies by visual inspection. However, single base pair deletions at 20% allele frequency are difficult to discern (Figure 6 and Figure 7).




Figure 6 Examples of small low frequency deletions in synthetic constructs.

The synthetic construct containing the wild-type NR-24 allele is displayed in brown. The synthetic construct containing the mutant NR-24 allele is displayed in blue.

- ① Unstable microsatellite with a 5 bp deletion at 20% allele frequency
- ② Unstable microsatellite with a 4 bp deletion at 20% allele frequency
- ③ Unstable microsatellite with a 3 bp deletion at 20% allele frequency
- ④ Unstable microsatellite with a 2 bp deletion at 20% allele frequency
- ⑤ Unstable microsatellite with a 1 bp deletion at 20% allele frequency

Automatic calling of low frequency small deletions in an endometrial tumor tissue sample

As an example of how the TrueMark™ MSI Analysis Software interprets low frequency small deletions, results from an endometrial tumor tissue sample with a relatively low fraction of tumor content (~25%) is shown in Figure 7.



Figure 7 Example of small low frequency deletions in a mismatch repair (MSH2 and MSH6) deficient endometrial tumor tissue sample.

The normal sample is displayed in brown. The tumor sample is displayed in blue. The tumor content of this specimen was 25%.

- | | |
|----------|-----------|
| ① BAT-25 | ⑧ ABI-19 |
| ② NR-24 | ⑨ ABI-20B |
| ③ NR-21 | ⑩ ABI-17 |
| ④ BAT-40 | ⑪ ABI-16 |
| ⑤ CAT-25 | ⑫ BAT-26 |
| ⑥ NR-22 | ⑬ ABI-20A |
| ⑦ NR-27 | |

The TrueMark™ MSI Analysis Software called 10 of 13 markers in the endometrial tumor tissue sample as unstable, with an overall call of **MSI-High** (Figure 8).

Overall Call: **MSI-High** | Unstable: **10 of 13** | Review Status: New

Marker Details | Peak Details | Review
















Dye	Marker	Auto Call	Review Flag	Include	Manual Call
	BAT-25	Unstable	✓	Yes ▾	▾
	NR-24	Unstable	✓	Yes ▾	▾
	NR-21	Unstable	✓	Yes ▾	▾
	TH01	Match	✓	Yes ▾	▾
	BAT-40	Unstable	✓	Yes ▾	▾
	CAT-25	Unstable	✓	Yes ▾	▾
	NR-22	Unstable	✓	Yes ▾	▾
	NR-27	Unstable	✓	Yes ▾	▾
	ABI-19	Unstable	✓	Yes ▾	▾
	ABI-20B	Unstable	✓	Yes ▾	▾
	PentaD	Match	✓	Yes ▾	▾
	ABI-17	Stable	✓	Yes ▾	▾
	ABI-16	Stable	✓	Yes ▾	▾
	BAT-26	Unstable	✓	Yes ▾	▾
	ABI-20A	Stable	✓	Yes ▾	▾

Figure 8 Marker calls in the TrueMark™ MSI Analysis Software

(Optional) Accept and approve a specimen

You must have permissions in the TrueMark™ MSI Analysis Software to accept (**Initial Review** and **Final Review**) and approve (**Approve Sample**) a specimen. For more information on set up and use of permissions, see the *TrueMark™ MSI Analysis Software User Guide* (Pub. No. MAN0018874).

1. In the **Specimen Data** tab, click the **Review** tab.
2. Click **Accept** to acknowledge that you have reviewed the specimen.
Note: The comment **Accepted** is displayed, and the **Review Status** changes to **Under Review**.
3. Select the **Ready for Approval** checkbox.
Note: The **Review Status** changes to **Ready for Approval**.
4. Open the **Batch Summary**.
5. Click the **Approve** button, then enter your user name and password.
Note: Users without **Approve** permissions will not be allowed to sign.

If a specimen is approved, the approval history of the specimen appears in the Specimen PDF report. See the *TrueMark™ MSI Analysis Software User Guide* (Pub. No. MAN0018874) for more information.

Generate a report

1. In the **Batch Summary** pane of the batch of interest, click **PDF Report**.
2. Select **Batch Summary** or **Specimen**.
A message is generated stating **Report generated successfully**. Click **Open folder location** to see where the report was saved, then click **OK**.

TrueMark™ MSI Analysis Software 1.0

Assay Tools Help

Technologist (Technologist) applied biosystems

Home Example7-20191... x

Batch ID: Example7-20191029091254 | Specimens: 7 | Approved: 0

Batch Summary Specimen Data Approvals Audit Records

Approve PDF Report Export Results

<input type="checkbox"/>	Specimen ID	Overall Call	Unstable	Stable	No Call	Total	Review Flag	Well Pair Status	Well	Review Status	Edited?	Comment
<input checked="" type="checkbox"/>	MSI1	MSS	0	13	0	13	⚠		B01	New	No	
<input type="checkbox"/>	MSI49	No Call	2	8	3	13	⚠	No Call	C02,D...	New	No	
<input type="checkbox"/>	MSI58	MSI-High	13	0	0	13	✓	Match	E07,F...	New	No	
<input type="checkbox"/>	NEG	MSS	0	13	0	13	✓		H10	New	No	
<input type="checkbox"/>	NTC	No Call	0	0	13	13	✓		A10	New	No	
<input type="checkbox"/>	S06-009596.2.1	MSI-High	13	0	0	13	✓	Match	H04,G...	New	No	
<input type="checkbox"/>	S07-004838.2.1	MSI-High	13	0	0	13	✓	Match	B02,A...	New	No	

Report Information

Report generated successfully.
[Open folder location](#)

OK

- ① PDF Report button
- ② Open folder location link
- ③ OK button

Export results

1. In the **Batch Summary** pane of the batch of interest, click **Export Results**.
2. Select the results format to export (**Batch Summary** or **Specimen**).
Batch Summary results is available in CSV format. Specimen results are available in CSV or VCF formats.
A message is generated stating **Export(s) generated successfully**. Click **Open folder location** to see where the results were saved, then click **OK**.

TrueMark™ MSI Analysis Software 1.0

Assay Tools Help Technologist (Technologist) applied biosystems

Home Example7-20191... x

Batch ID: Example7-20191029091254 | Specimens: 7 | Approved: 0

Batch Summary Specimen Data Approvals Audit Records

Approve PDF Report Export Results

<input type="checkbox"/>	Specimen ID	Overall Call	Unstable	Stable	No Call	Total	Review Flag	Well Pair Status	Well	Review Status	Edited?	Comment
<input checked="" type="checkbox"/>	MSI1	MSS	0	13	0	13	⚠		B01	New	No	
<input type="checkbox"/>	MSI49	No Call	2	8	3	13	⚠	No Call	C02,D...	New	No	
<input type="checkbox"/>	MSI58	MSI-High	13	0	0	13	✓	Match	E07,F...	New	No	
<input type="checkbox"/>	NEG	MSS	0	13	0	13	✓		H10	New	No	
<input type="checkbox"/>	NTC	No Call	0	0	13	13	✓		A10	New	No	
<input type="checkbox"/>	S06-009596.2.1	MSI-High	13	0	0	13	✓	Match	H04,G...	New	No	
<input type="checkbox"/>	S07-004838.2.1	MSI-High	13	0	0	13	✓	Match	B02,A...	New	No	

Export Information

Export(s) generated successfully.

[Open folder location](#)

OK

- ① Export Results button
- ② Open folder location link
- ③ OK button



Troubleshooting

Observation	Possible cause	Recommended action
Faint or no signal from both the TrueMark™ MSI Assay Amplification Control and the DNA test samples for all markers	The incorrect volume of Master Mix or Primer Set was used.	Use the correct volume of Master Mix or Primer Set.
	The Master Mix was not vortexed thoroughly before aliquoting.	Vortex the Master Mix thoroughly.
	The Primer Set was exposed to too much light.	Replace the Primer Set and store it protected from light.
	Evaporation.	Ensure that the adhesive seal or optical caps are properly sealed.
	The thermal cycler malfunctioned.	See the thermal cycler user guide and check the instrument calibration.
	Incorrect thermal cycler conditions were used.	Use correct thermal cycler conditions.
	The wrong PCR reaction tubes or plate were used.	Use MicroAmp™ Optical 96-Well Reaction Plate with Barcode with MicroAmp™ Clear Adhesive Film for the ProFlex™ Thermal Cycler.
	Insufficient PCR product was electrokinetically injected.	Use the correct genetic analyzer settings.
	Degraded formamide was used.	Check the storage of formamide. Do not thaw and refreeze multiple times. Use fresh Hi-Di™ Formamide.
Signal from the TrueMark™ MSI Assay Amplification Control but partial or no signal from DNA test samples	The test sample was diluted in the wrong buffer (for example, a TE buffer with an incorrect EDTA concentration).	Redilute the DNA using low-TE buffer (with 0.1 mM EDTA).
	Less than the recommended amount of DNA was added to the PCR reaction.	Quantitate the DNA sample using a Qubit™ dsDNA HS Assay Kit.
		Increase the injection time to boost the signal of the sample.
More than the expected number of peaks present for a sample identification marker (TH01 and PentaD)	Exogenous DNA is present in the sample.	Increase the number of PCR cycles to boost the signal of the sample.
		Use appropriate techniques to avoid introducing foreign DNA during laboratory handling.
		Include the final extension step of 60°C for 25 minutes in the PCR.
	Incomplete 3' A base addition (n-1 nt position) occurred.	Incubate the amplification products at 60°C for an additional 10 minutes.

Observation	Possible cause	Recommended action
More than the expected number of peaks present for a sample identification marker (TH01 and PentaD) <i>(continued)</i>	The signal exceeds the dynamic range of the instrument and is causing signal "pull-up" into adjacent channels.	Ensure the cycle number is optimized. User fewer PCR cycles or interpret the off-scale data according to your laboratory procedure.
		Decrease the injection time in the run module.
		Load less DNA in the PCR reaction.
	Poor spectral separation occurred.	Perform a spectral calibration.
		Confirm that Filter Set J6 modules are installed and used for analysis.
	The double-stranded DNA was not completely denatured.	Use the recommended amount of Hi-Di™ Formamide and heat the sample plate at 95°C for 3 minutes.
Some but not all markers visible on electropherogram of DNA Test Samples	Contamination was carried over from the DNA extraction.	Perform the DNA extraction again.
	MSI present in the sample is reflected in the tumor tissue sample identification markers.	Assess the level of instability in the overall sample to assist in determining whether the extra peak(s) are due to contamination or instability.
Marker profiles contain many off-scale peaks	The PCR reaction volume used is lower than the volume that is required for the amplification.	Use the correct PCR reaction volume: 10 µL
	There are fewer large DNA fragments due to FFPE fragmentation.	Use more PCR cycles, increase the injection time of the run module, or load more DNA into the PCR reaction to boost signal in the sample.
Marker profiles contain many off-scale peaks	The PCR cycle number used was too high.	Perform a sensitivity experiment to determine the optimal PCR cycle number based on the sample type.
		Decrease the injection time in the run module.

Extra peaks in the electropherogram

Causes of extra peaks

Peaks other than the target alleles may be detected on the electropherogram. Causes for the appearance of extra peaks include stutter products, incomplete 3' A nucleotide addition (at the n-1 position), dye artifacts, and non-specific amplification.

Stutter definition

Stutter is a well-characterized PCR artifact that refers to the appearance of a minor peak one repeat unit smaller than the target STR allele product (minus stutter), or less frequently, one repeat larger (plus stutter). Sequence analysis of stutter products at tetranucleotide STR loci has revealed that the minus stutter product is missing a single tetranucleotide core repeat unit relative to the main allele.

The level of stutter in this kit is normal and as expected for STR chemistries that are designed to overcome inhibitors and obtain robust performance with single source reference samples.

Artifact definition

Artifacts and anomalies are seen in all molecular biological systems. Artifacts are typically reproducible. Anomalies are non-reproducible, intermittent occurrences that are not observed consistently in a system (for example, spikes and baseline noise).

Artifact observation

Due to improvements in PCR primer manufacturing processes, the incidence of artifacts has been greatly reduced in the TrueMark™ MSI Assay. Electropherograms are free of reproducible dye artifacts in the kit read region of 80–215 nt for commonly used analytical thresholds.




Figure 9 Example of low baseline-level fluorescence in a typical no template control (NTC) PCR. Data produced on a 3500xL Genetic Analyzer (Y-axis scale 0 to 200 RFU).

Table 8 Marker-specific artifacts observed with the TrueMark™ MSI Assay

Artifact	Color	Size	Type
NR-22	Black	75–90 nt	Incomplete amplification
TH01	Green	80–100 nt	Incomplete amplification of BAT-40
SID208–213	Purple	208–213 nt	Non-specific amplification. This set of non-specific peaks is much larger than the nearest MSI allele, ABI-20A, and should not be ignored.

Marker-specific artifacts commonly observed in the TrueMark™ MSI Assay are shown in Figure 10 to Figure 12.




Figure 10 NR-22




Figure 11 TH01




Figure 12 Non-specific amplification in SID channel

Additional reproducible DNA-dependent artifacts have been characterized and documented. It is important to consider noise and other amplification-related artifacts when interpreting data.



Analyze the data with GeneMapper™ Software

Overview of the GeneMapper™ Software

GeneMapper™ Software is a flexible genotyping software package that provides DNA sizing and quality allele calls for genotyping systems. This software specializes in multi-application functionality, including amplified fragment length polymorphism (AFLP™), loss of heterozygosity (LOH), microsatellite and SNP genotyping analysis. GeneMapper™ Software uses Process Quality Values (PQVs) for automated identification that reduces data review time for high throughput genotyping.


Set up the GeneMapper™ Software (v5.0 or later) for analysis (before first use of the kit)

Define custom display and plot settings

1. Navigate to the **Plot Settings Editor**.
2. Select **Display Settings**, then enter the following display settings.
For Sample and Genotype plots
 - 5 panes
 - Select the graphs highlighted in blue
 - X-Axis—Basepairs
 - Y-Axis—Scale individually
 - Select **Toolbar**
For Sample plot
 - Select the first five colors, then the graphs highlighted in blue
 - Start Range—68.0 bp
 - End Range—225.0 bp

For Genotype plot

- Marker Margin—5 bp




- ① Display Settings tab
- ② For both Sample and Genotype plots pane
- ③ For Sample plot only pane
- ④ For Genotype plot only pane
- ⑤ OK button

3. Click **OK**.

Define the Size Standard

1. Navigate to the **Size Standard Editor**, then click **New**.
2. (Optional) Enter a description of the Size Standard.
3. Specify settings in the Size Standard Editor.
 - a. Enter a name as shown in the following figure [GS600(80-400)] or enter a new name.
 - b. In the **Security Group** field, select the **Security Group** appropriate for your software configuration.
 - c. In the **Size Standard Dye** field, select **Orange**.

- d. In the **Size Standard Table**, use sizes between **80.0** and **400.0** bp (copied from the GS600LIZ default size standard).



- ① Size Standard Table
- ② OK button

4. Click OK.

Create an analysis method

Create an analysis method

IMPORTANT! Analysis methods are version-specific, so you must create an analysis method for each version of the software.

1. Select **Tools ▸ GeneMapper Manager** to open the **GeneMapper Manager**.
2. Click the **Analysis Methods** tab, then click **New** to open the **Analysis Method Editor** with the **General** tab selected.

3. Enter the settings shown in the figures on the following pages.

Note: The **Analysis Method Editor** closes when you save your settings. To complete this step quickly, do not save the analysis method until you finish entering settings in all of the tabs.

4. After you enter the settings on all tabs, click **Save**.

Analysis method settings

General tab settings

The screenshot shows the 'Analysis Method Editor - Microsatellite' dialog box with the 'General' tab selected. The 'Analysis Method Description' section contains the following fields:

- Name:** ABI_MSI
- Description:** Analysis method for Applied Biosystems MSI Assay
- Instrument:** (empty text box)
- Analysis Type:** Microsatellite

At the bottom of the dialog are 'OK' and 'Cancel' buttons.

In the **Name** field, either type the name as shown or enter a name. In the **Security Group** field, select the security group appropriate to your software configuration from the list. The **Description** and **Instrument** fields are optional.

Allele tab settings

Analysis Method Editor - Microsatellite

General **Allele** Peak Detector Peak Quality Quality Flags

Bin Set: None

Marker Repeat Type

☐ Use marker-specific stutter ratio if available

Values for dinucleotide repeats are calculated automatically.

		Mono	Tri	Tetra	Penta	Hexa
Cut-off value		0.25	0.2	0.25	0.25	0.25
PlusA ratio		0.0	0.95	0.95	0.95	0.95
PlusA distance		0.0	1.6	1.6	1.6	1.6
Stutter ratio		0.0	0.95	0.15	0.15	0.15
Stutter distance	From	0.0	0.0	0.0	0.0	0.0
	To	0.0	3.5	4.5	5.5	6.5

Range Filter... Factory Defaults

OK Cancel

The following settings were used during development of the TrueMark™ MSI Assay:

- In the **Bin Set** field, select **None**.
- GeneMapper™ Software allows you to specify 4 types of marker repeat motifs: tri, tetra, penta, and hexa. You can enter parameter values for each type of repeat in the appropriate column.
- Deselect the **Use marker specific stutter ratio if available** checkbox (selected by default).

Perform appropriate internal studies to determine the appropriate filter setting to use.

Peak Detector tab settings

The screenshot shows the 'Analysis Method Editor - Microsatellite' dialog box with the 'Peak Detector' tab selected. The settings are as follows:

- Peak Detection Algorithm:** Advanced
- Ranges:**
 - Analysis:** Full Range
 - Sizing:** All Sizes
 - Start Pt:** 0
 - Stop Pt:** 10000
 - Start Size:** 0
 - Stop Size:** .000
- Smoothing and Baseline:**
 - Smoothing:** None (selected), Light, Heavy
 - Baseline Window:** 51 pts
- Size Calling Method:**
 - ☐ 2nd Order Least Squares
 - ☐ 3rd Order Least Squares
 - ☐ Cubic Spline Interpolation
 - ☒ Local Southern Method
 - ☐ Global Southern Method
- Peak Detection:**
 - Peak Amplitude Thresholds:**
 - B:** 200
 - G:** 200
 - Y:** 200
 - R:** 200
 - P:** 200
 - O:** 200
 - Min. Peak Half Width:** 2 pts
 - Polynomial Degree:** 3
 - Peak Window Size:** 5 pts
 - Slope Threshold:**
 - Peak Start:** 0.0
 - Peak End:** 0.0
 - Size Standard Normalization:**
 - ☐ Enable Normalization
 - Note: For 35XX series data collection normalization only.

Buttons: Factory Defaults, OK, Cancel

IMPORTANT! Perform the appropriate internal studies to determine the appropriate peak amplitude thresholds for interpretation of data.

Fill in the field information:

- **Peak Detection Algorithm**—Basic.
- **Minimum Peak Height**—Automatic.

Peak Quality tab settings

Analysis Method Editor - Microsatellite X

General
Allele
Peak Detector
Peak Quality
Quality Flags

Signal level

Homozygous min peak height

Heterozygous min peak height

Heterozygote balance

Min peak height ratio

Peak morphology

Max peak width (basepairs)

Pull-up peak

Pull-up ratio

Pull-up scan

Allele number

Max expected alleles

Cross-talk peak

Cross-talk ratio

Factory Defaults

OK
Cancel

IMPORTANT! Perform the appropriate internal studies to determine the heterozygous and homozygous minimum peak height thresholds, maximum peak height threshold, and the minimum peak height ratio threshold for interpretation of data.

Quality Flags tab settings

Analysis Method Editor - Microsatellite

General Allele Peak Detector Peak Quality **Quality Flags**

Quality weights are between 0 and 1.

Quality Flag Settings

Spectral Pull-up (SPU)	0.5	Control Concordance (CC)	0.5
Broad Peak (BD)	0.5	Low Peak Height (LPH)	0.5
Single Peak Artifact (SPA)	0.5	Off-scale (OS)	0.5
Sharp Peak (SHP)	0.5	Peak Height Ratio (PHR)	0.5
Cross Talk (XTLK)	0.5	One Basepair Allele (OBA)	0.5
Out of Bin Allele (BIN)	0.8	Split Peak (SP)	0.5

PQV Thresholds

	Pass Range:	Low Quality Range:
Sizing Quality:	From 0.75 to 1.0	From 0.0 to 0.5
Genotype Quality:	From 0.75 to 1.0	From 0.0 to 0.25

Assume Linearity From (bp): 0 To (bp): 800

Factory Defaults

OK Cancel

IMPORTANT! The values that are shown are the values used during assay development. Perform appropriate internal studies to determine the appropriate values to use.

Analyze and edit sample files with GeneMapper™ Software

1. In the **Project** window, select **Edit ► Add Samples to Project**, then navigate to the disk or directory that contains the sample files.
2. Apply analysis settings to the samples in the project.

Parameter	Settings
Sample Type	Select the sample type.
Analysis Method	Select the ABI MSI Assay analysis method.
Panel	Select the ABI MSI Panel .
Size Standard	Use GS600(80-400) with a size range of 80–400 bp for the Local Southern size calling method.

3. Click **Analyze**, enter a name for the project (in the **Save Project** dialog box), then click **OK** to start analysis.
 - The status bar displays the progress of analysis as a completion bar.
 - The table displays the row of the sample currently being analyzed in green (or red if analysis failed for the sample).



Safety



WARNING! GENERAL SAFETY. Using this product in a manner not specified in the user documentation may result in personal injury or damage to the instrument or device. Ensure that anyone using this product has received instructions in general safety practices for laboratories and the safety information provided in this document.

- Before using an instrument or device, read and understand the safety information provided in the user documentation provided by the manufacturer of the instrument or device.
 - Before handling chemicals, read and understand all applicable Safety Data Sheets (SDSs) and use appropriate personal protective equipment (gloves, gowns, eye protection, and so on). To obtain SDSs, see the “Documentation and Support” section in this document.
-

Chemical safety



WARNING! GENERAL CHEMICAL HANDLING. To minimize hazards, ensure laboratory personnel read and practice the general safety guidelines for chemical usage, storage, and waste provided below. Consult the relevant SDS for specific precautions and instructions:

- Read and understand the Safety Data Sheets (SDSs) provided by the chemical manufacturer before you store, handle, or work with any chemicals or hazardous materials. To obtain SDSs, see the "Documentation and Support" section in this document.
- Minimize contact with chemicals. Wear appropriate personal protective equipment when handling chemicals (for example, safety glasses, gloves, or protective clothing).
- Minimize the inhalation of chemicals. Do not leave chemical containers open. Use only with sufficient ventilation (for example, fume hood).
- Check regularly for chemical leaks or spills. If a leak or spill occurs, follow the manufacturer cleanup procedures as recommended in the SDS.
- Handle chemical wastes in a fume hood.
- Ensure use of primary and secondary waste containers. (A primary waste container holds the immediate waste. A secondary container contains spills or leaks from the primary container. Both containers must be compatible with the waste material and meet federal, state, and local requirements for container storage.)
- After emptying a waste container, seal it with the cap provided.
- Characterize (by analysis if needed) the waste generated by the particular applications, reagents, and substrates used in your laboratory.
- Ensure that the waste is stored, transferred, transported, and disposed of according to all local, state/provincial, and/or national regulations.
- **IMPORTANT!** Radioactive or biohazardous materials may require special handling, and disposal limitations may apply.

Biological hazard safety



WARNING! Potential Biohazard. Depending on the samples used on this instrument, the surface may be considered a biohazard. Use appropriate decontamination methods when working with biohazards.



WARNING! BIOHAZARD. Biological samples such as tissues, body fluids, infectious agents, and blood of humans and other animals have the potential to transmit infectious diseases. Conduct all work in properly equipped facilities with the appropriate safety equipment (for example, physical containment devices). Safety equipment can also include items for personal protection, such as gloves, coats, gowns, shoe covers, boots, respirators, face shields, safety glasses, or goggles. Individuals should be trained according to applicable regulatory and company/ institution requirements before working with potentially biohazardous materials. Follow all applicable local, state/provincial, and/or national regulations. The following references provide general guidelines when handling biological samples in laboratory environment.

- U.S. Department of Health and Human Services, *Biosafety in Microbiological and Biomedical Laboratories (BMBL)*, 5th Edition, HHS Publication No. (CDC) 21-1112, Revised December 2009; found at:
<https://www.cdc.gov/labs/pdf/CDC-BiosafetymicrobiologicalBiomedicalLaboratories-2009-P.pdf>
 - World Health Organization, *Laboratory Biosafety Manual*, 3rd Edition, WHO/CDS/CSR/LYO/2004.11; found at:
www.who.int/csr/resources/publications/biosafety/Biosafety7.pdf
-



Documentation and support

Related documentation

Document	Publication number
<i>TrueMark™ MSI Assay Quick Reference Guide</i>	MAN0018869
<i>TrueMark™ MSI Analysis Software User Guide</i>	MAN0018874
<i>3500/3500xL Genetic Analyzer with 3500 Series Data Collection Software 3.1 User Guide</i>	100031809
<i>3500/3500xL Genetic Analyzer with 3500 Series Data Collection Software v3.3 User Guide</i>	100079380
<i>SeqStudio™ Genetic Analyzer Instrument and Software User Guide</i>	MAN0016138
<i>SeqStudio™ Genetic Analyzer Instrument and Software User Guide</i>	MAN0018646
<i>GeneMapper™ Software v4.1 Quick Reference Guide</i>	4403615

Customer and technical support

Visit thermofisher.com/support for the latest service and support information.

- Worldwide contact telephone numbers
- Product support information
 - Product FAQs
 - Software, patches, and updates
 - Training for many applications and instruments
- Order and web support
- Product documentation
 - User guides, manuals, and protocols
 - Certificates of Analysis
 - Safety Data Sheets (SDSs; also known as MSDSs)

Note: For SDSs for reagents and chemicals from other manufacturers, contact the manufacturer.



Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale at www.thermofisher.com/us/en/home/global/terms-and-conditions.html. If you have any questions, please contact Life Technologies at www.thermofisher.com/support.

