
3

**Tristel**™

Rinse Wipe

## USER INSTRUCTIONS

- Disinfect hands and wear gloves when handling disinfectants and medical devices.
- Wear appropriate Personal Protective Equipment (PPE).
- Store between 10°C and 35°C in the original packaging. Do not freeze.
- Do not use if the Wipe is damaged, discoloured, or dry.
- Keep the empty Wipe sachet(s) to complete the Tristel Quality Audit Trail Record Book at the end of the decontamination procedure.
- Continuous wiping of the surface is not necessary and should be avoided to prevent shedding of the Wipe.
- For professional use only.
- These Instructions for Use (IFU) should be used in conjunction with the product label, Safety Data Sheet (SDS) and device manufacturer's instructions.
- Do not use if the Wipe sachet is damaged.
- Do not use past the expiry date given on the sachet. Discard expired Wipes according to local regulations.
- For batch number and expiry date see sachet.

### Intended Use

The third and final step in the decontamination procedure is rinsing of the medical device. The Rinse Wipe is impregnated with de-ionised water and a low-level of antioxidant which removes chemical residues from a surface. Each Rinse Wipe sachet is packed and then sterilised by gamma-irradiation.

#### STEP 1

Take one Rinse Wipe sachet.

#### STEP 2

Remove the Wipe from its sachet and lay it out in the palm of your hand.

#### STEP 3

Wipe the surface of the device that has been decontaminated to remove excess foam.

#### STEP 4

Discard the used Wipe and gloves in accordance with local regulations. Do not reuse the Wipe. Keep the empty Wipe sachet for traceability.

**Note:** Upon completion of the decontamination cycle the device should be left to air dry. Store the device in accordance with hospital protocols to prevent damage or recontamination.

**Please complete the Tristel Quality Audit Trail Record Book to keep record of the decontamination procedure.**

If you believe that this device has contributed to a serious incident (i.e. a death, the temporary or permanent serious deterioration of a person's state of health, or a serious public health threat), contact the device manufacturer and your national competent authority immediately.

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# IMPORTERS

**EUROPEAN UNION (EU) / EUROPEAN ECONOMIC AREA (EEA)  
/ NORTHERN IRELAND / SWITZERLAND / TURKEY**

**AUSTRIA, BULGARIA, CROATIA, CZECH REPUBLIC,  
ESTONIA, GERMANY, GREECE, HUNGARY, LATVIA,  
LITHUANIA, POLAND, SLOVAKIA, SLOVENIA**

**Tristel GmbH**

Karl-Marx-Allee 90A, 10243 Berlin, Germany

T +49 (0)30 54844226 E deutschland@tristel.com

W www.tristel.de

**BELGIUM, FRANCE, LUXEMBOURG, NETHERLANDS**

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**NORTHERN IRELAND**

**BRENNAN & CO NI**

Unit 6 & 7, Carrowreagh Business Park, Dundonald, Belfast, Northern Ireland, BT16 1QQ

T +44 2890 410 136 E info@brennanconi.com

W www.brennanconi.com

**NORWAY**

**VINGMED AS**

Solbråveien 13, N-1383 Asker, Norway

T +47 67 58 06 80 E info@vingmed-as.no

W www.vingmed-as.no

**SPAIN**

**VESISMIN, S.L.**

C/ Aribau 230-240, 6ª planta, L-O, 08006 Barcelona, Spain

T +34 934 09 53 01 E vesismin@vesismin.com

W www.vesismin.com

**SWITZERLAND**

**TRISTEL AG**

Sandgrube 29, 9050 Appenzell, Switzerland

T +41 (0)71 5670658 E schweiz@tristel.com

W www.tristel.com

**TURKEY**

**YEFE DIS TICARET LTD STI.**

Cennetoğlu Mah. Yeşillik Cad. 228 Karabağlar, İzmir 35400, Turkey

T +90 532 275 1133 E efeayvali@pelisLtdcom

W www.yefe.com.tr

For any further information, please contact us at mail@tristel.com.

# Tristel™

## WE HAVE CHEMISTRY.

**Manufactured by:** Tristel Solutions Limited  
Lynx Business Park, Fordham Road, Snailwell, Cambridgeshire CB8 7NY  
T +44 (0)1638 721500 - F +44 (0)1638 721911 - E [mail@tristel.com](mailto:mail@tristel.com) - W [www.tristel.com](http://www.tristel.com)

**Australian sponsor:** Tristel Pty Ltd. 40/328 Reserve Road, Cheltenham, VIC 3192, Australia  
T 1300 680 898 - E [mail-au@tristel.com](mailto:mail-au@tristel.com)

Tristel Rinse Wipes is entered into the Australian Register of Therapeutic Goods as a Class I medical device, ARTG number: 173046.

For Tristel patent information please visit [www.our-patents.info/tristel](http://www.our-patents.info/tristel)

Contact Tristel, your local distributor or visit [www.tristel.com](http://www.tristel.com) for supporting documents such as safety data sheets and reports.

Where necessary, external translations have been sourced from ISO 17100:2015 approved supplier: [www.tolingo.com](http://www.tolingo.com).

**EC REP** Tristel NV, Smallandlaan 14B,  
2660 Antwerp, Belgium

**CH REP** Tristel AG, Sandgrube 29,  
9050 Appenzell, Switzerland

RIN/098/5  
14-JUN-2022

**STERILE R**