
AQT90 FLEX

β hCG Test Kit

Contents

en	βhCG Test Kit.....	2
bg	Набор тестове за βhCG.....	7
cs	Souprava testů βhCG.....	12
da	βhCG-testkit.....	17
de	βhCG-Testkit.....	22
el	Κιτ εξέτασης βhCG.....	27
es	Kit de Test de βhCG.....	32
et	βhCG mõõtekomplekt.....	37
fi	βhCG-testaussarja.....	42
fr	Kit de tests βhCG.....	47
hr	Testni komplet βhCG.....	52
hu	βhCG tesztkészlet.....	57
it	Kit per il test βhCG.....	62
lt	βhCG testo rinkinys.....	67
lv	βhCG testēšanas komplekts.....	72
nl	βhCG-testkit.....	77
no	βhCG-testkit.....	82
pl	Zestaw testowy do oznaczania βhCG.....	86
pt	Kit de Testes de βhCG.....	91
ro	Kit de testare pentru βhCG.....	96
ru	Набор тестов на бета-ХГЧ.....	101
sk	Súprava na testovanie βhCG.....	107
sv	βhCG-testkit.....	112
tr	βhCG Test Kiti.....	117
zh	βhCG 检测试剂盒.....	121
References		126

en

βhCG Test Kit

REF 942-918

Intended use

For *in vitro* diagnostic use.

The βhCG Test is an *in vitro* diagnostic assay for the quantitative determination of human chorionic gonadotropin plus beta subunit human chorionic gonadotropin in EDTA or lithium-heparin whole blood or plasma specimens on the AQT90 FLEX analyzer in point-of-care and laboratory settings. It is indicated for use as an aid in the early detection of pregnancy. It is not indicated as a surrogate marker in the diagnosis or monitoring of cancer patients.

WARNING – Risk of incorrect results

The βhCG Test is not indicated as a surrogate marker in the diagnosis or monitoring of cancer patients.

The βhCG CAL cartridges are *in vitro* diagnostic reagents intended for the calibration adjustment of the βhCG Test on the AQT90 FLEX analyzer by establishing points of reference to estimate total βhCG values.

NOTE: The cartridges are ready for use.

Test

Human chorionic gonadotropin (hCG) is a glycoprotein hormone. It is secreted during pregnancy by the trophoblastic cells of the placenta, shortly after the implantation of the fertilized ovum into the uterine wall [1,2]. Physiologically, hCG appears to maintain the corpus luteum and its synthesis of progesterone and estrogens that support the endometrium until the 4th month of pregnancy, by which time the placenta takes over the production of the hormones.

hCG appears in the blood stream 6–8 days after conception. The plasma concentrations of biologically active hCG rise exponentially in the first trimester of pregnancy. It doubles every 48 hours and reaches a peak at about 10 weeks after the last menstrual period.

Concentrations decrease from the 10th to the 16th week of gestation, reaching approximately one-fifth of the peak. It remains around this concentration until the delivery [3], although there can be large individual variations in hCG concentrations. After normal delivery, the hCG concentrations decline to those observed in non-pregnant women in 1–3 weeks [4,5]. The rapid rise in hCG levels after conception makes it an excellent marker for early confirmation. hCG is also secreted by certain cancers [3]. However, this assay is not for use in the diagnosis or monitoring of cancer patients.

Human chorionic gonadotropin is composed of two highly glycosylated dissimilar subunits, α and β, which are joined non-covalently. The subunits are rapidly disassociated and inactivated in the blood stream. Thus, both free α and β subunits together with intact

hCG are found in circulation. Only intact hCG and β subunits have immunological specificity as the α subunit of hCG is similar to that of the α subunit of the pituitary glycoprotein hormones, luteinizing hormone (LH), follicle stimulating hormone (FSH) and thyroid stimulating hormone (TSH). The β subunit, in turn, is unique and distinguishes hCG from the other glycoprotein hormones [6].

The antibodies used in the βhCG assay recognize both intact hCG and free β subunits.

Calibration adjustment

The assay has been factory-calibrated using a 8-point calibration. The factory-defined calibration data is provided in the form of a barcode. This data is used as a lot-specific reference curve, which is adjusted to the correct analyzer-specific signal level by using a CAL Cartridge with the same lot number. The results obtained using the CAL Cartridge are combined with the factory-defined calibration data to produce a lot- and analyzer-specific calibration curve for the assay. When a CAL Cartridge with a new lot number is inserted, the analyzer requests the factory-defined calibration data so that it can do the calibration adjustment. If the user does not insert a CAL Cartridge when the Test Cartridge lot changes, the analyzer will request it as soon as it detects a Test Cartridge with a new lot number. Test Cartridges with a new lot number can only be used after a calibration adjustment is done with a CAL Cartridge with the same lot number.

See the AQT90 FLEX *instructions for use* document for details.

Principles of operation

The AQT90 FLEX analyzer is a fully automatic continuous-access analyzer utilizing immunoassay technology and time-resolved fluorometric detection. See the AQT90 FLEX *instructions for use* document for details.

The AQT90 FLEX analyzer is operated via a touch screen. The analyzer aspirates samples from capped sample tubes that are inserted manually.

The analyzer does all assay steps automatically and reports quantitative results. See the AQT90 FLEX *instructions for use* document for details.

Contents

- 10 βhCG Test Cartridges for 160 tests
- One βhCG CAL Cartridge for one calibration adjustment
- A factory-defined calibration data sheet with barcode

Components

- Each Test Cartridge contains 16 analyte-specific assay cups
- The CAL Cartridge contains 8 analyte-specific background cups and 8 cups with added antigen

Storage

Store unopened at 2-8 °C (36-46 °F).

Shelf life

If the protective pouch containing the cartridge remains sealed until just before the cartridge is inserted into the analyzer, it is stable until the end of the month stated in the expiry date shown on the pouch label.

On-board stability

The value (in days) is printed on the cartridge pouches alongside the following symbol:

Chemical composition

Reactive ingredients:

- Mouse monoclonal anti- β hCG capturing and tracing antibodies: approx. 600 ng/cup
- Native β hCG (only in the 8 CAL Cartridge cups with added antigen): approx. 0.01 IU/cup

Other ingredients:

- Bovine serum albumin
- Bovine γ -globulin
- Mouse IgG (blocker substances for heterophilic interferences)

Main buffer component: TRIS (tris hydroxymethylaminomethane)

Toxic ingredient: Sodium azide <0.1 %

Data in cartridge barcodes

- Checksum
- Barcode type
- Parameter code
- Lot no.
- Cartridge ID
- Expiry date
- Barcode version
- On-board stability

NOTE: The barcodes are lot-specific.

Product calibrator traceability

The calibration of the β hCG assay is traceable to WHO 75/589.

About using the cartridges

WARNING – Risk of incorrect results

Cartridges are sensitive to humidity and must therefore not be taken out of their sealed pouches until just before use. Only use cartridges that have been kept according to storage specifications and are not damaged.

WARNING – Risk of incorrect results

Do not use a cartridge that has already been used on another analyzer, otherwise you risk obtaining incorrect results.

WARNING – Risk of infection

This product contains materials from human source and should be handled and disposed of as potentially infectious. The plasma from human donors used in the manufacture of this product has been tested and found non-reactive to HbsAg, HIV-1, HIV-2 and HCV.

Sample collection and preparation

Sample tubes approved for use with the analyzer

See the *AQT90 FLEX instructions for use* document for details.

Anticoagulants

Use either EDTA or lithium heparin.

Sample size

The minimum sample volume is 2 mL (independent of the number of tests requested).

Sample type

Use whole-blood or plasma samples.

Sample collection

Collect blood samples by venipuncture.

WARNING – Risk of incorrect results

For sample collection, do not use tubes that contain a gel.

Sample handling and storage

Samples that cannot be analyzed immediately after collection must be handled and stored correctly before they are analyzed. The table provides an overview.

Sample type	Handling and storage	Analyze within...
Whole blood	1. Keep at an ambient storage temperature 18-25 °C (65-77 °F).	3 hours*
Plasma	1. Separate the plasma immediately after the sample is collected. 2. Keep at an ambient storage temperature 18-25 °C (65-77 °F).	3 hours*

Sample type	Handling and storage	Analyze within...
Plasma	1. Separate the plasma immediately after the sample is collected. 2. Refrigerate at 2-8 °C (36-46 °F).	24 hours*
Plasma	1. Separate the plasma immediately after the sample is collected. 2. Freeze at < -18 °C (< -1 °F).	>24 hours

* Sample stability.

NOTE: Never freeze whole-blood samples.

To thaw frozen plasma samples

1. Let the sample thaw at room temperature.
2. Gently invert the sample tube 5 times to mix the sample.
3. Centrifuge the thawed sample to remove possible micro clots.

NOTE: Do not freeze and thaw a sample more than once.

Maximum sample age

This is the period of time after a sample is collected within which whole-blood samples kept at an ambient temperature of 18-25 °C (65-77 °F) must be analyzed. See 'Sample handling and storage' for details.

Sample analysis

Samples must be analyzed immediately after they are collected. When this is not possible, samples must be handled and stored correctly, and analyzed within the time period given in 'Sample handling and storage'. See the *AQT90 FLEX instructions for use* document for detailed instructions.

NOTE: Keep the analyzer connected to the mains power supply so that it is always ready for use.

Results

Results are shown on the screen and can be printed as they become available. Results are saved in a data log. If the analyzer detects an error condition, a message is attached to the analysis and no result is given.

WARNING – Risk of incorrect results

Always interpret with caution any assay results that are inconsistent with other available clinical data and a patient's medical history. If samples contain heterophilic antibodies, immunoassays could give falsely elevated or depressed results. Blockers have been added to this assay to minimize potential interference from heterophilic antibodies.

Calibration adjustment

See the *AQT90 FLEX instructions for use* document for details.

Results are saved in a data log and can be printed. If a calibration adjustment fails, a message is shown.

NOTE: After every calibration adjustment, liquid quality control (LQC) measurements must be done.

Quality control

Liquid quality control (LQC) measurements must be done:

- after a calibration adjustment
- in accordance with local, national and/or international regulations or accreditation requirements

See the *AQT90 FLEX instructions for use* document for details.

LQC material

See the *AQT90 FLEX Ordering information* document for details.

Performance characteristics

Analytical specificity

The cross-reactivity with luteinizing hormone, thyroid stimulating hormone, follicle stimulating hormone, placental lactogen hormone and growth hormone at concentrations of 1000 IU/L, 175 IU/L, 1000 IU/L, 200 µg/L and 200 µg/L, respectively, was studied to determine the analytical specificity of the β hCG assay. The cross-reactivity was shown to be $\leq 1\%$ for thyroid stimulating hormone and $\leq 0.5\%$ for the other tested substances.

Analytical sensitivity and reportable range

Analytical sensitivity	β hCG, IU/L
Limit of blank (LoB)	<1
Limit of detection (LoD)	<2

The reportable range of the assay is 2-5,000 IU/L.

Reference values

EDTA whole blood and plasma were obtained from the following patient groups and analyzed with the β hCG assay.

Patient group	N*	β hCG, IU/L 95 % percentile
Male	129	<0.5
Female (all)	125	<6
Post-menopausal	38	<6
Pre-menopausal	87	<2

*N = number of patients

NOTE: These values should only be used as an example. Each laboratory should establish its own reference ranges.

Low level β hCG values ≥ 25 IU/L may be indicative of early pregnancy but the results should always be evaluated in the context of the clinical situation [7]. When borderline values are reported, a new patient sample should be collected 48 hours later.

Imprecision

Spiked plasma pools were analyzed over 20 days, twice a day, 4 replicates per run to determine within-run and within-lab imprecision.

Pool	β hCG*, IU/L	CV** (%)	
		Within-run	Within-lab
1	3.1	15.4	17.2
2	13.0	5.2	8.6
3	4173	2.3	3.3

* Mean concentration; ** Coefficient of variation.

Hook effect

No hook effect was found when β hCG concentrations up to 233,000 IU/L were measured.

Carry over

Carry over from a sample with high β hCG value ($\sim 100,000$ IU/L) to an adjacent negative sample was determined to be < 50 ppm (< 5 IU/L).

Interfering substances

Hemolytic, lipemic and icteric samples were shown not to interfere with the assay.

Samples with β hCG at two concentration levels (approximately 5 IU/L and 50 IU/L) were used in the interference tests for biotin. Samples with β hCG at a concentration level of approximately 13 IU/L were used in the interference tests for all other substances. The samples were spiked with interfering substances at concentrations above the upper therapeutic level.

The following interfering substances were found to have no notable effect on the β hCG assay (interference ≤ 20 %) at the concentrations listed in the table.

Substance	Concentration mg/L
Abciximab	200
Acetaminophen	200
Acetylcysteine	2500
Acetylsalicylic acid	2000
Allopurinol	100

Substance	Concentration mg/L
Ambroxol	75
Ampicillin	100
Ascorbic acid	100
Atenolol	10
Biotin	2.6 (2,600 ng/mL)
Caffeine	100
Captopril	50
Cefoxitin	1100
Cinnarizine	30
Cocaine	10
Cyclosporine	5
Diclofenac	50
Digoxin	0.010
Dopamine	130
Erythromycin	200
Ethanol	5 % (v/v)
Furosemide	150
Heparin low molecular weight	5000 IU/L
Heparin Sodium	8000 IU/L
Ibuprofen	500
Levodopa	20
Methyldopa	40
Metronidazole	200
Nicotine ($\pm$)	20
Nifedipine	60
Nitrofurantoin	10
Nitroglycerin	0.16
Nystatin	7
Oxytetracycline	100
Phenylbutazone	400
Phenytoin	100
Propranolol	5

Substance	Concentration mg/L
Quinidine	30
Rifampicin	150
Tetracycline	50
Theophylline	100
Trimethoprim	100
Verapamil	160
Warfarin	15

Recovery

Recovery throughout the reportable range was evaluated using a negative whole blood sample spiked with known amounts of the WHO 4th IS (75/589). The results are averages of 10 measurements.

Concentration IU/L		
On AQT90 FLEX analyzer	Expected	Bias (%)
12.67	12.63	-0.30
556.83	627.97	11.30
1123.57	1248.26	9.99
1802.64	1860.90	3.13
2444.87	2465.86	0.85
3125.17	3063.16	-2.00
3813.34	3652.78	-4.40
4366.27	4234.75	-3.10

Method comparison

The β hCG assay (y) was compared to the commercially available β hCG assay for Siemens ADVIA Centaur immunoassay system (x). Serum samples in the range 2-5,000 IU/L (as measured on the AQT90 FLEX analyzer) were used. A bridging study showed that serum samples gave the same results as plasma samples on the AQT90 FLEX analyzer.

The linear regression line and the correlation coefficient were calculated for the ranges 2-1,000 IU/L and 2-5,000 IU/L.

β hCG range, IU/L	Correlation equation	R ²	n
2-1,000	$y = 1.181x + 9.5$	0.991	79
2-5,000	$y = 1.228x - 20.7$	0.986	116

Patents

Radiometer products may be covered by one or more patents and patent applications, see www.radiometer.com/en/legal/patents.

Explanation of graphical symbols

See the AQT90 FLEX instructions for use document.

References

See 'References' on the last page of this insert.

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References

1. Braunstein GD, Rasor J, Adler D, Danzer H, Wade, ME. Serum human chorionic gonadotropin levels throughout normal pregnancy. *Am J Obstet Gynecol* 1976; 126: 678-81.
2. Brody S, Carlström G. Immunoassay of human chorionic gonadotropin in normal and pathologic pregnancy. *J Clin Endocrinol Metab* 1962; 22: 564-74.
3. Cole LA, Kardana A, Park S-Y, Braunstein GD. The deactivation of hCG by nicking and disassociation. *J Clin Endocrinol Metab* 1993; 76: 704-13.
4. Pastorfide GB, Goldstein DP, Kosasa TS, Levesue L. Serum chorionic gonadotropin activity after molar pregnancy, therapeutic abortion, and term delivery. *Am J Obstet Gynecol* 1974; 118: 293-94.
5. Reyes FI, Winter JSD, Faiman C. Postpartum disappearance of chorionic gonadotropin from the maternal and neonatal circulations. *Am J Obstet Gynecol* 1985; 153: 486-89.
6. Cole LA. Immunoassay of human chorionic gonadotropin, its free subunits, and metabolites. *Clin Chem* 1997; 43: 12: 223-43.
7. Burtis CA, Ashwood ER, Bruns DE. Clinical Chemistry of Pregnancy. In: Tietz textbook of clinical chemistry and molecular diagnostics, fourth edition. St. Louis: Elsevier Saunders, 2006.

Made in Finland for:
Radiometer Medical ApS
Åkandevej 21
2700 Brønshøj
Denmark
www.radiometer.com

Radiometer, the Radiometer logo, ABL, AQT, TCM, RADIANCE, AQUIRE, PICO, CLINITUBES and QUALICHECK are trademarks of or used under license by Radiometer Medical ApS.

© 2018 Radiometer Medical ApS. All rights reserved.

Code number: 994-700

Version: 201808P