
AQT90 FLEX

CRP Test Kit

Contents

en	CRP Test Kit.....	2
bg	Набор тестове за CRP.....	6
cs	Souprava testů CRP.....	11
da	Testkit f or CRP.....	16
de	CRP-Testkit.....	20
el	Κιτ εξέτασης CRP.....	25
es	Kit de Test de CRP.....	29
et	CRP mõõtekomplekt.....	34
fi	CRP-testaussarja.....	38
fr	Kit de tests CRP.....	42
hr	Testni komplet CRP.....	47
hu	CRP tesztkészlet.....	51
it	Kit per il test CRP.....	56
lt	CRP testo rinkinys.....	60
lv	CRP testēšanas komplekts.....	65
nl	CRP-testkit	69
no	CRP-testkit.....	74
pl	Zestaw testowy do oznaczania CRP.....	78
pt	Kit de Testes de CRP.....	83
ro	Kit de testare pentru CRP.....	88
ru	Набор тестов на CRP.....	92
sk	Súprava na testovanie CRP.....	97
sv	CRP-testkit.....	102
tr	CRP Test Kiti.....	106
zh	CRP 检测试剂盒.....	110
References		114

en

CRP Test Kit

REF 942-936

Intended use

For *in vitro* diagnostic use.

The CRP Test is an *in vitro* diagnostic assay for the quantitative determination of C-Reactive Protein in EDTA or lithium-heparin whole-blood or plasma specimens on the AQT90 FLEX analyzer in point-of-care and laboratory settings. It is intended for use as an aid in the detection and evaluation of infection, tissue injury, and inflammatory disorders and associated diseases.

The CRP CAL Cartridge is an *in vitro* diagnostic reagent intended for the calibration adjustment of the CRP Test on the AQT90 FLEX analyzer by establishing points of reference to estimate CRP values.

NOTE: The cartridges are ready for use.

Test

C-Reactive Protein (CRP) is an acute phase protein consisting of five identical polypeptide chains that form a five-membered ring with a molecular weight of 120,000 Daltons. CRP belongs to the pentraxin family of proteins and it is synthesized in the liver by hepatocytes. CRP has been shown to interact with the humoral and cellular effector system and to have an important role, for example in the elimination of foreign pathogens and potentially toxic substances produced as a result of tissue damage [1,2].

In healthy people, serum and plasma concentrations of CRP are normally below 5 mg/L and elevated levels are attributed to various tissue-damaging processes. Although the increase of CRP is non-specific, CRP measurements have proven to be useful in the detection and monitoring of for example, inflammatory conditions, infections and malignant diseases. CRP levels rise within 6-8 hours after an acute stimulus and may double approximately every 8 hours. The peak level, which in severe cases can be 100-1000 times the normal value, is usually reached within 50 hours [3,4,5].

Additionally, monitoring of CRP concentrations using ultrasensitive methods (hsCRP) has been reported to be prognostic in many conditions, including cardiovascular diseases [6].

Calibration adjustment

The assay has been factory-calibrated using a 8-point calibration. The factory-defined calibration data is provided in the form of a barcode. This data is used as a lot-specific reference curve, which is adjusted to the correct analyzer-specific signal level by using a CAL Cartridge with the same lot number. The results obtained using the CAL Cartridge are combined with the factory-defined calibration data to produce a lot- and analyzer-specific calibration curve for the assay.

When a CAL Cartridge with a new lot number is inserted, the analyzer requests the factory-defined calibration data so that it can do the calibration adjustment. If the user does not insert a CAL Cartridge when the Test Cartridge lot changes, the analyzer will request it as soon as it detects a Test Cartridge with a new lot number. Test Cartridges with a new lot number can only be used after a calibration adjustment is done with a CAL Cartridge with the same lot number.

See the AQT90 FLEX *instructions for use* document for details.

Principles of operation

The AQT90 FLEX analyzer is a fully automatic continuous-access analyzer utilizing immunoassay technology and time-resolved fluorometric detection. See the AQT90 FLEX *instructions for use* document for details.

The AQT90 FLEX analyzer is operated via a touch screen. The analyzer aspirates samples from capped sample tubes that are inserted manually.

The analyzer does all assay steps automatically and reports quantitative results. See the AQT90 FLEX *instructions for use* document for details.

Contents

- 10 CRP Test Cartridges for 80 tests
- One CRP CAL Cartridge for one calibration adjustment
- A factory-defined calibration data sheet with barcode

Components

- Each Test Cartridge contains 8 analyte-specific assay cups and 8 dilution cups.
- The CAL Cartridge contains 8 analyte-specific background cups and 8 cups with added antigen.

Storage

Store unopened at 2-8 °C (36-46 °F).

Shelf life

If the protective pouch containing the cartridge remains sealed until just before the cartridge is inserted into the analyzer, it is stable until the end of the month stated in the expiry date shown on the pouch label.

On-board stability

The value (in days) is printed on the cartridge pouches alongside the following symbol:

Chemical composition

Reactive ingredients:

- Mouse monoclonal CRP capturing and tracing antibodies: 800 ng/cup
- Native CRP (only in the 8 CAL Cartridge cups with added antigen): 5 ng/cup

Other ingredients:

- Bovine serum albumin
- Bovine γ -globulin
- Mouse IgG (blocker substances for heterophilic interferences)

Main buffer component: TRIS (Tris hydroxymethylaminomethane)

Toxic ingredient: Sodium azide <0.1 %

Data in cartridge barcodes

- Checksum
- Barcode type
- Parameter code
- Lot no.
- Cartridge ID
- Expiry date
- Barcode version
- On-board stability

NOTE: The barcodes are lot-specific.

About using the cartridges

WARNING – Risk of incorrect results

Cartridges are sensitive to humidity and must therefore not be taken out of their sealed pouches until just before use. Only use cartridges that have been kept according to storage specifications and are not damaged.

WARNING – Risk of incorrect results

Do not use a cartridge that has already been used on another analyzer, otherwise you risk obtaining incorrect results.

WARNING – Risk of infection

This product contains materials from human source and should be handled and disposed of as potentially infectious. The plasma from human donors used in the manufacture of this product has been tested and found non-reactive to HbsAg, HIV-1, HIV-2 and HCV.

Product calibrator traceability

Calibration of the CRP assay is traceable to ERM - DA472/IFCC.

Sample collection and preparation

Sample tubes approved for use with the analyzer

See the *AQT90 FLEX instructions for use* document for details.

Anticoagulants

Use either EDTA or lithium heparin.

Sample size

The minimum sample volume is 2 mL (independent of the number of tests requested).

Sample type

Use whole-blood or plasma samples.

Sample collection

Collect blood samples by venipuncture.

WARNING – Risk of incorrect results

For sample collection, do not use tubes that contain a gel.

Sample handling and storage

Samples that cannot be analyzed immediately after collection must be handled and stored correctly before they are analyzed. The table provides an overview.

Sample type	Handling and storage	Analyze within...
Whole blood	1. Keep at an ambient storage temperature 18-25 °C (65-77 °F).	3 hours ^a
Plasma	1. Separate the plasma immediately after the sample is collected. 2. Keep at an ambient storage temperature 18-25 °C (65-77 °F).	3 hours ^a
Plasma	1. Separate the plasma immediately after the sample is collected. 2. Refrigerate at 2-8 °C (36-46 °F).	24 hours ^a
Plasma	1. Separate the plasma immediately after the sample is collected. 2. Freeze at -18 °C (<-1 °F).	>24 hours

^a Sample stability.

NOTE: Never freeze whole-blood samples.

To thaw frozen plasma samples

1. Let the sample thaw at room temperature.
2. Gently invert the sample tube 5 times to mix the sample.
3. Centrifuge the thawed sample to remove possible micro clots.

NOTE: Do not freeze and thaw a sample more than once.

Maximum sample age

This is the period of time after a sample is collected within which whole-blood samples kept at an ambient temperature of 18-25 °C (65-77 °F) must be analyzed.

See 'Sample handling and storage' for details.

Sample analysis

Samples must be analyzed immediately after they are collected. When this is not possible, samples must be handled and stored correctly, and analyzed within the time period given in 'Sample handling and storage'.

See the *AQT90 FLEX instructions for use* document for detailed instructions.

NOTE: Keep the analyzer connected to the mains power supply so that it is always ready for use.

Results

Results are shown on the screen and can be printed as they become available. Results are saved in a data log. If the analyzer detects an error condition, a message is attached to the analysis and no result is given.

WARNING – Risk of incorrect results

Always interpret with caution any assay results that are inconsistent with other available clinical data and a patient's medical history. If samples contain heterophilic antibodies, immunoassays could give falsely elevated or depressed results. Blockers have been added to this assay to minimize potential interference from heterophilic antibodies.

Calibration adjustment

See the *AQT90 FLEX instructions for use* document for details.

Results are saved in a data log and can be printed. If a calibration adjustment fails, a message is shown.

NOTE: After every calibration adjustment, liquid quality control (LQC) measurements must be done.

Quality control

Liquid quality control (LQC) measurements must be done:

- after a calibration adjustment
- in accordance with local, national and/or international regulations or accreditation requirements

See the *AQT90 FLEX instructions for use* document for details.

LQC material

See the *AQT90 FLEX Ordering information* document for details.

Performance characteristics

Analytical sensitivity and reportable range

Analytical sensitivity	CRP mg/L
Limit of blank (LoB)	<0.5
Limit of detection (LoD)	<1

The reportable range of the assay is 5-500 mg/L.

Reference values

EDTA whole blood and plasma were obtained from 272 apparently healthy individuals (143 women and 129 men) and analyzed using the CRP assay. The 95th percentile for both whole-blood and plasma samples was determined to be <5 mg/L.

NOTE: This value should only be used as an example. Each laboratory should establish its own reference ranges.

Imprecision

Three levels of antigen-spiked plasma pools were analyzed to determine the within-run and total imprecision on plasma. Each level was analyzed in duplicate in 2 separate runs for 20 test days. The total precision is a combination of repeatability, between-run and between-day precision.

Plasma	CRP ^a mg/L	CV ^b (%)	
		Within-run	Total ^c
Pool 1	7.4	4.9	5.3
Pool 2	111	5.5	5.6
Pool 3	343	4.5	4.8

^a Mean concentration; ^b Coefficient of variation; ^c Includes within-run, between-run and between-day contributions

A total of 150-180 spiked whole-blood samples were analyzed to determine the within-run and total imprecision on whole-blood. Three Test Kit lots and 5-6 analyzers were used.

Whole blood	CRP ^a mg/L	CV ^b (%)	
		Within-run	Total ^c
Pool 1	6.5	8.0	12.8
Pool 2	99	7.3	11.1
Pool 3	439	7.0	12.8

^a Mean concentration; ^b Coefficient of variation; ^c Includes within-run, between-run, between-batch and between-day (found in plasma precision study) contributions

Hook effect

No hook effect was found when CRP concentrations up to 2,000 mg/L were measured.

Carry over

Carry over from a sample with a high CRP value (~2,000 mg/L) to an adjacent very low sample value (<1 mg/L) was determined to be <100 ppm (<0.2 mg/L).

Interfering substances

Blockers have been added to this assay to minimize potential interference from heterophilic antibodies, for example HAMA (human anti-mouse antibodies) and rheumatoid factor (RF). Tests showed that the blockers had the desired effect.

Hemolytic, lipemic and icteric samples were shown not to interfere with the assay.

Samples with CRP at two concentration levels (approximately 9 mg/L and 37 mg/L) were used in the interference tests for biotin. Samples with CRP at a concentration level of approximately 85 mg/L were used in the interference tests for all other substances.

The following interfering substances were found to have no notable effect on the CRP assay (interference $\leq 20\%$) at the concentrations listed in the table.

Substance	Concentration mg/L
Abciximab	200
Acetaminophen	200
Acetylcysteine	2,500
Acetylsalicylic acid	2,000
Allopurinol	100
Ambroxol	75
Ampicillin	100
Ascorbic acid	100
Atenolol	10
Biotin	2.6 (2,600 ng/mL)
Caffeine	100
Captopril	50
Cefoxitin	1,100
Cinnarizine	30
Cocaine	10
Cyclosporine	5
Diclofenac	50

Substance	Concentration mg/L
Digoxin	0.010
Dopamine	130
Erythromycin	200
Ethanol	5 % v/v
Furosemide	150
Heparin low molecular weight	5,000 IU/L
Heparin Sodium	8,000 IU/L
Ibuprofen	500
Levodopa	20
Methyldopa	40
Metronidazole	200
Nicotine ($\pm$)	20
Nifedipine	60
Nitrofurantoin	10
Nitroglycerin	0.16
Nystatin	7
Oxytetracycline	100
Phenylbutazone	400
Phenytoin	100
Propranolol	5
Quinidine	30
Rifampicin	150
Tetracycline	50
Theophylline	100
Trimethoprim	100
Verapamil	160
Warfarin	15

Method comparison

The CRP assay (y) was compared to a commercially available CRP assay for the c501 module of the Roche Cobas 6000 (x). EDTA whole-blood samples were analyzed on the AQT90 FLEX analyzer and lithium-heparin plasma samples from the same sample collections were analyzed on the Roche Cobas 6000.

The concentrations measured on the AQT90 FLEX analyzer were in the range 5-500 mg/L.

The linear regression line and correlation coefficient were found to be:

$$y = 0.88x + 11.6; R^2 = 0.98 (n = 110)$$

Patents

Radiometer products may be covered by one or more patents and patent applications, see www.radiometer.com/en/legal/patents.

Explanation of graphical symbols

See the AQT90 FLEX instructions for use document.

References

See 'References' on the last page of this insert.

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Test

C-Reactive Protein (CRP) is an acute phase protein consisting of five identical polypeptide chains that form a five-membered ring with a molecular weight of 120 000 Daltons. CRP belongs to the pentraxin family of proteins and it is synthesized in the liver by hepatocytes. CRP has been shown to interact with the humoral and cellular effector system and to have an important role, for example in the elimination of foreign pathogens and potentially toxic substances produced as a result of tissue damage [1,2].

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References

1. Kilpatrick JM, Volanakis JE. Molecular genetics, structure, and function of C-reactive protein. *Immunol Res* 1991; 10: 43-53.
2. Kolb-Bachofen V. A review on the biological properties of C-reactive protein. *Immunobiol* 1991; 183: 133-45.
3. Dati F, Schumann G, Thomas L et al. Consensus of a group of professional societies and diagnostic companies on guidelines for interim reference ranges for 14 proteins in serum based on the standardization against the IFCC/BCR/CAP Reference Material (CRM 470). *Eur J Clin Chem Clin Biochem* 1996; 34(6): 517-20.
4. Pepys MB, Baltz ML. Acute phase proteins with special reference to C-reactive protein and related proteins (pentaxins) and serum amyloid A protein. *Adv Immunol* 1983; 34: 141-212.
5. Clyne B, Olshaker JS. The C-reactive protein. *J Emerg Med* 1999; 17(6): 1019-25.
6. Bassuk SS, Rifai N, Ridker PM: High-sensitivity C-reactive protein: clinical importance. *Curr Probl Cardiol* 2004; 29(8): 439-93.

Made in Finland for:
Radiometer Medical ApS
Åkandevej 21
2700 Brønshøj
Denmark
www.radiometer.com

Radiometer, the Radiometer logo, ABL, AQT, TCM, RADIANCE, AQUIRE, PICO, CLINITUBES and QUALICHECK are trademarks of or used under license by Radiometer Medical ApS.

© 2018 Radiometer Medical ApS. All rights reserved.

Code number: 994-745

Version: 201808P