	REF	C €	IVD
Level 1 Kit	944-198	CC9111A	6 × 0,65 mL
Level 2 Kit	944-199	CC912A	6 × 0,65 mL

EN INTENDED USE
For *in vitro* diagnostic use.
Multi-CHECK Combi is intended for use as an assayed quality control material to monitor the precision of laboratory testing procedures for the analytes listed in the specifications insert.

PARAMETERS
C-Reactive Protein, D-Dimer, Human Chorionic Gonadotropin, Procalcitonin
Contents: ● 6 tubes with caps ● Specifications insert
Components: Each tube contains 0.65 mL of control solution, enough for 5 tests.
Composition: This product is prepared from human serum to which purified biochemical material (extracts of human origin), chemicals, drugs/ preservatives, and stabilizers have been added.
Storage: Store the unopened product at –20 °C (–4 °F) or colder.
Shelf life: When stored correctly the unopened product is stable until the expiry date shown on the label.
● The thawed product is stable for 2 hours (maximum), when kept unused at ambient temperature; maximum 30°C (86 °F)
● Analytes are stable for 4 days under these conditions with the exception of PCT (Procalcitonin) which will be stable for 1 day when kept unused at 2–8 °C (36–46 °F)
NOTE: Do not freezeze.

Limitations of use:
● Product stability is dependent on its storage temperature, detailed specifications are provided for: “Storage”, “Shelf Life” and “Shelf life of the thawed product”. Do not use the product if it has not been kept according to these specifications, or if it has been damaged.
● If there is evidence of microbial contamination or excessive turbidity in the product, discard the tube.
● This product is not intended for use as a standard or calibrator.
Data in both barcodes on the specifications insert contains:
● Checksum ● Level name ● Barcode type
● Lot no. ● Area code ● Expiry date
Extra data in the barcodes:
Upper and lower limits of the assigned control range.
Assigned value
The assigned value is the center value of the assigned control range. The values are lot-specific. Values are provided as guidelines, each laboratory should establish its own statistical limits.
Quality control
NOTE: LQC material must be registered on the analyzer before use by using the barcodes provided on the Specifications insert.
LQC measurements must be performed:
● after a calibration adjustment
● in accordance with local, national and/or international regulations or accreditation requirements
See the *AQT90 FLEX instructions for use* document.
This product is intended for near-patient use. It is to be used by trained professionals.

Preparation before use
1. Remove the LQC solution from the freezer.
2. Let material thaw completely at ambient temperature of 18–25 °C (65–77 °F).
3. Analyze it within the time period specified in the “Shelf life of the thawed product” section above.
NOTE: Do not freezeze.
WARNING – Risk of infection
Handle and discard used LQC tubes as biohazardous waste. Failure to do so may result in exposure to potentially infected material. Disposal of all waste should be in accordance with local regulations. Material Safety Data Sheet (MSDS) is available on request.
WARNING – Potentially infectious material
Each human donor unit used to manufacture this LQC was tested using FDA-accepted methods and found non-reactive for Hepatitis B Surface Antigen (HBsAg), antibody to Hepatitis C (HCV) and antibody to HIV-1/ HIV-2. This product may also contain other human source material for which there are no approved tests. In accordance with good laboratory practice, all human source material should be considered potentially infectious and handled with the same precautions used with patient specimens.

!

WARNING – Risk of infection

Handle and discard used LQC tubes as biohazardous waste. Failure to do so may result in exposure to potentially infected material. Disposal of all waste should be in accordance with local regulations. Material Safety Data Sheet (MSDS) is available on request.

WARNING – Potentially infectious material

Each human donor unit used to manufacture this LQC was tested using FDA-accepted methods and found non-reactive for Hepatitis B Surface Antigen (HBsAg), antibody to Hepatitis C (HCV) and antibody to HIV-1/ HIV-2. This product may also contain other human source material for which there are no approved tests. In accordance with good laboratory practice, all human source material should be considered potentially infectious and handled with the same precautions used with patient specimens.

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