
AQT90 FLEX

NT-proBNP CAL Cartridge

Contents

en	NT-proBNP CAL Cartridge.....	2
bg	Калибрираща касета NT-proBNP.....	5
cs	Zásobník testů pro kalibraci NT-proBNP.....	8
da	NT-proBNP KAL-kassette.....	11
de	NT-proBNP-KAL-Kassette.....	14
el	Κασέτα βαθμονόμησης NT-proBNP.....	17
es	Cartucho CAL de NT-proBNP.....	20
et	NT-proBNP CAL-kassett.....	23
fi	NT-proBNP-CAL-kasetti.....	26
fr	Cartouche de CAL NT-proBNP.....	29
hr	NT-proBNP kalibracijska kaseta.....	32
hu	NT-proBNP CAL kazetta.....	35
it	Cartuccia CAL NT-proBNP.....	38
ja	NT-proBNP CAL カートリッジ.....	41
lt	NT-proBNP CAL kasetė.....	44
lv	NT-proBNP CAL kārtidžs.....	47
nl	NT-proBNP CAL-cassette.....	50
no	NT-proBNP KAL-kassett.....	53
pl	Kaseta kalibracyjna NT-proBNP CAL.....	56
pt	Cartucho CAL NT-proBNP.....	59
ro	Cartușul CAL pentru NT-proBNP.....	62
ru	CAL-картридж для теста на NT-proBNP.....	65
sk	Kalibračná kazeta NT-proBNP.....	68
tr	NT-proBNP Kalibrasyon Kartuşu.....	71
sv	CAL-kassett för NT-proBNP.....	74
zh	NT-proBNP 定标卡.....	77

en

NT-proBNP CAL Cartridge

REF 944-258

Intended use

For *in vitro* diagnostic use.

The NT-proBNP CAL Cartridge is an *in vitro* diagnostic reagent intended for the calibration adjustment of the NT-proBNP Test on the AQT90 FLEX analyzer by establishing points of reference to estimate NT-proBNP values.

NOTE: The cartridge is ready for use.

Calibration adjustment

The assay has been factory-calibrated using a 9-point calibration. The factory-defined calibration data is provided in the form of a barcode. This data is used as a lot-specific reference curve, which is adjusted to the correct analyzer-specific signal level by using a CAL Cartridge with the same lot number. The results obtained using the CAL Cartridge are combined with the factory-defined calibration data to produce a lot- and analyzer-specific calibration curve for the assay. When a CAL Cartridge with a new lot number is inserted, the analyzer requests the factory-defined calibration data so that it can do the calibration adjustment. If the user does not insert a CAL Cartridge when the Test Cartridge lot changes, the analyzer will request it as soon as it detects a Test Cartridge with a new lot number. Test Cartridges with a new lot number can only be used after a calibration adjustment is done with a CAL Cartridge with the same lot number.

See the AQT90 FLEX *instructions for use* document for details.

Contents

- One NT-proBNP CAL Cartridge for one calibration adjustment
- A factory-defined calibration data sheet with barcode

Components

The CAL Cartridge contains 8 analyte-specific background cups and 8 cups with added antigen.

Storage

Store unopened at 2-8 °C (36-46 °F).

Shelf life

If the protective pouch containing the cartridge remains sealed until just before the cartridge is inserted into the analyzer, it is stable until the end of the month stated in the expiry date shown on the pouch label.

On-board stability

The value (in days) is printed on the cartridge pouches alongside the following symbol:

Chemical composition

Reactive ingredients:

- Sheep polyclonal anti-NT-proBNP capturing and tracing antibodies: approx. 150 ng/cup
- Synthetic NT-proBNP (only in the 8 CAL Cartridge cups with added antigen): approx. 50 pg/cup

Other ingredients:

- Bovine serum albumin
- Bovine γ -globulin
- Mouse IgG and sheep IgG (blocker substances for heterophilic interferences)
- Human serum

Main buffer component: TRIS (Tris hydroxymethylaminomethane)

Toxic ingredient: Sodium azide <0.1 %

Data in the cartridge barcode

- Checksum
- Barcode type
- Parameter code
- Lot no.
- Cartridge ID
- Expiry date
- Barcode version
- On-board stability

NOTE: The barcode is lot-specific.

About using the cartridge

WARNING – Risk of incorrect results

The cartridge is sensitive to humidity and must therefore not be taken out of its sealed pouch until just before use. Only use cartridges that have been kept according to storage specifications and are not damaged.

WARNING – Risk of infection

This product contains materials from human source and should be handled and disposed of as potentially infectious. The plasma from human donors used in the manufacture of this product has been tested and found non-reactive to HbsAg, HIV-1, HIV-2 and HCV.

Calibration adjustment

See the *AQT90 FLEX instructions for use* document for details.

Results are saved in a data log and can be printed. If a calibration adjustment fails, a message is shown.

NOTE: After every calibration adjustment, liquid quality control (LQC) measurements must be done.

Patents

Manufactured under licence from Roche Diagnostics. U.S. Patent Nos.US 5,786,163. Additional U.S. Patent pending.

Radiometer products may be covered by one or more patents and patent applications, see www.radiometer.com/en/legal/patents.

Explanation of graphical symbols

See the *AQT90 FLEX instructions for use* document.

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Radiometer, the Radiometer logo, ABL, AQT, TCM, RADIANCE, AQUIRE, PICO, CLINITUBES and QUALICHECK are trademarks of or used under license by Radiometer Medical ApS.
© 2018 Radiometer Medical ApS. All rights reserved.