	[REF]	CC CE	[VD]
Trilevel Kit	944-512	CC510A	3 × 2 × 0.65 mL
Level 1 Kit	944-513	CC511A	6 × 0.65 mL
Level 2 Kit	944-514	CC512A	6 × 0.65 mL
Level 3 Kit	944-515	CC513A	6 × 0.65 mL

EN INTENDED USE

For *in vitro* diagnostic use. For use with the AQT90 FLEX analyzer as an assayed liquid quality control (LQC) serum to monitor the precision of laboratory testing procedures for the parameters listed in the Specifications insert.

PARAMETERS
CKMB, Myoglobin, Troponin I, Troponin T, NT-proBNP
Contents:

- 6 tubes with caps
- Specifications insert

Components:
Each tube contains 0.65 mL of control solution, enough for 5 tests.
Composition:
This product is prepared from human serum to which purified biochemical material (extracts of human origin), chemicals, drugs/preservatives, and stabilizers have been added.

Storage: Store the unopened product at –20 °C (–4 °F) or colder.
Shelf life: When stored correctly the unopened product is stable until the expiry date shown on the label.

Shelf life of the thawed product:

- The thawed product is stable for 2 hours (maximum), when kept unused at ambient temperature; maximum 32 °C (89 °F)
- The thawed product is stable for 4 days (maximum), when kept unused at 2–8 °C (36–46 °F)

NOTE: Do not refreeze.

Limitations of use:

- Product stability is dependent on its storage temperature, detailed specifications are provided for: “Storage”, “Shelf Life” and “Shelf life of the thawed product”. Do not use the product if it has not been kept according to these specifications, or if it has been damaged.
- If there is evidence of microbial contamination or excessive turbidity in the product, discard the tube.
- This product is not intended for use as a standard or calibrator.

Data in both barcodes on the specifications insert contains:

- Checksum
- Level name
- Barcode type
- Lot no.
- Area code
- Expiry date

Extra data in the barcodes:
Upper and lower limits of the assigned control range.

Assigned value

The assigned value is the center value of the assigned control range. The values are lot-specific. Values are provided as guidelines, each laboratory should establish its own statistical limits.

Quality control

NOTE: LQC material must be registered on the analyzer before use by using the barcodes provided on the Specifications insert.

LQC measurements must be performed:

- after a calibration adjustment
- in accordance with local, national and/or international regulations or accreditation requirements

See the AQT90 FLEX *instructions for use* document.

Preparation before use

- Remove the LQC solution from the freezer.
- Let material thaw completely at ambient temperature of 18–25 °C (65–77 °F). Thawing takes up to 30 minutes (maximum).
- Analyze it within the time period specified in the “Shelf life of the thawed product” section above.

NOTE: Do not refreeze.

WARNING – Risk of infection
Handle and discard used LQC tubes as biohazardous waste. Failure to do so may result in exposure to potentially infected material. Disposal of all waste should be in accordance with local regulations. Material Safety Data Sheet (MSDS) is available on request.

WARNING – Potentially infectious material
Each human donor unit used to manufacture this LQC was tested using FDA-accepted methods and found non-reactive for Hepatitis B Surface Antigen (HBsAg), antibody to Hepatitis C (HCV) and antibody to HIV-1/HIV-2. This product may also contain other human source material for which there are no approved tests. In accordance with good laboratory practice, all human source material should be considered potentially infectious and handled with the same precautions used with patient specimens.

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