


Instructions for use

Ambu® Oval Plus Silicone Resuscitator

5

O ₂ flow (l/min)	V _T (ml) x f (pr. min.), I:E ratio = 1:2			
	250 x 12	600 x 12	750 x 12	1000 x 12
2	74	43	38	34
5	100	76	65	54
10	100	100	100	87
15	100	100	100	100

5.1

O ₂ flow (l/min)	V _T (ml) x f (pr. min.), I:E ratio = 1:2			
	40 x 40	100 x 20	200 x 20	400 x 15
1	70	60	40	34
2	100	100	60	47
4	100	100	100	73
6	100	100	100	100


5.2

O ₂ flow (l/min)	V _T (ml) x f (pr. min.), I:E ratio = 1:2				
	87%	70%	61%	42%	47%
1	100 %	100 %	100 %	63%	74%
2	100 %	100 %	100 %	100 %	100 %
4	100 %	100 %	100 %	100 %	100 %
6	100 %	100 %	100 %	100 %	100 %
8	100 %	100 %	100 %	100 %	100 %
15	100 %	100 %	100 %	100 %	100 %


5.3

6


6.1

6.2


6.3


7

7.1


7.2


8

8.1

8.2


8.3

8.4

8.5

8.6

8.7

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0086

CE mark. The product complies with the EU Council directive concerning Medical Devices 93/42/EEC

Ambu is a registered trademark of Ambu A/S, Denmark.
Ambu is certified according to ISO 13485

English Instructions for use8-15	Français Mode d'emploi.....80-87	Polski Instrukcja obsługi.....152-159
Български Указания за ползване16-23	Hrvatski Upute za upotrebu.....88-95	Português Manual de instruções.....160-167
Česky Návod k použití24-31	Magyar Használati útmutató.....96-103	Română Instrucțiuni de utilizare168-175
Dansk Brugsanvisning.....32-39	Italiano Manuale d'uso104-111	Русский Инструкция по применению.....176-183
Deutsch Bedienungsanleitung40-47	日本語 使用法112-119	Slovenčina Návod na použitie184-191
Ελληνικά δηγίες Χρήσεως48-55	Lietuviškai Naudojimo instrukcijos120-127	Slovenščina Navodila za uporabo192-199
Español Manual de instrucciones.....56-63	Latviski Lietošanas instrukcija.....128-135	Svenska Instruktionshandbok.....200-207
Eesti Kasutusjuhised64-71	Nederlands Gebruiksaanwijzing.....136-143	Türkçe Kullanım talimatları.....208-215
Suomi Käyttöohje.....72-79	Norsk Brukerveiledning144-151	中文 使用说明215-221

1. Intended use

The Ambu® Oval Plus Silicone resuscitator is a reusable resuscitator intended for pulmonary resuscitation.

The range of application for each version is:

- Adult: Adults and children with a body weight more than 30 kg (66 lbs).
- Paediatric: Infants and Children with a body weights up to 10 ~ 30 kg (22 ~ 66lbs).
- Neonate: Neonates and infant with a body weight up to 10 kg (22 lbs).

2. Warning and caution statements

Failure to observe these precautions may result in inefficient ventilation of the patient or damage to the equipment.

WARNING

Oil or grease should not be used in close proximity to oxygen equipment.
Do not smoke or use open flame when oxygen is in use – fire may result.

Never override the pressure-limiting valve (if available) unless medical and professional assessment indicates the necessity. High ventilation pressures may cause lung rupture to certain patients. If the pressure-limiting valve is overridden in patients with a bodyweight less than 10 kg (22 lbs.), a manometer must be used to monitor the ventilation pressure to avoid the possibility of a lung rupture.

By adding accessories, it may increase inspiratory and/or expiratory resistance.
Do not attach accessories if increased breathing resistance would be detrimental for the patient.

CAUTION

US federal law restricts this device to sale by or on the order of a physician (USA and Canada only)

For use by trained personnel only. The proper application of a facemask to obtain tight seal should be trained in particular. Make sure that the personnel are made familiar with the content of this manual.

Always inspect the resuscitator and perform a functional test after unpacking, cleaning, assembly and prior to use.

If the resuscitator with attachments is placed on standby for emergency use, it should be inspected at regular intervals to assure integrity.

Always watch the movement of the chest and listen for the expiratory flow from the valve in order to check the ventilation efficiency. Switch immediately to mouth-to-mouth ventilation if efficient ventilation cannot be obtained.

Insufficient, reduced, or no airflow may result in brain damage to the patient being ventilated.

Do not use the resuscitator in toxic or hazardous atmosphere.

Never store the resuscitator in a deformed state other than as folded when delivered by the manufacturer, otherwise permanent distortion of the bag will occur which may reduce the ventilation efficiency. The folding zone is clearly visible on the bag (only Adult and Paediatric versions may be folded).

3. Specifications

The Ambu Oval Plus Silicone resuscitator is in conformity with the product specific standard EN ISO 10651-4. The Ambu Oval Plus Silicone is in conformity with Council Directive 93/42/EEC concerning Medical Devices.

	Neonate	Paediatric	Adult
Resuscitator volume	approx. 260 ml	approx. 660 ml	approx. 1546 ml
Delivered volume one hand**	150 ml	450 ml	600 ml
Delivered volume two hands**	-	-	1100 ml
Dimensions without Reservoir (length x diameter)	approx. 198 x 75 mm	approx. 240 x 100 mm	approx. 269 x 130 mm
Weight without Reservoir and Mask	approx. 128 g	approx. 205 g	approx. 338 g
Pressure-limiting valve*	4.0 kPa (40 cm H ₂ O)	4.0 kPa (40 cm H ₂ O)	4.0 kPa (40 cm H ₂ O)
Dead space	≤ 5 ml + 10 % of the delivered volume	≤ 5 ml + 10 % of the delivered volume	≤ 5 ml + 10 % of the delivered volume
Inspiratory resistance	max 0,5 kPa (5,0 cm H ₂ O) at 5 l/min	max 0,5 kPa (5,0 cm H ₂ O) at 50 l/min	max 0,5 kPa (5,0 cm H ₂ O) at 50 l/min
Expiratory resistance	max 0,5 kPa (5,0 cm H ₂ O) at 5 l/min	max 0,5 kPa (5,0 cm H ₂ O) at 50 l/min	max 0,5 kPa (5,0 cm H ₂ O) at 50 l/min
Reservoir volume	approx. 1500 ml (bag) approx. 100 ml (tube)	approx. 1500 ml (bag)	approx. 1500 ml (bag)
Patient connector	Outside 22 mm male (ISO 5356-1) Inside 15 mm female (ISO 5356-1)		
Expiration connector (for PEEP valve attachment)	30 mm male (ISO 5356-1)		
Manometer Port connector	Ø 4,2 +/- 0,1 mm		
Bag inlet connector	-	Inside 32mm female (ISO10651-4)	Inside 32mm female (ISO10651-4)
Forward and backward leakage	Not measurable		
Operation temperature	-18 °C to +50 °C (-0,4 °F to +122 °F)		
Storage temperature	Tested at - 40 °C (-40 °F) and + 60 °C (+140 °F) according to EN ISO10651-4		
O ₂ inlet connector	according to EN 13544-2		
Long term storage	For long term storage the resuscitator should be kept in closed packaging in a cool place away from sunlight.		

* Higher delivery pressure can be obtained by overriding the pressure-limiting valve.

The Ambu Oval Plus Adult is also available without pressure-limiting valve.

** Tested according to EN ISO 10651-4.

4. Principle in operation ①

The illustration (1) shows the ventilation gas flow mixtures into the bag and to and from the patient during manual operation of the resuscitator. (a) Adult and paediatric resuscitator, (b) infant resuscitator with closed reservoir, (c) infant resuscitator with open reservoir.

The gas flow is similar when the patient is breathing spontaneously through the device. The O₂ reservoir assembly is fitted with two valves, one allowing ambient air to be drawn in when the reservoir is empty and one spilling out surplus oxygen when the reservoir bag is full.

1.1 Excess oxygen, 1.2 Air, 1.3 Oxygen inlet, 1.4 Patient, 1.5 Expiration, 1.6 Manometer port, 1.7 Pressure limiting valve,

5. Instructions for use

Preparation

- If the resuscitator is packed in a compressed state, unfold by pulling on the patient valve and the inlet valve.
- Fit the facemask and place all items in the plastic bag supplied with the resuscitator.
- The integrity of the kits issued for storage ready for use should be inspected at the interval established in the local protocol.
- Before use on the patient make a brief functional test as described in section 9.
- If connecting external devices to the resuscitator, make sure to test for functionality and consult the instructions for use accompanying the external device.

Patient use

- Clear the patient's mouth and airway using recommended techniques. Use recommended techniques to position the patient correctly to open the airway and to hold the mask firmly against the face. (2.1)
- Slide you hand (Adult Version) or ring and middle finger (Paediatric version) under the support strap. The infant version does not have a support strap. Ventilation without using the support strap can be achieved by turning the bag. (2.2)
- Ventilate the patient. During insufflation observe the rise of the patient's chest. Release the bag abruptly and listen for the expiratory flow from the patient valve and observe lowering of the chest.
- If continued resistance to insufflation is encountered, check the airway for obstruction or correct the backward tilt of the head.
- If the patient vomits during mask ventilation, immediately clear the patient's airway and then freely compress the bag a few times before resuming ventilation.
- The patient valve may be disassembles and cleaned if excessive amount of vomitus are obstructing the air flow. For details review the cleaning section of this manual.

5.2 Manometer port ③

WARNING

Use only for monitoring pressure. The cap must always be put on the connector when pressure is not being monitored.

A pressure gauge can be connected to the manometer port on the top of the patient valve. (This only applies to the version with manometer port). Remove the cap (3.1) and connect pressure manometer or the tube for the pressure gauge (3.2).

5.3 Pressure limitation system ④

WARNING

Never override the pressure-limiting valve (if available) unless medical and professional assessment indicates the necessity. High ventilation pressures may cause lung rupture to certain patients. If the pressure-limiting valve is overridden in patients with a bodyweight less than 10 kg (22 lbs.), a manometer must be used to monitor the ventilation pressure to avoid the possibility of a lung rupture.

If the resuscitator is equipped with a pressure limiting valve, the valve is set to open at 40 cm H₂O (4.0kPa) (4.1).

If medical and professional assessment indicates a pressure above 40 cm H₂O is required pressure limiting valve can be overridden by pressing the override cap onto the valve (4.2). Alternatively the pressure limiting valve can be overwritten by placing the index finger on the blue button while squeezing the bag.

6. Oxygen administration

Administer Oxygen according to medical indication.

Examples of O₂ percentages, which can be obtained with different volumes and frequencies, can be seen in ⑥. Adult (5.1), Pediatric (5.2), Infant (5.3).

VT: Ventilation volume, f: Frequency

Note: If high ventilation pressure is used, a higher O₂ flow settings is needed because part of the stroke volume is vented from the pressure-limiting valve. For infant version, use of supplementary Oxygen without reservoir attached will limit the oxygen concentration to 60-80 % at 15 L/min.

Oxygen enrichment of the inspiratory air in case of spontaneous respiration

- With oxygen reservoir bag or without an oxygen reservoir:
Attach a mask to the resuscitator, set the oxygen flow to between 2 and 4 l/min, and place the mask over the patient's mouth and nose, creating an airtight seal. If the spontaneous respiration is sufficient, an oxygen concentration close to **100 %** will be obtained.
- With an oxygen reservoir tube:
The end of the resuscitator with the oxygen reservoir that normally faces away from the patient can also be used for oxygen enrichment of the inspiratory air in case of spontaneous respiration. Set the oxygen flow should be set to between 2 and 4 l/min, and the free end of the oxygen reservoir should be held as tightly as possible against the patient's head.

The oxygen concentration will be very close to **100 %** at the end of the tube. Depending on the distance between the oxygen reservoir and the patient's head, there will be significantly more or less inhalation of atmospheric air, and thus lower oxygen concentrations should be expected.

7. Parts/Materials

Bag	Silicone rubber
Patient valve housing	Polysulphone
O-Ring (only patient valves with swivel)	Silicone rubber
Patient connector	Polysulphone
Expiratory connector	Polypropylene, reinforced
Outlet cap (option)	Polysulphone
Valve discs	Silicone rubber
Reservoir bag	Polyethersulfone/Polyurethane

Flange nut	Polypropylene, reinforced
Pressure limiting valve	Polysulphone/stainless steel
Override cap	Silicone rubber
Bag support connector	Polysulphone
Inlet valve cover	Polyoxymethylene
Inlet valve housing	Polyoxymethylene
Cap for manometer port	Polypropylene/EPDM rubber
Inlet valve cover (Neonate)	Polysulphone
Inlet valve housing (Neonate)	Polysulphone
Window with snap faster	Polysulphone
Reservoir housing	Polyoxymethylene/Polypropylene
Adapter for reservoir bag	Silicone rubber
Reservoir tube	Polyethylene

8. Cleaning-disinfection-sterilisation

Clean – disinfect – sterilize the complete resuscitator thoroughly after each patient use or if used in a patient/environment with infectious diseases.

8.1 Disassembly ⑥

The resuscitator must always be disassembled for cleaning as shown.
Adult (6.1), Paediatric (6.2), Neonate (6.3)

CAUTION

Do not try to disconnect the patient connector from the patient valve (permanent snap fit).
Do not disassemble parts further than shown.
The bag connector should not be detached from the bag during the cleaning and sterilisation procedure.

8.2 Cleaning procedure

The following cleaning method has been validated according to the AAMI TIR 12 standard. Ambu makes no warranties concerning other cleaning and sterilisation methods than the validated one described here. The validated method specifies cleaning by washing machine using detergents.

Disassemble the devices prior to placing in the washer and program the washer for the following parameters. Place devices into washer to allow for drainage:

Phase	Recirculation time (minutes)	Water temperature	Detergent type and concentration (if applicable)
Pre-Wash 1	02:00	Cold tap water	NA
Enzyme wash	02:00	Hot tap water	Enzymatic detergent
Wash 1	02:00	60.0 °C (140 °F) (Set point)	Enzymatic detergent
Rinse 1	0:15	Hot tap water (60.0 °C) (140 °F)	NA

Dry the devices appropriately.

Sterilisation by gravity steam autoclave running a full cycle at 134 degrees Celsius for 10 minutes at a standard pressure and then dry 30 minutes.

8.3 Parts that can be cleaned and sterilised

Applicable methods		
X Applicable O Not applicable	Machine Washing	Disinfection & sterilisation Autoclave 134°
Patient valve	X	X
Bag	X	X
Inlet valve	X	X
Extension tube	X	X
Silicone Facemask	X	X
Silicone Facemask size OA	X	O
Oxygen reservoir bag	X	X
Oxygen reservoir tube	X	O

8.4 Detergents and chemical disinfectants

Use only brands that are compatible with the resuscitator materials to avoid reduction in the lifetime of the materials. Follow the instructions of the manufacturer of the detergent or chemical disinfectant for dilution and exposure time.

CAUTION

Avoid using substances containing phenol to clean the product. Phenol will cause premature wear and degradation of the materials or reduced product life span. Promptly remove all residue of cleaning materials from the resuscitator. Residues may cause premature wear or reduced product life span. The following steps are generally recommended. Select proper methods for the resuscitator parts in question according to the table above.

8.5 Disinfecting and/or sterilising

Select heated or chemical disinfection according to local standards for disinfection and the table of applicable methods.

Washing machine - heat-disinfecting program:

Automatic washing machines designed for medical equipment will normally include program cycles for heat disinfection.

Autoclaving: Use standard autoclaving equipment adjusted for 134 °C (273 °F) for resuscitator parts according to the table above.

Sterilisation by gravity steam autoclave running a full cycle at 134 degrees Celsius for 10 minutes at a standard pressure and then dry 30 minutes.

For drying and cooling leave the parts to dry and/or cool completely before reassembling the resuscitator.

Ambu product testing has shown that the resuscitator is still fully functional following autoclaving 30 times. The actual number of times that the product can be autoclaved while maintaining full functionality will vary, and can be both higher and lower than 30 times, depending on product use, storage, and wear and tear. Always perform a test of function prior to each use.

Note: that it is the user's responsibility to qualify any deviations from the recommended method of processing

8.6 Inspection of parts

After cleaning, disinfecting and/or sterilising carefully inspect all parts for damage and residuals or excessive wear and replace if necessary. Some methods may cause discoloration of rubber parts without impact on their lifetime. In case of material deterioration e.g. cracking, the parts should be replaced.

When inserting the valve housing of the inlet valve: Make sure that the bag opening seats smoothly against the flange.

When inserting the patient valve disc: Make sure that the butt of the stem is pushed through the hole in the middle of the patient valve housing.

When mounting the outlet cap (option): Note that the slot of the outlet cap should face the patient connector.

8.7 Reassembly

Reassemble the parts of the resuscitator as shown in the illustration © Adult (6.1), Paediatric (6.2), Neonate (6.3)

CAUTION

For neonate version, the adapter for reservoir bag shall be mounted on the inlet connector, and the cap on the adapter must be put on the oxygen nipple before the reservoir bag can be connected to the resuscitator.

9. Test of function ⑦

Resuscitator

Close the pressure-limiting valve with the override cap (This only applies to the version with pressure limiting valve) and close the patient connector with the thumb (7.1). Briskly squeeze the bag. The resuscitator shall offer resistance to the squeeze.

Open the pressure-limiting valve by opening the override cap or by removing the finger and repeating the procedure. The pressure limiting valve should now be activated and it should be possible to hear the expiratory flow from the valve.

Squeeze and release the resuscitator a few times to ensure that air is moving through the valve system and out of the patient valve. (7.2)
Note: As the valve plates are moving during functional test or during ventilation a slight sound may appear. This does not compromise the functionality of the resuscitator.

Oxygen reservoir bag

Supply a gas flow of 5 l/min to the oxygen bag. Check that the reservoir fills. If not, check the integrity of the two valve shutters or for a torn reservoir.

Oxygen reservoir tube

Supply a gas flow of 10 l/min to the oxygen tube. Check that the oxygen flows out at the end of the reservoir tube. If not, check for a blocked oxygen tube.

10. Service

The resuscitator requires no scheduled maintenance other than regular cleaning, inspection and testing.

11. Accessories and Spare parts

Accessories ®

For further information please refer to the directions for use of the specific accessory

- Ambu reusable PEEP valve 10 (8.1)	A000 137 000
- Ambu reusable PEEP valve 20 (8.1)	A000 213 000
- Ambu disposable Pressure Manometer (8.2)	322 003 000
- Splash guard – Only Adult resuscitator (8.4)	245 000 707
- Extension tube – 30 cm (8.5)	209 000 701
- Extension tube – 90 cm (8.5)	209 000 702
- Extension tube – 105 cm (8.5)	209 000 703
- Adapter 28 mm (8.6)	209 000 506
- Adapter 24 mm (8.7)	209 000 507

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