22mm Breathing System

en	22MM BREATHING SYSTEM
fr	22MM CIRCUIT DE RESPIRATION
de	22MM BEATMUNGS-SYSTEM
es	22MM CIRCUITO RESPIRATORIO
pt	22MM SISTEMA RESPRATÓRIO
it	22MM CIRCUITO RESPIRATORIO
nl	22MM BEADEMINGSSYSTEEM
no	22MM RESPIRATOR SLANGESYSTEM
fi	22MM HENGITYSLETKUSTO
sv	22MM ANDNINGSSYSTEM
da	22MM RESPIRATIONSSYSTEM
el	22MM ΣΥΣΤΗΜΑ ΑΝΑΠΝΟΗΣ
lt	22MM KŲĖPAVIMO SISTEMA
pl	22MM SYSTEM ODDECHOWY
ru	22MM ДЫХАТЕЛЬНЫЙ КОНТУР
cs	22MM DÝCHACÍ SYSTÉM
hu	22MM LÉLEGEZTETŐ RENDSZER
sl	22MM DIHALNI SISTEM
lv	22MM ELPOŠANAS SISTĒMA
et	22MM HINGAMISAPARAAT
bg	22 MM РЕСПИРАТОРНА СИСТЕМА
sr	22MM SISTEM ZA DISANJE
ro	22MM SISTEM DE RESPIRAȚIE
sk	22MM D/SYSTÉM
hr	22MM SUSTAV ZA DISANJE
tr	22MM SOLUNUM DEVRESİ


Figure 1. / Figure 1. / Abbildung 1. / Figura 1. / Figura 1. / Figura 1. / Figur 1. / Figur 1. / Kuva 1. / Figur 1. / Figur 1. / Εικόνα 2. / 1 pav. / Rysunek 1. / Рисунок 1. / Obr. 1. / 1. ábra. / 1. attēls. / Joonis 1. / Фигура 1. / Slika 1. / Figura 1. / Obrázok 1. / Slika 1. / Şekil 1.


Figure 2. / Figure 2. / Abbildung 2. / Figura 2. / Figura 2. / Figura 2. / Figur 2. / Figur 2. / Kuva 2. / Figur 2. / Figur 2. / Εικόνα 2. / 2 pav. / Rysunek 2. / Рисунок 2. / Obr. 2. / 2. ábra. / 2. attēls. / Joonis 2. / Фигура 2. / Slika 2. / Figura 2. / Obrázok 2. / Slika 2. / Şekil 2.

C- Compliance R_I - Inspiratory limb resistance to flow R_E - Expiratory limb resistance to flow / C- Conformité R_I - Résistance de la branche inspiratoire à l'écoulement R_E - Résistance de la branche / C- Compliance R_I - Strömungswiderstand des Inspirationsschlauches R_E - Strömungswiderstand des / C- cumplimiento normativo; R_I - resistencia de la extremidad inspiratoria al flujo; R_E - resistencia al flujo / C- Complacência R_I - Resistência ao fluxo do ramo inspiratório R_E - Resistência ao fluxo do ramo / C- Conformità. R_I - resistenza della diramazione inspiratoria al flusso. R_E - resistenza della diramazione espiratoria al flusso / C- conformiteit R_I - Weerstand van de inspiratie-arm tegen doorstroming R_E - Weerstand van de expiratie-arm tegen doorstroming / C- Ettergivighet R_I - Flowmotstand i inspiratorisk gren R_E - Flowmotstand i expiratorisk gren / C- kompianssi R_I - sisäänhengityshaaran virtausvastus R_E - uloshengityshaaran virtausvastus / C- Efterlevnad R_I - Inandningsrörets flödesmotstånd R_E - Utandningsrörets flödesmotstånd / C- Overholdelse af R_I - Inspiratorisk slange flow modstand R_E - Ekspiratorisk slange flow modstand / C- Συμμόρφωση R_I - Αντίσταση εισπνευστικού σκέλους στη ροή R_E - Αντίσταση εκπνευστικού σκέλους στη ροή / C- atitiktis R_I - įkvėpimo galo atsparumas srautui R_E - iškvėpimo galo atsparumas srautui / C- podatnosť R_I - opór odnogi vdechowej na przepływ R_E - opór odnogi wydechowej na przepływ / C- Соответствие R_I - Сопротивление сегмента вдоха потока R_E - Сопротивление сегмента выдоха потока / C- soulad R_I - odpor proudění vdechovací části R_E - odpor proudění výdechové části / C- tágulékonyág, R_I - belégzési cső áramlási ellenállása, R_E - kilégzési cső áramlási ellenállása / Č- podajnosť R_I - upornost na pretok vdihovalnega kraka R_E - upornost na pretok izdihovalnega kraka / C- atbilstība R_I - ielpas daļas plūsmas pretestība R_E - izelpas daļas plūsmas pretestība / C- vastavus R_I - inspiratsiooniosa resistentsus voolule R_E - ekspirsatsiooniosa resistentsus voolule / C- Соответствие с R_I - Съпротивление на потока на разклонението за инспирация / R_E - Съпротивление на потока на разклонението за експирация / C- Usaglašenost R_I - otpornost na protok inspiratorne grane R_E - otpornost na protok ekspiratorne grane / C- Conform R_I - Rezistență la debit pentru tubul de inspirare R_E - Rezistență la debit pentru tubul de expirare / C- podajnosť R_I - prietokový odpor inspiračnej vetvy R_E - prietokový odpor expiračnej vetvy / C- sukladnost R_I - otpornost protoka u kraku za udisanje R_E - otpornost protoka u kraku za izdisanje / C- Uygunluk R_I - Nefes alma uzantisinin akış direnci R_E - Nefes verme uzantisinin akış direnci

Patient category ¹	Intended delivered volume ² ml	Flow resistance limit ³ cmH ₂ O/l/min	At flow ⁴ l/min
Adult	≥ 300 ml	0.06	30

Patient category ¹	Intended delivered volume ² ml	Compliance ⁵ ml/cmH ₂ O	At Pressure ⁶ cmH ₂ O
Adult	≥ 300 ml	5	60 ± 3

1 Patient category/ Catégorie de patients/ Patientenkategorie/ Categoría de paciente/ Categoría de doentes/ Categorie di pazienti/ Patiëntencategorie/ Pasientkategori/ Potilastyypit/ Patientkategori/ Patientkategori/ Κατηγορία ασθενών/ Pacientų kategorija/ Kategorija pacjentów/ Категорія пацієнтів/ Kategorije pacientů/ Betegek köre/ Kategorija bolnika/ Pacientu kategorija/ Patsiendi kategooria/ Категорія на пацієнта/ Kategorija pacijenata/ Categorie de pacienți/ Kategoría pacientov/ Kategorija pacijenata/ Hasta kategorisi

2 Intended delivered volume ml/ Volume administré prévu ml/ Zuzuführendes Volumen ml/ Volumen de suministro previsto ml/ Volume administrado previsto ml/ Volume erogato previsto ml/ Beoogd volume ml/ Tiltentkt leveret volum ml/ Aiottu antotilavuus ml/ Avsedd levererad volym ml/ Tilsigter leveret volumen ml/ Προβλεπόμενος χορηγούμενος όγκος ml/ Numatomas pristatomas tūris ml/ Planowana objętość dostarczana ml/ Планируемый дыхательный объем ml/ Předpokládaný přijatý objem ml/ Leadandó légzési térfogat ml/ Predvidena dovedena prostornina ml/ Paredzētais piegādes tilpums ml/ / Kavandatud manustatav maht ml/ Планиран осигуряван капацитет ml/ Planirana isporučena zapremina ml/ Volumul estimat livrat ml/ Určený podávaný objem ml/ Predvideni isporučenii volumen ml/ Amaçlanan uygulama hacmi ml

3 Flow resistance limit cmH₂O/l/min/ Limite de résistance à l'écoulement cmH₂O/l/min/ Grenzwert des Strömungswiderstands cmH₂O/l/min/ Limite de resistencia de flujo cmH₂O/l/min/ Limite de resistência ao fluxo cmH₂O/l/min/ Limite della resistenza al flusso cmH₂O/l/min/ Limiet stromingsweerstand cmH₂O/l/min/ Grense for flow-motstand cmH₂O/l/min/ Virtausvastuksen raja cmH₂O/l/min/ Gräns för flödesmotstånd cmH₂O/l/min/ Grænse for flow-modstand cmH₂O/l/min./ Όριο αντίστασης στη ροή cm H₂O/l/λεπτό/ Pasipriešinimo srautui riba cmH₂O/l/min/ Granica oporu przepływu cmH₂O/l/min/ Предел сопротивления потоку см H₂O/л/мин/ Limit odporu proudění cmH₂O/l/min/ Aramlási ellenállás határértéke cmH₂O/l/perc/ Омејитев одржности за pretok cmH₂O/l/min/ Plūsmas pretestības robeža cmH₂O/l/min/ Voolutakistuse piirväärtus cmH₂O/l/min/ Граница на съпротивление на потока cmH₂O/l/min/ Granična vrednost otpornosti protoku cmH₂O/l/min/ Limită de rezistență la curgere cmH₂O/l/min/ Limit odporu voči prúdeniu cmH₂O/l/min/ Granica otpora na protok cmH₂O/l/min/ Akış direnci sınırı cmH₂O/dakika

4 At flow l/min/ Débit l/min/ Bei Durchfluss l/min/ Flujo l/min/ Com o fluxo l/min/ Al flusso l/min/ Bij het stromen l/min/ Ved flow l/min/ Virtauksella l/min/ Vid flödet l/min/ Ved flow L/min./ Σε ροή l/λεπτό/ Esant srautui l/min/ Przy przepływie l/min/ В потоке л/мин/ Při proudění l/min/ Áramlási sebesség l/perc/ Pri pretoku l/min/ Pie aktīvas plūsmas l/min/ Voolukirusel l/min/ При поток l/min/ Pri protoku od l/min/ La curgere L/min/ Pri prúdeni l/min/ Pri protoku l/min/ Akışta l/ dakika/

5 Compliance ml/cmH₂O/ Tolérance ml/cmH₂O/ Compliance ml/cmH₂O/ Compliance ml/cmH₂O/ Complacência ml/cmH₂O/ Compliance ml/cmH₂O/ Naleving ml/cmH₂O/ Compliance ml/cmH₂O/ Kompiianssi ml/cmH₂O/ Compliance ml/cmH₂O/ Overensstemmelse ml/cmH₂O/ Συμμόρφωση ml/cm H₂O/ Atitiktis ml/cmH₂O/ Podatność ml/cmH₂O/ Растяжимость мл/см H₂O/ Shoda ml/cmH₂O/ Megfelelőség ml/cmH₂O/ Kompiancia ml/cmH₂O/ Atbilstība ml/cmH₂O/ Vastavus ml/cmH₂O/ Еластичность ml/cmH₂O/ Kompijansa ml/cmH₂O/ Compliantă ml/cmH₂O/ Poddajność ml/cmH₂O/ Rastezljivost ml/cmH₂O/ Uygunluk ml/cmH₂O/

6 At Pressure cmH₂O/ Pression cmH₂O/ Bei Druck cmH₂O/ Presión cmH₂O/ Com a pressão cmH₂O/ Alla pressione cmH₂O/ Onder druk cmH₂O/ Ved trykk cmH₂O/ Paineella cmH₂O/ Vid tryck cmH₂O/ Ved tryk cmH₂O/ Σε πίεση cm H₂O/ Esant slėgiui cmH₂O/ Przy ciśnieniu cmH₂O/ Под давлением см H₂O/ Při tlaku cmH₂O/ Nyomás cmH₂O/ Pri tlaku cmH₂O/ Pie spiediena cmH₂O/ Rõhul cmH₂O/ При налягане cmH₂O/ Pri pritisku od cmH₂O/ La presiune cmH₂O/ Pri tlaku cmH₂O/ Pri tlaku cmH₂O/ Basiñçta cmH₂O/

7 Adult/ Adulte/ Für Erwachsene/ Adulto/ Adulto/ Adulto/ Volwassene/ Voksne/ Aikuinen/ Vuxen/ Voksen/ Ενήλικες/ Su-augusieji/ Dorosły/ Для взрослых/ Dospělý/ Felnőtt/ Odrasli/ Pieaugušajiem/ Tāiskasvanu/ Взрaстни/ Odrasli/ Adultji/ Dospelí/ Odrasla osoba/ Yetişkin/

en 22mm Adult Breathing System

The instruction for use (IFU) is not supplied individually. Only 1 copy of the instructions for use is provided per box and should therefore be kept in an accessible location for all users. More copies are available upon request.

NOTE: Distribute instructions to all product locations and users. These instructions contain important information for the safe use of the product. Read these instructions for use in their entirety, including warnings and cautions before using this product. Failure to properly follow warnings, cautions and instructions could result in serious injury or death to the patient. Serious incidents to be reported as per the local regulation.

WARNING: A WARNING statement provides important information about a potentially hazardous situation, which, if not avoided could result in death or serious injury.

CAUTION: A CAUTION statement provides important information about a potentially hazardous situation, which if not avoided could result in minor or moderate injury.

Intended Use

To deliver and remove respiratory gases to and from a patient via a breathing system comprised of tubing and connectors.

Intended user

For use by qualified medical personnel.

Operational notes:

The breathing system should be connected between the ventilator or the anaesthetic machine and the patient airway connection. The product should only be used on a single patient up to a maximum of 7 days. The start of use is defined from the first connection to the equipment. The start of use is defined from the first connection to the equipment. The breathing system may be provided in Flex-tube™, or Smoothbore tubing.

Limitations of use:

Do not use with equipment that does not comply to international standards.

Product Use:

The product is intended for a maximum of 7 days use with a variety of anaesthetic and ventilator equipment meeting international standards. Use in the acute and non-acute environment under the supervision of suitably or appropriately qualified medical personnel. The breathing system is safe for use with atmospheric gases: N₂O, O₂ and CO₂. Anaesthetic agents: isoflurane, sevoflurane, desflurane and other comparable volatile agents.

Patient Category

The range of 22mm breathing systems are suitable for use on all adult patients with intended delivered tidal volumes of ≥ 300ml.

Instruction for use

Connect the breathing system

1. Attach the breathing system limbs to the inspiratory and expiratory ports of the ventilator or anaesthetic machine using a push and twist action (fig.1).
2. Remove any red dust cap from the Y-Piece of the breathing system.
3. Perform pre-use checks as recommended by the equipment manufacturer.
4. Attach the Y-Piece to the patient interface.

Accessories:

Where the 22mm breathing system is supplied with additional components the clinician should assess any potential impact upon the patient's respiratory function, and be aware of the subsequent increased resistance to flow, compressible volume and weight of the product.

When the ventilator and anaesthetic machine parameters are set, any impact that the internal volume and resistance to flow of the device may have upon the effectiveness of the prescribed ventilation, should be assessed on an individual patient basis by the clinician. Product configuration can be identified on the individual product labelling.

Water traps and breathing filters, HMEs, HMEFs (if included in the configuration)

Follow the instructions in individual IFUs supplied.

Monitoring lines (if included in the configuration)

Monitoring lines provide a leak free connection between the breathing system and a respiratory/anaesthetic gas monitor via a luer lock fitting. The connection should be checked prior to use. The line should be checked to confirm it is kink free and a patent pathway exists.

Pre-Use Check

1. Prior to attachment, check that all breathing system components are free of any obstruction or foreign bodies and that there is a patent airway.
2. Remove any red dust cap before connecting the breathing system between the equipment and the patient.
3. Connect/disconnect the tapered connections of the breathing system to the tapered connection of the ventilator and/or respiratory equipment using a push and twist action (fig.1).
4. Following attachment the breathing system and all accessories must be checked for leaks and occlusions prior to use and that all connections are secure.

5. When setting ventilator or anaesthetic machine parameters consider any impact that the internal volume of the breathing system and accessories may have upon the effectiveness of the prescribed ventilation. The clinician should assess this on an individual patient basis.

6. Check that any port cap or accessory is securely attached to prevent leakage.

In-Use Check

WARNING Monitor the clinical parameters of the patient during the treatment.

Warning

1. This system is not suitable for use with flammable anaesthetic agents.

2. Monitor the breathing system throughout use and replace if increased resistance or contamination occurs.

3. Refer to equipment manufacturers guidelines for machine settings and pre-use check requirements.

4. Perform a self-test of the equipment after full assembly of the breathing system and associated accessories before use on each patient.

5. Prior to use, consult the respective documentation and instructions for use of all connected devices, or combination of devices.

6. Do not stretch, cut or reconfigure the tubing.

7. Cleaning can lead to device malfunction and breakage.

8. Ventilator accuracy can be affected by use of a gas driven nebuliser.

9. Adding attachments or other devices to the breathing system can change the pressure gradient across the breathing system and adversely affect the ventilator performance. Perform a self-test of the ventilator after full assembly of the breathing system and all connected devices before use on each patient.

10. Use of this breathing system with equipment not compliant with appropriate International Standards can result in degraded performance. Prior to use, consult the respective documentation and instructions for use of all connected devices, or combination of devices to ensure compatibility.

11. Consult individual product labelling to determine if device is MR safe or unsafe.

12. Consult individual product labelling for presence or absence of phthalates.

13. If a luer port connector is included in the system, it is for the exclusive use of monitoring respiratory gases. If applicable, check that the cap is securely attached to prevent inadvertent leakage when the port is not in use. Do not introduce gases or liquids via the misconnection of the luer lock port as this can result in serious injury to the patient.

14. Do not use product past the expiry date.

Caution

1. The breathing system is only to be used under medical supervision by suitably or appropriately qualified medical personnel.

2. To avoid contamination, the product should remain packaged until ready for use. Do not use the product if the packaging is damaged.

3. When connecting or disconnecting the system, hold connector and use a push and twist action.

FOR USA: Rx Only: Federal law restricts this device to sale by or on the order of a physician.

Working pressure: <60 cmH₂O

Storage Conditions

Recommended storage at room temperature over the indicated product shelf-life.

Disposal

Following use, the product must be disposed of in accordance with local hospital, infection control and waste disposal regulations.

Note

The length referenced in the product description on the product label denotes the minimum working length of the breathing system. The full breathing system length is shown in the technical data section of the product label.

Technical and performance data

See product label schematic for full breathing system length (l), compliance (C) and resistance to flow (R) data.

WARNING Ensure that the breathing system technical and performance data are appropriate to the intended patient and treatment. Patient ventilation and vital signs should be monitored. Please refer to product label for full information.

Appropriate standards

This product complies with relevant requirements of the following standards:

- BS EN ISO 5367: Anaesthetic and respiratory equipment. Breathing tube and sets
- BS EN ISO 5356-1: Anaesthetic and respiratory equipment. Conical connectors. Cones and sockets

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