22mm adult heated wire breathing system

en	22mm adult heated wire breathing system
fr	Circuit respiratoire 22 mm adultes avec fil chauffant
de	Beheiztes Beatmungssystem mit 22-mm-Schlauch für Erwachsene
es	Sistema de respiración con cable caliente de 22 mm para adultos
pt	Sistema de respiração de fio aquecido de 22 mm para doentes adultos
it	Sistema di respirazione con cavo riscaldato per adulti da 22 mm
nl	22 mm ademhalingssysteem met verwarmde draad voor volwassenen
no	22 mm pustesystem med varmetråd for voksne
fi	22 mm:n lämmitysjohtimellinen hengitysjärjestelmä aikuisille
sv	22 mm andningssystem med värmeslinga för vuxna
da	22 mm voksen respirations-system med varmetråd
el	Σύστημα αναπνοής με θερμαινόμενο σύρμα 22 mm για ενήλικες
lt	22 mm suaugusiesiems skirta kvėpavimo sistema su šildoma viela
pl	System oddychania o średnicy 22 mm z drutem grzewczym dla dorosłych
ru	Дыхательная система для взрослых 22 мм с подогревом
cs	22mm dýchací systém s vyhřívacím drátem pro dospělé
hu	22 mm-es fűtőcsöves felhőtt lélegeztető rendszer
sl	Dihalni sistem za odrasle z grelno žico, 22 mm
lv	22 mm apsildāmā vada elpošanas sistēma pieaugušiem pacientiem
et	22 mm täiskasvanute hingamissüsteem soojendusjuhtmega
bg	Система за дишане с нагревател за възрастни 22 mm
sr	Sistem za disanje od 22 mm sa grejnim navojem za odrasle
ro	Circuit de ventilație pentru adulți, cu fir încălzit, 22mm
sk	22 mm dýchací systém s ohrievaným drôtom pre dospelých
hr	Sustav za disanje za odrasle sa zagrijanom žicom od 22 mm
tr	22 mm ısıtma telli yetişkin solunum sistemi
ja	22mm成人用熱線入り呼吸回路
zh-s	22 mm 成人加热丝呼吸系统
zh-t	22 mm 成人加熱線呼吸系統
ar	نظام تنفس مزود سخان سلكي للبالغين مقاس 22 مم


Figure 1. System set-up / Figure 1. Configuration du système / Abbildung 1. Systemaufbau / Figura 1. Configuración del sistema / Figura 1. Configuração do sistema / Figura 1. Configurazione del sistema / Figure 1. Systeem instellen / Figure 1. Systemoppsett / Kuva 1. Järjestelmän määrittäminen / Figur 1. Systeminställning / Figur 1. Opsætning af systemet / Εικόνα 1. Πύθμιση συστήματος / 1 pav. Sistemos nustatymas / Rysunek 1. Konfiguracja systemu / Рисунок 1. Настройка системы / Obr. 1. Sestavení systému / 1. ábra. A rendszer felépítése / 1. attēls. Sistēmas iestatīšana / Joonis 1. Süsteemi seadistamine / Фигура 1. Конфигурация на системата / Slika 1. Postavljanje sistema / Figura 1. Configurarea sistemului / Obrázok 1. Zostavenie systému / Slika 1. Postavljanje sustava / Şekil 1. Sistemin kurulumu / 図1. 回路のセットアップ / 图1. 系统设置 / 圖1. 系統設定 / 1. ماطنل داداغ / 1. لفشلال


Figure 2. Correct distal temperature probe insertion / Figure 2. Insertion correcte de l'extrémité distale de la sonde de température / Abbildung 2. Korrektes Einführen des distalen Temperaturfühlers / Figura 2. Inserción correcta de la sonda térmica distal / Figura 2. Inserção correta da sonda de temperatura distal / Figura 2. Inserimento corretto della sonda di temperatura distale / Figure 2. Correcte plaatsing van de distale temperatuursonde / Figur 2. Riktig innsetning av distal temperatuursonde / Kuva 2. Distaalisen lämpötila-anturin oikea asennus / Figur 2. Korrekt införing av distal temperaturgivare / Figure 2. Korrekt montering af distal temperatur probe / Εικόνα 2. Σωστή εισαγωγή του περιφερικού άκρου του αισθητήρα θερμοκρασίας / 2 pav. Tinkamas toliau esančio temperatūros zondo įvedimas / Rysunek 2. Poprawne wprowadzenie dystalnej sondy temperatury / Рисунок 2. Правильная установка дальнего датчика температуры / Obr. 2. Správne zasunutí distálného konce teplotní sondy / 2. ábra. A distális hőmérséklet-érzékelő megfelelő behelyezése / Slika 2. Pravilno vstavljanje distalne temperature sonde / 2. attēls. Pareiza distālās temperatūras zondes ievietošana / Joonis 2. Distaalse temperatuurianduri õige sisestamine / Фигура 2. Правильно поставлена дистална термосонда / Slika 2. Pravilno umetanje distalne temperature sonde / Figura 2. Introducerea corectă a sondei de temperatură distale / Obrázok 2. Správne zasunutie distálnej teplotnej sondy / Slika 2. Ispravno umetanje sonde za distalnu temperaturu / Şekil 2. Distal sıcaklık probunun doğru bir biçimde yerleştirilmesi / 図2. 温度プローブの正しい挿入 / 图2. 远端温度探头的正确插入 / 圖2. 正確的遠端溫度探針插入方式 / 2. ڊي جيبل ڤرارجل ۽ جرد سچ ۾ ڊيڊل ۽ جي ڇيصل ڦي رظلا / 2. لفشلال

en

The instruction for use is not supplied individually. Only 1 copy of the instructions for use is provided per box and should therefore be kept in an accessible location for all users. More copies are available upon request. **NOTE:** Distribute instructions to all product locations and users.

These instructions contain important information for the safe use of the product. Read these instructions for use in their entirety, including warnings and cautions before using this product. Failure to properly follow warnings, cautions and instructions could result in serious injury or death to the patient. Serious incidents to be reported as per the local regulation.

WARNING:

A WARNING statement provides important information about a potentially hazardous situation which, if not avoided could result in death or serious injury.

CAUTION:

A CAUTION statement provides important information about a potentially hazardous situation which, if not avoided could result in minor or moderate injury.

Description

22mm adult heated wire breathing system

Intended use

To deliver and remove humidified respiratory gases from a patient via a system of tubing and connectors.

Intended user

For use by or under supervision of appropriately qualified clinician or healthcare professional.

Indications

For use in a hospital environment.

For use with the Intersurgical PMH7000, Fisher & Paykel® MR850 and MR730 humidifiers.

For use with ventilators compliant to ISO 80601-2-12 Medical electrical equipment. Particular requirements for basic safety and essential performance of critical care ventilators.

For use with oxygen rich respiratory gases.

The breathing system may be provided in Flextube™, or Smoothbore tubing. Silver Knight™ breathing systems contain anti-microbial additive in the Flextube.

* Fisher and Paykel is a registered trademark of Fisher & Paykel®

Healthcare Limited.

Contraindications

CAUTION This product is not suitable for use in an MRI environment.

WARNING To ensure compatibility with this breathing system, only use ventilators that are compliant to ISO 80601-2-12 or ISO 80601-2-72. Prior to use consult the respective documentation and instructions for use of all connected devices, or combination of devices and perform system tests as required.

Patient category

Intended for use with adult patients.

Instructions for use

Preparation

1. Fit the chamber

- Slide the humidification chamber on to base unit top plate and remove cassette from chamber ports.

2. Hang water bag from pole

- Unwind the fill set from the cassette and remove the red cap from spike, then spike the bag. Hang bag on pole above the humidification chamber.

See humidification chamber IFU for details

- Ensure that the feed line is not kinked and chamber begins to fill.

3. Connect the breathing system (see figure 1)

- Attach dry limb (A) to ventilator inspiratory port and to chamber port.
- Connect the elbow on the green inspiratory limb (B) to other port on chamber.
- Connect expiratory limb (C) to expiratory port on ventilator.

4. Connect temperature probe

- Connect the temperature probe connection to the corresponding socket on humidifier base unit.
- Insert the distal temperature probe into the port on the elbow at the chamber. See figure 2.
- Insert the proximal temperature probe into the port at the patient end of the green inspiratory limb.

WARNING Ensure proximal temperature probe is located in the inspiratory limb and not in the y-piece.

5. Connect heater wire electrical adapter

- Connect the heated wire adapter into the socket on humidifier base unit.
- Connect the clover leaf end to the socket on the breathing system elbow above chamber.

Dual heated wires breathing systems

- Connect the oval end to socket on elbow at ventilator expiratory port.

6. Turn on humidifier base unit

- Ensure correct mode is activated.

7. Connect breathing system

- Remove red dust cap (D) from y-piece and attach y-piece to patient interface.

Pre-use Checks

WARNING

1. Every component of the breathing system must be visually inspected and checked for function, leakage and occlusions, immediately before use on each patient. Ensure that connections are secure.

2. To prevent disconnection of the tubing system during use, use only with devices, and/or accessories in compliance with the applicable international standards.

3. Perform a self-test of the ventilator after full assembly of the breathing

system and associated accessories before use on each patient.

4. Prior to use, consult the respective documentation and instructions for use of all connected devices, or combination of devices.

5. Check coiled heated wire. Do not use heated wire breathing system if wire is bunched within the tubes.

6. Ensure that the humidifier is mounted at a lower level than the patient and ventilator. See figure 1.

CAUTION

1. Prior to installation, check that all breathing system components are free of any obstruction or foreign bodies and that there is a patent airway.

2. When connecting or disconnecting the system, hold connector and use a push and twist action.

3. To avoid contamination, the product should remain packaged until ready for use. Do not use the product if the packaging is damaged

4. Use the red dust cap to occlude the breathing system patient connection port as and when required in the ventilator self-test protocol. Ensure red dust cap is held securely to prevent unwanted disconnection during self-test.

In-use checks

WARNING

1. Do not cover the breathing system with any thing that may insulate the air flow around the tube as this may impede normal heat transfer and cause damage to the system or affect gas delivery to the patient.

2. The 25cm of breathing system from the patient connection is an applied part, avoid extended contact with patient's skin for the remainder of the system.

3. Do not stretch the tubing.

4. Monitor the clinical parameters of the patient during the treatment at regular intervals.

CAUTION

1. Check humidification chamber water level before starting and during treatment. See chamber IFU for more details.

2. Check humidifier status before starting and during treatment. See humidifier IFU for more details.

Duration of use

Single use product. Maximum duration of use 7 days. The start of use is defined from the first connection to the equipment. Do not reuse. Reuse poses a threat to patient safety in terms of cross infection and device malfunction and breakage. Cleaning can lead to device malfunction and breakage.

WARNING

1. Disconnect breathing system from patient during defibrillation.

2. The materials used in the airways and related equipment may not be compatible with anaesthetic or respirable gases, solutions/suspensions/emulsions that have not been evaluated.

3. The use of nebulised drugs and humidification may increase the expiratory filter resistance to flow (where used). Frequently monitor resistance and replace as required.

4. Ventilator accuracy can be affected by use of a gas driven nebuliser.

5. Adding attachments or other devices to the breathing system can change the pressure gradient across the breathing system and adversely affect the ventilator performance. Perform a self-test of the ventilator after full assembly of the breathing system and all connected devices before use on each patient.

6. No modification of this equipment is allowed.

7. Monitor breathing system condensation and clear any mobile condensate that may build up in the tubes during treatment. Ensure condensate is drained towards the machine and away from the patient.

CAUTION

1. This product has been verified for a maximum duration of 7 days use.

2. Following use, the product must be disposed of in line with the local healthcare, hygiene and waste disposal protocols and regulations.

3. Interference can occur between humidifier electrical devices in close proximity. Ensure normal operation of devices prior to use.

Note

WARNING Do not use outside the specified environmental conditions. Using the device outside this range can affect the quality of the therapy or injure the patient.

Ensure compatibility of all parts of the system when connecting accessories and breathing systems to a ventilator and/ or humidifier.

Patient ventilation and vital signs should be monitored. Please refer to product label for full information.

Humidity output:

- **Invasive mode:** > 33mg/L* @ 5 - 50L/min

- **Non-invasive mode:** > 12mg/L* @ 5 - 50L/min

Working pressure: <125 cmH₂O

Maximum operating temperature: 26 °C

Minimum operating temperature: 18 °C

Breathing system leakage: < 70mL/min, additional limbs < 25mL/min at 60cmH₂O

* At sea level when used with Intersurgical humidification chamber.

Recommended storage at room temperature over the indicated product shelf-life.

For further information contact Intersurgical - info@intersurgical.com.

Standards

This product complies with the following relevant standards:

- BS EN ISO 5367: Anaesthetic and respiratory equipment. Breathing tube and sets

- BS EN ISO 5356-1: Anaesthetic and respiratory equipment. Conical connectors. Cones and sockets

- ISO 80601-2-74: Medical electrical equipment - Particular requirements for basic safety and essential performance of respiratory humidifying equipment

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INTERSURGICAL®
COMPLETE RESPIRATORY SYSTEMS

Intersurgical Ltd, Crane House, Molly Millars Lane,
Wokingham, Berkshire, RG41 2RZ, UK
T: +44 (0)118 9656 300 info@intersurgical.com
www.intersurgical.com

Ireland

T: +353 (0)1 697 8540
info@intersurgical.ie

France

T: +33 (0)1 48 76 72 30
info@intersurgical.fr

Deutschland

T: +49 (0)2241 25690
info@intersurgical.de

España

T: +34 91 665 73 15
info@intersurgical-es.com

Portugal

T: +351 219 108 550
info@intersurgical.pt

Italia

T: +39 0535 20836
info@intersurgical.it

Nederland

T: +31 (0)413 243 860
info@intersurgical.nl

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T: +32 5 2456991
info@intersurgical.be

Luxembourg

T: +32 5 2456991
info@intersurgical.be

Sverige

T: +46 08 514 30 600
info@intersurgical.se

Danmark

T: +45 28 69 80 00
info@intersurgical.dk

Lietuva

T: +370 387 66611
info@intersurgical.lt

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T: +7 495 7716809
info@intersurgical.ru

Česká Republika

T: +420 272 911 814
info@intersurgical.cz

Türkiye

T: +90 216 468 88 28
info@intersurgical.com.tr

South Africa

T: +27 (0)11 444 7968
info@intersurgical.co.za

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T: +86 519 69817000
info@intersurgical-cn.com

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T: 81 (0)3 6863 4359
info@intersurgical.co.jp

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T: +886 4 2380 5430
info@intersurgical.com.tw

Pilipinas

T: +63 2820 4124
info@intersurgical.ph

USA

T: +1 800 828 9633
support@intersurgicalinc.com

Canada

T: +1 905-319-6500
info@intersurgical.ca

Colombia

T: +57 1 742 5161
info@intersurgical.com

Australia

T: +61 (0)2 8048 3300
info@intersurgical.com.au

patent: <http://us.intersurgical.com/support>

Distributed in the USA by Intersurgical Incorporated, 6757 Kinne St., East Syracuse, NY 13057
T: (800) 828-9633 F: (315) 451-3696 support@intersurgicalinc.com www.intersurgical.com

2490A Issue 7 RO 08.22

Dual float auto fill humidification chamber

en	DUAL FLOAT AUTO FILL HUMIDIFICATION CHAMBER
fr	CHAMBRE D'HUMIDIFICATION À REMPLISSAGE AUTOMATIQUE AVEC DEUX FLOTTEURS
de	BEFEUCHTERKAMMER MIT AUTOMATISCHER FÜLLUNG UND DOPPELSCHWIMMER
es	CÁMARA DE HUMIDIFICACIÓN DE LLENADO AUTOMÁTICO CON DOBLE FLOTADOR
pt	CÂMARA DE HUMIDIFICAÇÃO DE ENCHIMENTO AUTOMÁTICO COM FLUTUADOR DUPLO
it	CAMERA DI UMIDIFICAZIONE A RIEMPIMENTO AUTOMATICO A DOPPIO GALLEGGIANTE
nl	AUTOMATISCHE VUL BEVOCHTIGINGSKAMER MET DUBBELE VLOTTER
no	FUKTINGSKAMMER MED DOBBEL FLOTTØR OG AUTOMATISK FYLING
fi	AUTOMAATTISESTI TÄYTETTÄVÄ KAKSOISVENTTIILILLÄ VARUSTETTU KOSTUTUSKAMMIO
sv	BEFUKTNINGSKAMMARE MED DUBBELT FLÖDE OCH AUTOMATISK Fyllning
da	BEFUGTNINGSKAMMER MED TO FLYDERE, DER PÅFYLDES AUTOMATISK
el	ΘΑΛΑΜΟΣ ΥΓΡΑΝΣΗΣ ΑΥΤΟΜΑΤΗΣ ΠΛΗΡΩΣΗΣ ΜΕ ΔΙΠΛΟ ΠΛΩΤΗΡΑ
lt	DVIAUKŠTĖ AUTOMATINIO UŽPILDYMO DRĖKINIMO KAMERA
pl	AUTOMATYCZNE NAPEŁNIANA KOMORA NAWILŻANIA Z DWOMA PŁYWAKAMI
ru	САМОЗАПОЛНЯЮЩАЯСЯ КАМЕРА УВЛАЖНЕНИЯ С ДВОЙНЫМ ПОПЛАВКОМ
cs	AUTOMATICKY PLNĚNÁ ZVLHČOVACÍ KOMORA SE SYSTÉMEM DVOU PLOVÁKŮ
hu	KETTŐS, ÚSZÓVAL ELLÁTOTT AUTOMATIKUS FELTÖLTÉSŰ PÁRÁSÍTÓ KAMRA
sl	VLAŽILNA KOMORA Z DVOJNIM PLOVCEM IN SAMODEJNIM POLNJENJEM
lv	AUTOMĀTISKI UZPILDĀMA MITRINĀŠANAS KAMERA AR DIVIEM PLUDIŅIEM
et	KAHE UJUKIGA AUTOMAATSE TÄITUMISEGA NIISUTIKAMBER
bg	ОВЛАЖНИТЕЛНА КАМЕРА С АВТОМАТИЧНО ПЪЛНЕНЕ С ДВОЕН ПОПЛАВЪК
sr	KOMORA ZA OVLAŽIVANJE SA DVA VENTILA I AUTOMATSKIM PUNJENJEM
ro	CAMERĂ DE UMIDIFICARE CU UMPLERE AUTOMATĂ CU FLOTOARE DUALE
sk	ZVLHČOVACIA KOMORA S DVOJPLAVÁKOVÝM SYSTÉMOM AUTOMATICKÉHO PLNENIA
hr	KOMORA ZA OVLAŽIVANJE S AUTOMATSKIM PUNJENJEM I DVA PLOVKA
tr	ÇİFT ŞAMANDIRALI, OTOMATİK DOLUMLU NEMLENDİRME HAZNESİ
ja	デュアルフロート自動充填加湿器チャンバ
zh-s	双浮子自动加水式湿化罐
zh-t	雙浮式自動填充加湿瓶
ar	غرفة ترطيب تلقائية الملء مزدوجة العوامة

The instruction for use is not supplied individually. Only 1 copy of the instructions for use is provided per box and should therefore be kept in an accessible location for all users. More copies are available upon request.

NOTE: Distribute instructions to all product locations and users.

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WARNING: A WARNING statement provides important information about a potentially hazardous situation which, if not avoided could result in death or serious injury.

CAUTION: A CAUTION statement provides important information about a potentially hazardous situation which, if not avoided could result in minor or moderate injury.

Description. Auto fill humidification chamber with fill set. **Intended use.** To humidify the respiratory gases being delivered to a patient.

Intended user. For use by qualified medical personnel. This product may be used in a homecare setting. Appropriately qualified medical personnel to instruct homecare users on the correct usage.

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* Fisher and Paykel is a registered trademark of Fisher & Paykel® Healthcare Limited.

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Patient category. Intended for use with adult, paediatric and neonatal patients.

INSTRUCTIONS FOR USE

Pre-use Checks

CAUTION When connecting or disconnecting the system, use a push and twist action.

WARNING To prevent disconnection of the tubing system during use, use only tubes in compliance with applicable ISO standards.

WARNING Every component of the breathing system must be visually inspected and checked for function, leakage and occlusions, immediately before use on each patient. Ensure that connections are secure.

WARNING Perform a self-test of the ventilator after full assembly of the breathing system and associated accessories before use on each patient.

WARNING Prior to use, consult the respective documentation and instructions for use of all connected devices, or combination of devices.

WARNING Ensure humidifier is level, maximum tilt angle < 10 degrees. See figure a.

WARNING Only use the chamber with sterile water suitable for inhalation or equivalent.

Preparation

1. Fit the chamber

- Slide the humidification chamber on to base unit top plate and remove the green cassette from the humidification ports. See figure a - 3.

2. Hang water bag from pole

- Unwind the fill set from the cassette and remove the red cap from spike (see figure a - 1), then spike the bag (see figure a - 2). Hang bag so that the height (H)

above the chamber is greater than maximum ventilated pressure (Pmax) in cmH₂O. See figure a.

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- Open spike vent when using rigid water reservoir. See figure a - 4.

CAUTION: Ensure that the feed line is free from trapped air bubbles before commencing treatment.

3. Turn on humidifier base unit

- Ensure correct mode is activated.

In-use checks

CAUTION Check humidification chamber water level before starting and during treatment. See figure a.

CAUTION Check humidifier status before starting and during treatment. See humidifier IFU for more details.

WARNING Replace the humidifier if water level exceeds maximum fill line. See figure a.

CAUTION When using Infant nasal CPAP please use a pressure infuser to ensure that reservoir bag pressure exceeds 140 cmH₂O. Check fluid levels after 2 minutes and while in use.

WARNING Metal surfaces may be hot when exposed. See figure a.

WARNING Monitor the clinical parameters of the patient, i.e O₂ concentration measured at the point of delivery to the patient, during the treatment.

Duration of use

Single use product. Maximum duration of use 7 days. Do not reuse. Reuse poses a threat to patient safety in terms of cross infection and device malfunction and breakage.

Cleaning can lead to device malfunction and breakage.

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Environmental Operating Conditions: Performance verified at room/ ambient temperature at sea level.

Resistance to flow

	Resistance to Flow (cm H ₂ O)					
	2.5 L/min	15 L/min	30 L/min	60 L/min	90 L/min	120 L/min
Max Fill Chamber	0.0	0.0	0.0	0.4	0.9	1.5

Compliance

	Compliance at 60cmH ₂ O
Min Fill Chamber	0.5
Max Fill Chamber	0.3

Min Humidification output: 33 mg/L @ BTPS

WARNING Please refer to breathing system IFU for humidification output data. Ensure technical and performance data are appropriate to the intended patient and treatment.

Patient ventilation and vital signs should be monitored.

For further information contact Intersurgical – info@intersurgical.com.

CAUTION Only to be used by qualified medical personnel or a user suitably trained to use the device by qualified medical personnel.

CAUTION Do not stretch or cut the fill set.

CAUTION This product has been verified for a maximum duration of 7 days use.

CAUTION Following use, the product must be disposed of in line with the local healthcare, hygiene and waste disposal protocols and regulations.

WARNING The materials used in the airways and related equipment may not be compatible with anaesthetic or respirable gases, solutions/suspensions/emulsions that have not been evaluated.

WARNING The use of nebulised drugs may cause discolouration of the chamber and / or water.

WARNING Do not use sterile products if the packaging is damaged, refer to the product label for information.

WARNING Use of this breathing system with systems and accessories not compliant with applicable ISO standards can result in degraded performance. Prior to use, consult the respective documentation and instructions for use of all connected devices, or combination of devices to ensure

compatibility.

WARNING Do not use outside the specified environmental conditions, as per the Humidifier IFU. Using the device outside this range can affect the quality of the therapy or injure the patient.

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- BS EN ISO 5356-1: Anaesthetic and respiratory equipment. Conical connectors. Cones and sockets.
- ISO 80601-2-74: Medical electrical equipment – Particular requirements for basic safety and essential performance of respiratory humidifying equipment.
- BS EN ISO 18190: Anaesthetic and respiratory equipment. General requirements for airways and related devices.

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INTERSURGICAL[®]
COMPLETE RESPIRATORY SYSTEMS

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Česká Republika
T: +420 272 911 814
info@intersurgical.cz

Türkiye
T: +90 216 468 88 28
info@intersurgical.com.tr

South Africa
T: +27 (0)11 444 7968
info@intersurgical.co.za

中国
T: +86 519 69817000
info@intersurgical-cn.com

日本
T: 81 (0)3 6863 4359
info@intersurgical.co.jp

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Pilipinas
T: +63 2820 4124
info@intersurgical.ph

USA
T: +1 800 828 9633
support@intersurgicalinc.com

Canada
T: +1 905-319-6500
info@intersurgical.ca

Colombia
T: +57 1 742 5161
info@intersurgical.com



Australia
T: +61 (0)2 8048 3300
info@intersurgical.com.au

patent: <http://us.intersurgical.com/support>
MADE BY INTERSURGICAL IN LITHUANIA

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High Volume Self-sealing water trap

CE
1639 Rx ONLY

Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK
T: +44 (0)118 9656 300 info@intersurgical.com www.intersurgical.com

EC REP

Ireland
T: +353 (0)1 697 8540
info@intersurgical.ie

España
T: +34 91 665 73 15
info@intersurgical-es.com

Nederland
T: +31 (0)413 243 860
info@intersurgical.nl

Sverige
T: +46 08 514 30 600
info@intersurgical.se

Россия
T: +7 495 7716809
info@intersurgical.ru

South Africa
T: +27 (0)11 444 7968
info@intersurgical.co.za

中国
T: +86 519 69817000
info@intersurgical-cn.com

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Canada
T: +1 905-319-6500
info@intersurgical.ca

France
T: +33 (0)1 48 76 72 30
info@intersurgical.fr

Portugal
T: +351 219 108 550
info@intersurgical.pt

België/Belgique
T: +32 5 2456991
info@intersurgical.be

Danmark
T: +45 28 69 80 00
info@intersurgical.dk

Česká Republika
T: +420 272 911 814
info@intersurgical.cz

日本
T: 81 (0)3 6863 4359
info@intersurgical.co.jp

Philippines
T: +63 2820 4124
info@intersurgical.ph

USA
T: +1 800 828 9633
support@intersurgicalinc.com

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MADE BY INTERSURGICAL IN LITHUANIA
9665 Issue 8 RO 08.21