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
en Breathing System

Instructions for use

1. Every component of the breathing system must be visually inspected and checked for function, leakage, and occlusions,

immediately before use on each patient. Ensure that connections are secure.

2. Remove the safety cap before connecting the system between the equipment and the patient.

3. When connecting or disconnecting the system, use a push and twist action.

4. If the system contains an Adjustable Pressure Limiting (APL) valve, set and adjust it as clinically required throughout the procedure. System pressure is maximised when the APL valve green cap is turned fully clockwise. The 60 cm H₂O pressure release cannot be overridden. It may be necessary to squeeze the reservoir bag to achieve 60 mbar system pressure.

5. If this system is used without a respiratory filter at the patient connection, then it must only be used on a single patient.

6. If this system is protected by an Intersurgical filter, or other appropriately validated respiratory filter, at the patient connection, and the filter and patient connection are changed for each patient, then the system may be used on more than one patient.

7. The system should be changed in accordance with the hospital department cleaning and disinfection procedures, or when it is no longer fit for its intended purpose.

8. When the system is not in use, the patient connection should be protected by the red cap supplied to prevent

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Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK

T: +44 (0)118 9656 300 info@intersurgical.com www.intersurgical.com

Distributed in the USA by Intersurgical Incorporated, 6751 Kinne Street, East Syracuse, NY 13057

T: (800) 828-9633 F: (315) 451-3696 support@intersurgicalinc.com www.intersurgical.com

Ireland T: +353 (0)1 697 8540 info@intersurgical.ie

France T: +33 (0)1 48 76 72 30 info@intersurgical.fr

Germany T: +49 (0)2241 25690 info@intersurgical.de

Spain T: +34 (0)11 665 73 15 info@intersurgical.es

Portugal T: +351 219 108 550 info@intersurgical.pt

Italy T: +39 0535 20836 info@intersurgical.it

Nederland T: +31 (0)143 243 860 info@intersurgical.nl

België/Belgique T: +32 5 2456991 info@intersurgical.be

Luxembourg T: +32 5 2456991 info@intersurgical.be

Sweden T: +46 08 514 30 600 info@intersurgical.se

Denmark T: +45 28 69 80 00 info@intersurgical.dk

Lietuva T: +370 387 66611 info@intersurgical.lt

South Africa T: +27 (0)11 444 7968 info@intersurgical.co.za

Ceská Republika T: +420 272 911 814 info@intersurgical.cz

Japan T: +81 (0)3 6863 4359 info@intersurgical.co.jp

Turkiye T: +90 216 468 88 28 info@intersurgical.com.tr

USA T: +1 800 828 9633 support@intersurgicalinc.com

Canada T: +1 905 319 6500 info@intersurgical.ca

Columbia T: +51 1 742 5161 info@intersurgical.pe

Australia T: +61 (0)2 8048 3300 info@intersurgical.com.au

Taiwan T: +886 4 2890 5410 info@intersurgical.com.tw

Philipinas T: +63 2820 4124 info@intersurgical.ph

USA T: +1 800 828 9633 support@intersurgicalinc.com

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