Mapleson D Anaesthetic Breathing System

en	MAPLESON D ANAESTHETIC BREATHING SYSTEM
fr	CIRCUIT RESPIRATOIRE D'ANESTHÉSIE MAPLESON D
de	MAPLESON D ANÄSTHESIE BEATMUNGSSYSTEM
es	CIRCUITO RESPIRATORIO PARA ANESTESIA MAPLESON D
pt	MAPLESON D SISTEMA RESPIRATÓRIO ANESTÉSICO
it	UNITÀ RESPIRATORIA MANUALE PER ANESTESIA MAPLESON D
nl	MAPLESON D BEADEMINGSSYSTEEM
no	MAPLESON D ANESTETISK PUSTESYSTEM
fi	MAPLESON D HENGITYSLETKUSTO
sv	MAPELSON D ANDNINGSSYSTEM FÖR ANESTESI
da	MAPLESON D ANÆSTESI RESPIRATIONSSYSTEM
el	MAPLESON D ΑΝΑΠΝΕΥΣΤΙΚΟ ΣΥΣΤΗΜΑ ΑΝΑΙΣΘΗΣΙΑΣ
lt	MAPLESON D ANESTEZINĖ KVĖPAVIMO SISTEMA
pl	MAPLESON D ANESTETYCZNY UKŁAD ODDECHOWY
ru	АНЕСТЕЗИОЛОГИЧЕСКАЯ ДЫХАТЕЛЬНАЯ СИСТЕМА MAPLESON D
cs	ANESTEZIOLOGICKÝ DÝCHACÍ SYSTÉM MAPLESON D
hu	MAPLESON D ANESZTÉZIÁS LÉGZŐ RENDSZER
sl	MAPLESON D ANESTEZIJSKI DIHALNI SISTEM
lv	MAPLESON D ANESTĒZIJAS ELPOŠANAS SISTĒMA
et	ANESTESIA HINGAMISSÜSTEEM MAPLESON D
bg	MAPLESON D АНЕСТЕТИЧНА СИСТЕМА ЗА ДИШАНЕ
sr	MAPLESON D SISTEM ZA DISANJE, ZA ANESTEZIJU
ro	CIRCUIT PENTRU ANESTEZIE MAPLESON D
sk	ANESTETICKÝ DÝCHACÍ SYSTÉM MAPLESON D
hr	MAPLESON D SISTEM ZA DISANJE ZA ANESTEZIJU
tr	MAPLESON D ANESTEZİ SOLUNUM SİSTEMİ


Fig. 1 / Figure 1. / Abbildung 1. / Figura 1. / Figur 1. / Kuva 1. / Εικόνα 1. / 1 pav. / Rysunek 1. / Рис. 1. /
 Obrázek 1. / 1. ábra / Slika 1 / 1. attēls. / Joonis 1. / Фигура 1. / Slika 1. / Obrázok 1. / Şekil 1

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en Mapleson D Anaesthetic Breathing System

The instruction for use (IFU) is not supplied individually. Only 1 copy of the instructions for use is provided per box and should therefore be kept in an accessible location for all users. More copies are available upon request.

NOTE: Distribute instructions to all product locations and users.

These instructions contain important information for the safe use of the product. Read these instructions for use in their entirety, including warnings and cautions before using this product. Failure to properly follow warnings, cautions and instructions could result in serious injury or death to the patient. Serious incidents to be reported as per the local regulation.

WARNING: A WARNING statement provides important information about a potentially hazardous situation, which, if not avoided could result in death or serious injury.

CAUTION: A CAUTION statement provides important information about a potentially hazardous situation, which if not avoided could result in minor or moderate injury.

Intended Use

To deliver and remove anaesthetic and respiratory gases to and from a patient via a breathing system comprised of tubing and connectors.

Intended user

For use by qualified medical personnel.

Indications

The breathing system should be connected to the anaesthetic machine and the patient airway connection. The product can be used repeatedly on several patients up to 7 days, if a suitably validated single-use respiratory filter is used between the product and each patient. The breathing system is provided in Flextube™ tubing. Silver Knight™ breathing systems contain anti-microbial additive in the Flextube.

Contraindications

No known contraindications.

Product Use:

The product is intended for a maximum of 7 days use with a variety of anaesthetic machines, meeting international standards. The start of use is defined from the first connection to the equipment. Use in the acute and non-acute environment under the supervision of suitably or appropriately qualified medical personnel. The breathing system is safe for use with atmospheric gases: N₂O, O₂ and CO₂; anaesthetic agents: isoflurane, sevoflurane, desflurane and other comparable volatile agents.

Patient Category

The range of Mapleson D anaesthetic breathing systems are suitable for use on all patients with intended delivered tidal volumes of ≥300 ml.

Connect the breathing system

Attach the breathing system to the fresh gas outlet of the anaesthetic machine using a push and twist action (Fig. 1).

Remove any red dust cap from the patient connection of the breathing system.

Perform pre-use checks as recommended by the equipment manufacturer.

Attach the patient connection to the patient airway.

Accessories:

Where the Mapleson anaesthetic breathing system is supplied with additional components the clinician should assess any potential impact upon the patient's respiratory function and be aware of the subsequent increased resistance to flow, compressible volume and weight of the product.

When the anaesthetic machine parameters are set, any impact that the system length, internal volume and resistance to flow of the device may have upon the effectiveness of the prescribed ventilation, should be assessed on an individual patient basis by the clinician.

Water traps and breathing filters, HMEs, HMEFs if included in the configuration

Follow the instructions in individual IFUs supplied.

Reservoir bags (if included in the system)

1. Ensure that the reservoir bag is securely attached to the breathing system via a push and twist action.

2. Consider the impact of increased compliance of the reservoir bag upon prescribed ventilation.

3. The reservoir bag allows manual ventilation, and gives a visual indication of rate and depth of breathing.

Monitoring lines (if included in the system)

Monitoring lines provide a leak free connection between the breathing system and a respiratory/anaesthetic gas monitor via a luer lock fitting. The connection should be checked prior to use. The line should be checked to confirm it is kink free and a patent pathway exists.

Pre-Use Check

1. Prior to installation, and each use, check that all breathing system components are free of any obstruction or foreign bodies and that there is a patent airway.

2. Remove any red dust cap before connecting the breathing system between the equipment and the patient.

3. Test the integrity of the green inner lumen tubing of the breathing system. Please follow steps indicated in figure 2. Further information on how to correctly perform this leak test can be found on the Intersurgical YouTube channel by searching *Using the Uniflow/Bain Leak Tester*.

4. Connect/disconnect the tapered connection of the breathing system to the tapered connection of the anaesthesia machine using a push and twist action.

5. Check function of the Adjustable Pressure Limiting (APL) valve by positioning the valve in the open position and set the fresh gas flow. Occlude the patient connection end to demonstrate gas flow through the APL valve. Adjust the APL valve to pressurise the system and the reservoir bag should be observed to inflate. Position APL valve in the closed position, squeeze firmly the inflated reservoir bag and listen for the APL

valve to open to confirm the safety pressure relief function of the valve.

6. The breathing system should be protected from cross contamination by an Intersurgical breathing filter, or other appropriately validated respiratory filter, at the patient connection. The filter and patient connection should be changed for each patient to minimise the risk of cross contamination.

7. Following installation and each use the breathing system and all accessories must be checked for leaks and occlusions and that all connections are secure. To complete a pre-use leak test, the APL valve must be in the closed position.

8. When setting prescribed fresh gas flow, consider any impact that the internal volume of the device may have upon the effectiveness of the patient's respiratory function. Ensure that sufficient fresh gas flow is set to avoid re-breathing of alveolar gas. The clinician should assess this on an individual patient basis.

9. Adjust the APL valve in accordance with patient's required respiratory function. Please refer to diagram below for further guidance. (Fig. 3)

10. Check that any port cap or accessory is securely attached to prevent leakage.

11. The breathing system should be changed in accordance with the hospital policy (max 7 days), or earlier if it is no longer fit for its intended purpose; it may pose a threat to patient safety in terms of cross infection, device malfunction and breakage.

In-Use Check

Monitor the clinical parameters of the patient during the treatment.

Warning

1. This system is not suitable for use with flammable anaesthetic agents.

2. Monitor the breathing system throughout use and replace if increased resistance or contamination occurs.

3. Refer to equipment manufacturers guidelines for machine settings and pre-use check requirements.

4. Perform a self-test of the equipment after full assembly of the breathing system and associated accessories before use on each patient.

5. Prior to use, consult the respective documentation and instructions for use of all connected devices, or combination of devices.

6. Do not stretch, cut or reconfigure the tubing.

7. Cleaning can lead to device malfunction and breakage.

8. Adding attachments or other devices to the breathing system can change the pressure gradient across the breathing system and adversely affect the ventilator performance. Perform a self-test of the ventilator (or anaesthesia machine) after full assembly of the breathing system and all connected devices before use on each patient.

9. Reuse of a breathing system in anaesthesia, that hasn't used a suitably validated single-use respiratory filter between the product and each patient, poses a threat to patient safety in terms of cross infection, device malfunction and breakage.

10. Breathing systems, their parts and accessories are validated for use with specific ventilators and anaesthesia machines. Incompatible parts can result in degraded performance. The responsible organisation is responsible for ensuring the compatibility of the ventilator or anaesthesia machine and all of the parts used to connect to the patient before use.

11. Use of this breathing system with equipment not compliant with appropriate International Standards can result in degraded performance. Prior to use, consult the respective documentation and instructions for use of all connected devices, or combination of devices to ensure compatibility.

12. If the system includes an APL valve, it is not MRI safe due to metal components in the valve. Please check product labelling.

13. Do not use with equipment that does not comply to international standards.

14. Do not use on patients with a tidal volume < 300ml.

Caution

1. The breathing system is only to be used under medical supervision by suitably or appropriately qualified medical personnel.

2. To avoid contamination, the product should remain packaged until ready for use. Do not use the product if the packaging is damaged.

3. The breathing system should not be re-used without an Intersurgical filter or other appropriately validated respiratory filter at the patient connection.

4. When connecting or disconnecting the system, use a push and twist action. (Fig. 1)

5. When the breathing system is not in use, the patient connection which is changed for each patient should be protected by the red dust cap supplied to prevent contamination.

FOR USA: Rx Only: Federal law restricts this device to sale by or on the order of a physician.

Working pressure: <60 cmH₂O

Storage Conditions

Recommended storage at room temperature over the indicated product shelf-life.

Disposal

Following use, the product must be disposed of in accordance with local hospital, infection control and waste disposal regulations.

Note

The length referenced in the product description on the product label denotes the minimum working length of the breathing system. The full breathing system length is shown in the technical data section of the product label.

Technical and performance data

See product label schematic for full breathing system length (l), compliance (C) and resistance to flow (R) data.

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