22mm Anaesthetic Breathing System

en	22MM ANAESTHETIC BREATHING SYSTEM
fr	22MM CIRCUIT D'ANESTHESIE
de	22MM ANAESTHESIE-BEATMUNGSSYSTEM
es	22MM CIRCUITO RESPIRATORIO PARA ANESTESIA
pt	22MM SISTEMA RESPIRATÓRIO ANESTÉSICO
it	22MM CIRCUITO PER ANESTESIA
nl	22MM BEADEMINGSSYSTEEM ANESTHESIE
no	22MM ANESTESISLANGESYSTEM
fi	22MM ANESTESIALETKUSTO
sv	22MM NARKOSSYSTEM
da	22MM ANÆSTESI RESPIRATIONSSYSTEM
el	22MM ΣΥΣΤΗΜΑ ΑΝΑΙΣΘΗΣΙΑΣ
lt	22MM ANESTEZINĖ KVĖPAVIMO SISTEMA
pl	22MM ANESTETYCZNY UKŁAD ODDECHOWY
ru	22MM КОМПОНАТ АНЕСТЕЗИОЛОГИЧЕСКИЙ
cs	22MM ANESTEZOLOGICKÝ DÝCHACÍ SYSTÉM
hu	22MM ANESZTEZIOLÓGIAI LÉLEGEZTETŐ RENDSZER
sl	22MM DIHALNI SISTEM ZA ANASTEZIJU
lv	22MM ANESTĒZIJAS KONTURS
et	22MM NARKOOSI HINGAMISSÜSTEEM
bg	22 MM АНЕСТЕТИЧНА СИСТЕМА ЗА ДИШАНЕ
sr	22MM SISTEM ZA DISANJE, ZA ANESTEZIJU
ro	22MM SISTEM DE RESPIRAȚIE ANESTEZIC
sk	22MM ANESTETICKÝ DÝCHACÍ SYSTÉM
hr	22MM ANESTETIČKI SUSTAV ZA DISANJE
tr	22MM ANESTEZİK SOLUNUM SİSTEMİ


Fig. 1 / Figure 1. / Abbildung 1. / Figura 1. / Figur 1. / Kuva 1. / Εικόνα 1. / 1 pav. / Rysunek 1. / Рис. 1. /Obrázek 1. / 1. ábra / Slika 1 / 1. attēls. / Joonis 1. / Фигура 1. / Slika 1. / Obrázok 1. / Šekil 1

Patient category ¹	Intended delivered volume ² ml	Flow resistance limit ³ cmH ₂ O/l/min	At flow ⁴ l/min
Adult ⁷	≥ 300 ml	0.06	30

Patient category	Intended delivered volume ml	Compliance ⁵ ml/cmH ₂ O	At Pressure ⁶ cmH ₂ O
Adult	≥ 300 ml	5	60 ± 3

1 Patient category/ Catégorie de patients/ Patientenkatēgorie/ Categoría de paciente/ Categoria de doentes/ Categorie di pazienti/ Παθιέντεκατēgorie/ Pasientekategorii/ Potilastyyppi/ Patientkatēgorii/ Patientkatēgorii/ Κατēgoria ασθενών/ Pacientų kategorija/ Kategoria pacjentów/ Категория пациентов/ Kategorie pacientů/ Betegek köre/ Kategorija bolnika/ Pacientu kategorija/ Patsiendi kategooria/ Категория на пациента/ Kategorija pacijenata/ Categorie de pacienți/ Kategoría pacientov/ Kategorija pacijenata/ Hasta kategorisi

2 Intended delivered volume ml/ Volume administré prévu ml/ Zuzuführendes Volumen ml/ Volumen de suministro previsto ml/ Volume administrado previsto ml/ Volume erogato previsto ml/ Beoogd volume ml/ Tiltent levert volum ml/ Aiottu antotilavuus ml/ Avsedd levererad volym ml/ Tilsigtet leveret volumen ml/ Προβλεπόμενος χορηγούμενος όγκος ml/ Numatomas pristatomas tūris ml/ Planowana objętość dostarczana ml/ Планируемый дыхательный объем ml/ Předpokládaný přijatý objem ml/ Leadandó légzési térfogat ml/ Predvidena dovedena prostornina ml/ Paredzētais piegādes tilpums ml/ / Kavandatud manustatav maht ml/ Планиран осигуряван капацитет ml/ Planirana isporučena zapremina ml/ Volumul estimat livrat ml/ Určený podávaný objem ml/ Predvideni isporučenii volumen ml/ Amaçlanan uygulama hacmi ml

3 Flow resistance limit cmH₂O/l/min/ Limite de résistance à l'écoulement cmH₂O/l/min/ Grenzwert des Strömungswiderstands cmH₂O/l/min/ Limite de resistenza de fluxo cmH₂O/l/min/ Limite de resistência ao fluxo cmH₂O/l/min/ Limite della resistenza al flusso cmH₂O/l/min/ Limiet stromingsweerstand cmH₂O/l/min/ Grense for flow-motstand cmH₂O/l/min/ Virtausvastuksen raja cmH₂O/l/min/ Gräns för flödesmotstånd cmH₂O/l/min/ Grænse for flow-motstand cmH₂O/l/min./ Όριο αντίστασης στη ροή cm H₂O/l/λεπτό/ Pasipriešinimo srautui riba cmH₂O/l/min/ Granica oporu przepływu cmH₂O/l/min/ Предел сопротивления потоку cm H₂O/l/мин/ Limit odporu proudění cmH₂O/l/min/ Aramlási ellenállás határértéke cmH₂O/l/perc/ Omejitvei odpornosti za pretok cmH₂O/l/min/ Plüsmas pretestības robeža cmH₂O/l/min/ Voolutakistuse piirväärtus cmH₂O/l/min/ Граница на съпротивление на потока cmH₂O/l/min/ Granična vrednost otpornosti protoku cmH₂O/l/min/ Limită de rezistență la curgere cmH₂O/L/min/ Limit odporu voči prúdeniu cmH₂O/l/min/ Granica otpora na protok cmH₂O/l/min/ Akış direnci sınırı cmH₂O/dakika

4 At flow l/min/ Débit l/min/ Bei Durchfluss l/min/ Flujo l/min/ Com o fluxo l/min/ Al flusso l/min/ Bij het stromen l/min/ Ved flow l/min/ Virtauksella l/min/ Vid flödet l/min/ Ved flow l/min./ Σε ροή l/λεπτό/ Esant srautui l/min/ Przy przepływie l/min/ В потоке л/мин/ Pri proudění l/min/ Aramlási sebesség l/perc/ Pri pretoku l/min/ Pie aktīvas plūsmas l/min/ Voolutiirusel l/min/ При поток l/min/ Pri protoku od l/min/ La curgere L/min/ Pri prúdení l/min/ Pri protoku l/min/ Akışta l/dakika/

5 Compliance ml/cmH₂O/ Tolérance ml/cmH₂O/ Compliance ml/cmH₂O/ Compliance ml/cmH₂O/ Complacência ml/cmH₂O/ Compliance ml/cmH₂O/ Naleving ml/cmH₂O/ Compliance ml/cmH₂O/ Komplianssi ml/cmH₂O/ Compliance ml/cmH₂O/ Overensstemmelse ml/cmH₂O/ Συμμόρφωση ml/cm H₂O/ Atitiklis ml/cmH₂O/ Podatnosť ml/cmH₂O/ Растяжимость мл/см H₂O/ Shoda ml/cmH₂O/ Megfelelőség ml/cmH₂O/ Komplianca ml/cmH₂O/ Atbilstība ml/cmH₂O/ Vastāvus ml/cmH₂O/ Еластичност ml/cmH₂O/ Komplijansa ml/cmH₂O/ Compliantă ml/cmH₂O/ Poddajnosť ml/cmH₂O/ Rastezljivost ml/cmH₂O/ Uygünlük ml/cmH₂O/

6 At Pressure cmH₂O/ Pression cmH₂O/ Bei Druck cmH₂O/ Presión cmH₂O/ Com a pressão cmH₂O/ Alla pressione cmH₂O/ Onder druk cmH₂O/ Ved trykk cmH₂O/ Painéella cmH₂O/ Vid tryck cmH₂O/ Ved tryk cmH₂O/ Σε πίεση cm H₂O/ Esant slėgiui cmH₂O/ Przy ciśnieniu cmH₂O/ Под давлением см H₂O/ Pri tlaku cmH₂O/ Nyomás cmH₂O/ Pri tlaku cmH₂O/ Pie spiediena cmH₂O/ Rõhul cmH₂O/ При налягане cmH₂O/ Pri pritisku od cmH₂O/ La presiune cmH₂O/ Pri tlaku cmH₂O/ Pri tlaku cmH₂O/ Basınçta cmH₂O/

7 Adult/Adulte/ Für Erwachsene/Adulto/Adulto/Adulto/ Volwassene/ Voksne/ Aikuisen/ Vuxen/ Voksen/ Ενήλικες/ Suaugusieji/ Dorosły/ Для взрослых/ Dospělý/ Felnőtt/ Odrasli/ Pieaugušajiem/ Täiskasvanu/ Взрастни/ Odrasli/ Adulti/ Dospeli/ Odrasla osoba/ Yetişkin/

en 22mm Adult Anaesthetic Breathing System

The instruction for use (IFU) is not supplied individually. Only 1 copy of the instructions for use is provided per box and should therefore be kept in an accessible location for all users. More copies are available upon request.

NOTE: Distribute instructions to all product locations and users. These instructions contain important information for the safe use of the product. Read these instructions for use in their entirety, including warnings and cautions before using this product. Failure to properly follow warnings, cautions and instructions could result in serious injury or death of the patient. Serious incidents to be reported as per the local regulation.

WARNING: A WARNING statement provides important information about a potentially hazardous situation, which, if not avoided could result in death or serious injury.

CAUTION: A CAUTION statement provides important information about a potentially hazardous situation, which if not avoided could result in minor or moderate injury.

Intended Use

To deliver and remove anaesthetic and respiratory gases to and from a patient via a breathing system comprised of tubing and connectors.

Intended user

For use by qualified medical personnel.

Indications

The anaesthetic breathing system should be connected between the ventilator or the anaesthetic machine, and the patient airway connection. The product can be used repeatedly on several patients up to 7 days, if a suitably validated single-use respiratory filter is used between the product and each patient. The breathing system may be provided in Flextube™, or Smoothbore tubing. Silver Knight™ breathing systems contain anti-microbial additive in the Flextube.

Contraindications

No known contraindications.

Product Use:

The product is intended for a maximum of 7 days use with a variety of anaesthetic and ventilator equipment meeting international taper standards. The start of use is defined from the first connection to the equipment. Use in the acute and non-acute environment under the supervision of suitably or appropriately qualified medical personnel. The breathing system is safe for use with atmospheric gases: N₂O, O₂, and CO₂; anaesthetic agents: isoflurane, sevoflurane, desflurane and other comparable volatile agents.

Patient Category

The range of 22mm breathing systems are suitable for use on all adult patients with intended delivered tidal volumes of ≥ 300ml.

Connect the breathing system

Attach the breathing system limb's to the inspiratory and expiratory ports of the ventilator or anaesthetic machine using a push and twist action (fig.1). Remove any red dust cap from the Y-Piece of the breathing system.

Perform pre-use checks as recommended by the equipment manufacturer.

Attach the Y-Piece to the patient interface.

Accessories:

Where the 22mm breathing system is supplied with additional components the clinician should assess any potential impact upon the patient's respiratory function, and be aware of the subsequent increased resistance to flow, compressible volume and weight of the product.

When the ventilator and anaesthetic machine parameters are set, any impact that the internal volume and resistance to flow of the device may have upon the effectiveness of the prescribed ventilation, should be assessed on an individual patient basis by the clinician.

Product configuration can be identified on the individual product labelling.

Water traps and breathing filters, HMEs, HMEFs (if included in the configuration)

Follow the instructions in individual IFUs supplied.

Reservoir bags (if included in the configuration)

1. Ensure that the reservoir bag is securely attached to the breathing system via a push and twist action.
2. Consider the impact of increased compliance of the reservoir bag upon prescribed ventilation.
3. The reservoir bag allows manual ventilation, and gives a visual indication of rate and depth of breathing.

Monitoring lines (if included in the configuration)

Monitoring lines provide a leak free connection between the breathing system and a respiratory/anaesthetic gas monitor via a luer lock fitting. The connection should be checked prior to use. The line should be checked to confirm it is kink free and a patent pathway exists.

Pre-Use Check

1. Prior to installation, and each use, check that all breathing system components are free of any obstruction or foreign bodies and that there is a patent airway.
2. Remove any red dust cap before connecting the breathing system between the equipment and the patient.
3. Connect/disconnect the tapered connection of the breathing system to the tapered connection of the anaesthesia machine using a push and twist action.
4. The breathing system should be protected from cross contamination by an Intersurgical breathing filter, or other appropriately validated respiratory filter, at the patient connection. The filter and patient connection should be changed for each patient to minimise the risk of cross contamination. If the product is visibly contaminated dispose and replace with a new product.
5. Following installation and each use the breathing system and all accessories must be checked for leaks and occlusions and that all connections are secure.
6. When setting ventilator parameters consider any impact that the internal volume of the device may have upon the effectiveness of the prescribed ventilation. The clinician should assess this on an individual patient basis.
7. Check that any port cap or accessory is securely attached to prevent leakage.
8. The breathing system should be changed in accordance with the hospital policy (max 7 days), or earlier if it is no longer fit for its intended purpose.

In-Use Check

WARNING Monitor the clinical parameters of the patient during the treatment.

Warning

1. This system is not suitable for use with flammable anaesthetic agents.
 2. Monitor the breathing system throughout use and replace if increased resistance or contamination occurs.
 3. Refer to equipment manufacturers guidelines for machine settings and pre-use check requirements.
 4. Perform a self-test of the equipment after full assembly of the breathing system and associated accessories before use on each patient.
 5. Prior to use, consult the respective documentation and instructions for use of all connected devices, or combination of devices.
 6. Do not stretch, cut or reconfigure the tubing.
 7. Cleaning can lead to device malfunction and breakage.
 8. Ventilator accuracy can be affected by use of a gas driven nebuliser.
 9. Adding attachments or other devices to the breathing system can change the pressure gradient across the breathing system and adversely affect the ventilator performance. Perform a self-test of the ventilator (or anaesthesia machine) after full assembly of the breathing system and all connected devices before use on each patient.
 10. Reuse of a breathing system in anaesthesia, that hasn't used a suitably validated single-use respiratory filter between the product and each patient, poses a threat to patient safety in terms of cross infection, device malfunction and breakage.
 11. Breathing systems, their parts and accessories are validated for use with specific ventilators and anaesthesia machines. Incompatible parts can result in degraded performance. The responsible organisation is responsible for ensuring the compatibility of the ventilator or anaesthesia machine and all of the parts used to connect to the patient before use.
 12. Use of this breathing system with equipment not compliant with appropriate International Standards can result in degraded performance. Prior to use, consult the respective documentation and instructions for use of all connected devices, or combination of devices to ensure compatibility.
 13. Consult individual product labelling to determine if device is MR safe or unsafe.
 14. Consult individual product labelling for presence or absence of phthalates.
 15. If a luer port connector is included in the system, it is for the exclusive use of monitoring respiratory and anaesthetic gases. If applicable, check that the cap is securely attached to prevent inadvertent leakage when the port is not in use. Do not introduce gases or liquids via the connection of the luer lock port as this can result in serious injury to the patient.
 16. Do not use with equipment that does not comply to international taper connection standards.
 17. Do not use on patients with a tidal volume < 300ml.
 18. Ensure that the breathing system technical and performance data are appropriate to the intended patient and treatment. Patient ventilation and vital signs should be monitored. Please refer to product label for full information.
- ### Caution
1. The breathing system is only to be used under medical supervision by suitably or appropriately qualified medical personnel.
 2. To avoid contamination, the product should remain packaged until ready for use. Do not use the product if the packaging is damaged.
 3. When used on more than one patient, the anaesthetic breathing system should not be re-used without an Intersurgical filter or other appropriately validated respiratory filter at the patient connection.
 4. When connecting or disconnecting the system, hold connector and use a push and twist action.
 5. When the breathing system is not in use, the patient connection which is changed for each patient when used in anaesthesia should be protected by the red dust cap supplied to prevent contamination.

FOR USA: Rx Only: Federal law restricts this device to sale by or on the order of a physician.

Working pressure: <60 cmH₂O

Storage Conditions

Recommended storage at room temperature over the indicated product shelf-life.

Disposal

Following use, the product must be disposed of in accordance with local hospital, infection control and waste disposal regulations.

Note

The length referenced in the product description on the product label denotes the minimum working length of the breathing system. The full breathing system length is shown in the technical data section of the product label.

Appropriate standards

This product complies with relevant requirements of the following standards:

- BS EN ISO 5367: Anaesthetic and respiratory equipment. Breathing tube and sets
- BS EN ISO 5356-1: Anaesthetic and respiratory equipment. Conical connectors. Cones and sockets
- ISO 80601-2-13: Particular requirements for basic safety and essential performance of an anaesthetic workstation

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* Explanation of packing symbols can be found on the Intersurgical website <http://www.intersurgical.com/support> or in standard ISO 15223-1:2016.

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CE
1639

Rx ONLY

INTERSURGICAL[®]
COMPLETE RESPIRATORY SYSTEMS

Intersurgical Ltd, Crane House, Molly Millars Lane,
Wokingham, Berkshire, RG41 2RZ, UK
T: +44 (0)118 9656 300 info@intersurgical.com
www.intersurgical.com

Ireland
T: +353 (0)1 697 8540
info@intersurgical.ie

France
T: +33 (0)1 48 76 72 30
info@intersurgical.fr

Deutschland
T: +49 (0)2241 25690
info@intersurgical.de

España
T: +34 91 665 73 15
info@intersurgical-es.com

Portugal
T: +351 219 108 550
info@intersurgical.pt

Italia
T: +39 0535 20836
info@intersurgical.it

Nederland
T: +31 (0)413 243 860
info@intersurgical.nl

België/Belgique
T: +32 5 2456991
info@intersurgical.be

Luxembourg
T: +32 5 2456991
info@intersurgical.be

Sverige
T: +46 08 514 30 600
info@intersurgical.se

Danmark
T: +45 28 69 80 00
info@intersurgical.dk

Lietuva
T: +370 387 66611
info@intersurgical.lt

Россия
T: +7 495 7716809
info@intersurgical.ru

Česká Republika
T: +420 272 911 814
info@intersurgical.cz

Türkiye
T: +90 216 468 88 28
info@intersurgical.com.tr

South Africa
T: +27 (0)11 444 7968
info@intersurgical.co.za

中国
T: +86 519 69817000
info@intersurgical-cn.com

日本
T: 81 (0)3 6863 4359
info@intersurgical.co.jp

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Pilipinas
T: +63 2820 4124
info@intersurgical.ph

USA
T: +1 800 828 9633
support@intersurgicalinc.com



Canada
T: +1 905-319-6500
info@intersurgical.ca

Colombia
T: +57 1 742 5161
info@intersurgical.com

Australia
T: +61 (0)2 8048 3300
info@intersurgical.com.au

WATER TRAP, SELF SEALING, MEDIUM VOLUME

CE
1639 Rx ONLY

INTERSURGICAL
COMPLETE RESPIRATORY SYSTEMS

Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK

T: +44 (0)118 9656 300 info@intersurgical.com www.intersurgical.com

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Česká Republika
T: +420 272 911 814
info@intersurgical.cz

Türkiye
T: +90 216 468 88 28
info@intersurgical.com.tr

South Africa
T: +27 (0)11 444 7968
info@intersurgical.co.za

中国
T: +86 519 6981 7000
info@intersurgical-cn.com

日本
T: 81 (0)3 6863 4359
info@intersurgical.co.jp

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Filipinas
T: +63 2820 4124
info@intersurgical.ph

USA
T: +1 800 828 9633
support@intersurgicalinc.com

Canada
T: +1 905-319-6500
info@intersurgical.ca

Colombia
T: +57 1 742 5161
info@intersurgical.com

Australia
T: +61 (0)2 8048 3300
info@intersurgical.com.au

Distributed in the USA by Intersurgical Incorporated, 6757 Kinne St., East Syracuse, NY 13057
T: (800) 828-9633 F: (315) 451-3696 support@intersurgicalinc.com www.intersurgical.com

7705 Issue 9 RO 08.21