15 mm Breathing System

en	15MM BREATHING SYSTEM
fr	15MM CIRCUIT DE RESPIRATION
de	15MM BEATMUNGS-SYSTEM
es	15MM CIRCUITO RESPIRATORIO
pt	15MM SISTEMA RESPRATÓRIO
it	15MM CIRCUITO RESPIRATORIO
nl	15MM BEADEMINGSSYSTEEM
no	15MM RESPIRATOR SLANGESYSTEM
fi	15MM HENGITYSLETKUSTO
sv	15MM ANDNINGSSYSTEM
da	15MM RESPIRATIONSSYSTEM
el	15MM ΣΥΣΤΗΜΑ ΑΝΑΠΝΟΗΣ
lt	15MM KVĖPAVIMO SISTEMA
pl	15MM SYSTEM ODDECHOWY
ru	15MM ДЫХАТЕЛЬНЫЙ КОНТУР
cs	15MM DÝCHACÍ SYSTÉM
hu	15MM LÉLEGEZTETŐ RENDSZER
sl	15MM DIHALNI SISTEM
lv	15MM ELPOŠANAS SISTĒMA
et	15MM HINGAMISAPARAAT
bg	15MM РЕСПИРАТОРНА СИСТЕМА
sr	15MM SISTEM ZA DISANJE
ro	15MM SISTEM DE RESPIRAȚIE
sk	15MM D/SYSTÉM
hr	15MM SUSTAV ZA DISANJE
tr	15MM SOLUNUM DEVRESİ
zh-s	15 MM 呼吸回路
zh-t	15MM 呼吸回路
ar	15MM ملء نظام تنفس



Figure 1. / Figure 1. / Abbildung 1. / Figura 1. / Figura 1. / Figuur 1. / Kuva 1. / Figur 1. / Figur 1. / Eikóna 1. / 1 pav. / Rysunek 1. / Рисунок 1. / Obr. 1. / 1. ábra. / 1. attēls. / Joonis 1. / Фигура 1. / Slika 1. / Figura 1. / Obrázok 1. / Slika 1. / Şekil 1. / 图 1 / 圖 1 / لفشال 1

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Patient category	Intended delivered volume ml	Compliance ⁵ ml/cmH ₂ O	At Pressure ⁶ cmH ₂ O
Adult	≥ 300 ml	5	60 ± 3
Paediatric	50 ml < 300 ml	4	60 ± 3
Neonatal	≤ 50 ml	1.5	60 ± 3

9 Neonatal/ Nouveau-nés/ Neonatalen/ Neonatal/ Recém-nascidos/ Neonatale/ Neonataal/ Nyföddte/ Vastasyntyneille/ Nyfödda/ Neonatal/ Neovvđ/ Naujagimiai/ Noworodki/ Для новорожденных/ Novorozenci/ Újszülött/ Neonatalni/ Jaundzimušiem/ Vastsündinute/ Неонатални/ Nou-născuți/ Новорождённые/ Novorođenčad/ Yenidoğan/ 新生儿/ دالوال نیو-ناسد

en 15mm Breathing System

The instruction for use (IFU) is not supplied individually. Only 1 copy of the instructions for use is provided per box and should therefore be kept in an accessible location for all users. More copies are available upon request.

NOTE: Distribute instructions to all product locations and users.

These instructions contain important information for the safe use of the product. Read these instructions for use in their entirety, including warnings and cautions before using this product. Failure to properly follow warnings, cautions and instructions could result in serious injury or death to the patient. Serious incidents to be reported as per the local regulation.

WARNING: A WARNING statement provides important information about a potentially hazardous situation, which, if not avoided could result in death or serious injury.

CAUTION: A CAUTION statement provides important information about a potentially hazardous situation, which if not avoided could result in minor or moderate injury.

Intended Use

When used in Critical Care: To deliver and remove respiratory gases to and from a patient via a breathing system comprised of tubing and connectors. When used in Anaesthesia: To deliver and remove anaesthetic and respiratory gases to and from a patient via a breathing system comprised of tubing and connectors.

Intended user

For use by qualified medical personnel.

Indications

The breathing system should be connected between the ventilator or the anaesthetic machine, and the patient airway connection. The product should only be used on a single patient up to a maximum of 7 days. The breathing system may be provided in Flextube™, Compact™ extendable and Smoothbore tubing.

Contraindications

Do not use with equipment that does not comply to international taper connection standards.

Product Use:

The product is intended for a maximum of 7 days use with a variety of anaesthetic and ventilator equipment meeting international taper standards. The start of use is defined from the first connection to the equipment. Use in the acute and non-acute environment under the supervision of suitably or appropriately qualified medical personnel. The breathing system is safe for use with atmospheric gases: N₂O, O₂ and CO₂; anaesthetic agents: isoflurane, sevoflurane, desflurane and other comparable volatile agents.

Patient Category

The range of 15mm breathing systems are suitable for use on all paediatric patients with intended delivered tidal volumes of 50ml < 300ml.

Connect the breathing system

Attach the breathing system limb's to the inspiratory and expiratory ports of the ventilator or anaesthetic machine using a push and twist action (fig.1).

Remove any red dust cap from the Y-Piece of the breathing system.

Perform pre-use checks as recommended by the equipment manufacturer. Attach the Y-Piece to the patient interface.

Accessories:

Where the 15mm breathing system is supplied with additional components the clinician should assess any potential impact upon the patient's respiratory function, and be aware of the subsequent increased resistance to flow, compressible volume and weight of the product. When the ventilator and anaesthetic machine parameters are set, any impact that the internal volume and resistance to flow of the device may have upon the effectiveness of the prescribed ventilation, should be assessed on an individual patient basis by the clinician.

Water traps and breathing filters, HMEs, HMEFs (if included in the configuration)

Follow the instructions in individual IFUs supplied.

Reservoir bags (if included in the configuration)

1. Ensure that the reservoir bag is securely attached to the breathing system via a push and twist action.

2. Consider the impact of increased compliance of the reservoir bag upon prescribed ventilation.

3. The reservoir bag allows manual ventilation, and gives a visual indication of rate and depth of breathing.

Monitoring lines (if included in the configuration)

Monitoring lines provide a leak free connection between the breathing system and a respiratory/anaesthetic gas monitor via a luer lock fitting. The connection should be checked prior to use. The line should be checked to confirm it is kink free and a patent pathway exists.

Pre-Use Check

1. Prior to attachment, check that all breathing system components are free of any obstruction or foreign bodies and that there is a patent airway.

2. Remove any red dust cap before connecting the breathing system between the equipment and the patient.

3. Connect/disconnect the tapered connections of the breathing system to the tapered connection of the ventilator and/or respiratory equipment using a push and twist action (fig.1).

4. Following attachment the breathing system must be checked for leaks and occlusions prior to use and that all connections are secure.

5. When setting ventilator parameters consider any impact that the internal volume of the breathing system and accessories may have upon the effectiveness of the prescribed ventilation. The clinician should assess this on an individual patient basis.

6. Check that any port cap or accessory is securely attached to prevent leakage.

In-Use Check

The system should be changed if visibly contaminated or if resistance to flow increases.

Monitor the clinical parameters of the patient during the treatment.

Warning

1. This system is not suitable for use with flammable anaesthetic agents.

2. Monitor the breathing system throughout use and replace if increased resistance or contamination occurs.

3. Refer to equipment manufacturers guidelines for machine settings and pre-use check requirements.

4. Following use, the product must be disposed of in accordance with local hospital, infection control in line with the local healthcare, hygiene and waste disposal protocols and regulations.

5. Perform a self-test of the equipment after full assembly of the breathing system and associated accessories before use on each patient.

6. Prior to use, consult the respective documentation and instructions for use of all connected devices, or combination of devices.

7. Do not stretch, cut or reconfigure the tubing.

8. Cleaning can lead to device malfunction and breakage.

9. Ventilator accuracy can be affected by use of a gas driven nebuliser.

10. Adding attachments or other devices to the breathing system can change the pressure gradient across the breathing system and adversely affect the ventilator performance. Perform a self-test of the ventilator (or anaesthesia machine) after full assembly of the breathing system and all connected devices before use on each patient.

11. Breathing systems, their parts and accessories are validated for use with specific ventilators and anaesthesia machines. Incompatible parts can result in degraded performance. The responsible organisation is responsible for ensuring the compatibility of the ventilator or anaesthesia machine and all of the parts used to connect to the patient before use.

12. Use of this breathing system with equipment not compliant with appropriate International Standards can result in degraded performance. Prior to use, consult the respective documentation and instructions for use of all connected devices, or combination of devices to ensure compatibility.

13. Presence or absence of Phthalates is indicated on the individual product labelling.

14. Consult individual product labelling for presence or absence of phthalates.

15. If a luer port connector is included in the system, it is for the exclusive use of monitoring respiratory and anaesthetic gases. If applicable, check that the cap is securely attached to prevent inadvertent leakage when the port is not in use. Do not introduce gases or liquids via the misconnection of the luer lock port as this can result in serious injury to the patient.

Caution

1. The breathing system is only to be used under medical supervision by suitably or appropriately qualified medical personnel.

2. To avoid contamination, the product should remain packaged until ready for use. Do not use the product if the packaging is damaged.

3. When connecting or disconnecting the system, hold connector and use a push and twist action.

4. When the breathing system is not in use, the patient connection which is changed for each patient should be protected by the red dust cap supplied to prevent contamination.

FOR USA: Rx Only: Federal law restricts this device to sale by or on the order of a physician.

Working pressure: <60 cmH₂O

Storage Conditions

Store out of direct sunlight.

Recommended storage at room temperature over the indicated product shelf-life.

Environmental Operating Conditions: Performance verified at room/ambient temperature at sea level.

Disposal

Following use, the product must be disposed of in accordance with local hospital, infection control and waste disposal regulations.

Note

The length referenced in the product description on the product label denotes the minimum working length of the breathing system. The full breathing system length is shown in the technical data section of the product label.

Technical and performance data

See product label schematic for full breathing system length (l), compliance (C) and resistance to flow (R) data.

WARNING Ensure that the breathing system technical and performance data are appropriate to the intended patient and treatment.

Patient ventilation and vital signs should be monitored. Please refer to product label for full information.

Appropriate standards

This product complies with relevant requirements of the following standards:

- BS EN ISO 5367: Anaesthetic and respiratory equipment. Breathing tube and sets
- BS EN ISO 5356-1: Anaesthetic and respiratory equipment. Conical connectors. Cones and sockets
- ISO 80601-2-13: Particular requirements for basic safety and essential performance of an anaesthetic workstation
- BS EN ISO 80601-2-12: MEE Critical Care Ventilators

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* Explanation of packing symbols can be found on the Intersurgical website <http://www.intersurgical.com/support> or in standard ISO 15223-1:2016.

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INTERSURGICAL®
COMPLETE RESPIRATORY SYSTEMS

Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK

T: +44 (0)118 9656 300 info@intersurgical.com www.intersurgical.com

Ireland
T: +353 (0)1 897 8540
info@intersurgical.ie

España
T: +34 91 865 73 15
info@intersurgical-es.com

Portugal
T: +351 219 108 550
info@intersurgical.pt

Deutschland
T: +49 (0)2241 25690
info@intersurgical.de

Italia
T: +39 0535 20836
info@intersurgical.it

Nederland
T: +31 (0)413 243 860
info@intersurgical.nl

België/Belgique
T: +32 5 2456991
info@intersurgical.be

Luxembourg
T: +32 5 2456991
info@intersurgical.be

Sverige
T: +46 08 514 30 600
info@intersurgical.se

Danmark
T: +45 28 69 80 00
info@intersurgical.dk

Lietuva
T: +370 387 66611
info@intersurgical.lt

Россия
T: +7 495 7716809
info@intersurgical.ru

Česká Republika
T: +420 272 911 814
info@intersurgical.cz

Türkiye
T: +90 216 468 88 28
info@intersurgical.com.tr

South Africa
T: +27 (0)11 444 7968
info@intersurgical.co.za

中国
T: +86 519 69817000
info@intersurgical-cn.com

日本
T: 81 (0)3 6863 4359
info@intersurgical.co.jp

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Philippines
T: +63 2820 4124
info@intersurgical.ph

USA
T: +1 800 828 9633
support@intersurgicalinc.com

Canada
T: +1 905-319-6500
info@intersurgical.ca



Colombia
T: +57 1 742 5161
info@intersurgical.com

Australia
T: +61 (0)2 8048 3300
info@intersurgical.com.au

Distributed in the USA by Intersurgical Incorporated, 6757 Kinne St., East Syracuse, NY 13057
T: (800) 828-9633 F: (315) 451-3696 support@intersurgicalinc.com www.intersurgical.com

2665 Issue 5 RO 02.22

WATER TRAP, SELF SEALING, MEDIUM VOLUME

CE
1639 Rx ONLY

Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK

T: +44 (0)118 9656 300 info@intersurgical.com www.intersurgical.com

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