22mm Adult Compact™ Breathing System

en	22MM ADULT COMPACT BREATHING SYSTEM
fr	CIRCUIT RESPIRATOIRE ADULTE 22MM COMPACT
de	22MM COMPACT BEATMUNGSSYSTEM, ERWACHSENE
es	CIRCUITO RESPIRATORIO COMPACT PARA ADULTO DE 22MM
pt	SISTEMA RESPIRATÓRIO COMPACT ADULTO 22MM
it	CIRCUITO PER VENTILAZIONE COMPACT 22MM, ADULTO
nl	22MM COMPACT BEADEMINGSSYSTEEM VOOR VOLWASSENEN
no	22MM COMPACT SLANGESYSTEM, VOKSEN
fi	22MM AIKUISTEN COMPACT HENGITYS LETKUSTO
sv	22MM COMPACT ANDNINGSSYSTEM VUXEN
da	22 MM COMPACT RESPIRATIONSSÆT, VOKSEN
el	COMPACT ΚΥΚΛΩΜΑ ΑΝΑΙΣΘΗΣΙΑΣ ΣΥΜΠΙΕΣΙΜΟ 22Μ ΕΝΗΛΙΚΩΝ
lt	22MM DIAMETRO COMPACT REGULIUOJAMO ILGIO KVĖPAVIMO SISTEMA SUAUGUSIEMS
pl	22MM UKŁAD ODDECHOWY COMPACT DLA DOROSŁYCH
ru	ДЫХАТЕЛЬНЫЙ КОНТУР АНЕСТЕЗИОЛОГИЧЕСКИЙ ДЛЯ ВЗРОСЛЫХ COMPACT
cs	22MM DÝCHACÍ SYSTÉM COMPACT PRO DOSPĚLÉ
hu	22 MM FELNŐTT COMPACT LÉGZŐKÖR
sl	COMPACT DIHALNI SISTEM ZA ODRASLE 22MM, RAZTEGLJIV
lv	ELPOŠANAS KONTŪRA PIEAUGUŠO 22MM COMPACT
et	22MM TÄISKSVANU COMPACT HINGAMISKONTUUR
bg	22MM, COMPACT, ДИХАТЕЛНА СИСТЕМА ЗА ВЪЗРАСТНИ
sr	22MM SISTEM ZA DISANJE ODRASLE COMPACT
ro	CIRCUIT DE VENTILATIE COMPACT PENTRU ADULTI, 22MM
sk	COMPACT DÝCHACÍ SYSTÉM PRE DOSPELÝCH 22 MM
hr	COMPACT SUSTAV ZA DISANJE PROMJERA 22 MM, ZA ODRASLE
tr	22MM YETİŞKİN COMPACT SOLUNUM SİSTEMİ


Figure 1. / Figure 1. / Abbildung 1. / Figura 1. / Figura 1. / Figura 1. / Figuur 1. / Figur 1. / Kuva 1. / Figur 1. / Figur 1. / Εικόνα 1. / 1 pav. / Rysunek 1. / Рисунок 1. / Obr. 1. / 1. ábra. / 1. attēls. / Joonis 1. / Фигура 1. / Slika 1. / Figura 1. / Obrázok 1. / Slika 1.


Figure 2. / Figure 2. / Abbildung 2. / Figura 2. / Figura 2. / Figura 2. / Figuur 2. / Figur 2. / Kuva 2. / Figur 2. / Figur 2. / Εικόνα 2. / 2 pav. / Rysunek 2. / Рисунок 2. / Obr. 2. / 2. ábra. / Slika 2: / 2. attēls. / Joonis 2. / Фигура 2. / Slika 2. / Figura 2. / Obrázok 2. / Slika 2. / Şekil 2. /

¹ l: Length/ l — Longueur/ l — Länge/ l: longitud/ l — Comprimento/ l — Comprimento/ l: lengte/ l — Lengde/ l — pituus/ l — Längd/ l-Længde/ l — μήκος/ l — ilgis/ l — długość/ l — Длина/ l — Délka/ l — hosszúság/ l — dolžina/ l — garums/ l — pikkus/ l — дължина/ l — Dužina/ l — Lungime/ l — dĺžka/ l = duljina/ l — Uzunluk/

²C- Compliance R_I - Inspiratory limb resistance to flow R_E - Expiratory limb resistance to flow / C- Conformité R_I - Résistance de la branche inspiratoire à l'écoulement R_E - Résistance de la branche / C- Compliance R_I - Strömungswiderstand des Inspirationsschlauches R_E - Strömungswiderstand des / C- cumplimiento normativo; R_I - resistencia de la extremidad inspiratoria al flujo; R_E - resistencia al flujo / C- Complacência R_I - Resistência ao fluxo do ramo inspiratório R_E - Resistência ao fluxo do ramo / C- Conformità. R_I - resistenza della diramazione inspiratoria al flusso. R_E - resistenza della diramazione espiratoria al flusso / C- conformiteit R_I - Weerstand van de inspiratie-arm tegen doorstroming R_E - Weerstand van de expiratie-arm tegen doorstroming / C- Ettergivelighet R_I - Flowmotstand i inspiratorisk gren R_E - Flowmotstand i eksploratorisk gren / C- komplianssi R_I - sisäänhengityshaaran virtausvastus R_E - uloshengityshaaran virtausvastus / C- Etterlevnad R_I - Inandningsrørets flødesmotstand R_E - Utandningsrørets flødesmotstand / C- Overholdelse af R_I - Inspiratorisk slange flow modstand R_E - Eksploratorisk slange flow modstand / C- Συμμόρφωση R_I - Αντίσταση εισπνευστικού σκέλους στη ροή R_E - Αντίσταση εκπνευστικού σκέλους στη ροή / C- atitiktis R_I - įkvėpimo galo atsparumas srautui R_E - iškvėpimo galo atsparumas srautui / C- podatność R_I - opór odnogi wdechowej na przepływ R_E - opór odnogi wydechowej na przepływ / C- Соотвѣствие R_I - Сопротивление сегмента вдоха потока R_E - Сопротивление сегмента выдоха потока / C- soulad R_I - odpor proudění vdechovací části R_E - odpor proudění vdechovací části / C- táglékonyág, R_I - belégzési cső áramlási ellenállása, R_E - kilégzési cső áramlási ellenállása / C- podajnost R_I - upornost na pretok vdihovalnega kraka R_E - upornost na pretok izdihovalnega kraka / C- atbilstība R_I - ieelpas daļas plūsmas pretestība R_E - izelpas daļas plūsmas pretestība / C- vastavus R_I - inspiratsiooniosa resistentsus voolule R_E - ekspiratsiooniosa resistentsus voolule / C- Съответствие с R_I - Съпротивление на потока на разклонението за инспирация / R_E - Съпротивление на потока на разклонението за експирация / C- Usaglašenost R_I - otpornost na protok inspiratorne grane R_E - otpornost na protok ekspiratorne grane / C- Conform R_I - Rezistență la debit pentru tubul de inspirare R_E - Rezistență la debit pentru tubul de expirare / C- podajnost R_I - prietokový odpor inspiračnej vetvy R_E - prietokový odpor expiračnej vetvy / C- sukladnost R_I - otpornost protoka u kraku za udisanje R_E - otpornost protoka u kraku za izdisanje / C- Uygunluk R_I - Nefes alma uzantisinin akış direnci R_E - Nefes verme uzantisinin akış direnci /

Patient category ¹	Intended delivered volume ² ml	Flow resistance limit ³ cmH ₂ O//min	At flow ⁴ l/min
Adult ⁷	≥ 300 ml	0.06	30
Paediatric ⁸	50 ml < 300 ml	0.12	15
Neonatal ⁹	≤ 50 ml	0.74	2.5

Patient category	Intended delivered volume ml	Compliance ⁵ ml/cmH ₂ O	At Pressure ⁶ cmH ₂ O
Adult	≥ 300 ml	5	60 ± 3
Paediatric	50 ml < 300 ml	4	60 ± 3
Neonatal	≤ 50 ml	1.5	60 ± 3

1 Patient category/ Catégorie de patients/ Patientenkategorie/ Categoría de paciente/ Categoria de doentes/ Categorie di pazienti/ Πατιέντκατєrie/ Pasientkategorii/ Potilastyypii/ Patientkategorii/ Κατєγορία ασθενών/ Pacientų kategorija/ Kategoria pacjentów/ Категория пациентов/ Kategorie pacientů/ Betegek köre/ Kategorija bolnika/ Pacientu kategorija/ Patsiendi kateooria/ Категория на пациента/ Kategorija pacijenata/ Categorie de pacienți/ Категорія pacientov/ Kategorija pacijenata/ Hasta kategorisi

2 Intended delivered volume ml/ Volume administré prévu ml/ Zuführendes Volumen ml/ Volumen de suministro previsto ml/ Volume administrado previsto ml/ Volume erogato previsto ml/ Beoogd volume ml/ Tiltentk levert volum ml/ Aiottu antotilavuus ml/ Avsedd levererd volym ml/ Tilsigtet leveret volumen ml/ Προβλεπόμενος χορηγούμενος όγκος ml/ Numatomas pristatoma tūris ml/ Planowana objętość dostarczana ml/ Планируемый дыхательный объем мл/ Předpokládaný přijatý objem ml/ Leadandó légzési térfogat ml/ Predvidena dovedena prostornina ml/ Paredzētais piegādes tilpums ml/ / Kavandatud manus-tatav maht ml/ Планиран осигуряван капацитет ml/ Planirana isporučena zapremina ml/ Volumul estimat livrat ml/ Určený podávaný objem ml/ Predviđeni isporučeni volumen ml/ Amaçlanan uygulama hacmi ml

3 Flow resistance limit cmH₂O//min/ Limite de résistance à l'écoulement cmH₂O//min/ Grenzwert des Strömungswiderstands cmH₂O//min/ Limite de resistencia de flujo cmH₂O//min/ Limite de resistência ao fluxo cmH₂O//min/ Limite della resistenza al flusso cmH₂O//min/ Limiet stromingsweerstand cmH₂O//min/ Grense for flow-motstand cmH₂O//min/ Virtausvastuksen raja cmH₂O//min/ Gräns för flödesmotstånd cmH₂O//min/ Grænse for flow-modstand cmH₂O//min./ Όριο αντίστασης στη ροή cm H₂O//λεπτό/ Pasipriešinio srautui riba cmH₂O//min/ Granica oporu przepływu cmH₂O//min/ Предел сопротивления потоку см H₂O//л/мин/ Limit odporu proudění cmH₂O//min/ Aramlási ellenállás határértéke cmH₂O//perc/ Omejitven odpor-nosti za pretok cmH₂O//l/min/ Plūsmas pretestības robeža cmH₂O//l/min/ Voolutakistuse piirväärtus cmH₂O//min/ Граница на съпротивление на потока cmH₂O//l/min/ Granična vrednost otpornosti protoku cmH₂O//l/min/ Limită de rezistență la curgere cmH₂O//l/min/ Limit odporu voči prúdeniu cmH₂O//l/min/ Granica otpora na protok cmH₂O//l/min/ Akış direnci sınırı cmH₂O// dakika

4 At flow l/min/ Débit l/min/ Bei Durchfluss l/min/ Flujo l/min/ Com o fluxo l/min/ Al flusso l/min/ Bij het stromen l/min/ Ved flow l/ min/ Virtauksella l/min/ Vid flödet l/min/ Ved flow L/min./ Σε ροή l/λεπτό/ Esant srautui l/min/ Przy przepływie l/min/ В потоке л/ мин/ Pri proudění l/min/ Aramlási sebesség l/perc/ Pri pretoku l/min/ Pie aktīvas plūsmas l/min/ Voolutkiirusel l/min/ При поток l/min/ Pri protoku od l/min/ La curgere L/min/ Pri prúdení l/min/ Pri protoku l/min/ Akışta l/dakika

5 Compliance ml/cmH₂O/ Tolérance ml/cmH₂O/ Compliance ml/cmH₂O/ Compliance ml/cmH₂O/ Complacência ml/cmH₂O/ Compliance ml/cmH₂O/ Naleving ml/cmH₂O/ Compliance ml/cmH₂O/ Complianssi ml/cmH₂O/ Compliance ml/cmH₂O/ Överens-stemmelse ml/cmH₂O/ Συμμόρφωση ml/cm H₂O/ Atitiktis ml/cmH₂O/ Podatność ml/cmH₂O/ Растяжимость мл/см H₂O/ Shoda ml/cmH₂O/ Megfelelőség ml/cmH₂O/ Kompliance ml/cmH₂O/ Atbilstība ml/cmH₂O/ Vastāvus ml/cmH₂O/ Elasticност ml/ cmH₂O/ Komplijansa ml/cmH₂O/ Complianță ml/cmH₂O/ Poddajnost ml/cmH₂O/ Rastezljivost ml/cmH₂O/ Uygunluk ml/cmH₂O

6 At Pressure cmH₂O/ Pression cmH₂O/ Bei Druck cmH₂O/ Presión cmH₂O/ Com a pressão cmH₂O/ Alla pressione cmH₂O/ Onder druk cmH₂O/ Ved trykk cmH₂O/ Paineella cmH₂O/ Vid tryck cmH₂O/ Ved tryk cmH₂O/ Σε πίεση cm H₂O/ Esant slėgiui cmH₂O/ Przy ciśnieniu cmH₂O/ Под давлением см H₂O/ Pri tlaku cmH₂O/ Nyomás cmH₂O/ Pri tlaku cmH₂O/ Pie spiediena cmH₂O/ Rõhul cmH₂O/ При налягане cmH₂O/ Pri pritisku od cmH₂O/ La presiune cmH₂O/ Pri tlaku cmH₂O/ Pri tlaku cmH₂O/ Basiñeta cmH₂O

7 Adult/ Adulte/ Für Erwachsene/ Adulto/ Adulto/ Adulto/ Volwassene/ Voksne/ Aikuinen/ Vuxen/ Voksen/ Ενήλικες/ Suau-gusiej/ Dorosły/ Для взрослых/ Dospělý/ Feinótt/ Odrasli/ Pieaugušajiem/ Tāiskasvanu/ Взрастни/ Odrasli/ Adulti/ Dospeli/ Odrasla osoba/ Yetişkin

8 Paediatric/ Pédiatriques/ Pädiastrischen/ Pediátrico/ Pediátricos/ Pediatrico/ Pediatricisch/ Barn/ Lapsille/ Barn/ Pædia-trisk/ Ταϊδιatrikούς/ Vaikai/ Dzieci/ Для педиатрических/ Pediatrické/ Gyermek/ Pediatricni/ Pediatricas/ Pediaatriline/ Педиатрична/ Pedijatrijske/ Copii/ Pediatrických/ Pedijatriju/ Pediatrik

9 Neonatal/ Nouveau-nés/ Neonatalen/ Neonatal/ Recém-nascidos/ Neonatale/ Neonataal/ Nyfødte/ Vastasyntyneille/ Nyföd-da/ Neonatal/ Νεογνά/ Naujagimiai/ Noworodki/ Для новорожденных/ Novorozenci/ Ujszülött/ Neonatalni/ Jaundzimušiem/ Vastsündinute/ Неонатална/ Neonatalne/ Nou-născuți/ Novorodenci/ Novorođenčad/ Yenidoğan

en 22mm Adult Compact™ Breathing System

The instruction for use (IFU) is not supplied individually. Only 1 copy of the instructions for use is provided per box and should therefore be kept in an accessible location for all users. More copies are available upon request.

NOTE: Distribute instructions to all product locations and users. These instructions contain important information for the safe use of the product. Read these instructions for use in their entirety, including warnings and cautions before using this product. Failure to properly follow warnings, cautions and instructions could result in serious injury or death to the patient. Serious incidents to be reported as per the local regulation.

WARNING: A WARNING statement provides important information about a potentially hazardous situation, which, if not avoided could result in death or serious injury.

CAUTION: A CAUTION statement provides important information about a potentially hazardous situation, which if not avoided could result in minor or moderate injury.

Intended Use

To deliver and remove anaesthetic and respiratory gases to and from a patient via a breathing system comprised of tubing and connectors.

Intended user

For use by qualified medical personnel.

Indications

The breathing system should be connected between the ventilator or the anaesthetic machine, and the patient airway connection. The product should only be used on a single patient up to a maximum of 7 days. The breathing system is provided in Compact expandable tubing.

Contraindications

Do not use with equipment that does not comply to international taper connection standards.

Product Use:

The product is intended for use on a single patient for a maximum of 7 days use with a variety of anaesthetic and ventilator equipment meeting international standards. The start of use is defined from the first connection to the equipment. Use in the acute and non-acute environment under the supervision of appropriately qualified medical personnel. The breathing system is safe for use with atmospheric gases: N₂O, O₂ and CO₂; anaesthetic agents: isoflurane, sevoflurane, desflurane and other comparable volatile agents.

Patient Category

The range of 22mm breathing systems are suitable for use on all adult patients with intended delivered tidal volumes of ≥ 300ml.

Connect the breathing system

Attach the breathing system limbs to the inspiratory and expiratory ports of the ventilator or anaesthetic machine using a push and twist action (fig. 1).

Where applicable remove any red dust cap from the Y-Piece of the breathing system.

Fully extend the breathing system and perform pre-use checks as recommended by the equipment manufacturer.

Adjust the length and position of the breathing system to meet the clinical requirements.

Attach the Y-Piece to the patient interface.

Accessories:

Where the 22mm breathing system is supplied with additional components the clinician should assess any potential impact upon the patient's respiratory function, and be aware of the subsequent increased resistance to flow, compressible volume and weight of the product.

When the ventilator and anaesthetic machine parameters are set, any impact that the internal volume and resistance to flow of the device may have upon the effectiveness of the prescribed ventilation, should be assessed on an individual patient basis by the clinician.

Water traps and breathing filters, HMEs, HMEFs if included in the configuration

Follow the instructions in individual IFUs supplied.

Reservoir bags

1. Ensure that the reservoir bag is securely attached to the breathing system via a push and twist action (fig. 1).

2. Consider the impact of increased compliance of the reservoir bag upon prescribed ventilation.

3. The reservoir bag allows manual ventilation, and gives a visual indication of rate and depth of breathing.

Monitoring lines

Monitoring lines provide a leak free connection between the breathing system and a respiratory/anaesthetic gas monitor via a luer lock fitting.

The connection should be checked prior to use. The line should be checked to confirm it is kink free and a patent pathway exists.

Pre-Use Check

1. Prior to installation, check that all breathing system components are free of any obstruction or foreign bodies and that there is a patent airway.

2. Remove any red dust cap before connecting the breathing system between the equipment and the patient.

3. Connect/disconnect the tapered connections of the breathing system to the tapered connection of the ventilator and/or respiratory equipment using a push and twist action (fig. 1).

4. Following installation the breathing system must be checked for leaks and occlusions prior to use and that all connections are secure.

5. When setting ventilator parameters consider any impact that the internal volume of the breathing system and accessories may have upon the effectiveness of the prescribed ventilation. The clinician should be aware that the compressible volume of the breathing system will vary depending upon the position and extension of the system. The clinician should assess this on an individual patient basis.

6. Check that any port cap or accessory is securely attached to prevent

leakage.

In-Use Check

The system should be changed if visibly contaminated or if resistance to flow increases.

Monitor the clinical parameters of the patient during the treatment.

Warning

1. This system is not suitable for use with flammable anaesthetic agents.

2. Monitor the breathing system throughout use and replace if increased resistance or contamination occurs.

3. Refer to equipment manufacturers guidelines for machine settings and pre-use check requirements.

4. Perform a self-test of the equipment after full assembly of the breathing system and associated accessories before use on each patient.

5. Before use, consult the respective documentation and instructions for use of all connected devices, or combination of devices.

6. Do not over extend, cut or reconfigure the tubing.

7. Reuse of a breathing system in anaesthesia, that hasn't used a suitably validated single-use respiratory filter between the product and each patient, poses a threat to patient safety in terms of cross infection, device malfunction and breakage.

8. Ventilator accuracy can be affected by use of a gas driven nebuliser.

9. Adding attachments or other devices to the breathing system can change the pressure gradient across the breathing system and adversely affect the ventilator performance. Adjustment of the length of the Compact breathing system may change compliance and compressible volume of the system and therefore may have an impact upon the prescribed ventilation. Perform a self-test of the ventilator (or anaesthesia machine) after full assembly of the breathing system and all connected devices before use on each patient.

10. Breathing systems, their parts and accessories are validated for use with specific ventilators and anaesthesia machines. Incompatible parts can result in degraded performance. The responsible organization is responsible for ensuring the compatibility of the ventilator or anaesthesia machine and all of the parts used to connect to the patient before use.

11. Use of this breathing system with equipment not compliant with appropriate International Standards can result in degraded performance. Prior to use, consult the respective documentation and instructions for use of all connected devices, or combination of devices to ensure compatibility.

12. The Luer Lock port connector is for the exclusive use of monitoring respiratory and anaesthetic gases. Check that the cap is securely attached to prevent inadvertent leakage when the port is not in use. Do not introduce gases or liquids via the misconnection of the luer lock port as this can result in serious injury to the patient.

13. Consult individual product labelling for presence or absence of latex.

14. Consult individual product labelling for presence or absence of phthalates.

15. If a luer port connector is included in the system, it is for the exclusive use of monitoring respiratory and anaesthetic gases. If applicable, check that the cap is securely attached to prevent inadvertent leakage when the port is not in use. Do not introduce gases or liquids via the misconnection of the luer lock port as this can result in serious injury to the patient.

16. Consult individual product labelling for presence or absence of phthalates.

17. Consult individual product labelling for presence or absence of phthalates.

18. Consult individual product labelling for presence or absence of phthalates.

19. Consult individual product labelling for presence or absence of phthalates.

20. Consult individual product labelling for presence or absence of phthalates.

21. Consult individual product labelling for presence or absence of phthalates.

22. Consult individual product labelling for presence or absence of phthalates.

23. Consult individual product labelling for presence or absence of phthalates.

24. Consult individual product labelling for presence or absence of phthalates.

25. Consult individual product labelling for presence or absence of phthalates.

26. Consult individual product labelling for presence or absence of phthalates.

27. Consult individual product labelling for presence or absence of phthalates.

28. Consult individual product labelling for presence or absence of phthalates.

29. Consult individual product labelling for presence or absence of phthalates.

30. Consult individual product labelling for presence or absence of phthalates.

31. Consult individual product labelling for presence or absence of phthalates.

32. Consult individual product labelling for presence or absence of phthalates.

33. Consult individual product labelling for presence or absence of phthalates.

34. Consult individual product labelling for presence or absence of phthalates.

35. Consult individual product labelling for presence or absence of phthalates.

36. Consult individual product labelling for presence or absence of phthalates.

37. Consult individual product labelling for presence or absence of phthalates.

38. Consult individual product labelling for presence or absence of phthalates.

39. Consult individual product labelling for presence or absence of phthalates.

40. Consult individual product labelling for presence or absence of phthalates.

41. Consult individual product labelling for presence or absence of phthalates.

42. Consult individual product labelling for presence or absence of phthalates.

43. Consult individual product labelling for presence or absence of phthalates.

44. Consult individual product labelling for presence or absence of phthalates.

45. Consult individual product labelling for presence or absence of phthalates.

46. Consult individual product labelling for presence or absence of phthalates.

47. Consult individual product labelling for presence or absence of phthalates.

48. Consult individual product labelling for presence or absence of phthalates.

49. Consult individual product labelling for presence or absence of phthalates.

50. Consult individual product labelling for presence or absence of phthalates.

51. Consult individual product labelling for presence or absence of phthalates.

52. Consult individual product labelling for presence or absence of phthalates.

53. Consult individual product labelling for presence or absence of phthalates.

54. Consult individual product labelling for presence or absence of phthalates.

55. Consult individual product labelling for presence or absence of phthalates.

56. Consult individual product labelling for presence or absence of phthalates.

57. Consult individual product labelling for presence or absence of phthalates.

58. Consult individual product labelling for presence or absence of phthalates.

59. Consult individual product labelling for presence or absence of phthalates.

60. Consult individual product labelling for presence or absence of phthalates.

61. Consult individual product labelling for presence or absence of phthalates.

62. Consult individual product labelling for presence or absence of phthalates.

63. Consult individual product labelling for presence or absence of phthalates.

64. Consult individual product labelling for presence or absence of phthalates.

65. Consult individual product labelling for presence or absence of phthalates.

66. Consult individual product labelling for presence or absence of phthalates.

67. Consult individual product labelling for presence or absence of phthalates.

68. Consult individual product labelling for presence or absence of phthalates.

69. Consult individual product labelling for presence or absence of phthalates.

70. Consult individual product labelling for presence or absence of phthalates.

71. Consult individual product labelling for presence or absence of phthalates.

72. Consult individual product labelling for presence or absence of phthalates.

73. Consult individual product labelling for presence or absence of phthalates.

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Rx ONLY

INTERSURGICAL
COMPLETE RESPIRATORY SYSTEMS

Intersurgical Ltd, Crane House, Molly Millars Lane,
Wokingham, Berkshire, RG41 2RZ, UK
T: +44 (0)118 9656 300 info@intersurgical.com
www.intersurgical.com

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T: +7 495 7716809
info@intersurgical.ru

Česká Republika

T: +420 272 911 814
info@intersurgical.cz

Türkiye

T: +90 216 468 88 28
info@intersurgical.com.tr

South Africa

T: +27 (0)11 444 7968
info@intersurgical.co.za

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T: +86 519 69817000
info@intersurgical-cn.com

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T: 81 (0)3 6863 4359
info@intersurgical.co.jp

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T: +886 4 2380 5430
info@intersurgical.com.tw

Pilipinas

T: +63 2820 4124
info@intersurgical.ph

USA

T: +1 800 828 9633
support@intersurgicalinc.com

Canada

T: +1 905-319-6500
info@intersurgical.ca

Colombia

T: +57 1 742 5161
info@intersurgical.com

Australia

T: +61 (0)2 8048 3300
info@intersurgical.com.au

Ireland

T: +353 (0)1 697 8540
info@intersurgical.ie

France

T: +33 (0)1 48 76 72 30
info@intersurgical.fr

Deutschland

T: +49 (0)2241 25690
info@intersurgical.de

España

T: +34 91 665 73 15
info@intersurgical-es.com

Portugal

T: +351 219 108 550
info@intersurgical.pt

Italia

T: +39 0535 20836
info@intersurgical.it

Nederland

T: +31 (0)413 243 860
info@intersurgical.nl

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T: +32 5 2456991
info@intersurgical.be

Luxembourg

T: +32 5 2456991
info@intersurgical.be

Sverige

T: +46 08 514 30 600
info@intersurgical.se

Danmark



T: +45 28 69 80 00
info@intersurgical.dk

Lietuva

T: +370 387 66611
info@intersurgical.lt

WATER TRAP, SELF SEALING, MEDIUM VOLUME

CE
1639 Rx ONLY

INTERSURGICAL
COMPLETE RESPIRATORY SYSTEMS

Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK

T: +44 (0)118 9656 300 info@intersurgical.com www.intersurgical.com

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Česká Republika
T: +420 272 911 814
info@intersurgical.cz

Türkiye
T: +90 216 468 88 28
info@intersurgical.com.tr

South Africa
T: +27 (0)11 444 7968
info@intersurgical.co.za

中国
T: +86 519 6981 7000
info@intersurgical-cn.com

日本
T: 81 (0)3 6863 4359
info@intersurgical.co.jp

台灣
T: +886 4 2380 5430
info@intersurgical.tw

Filipinas
T: +63 2820 4124
info@intersurgical.ph

USA
T: +1 800 828 9633
support@intersurgicalinc.com

Canada
T: +1 905-319-6500
info@intersurgical.ca

Colombia
T: +57 1 742 5161
info@intersurgical.com

Australia
T: +61 (0)2 8048 3300
info@intersurgical.com.au

Distributed in the USA by Intersurgical Incorporated, 6757 Kinne St., East Syracuse, NY 13057
T: (800) 828-9633 F: (315) 451-3696 support@intersurgicalinc.com www.intersurgical.com

7705 Issue 9 RO 08.21