
The Drum™

Instructions For Use

en Instruction for use: The Intersurgical Drum™ CO₂ absorbent cartridge

These instructions must be carefully read before using the Intersurgical Drum Cartridge. Also refer to the instructions of use provided by the anaesthetic machine manufacturer.

WARNING: A WARNING statement provides important information about a potentially hazardous situation which, if not avoided could result in death or serious injury.

CAUTION: A CAUTION statement provides important information about a potentially hazardous situation which, if not avoided could result in minor or moderate injury.

TECHNICAL DATA

Product Code	Absorbent Volume	Resistance to flow of canister at 60 lpm	Casing	Typical pH	Absorbent details
2186000	1.2 L / 1.0 kg / 2.2 lb	1.0 mbar	Acrylonitrile Butadiene Styrene (ABS)	Used below 10 Unused 12 – 14 depending on absorbent	For full details please contact Intersurgical
2187000	1.2 L / 1.0 kg / 2.2 lb				
2188000	1.2 L / 0.9 kg / 2.0 lb				

INTENDED USE

To remove carbon dioxide from anaesthetic and respiratory gases delivered to a patient. The Intersurgical Drum CO₂ absorbent cartridge, is for the absorption of carbon dioxide within an anaesthetic breathing system.

INTENDED USER

For use by qualified medical personnel.

The Intersurgical Drum cartridge is compatible with numerous anaesthetic machines and absorbers, single or double, that are designed to house this standard cartridge design. These include:

- GE Healthcare® Aestiva® 5 machine
- GMS™ absorber
- Spacelabs® Focus™ and Sirius™ machines and CAS™ absorber
- Penlon® A100™, SP100™, SP200™ machines and AP100™ absorber
- Mindray® A4, A5, A7, WATO EX65 machines
- Any other machine that uses an absorber designed to house this standard cartridge design

CHOICES AVAILABLE

Product codes:

2186000 The Drum, Spherasorb™ CO₂ absorbent cartridge, white to violet colour change, 1kg

2187000 The Drum, Spherasorb™ CO₂ absorbent cartridge, pink to white colour change, 1 kg

2188000 The Drum, LoFloSorb™ CO₂ absorbent cartridge, green to violet colour change, 1kg

For full details of Spherasorb and LoFloSorb please contact Intersurgical or refer to the Intersurgical Absorbents catalogue.

PRINCIPLE OF OPERATION

The Drum cartridge is designed to fit inside reusable absorbers that have been designed for this type of cartridge.

The absorbent within the Drum cartridge undergoes a chemical reaction with the CO₂ within the breathing system. The hydroxides within the absorbent are converted to carbonates and water in an exothermic process. The CO₂ is therefore, bound within the absorbent.

WARNING If the moisture content of the absorbent is outside that set during the manufacturing process, this will affect the absorbent performance.

For more information on this issue please see the Safety Information section or contact Intersurgical.

FITTING THE INTERSURGICAL DRUM CARTRIDGE INTO THE REUSABLE ABSORBER

WARNING Do not use the Drum cartridge if it is cracked or damaged, or if the bag it is contained in is unsealed or torn. Also do not use if there is any unusual discolouration of the granules.

TO FIT THE DRUM CARTRIDGE INTO PLACE

WARNINGS

1. Always insert the Drum cartridge into the absorber in the correct orientation, with the arrows pointing upwards.
2. Remove the Drum from the laminate bag completely before use.

CAUTIONS

1. Always insert the Drum cartridge following instructions of the anaesthetic machine manufacturer (double absorbers will contain one cartridge in each absorber chamber).
2. Before use, ensure that the granules inside the Drum are settled and level, this is to provide optimum performance. This can be achieved by a few gentle taps.
3. When the Drum cartridge is contained within the absorber, the absorber must complete and pass a pre-use leakage test before use. This must be done in accordance with the instruction for use of the anaesthetic machine being used.
4. The Drum cartridge must be correctly fitted inside the absorber, before starting these tests. If a double absorber is being used, a Drum cartridge must be inserted in each chamber.

REMOVING THE DRUM CARTRIDGE

CAUTIONS

1. When exhausted, the Drum cartridge must be replaced to ensure that CO₂ absorption in the breathing system is continued.
2. Always remove the Drum cartridge from the absorber, following the instructions for use of the anaesthetic machine manufacturer.

WHEN TO REPLACE THE DRUM CARTRIDGE

WARNING

- A new filter or HMEF should be used on the patient Y piece of the breathing system for every patient. When this is done, the Drum

cartridge is not restricted to single patient use and can be used until exhaustion is identified. It is advised that a filter or HMEF is also placed at the expiratory port of the machine.

CAUTIONS

1. If the Drum cartridge has been attached to the machine for one month, it is recommended to replace it regardless of extent of usage.
2. The degree of colour change at exhaustion depends on the conditions of use and is an indication only. The Drum cartridge must be used in conjunction with inspiratory FiCO_2 monitoring. The Drum cartridge needs replacing when the inspiratory FiCO_2 concentration exceeds 5mmHg (approximately 0.5% volume FiCO_2).

NOTE: With some absorbents, the colour indicator can sometimes regenerate slowly in between uses and the absorbent may revert back to its original colour. However, this is not a regeneration of the absorption capacity.

Therefore, the Intersurgical Drum cartridge should be discarded in the first instance that exhaustion has been determined. Always refer to and follow the Instructions For Use of the anaesthetic machine being used.

SAFETY INFORMATION

The Drum cartridge may be used with oxygen, nitrous oxide, Xenon, Enflurane, Halothane, Isoflurane, Desflurane and Sevoflurane.

CAUTIONS

1. Before using the Drum cartridge it is essential to become familiar with and follow the Instructions For Use of the anaesthetic machine on which it is to be used.
2. The Drum cartridge must be used in conjunction with monitoring of FiCO_2 inspiratory carbon dioxide.
3. Once exhausted, the Drum cartridge cannot be reused. Do not leave within the machine once exhaustion is identified.

WARNINGS

1. Do not allow contact with alcohol based cleaning products.
2. Once the Drum cartridge is inserted within the absorber, always conduct a pre-use leakage test, according to the Instructions For Use of the anaesthetic machine on which it is to be used.
3. Do not use the Drum cartridge if the machine pre-use leakage test fails.
4. Do not use with trichloroethylene or chloroform because toxic products may be produced by chemical reactions.
5. The carbon dioxide absorbent inside the Drum cartridge has a moisture content set in the manufacturing process. Do not flush the Drum cartridge with dry gas for long periods. This may result in loss of moisture and reduced CO_2 absorption efficiency. Excessive loss of moisture can lead to failure to absorb CO_2 and to clinical problems.
6. Do not sterilize, disinfect or clean the Drum cartridge before or during use. Cleaning can lead to device malfunction and breakage.
7. Do not use the Drum cartridge after expiry date.
8. Never allow pure carbon dioxide to pass through any Carbon Dioxide Absorbents or disposable absorbers, as this will generate high temperatures.
9. Do not use if the Drum cartridge becomes excessively hot.
10. Do not use if it is believed that dry gas has been left running during long periods of non-use at high flows.
11. Do not use if there is a delay in the increase of inspired agent concentration.
12. It is recommended to use water traps in the patient circuit during low flow anaesthesia (fresh gas flows of 1 L/min or less) and when excessive condensed moisture accumulates in the limbs of the breathing system. Use water traps in both the expiratory and inspiratory limbs to remove water condensation.

For more information please refer to the Intersurgical Absorbents Catalogue or contact Intersurgical.

RISK AND SAFETY PHRASES RELATING TO THE ABSORBENTS WITHIN THE DRUM CARTRIDGE

H315: Causes skin irritation.

H318 Causes serious eye damage.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P302/P352 IF ON SKIN: Wash with plenty of soap and water.

P305/351/338. IF IN EYES: Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do.

CONTINUE RINSING

P332/313: If skin irritation occurs: Get medical advice/attention.

STORAGE

The Drum cartridge must be stored, sealed, in a clean, dry environment away from direct sunlight and not under heavy load.

Recommended storage at room temperature over the indicated product shelf-life.

DISPOSAL

Intersurgical advises that after clinical use, the waste Drum cartridge is treated at source with minimal handling; by incineration within a mixture of other clinical waste. In any case, the Drum cartridge must be disposed of according to local regulations. Refer to the MSDS for further details.

DATA SHEET

The Material Safety Data Sheet is available on request.

GE Healthcare is a registered trademark of General Electric Company. Aestiva is a registered trademark of Datex-Ohmeda, Inc. GMS is a trademark of GE Healthcare. Spacelabs is a registered trademark of Spacelabs Healthcare, LLC. Focus, Sirius and CAS are trademarks of Spacelabs Healthcare, LLC. Penlon is a registered trademark of Penlon Limited. A100, SP100, SP200 and AP100 are trademarks of Penlon Limited. Mindray is a registered trademark of Shenzhen Mindray Bio-Medical Electronics Co., Ltd. Datex-Ohmeda, Inc. is a General Electric Company.

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