
IS Can™

Instructions For Use

en Instruction for use: The Intersurgical IS Can™ disposable CO₂ absorber.

These instructions must be carefully read before using the Intersurgical IS Can absorber. Also refer to the instructions of use provided by the anaesthetic machine manufacturer.

WARNING: A WARNING statement provides important information about a potentially hazardous situation which, if not avoided could result in death or serious injury.

CAUTION: A CAUTION statement provides important information about a potentially hazardous situation which, if not avoided could result in minor or moderate injury.

TECHNICAL DATA

Product Code	Absorbent Volume	Resistance to flow of canister at 60 L/Min	Casing	Dust filters	Typical pH	Absorbent details
2196000	1.2 L / 1.0 kg / 2.2 lb	2.0 cm H ₂ O	Acrylonitrile Butadiene Styrene (ABS)	Polyester	Used below 10 Unused 12 – 14 depending on absorbent	For full details please contact Intersurgical
2197000	1.2 L / 1.0 kg / 2.2 lb					
2198000	1.2 L / 0.9 kg / 2.0 lb					

INTENDED USE

To remove carbon dioxide from anaesthetic and respiratory gases delivered to a patient.

The IS Can absorber is compatible with the GE Healthcare® Aespire®Avance®Aisys® machine ranges and the GE Healthcare® Compact Block ADU™.

INTENDED USER

For use by qualified medical personnel.

CHOICES AVAILABLE

Product codes:

2196000 IS Can, Spherasorb™ disposable CO₂ absorber, white to violet colour change

2197000 (Not available in the US market) IS Can, Spherasorb™ disposable CO₂ absorber, pink to white colour change

2198000 IS Can, LoFloSorb™ disposable CO₂ absorber, green to violet colour change

For full details of Spherasorb and LoFloSorb please contact Intersurgical or refer to the Intersurgical catalogue.

PRINCIPLE OF OPERATION


The absorbent within the IS Can absorber undergoes a chemical reaction with the CO₂ within the breathing system. The hydroxides within the absorbent are converted to carbonates and water in an exothermic process. The CO₂ is therefore, bound within the absorbent.

WARNING If the moisture content of the absorbent is outside of the moisture set during the manufacturing process, this will affect CO₂ absorption performance. For more information on this issue please refer to warnings under the Safety Information section or contact Intersurgical.


FITTING THE IS CAN

Fitting to the Aespire®Avance®Aisys® machine ranges (Figure 1)

Fitting to the Compact Block ADU™ (Figure 2)

(Figure 1)

(Figure 2)


Before fitting ensure that the granules inside the IS Can are settled and level by gently tapping the canister. This is to ensure optimum performance.

WARNING: Do not use the IS Can if the twin seal cap is missing, if the canister is cracked or damaged or if the absorbent granules are showing any color change. **WARNING:** The IS Can must be properly fitted to the anesthetic machine prior to performing the pre-use leak test. Once fitted and prior to use, conduct a system leak test in accordance to the instructions provided by the anaesthetic machine manufacturer. **WARNING:** Do not use if this leak test fails. Please note that failure of the pre-use leakage test may be due to leakage elsewhere in the machine or breathing system. In the event of of repeated failure of the pre-use leak test, examine other possible causes for leakage within the anesthetic machine or breathing system.


REMOVING THE IS CAN

Removing the IS Can from the Aespire®Avance®Aisys® machine ranges (Figure 3)

Removing the IS Can from the Compact Block ADU™ (Figure 4)

(Figure 3)

(Figure 4)


CAUTION When exhausted, the IS Can must be removed and replaced to ensure that CO₂ removal from the breathing system is continued. See Figure 3 and 4 for disconnection.

- Remove the twin seal cap from a new Intersurgical IS Can and attach as shown in Figure 1. Always ensure that the pre-use leakage test is successfully completed in accordance with the anaesthetic machine instructions for use using the IS Can.

WARNINGS

1. DO NOT use the IS Can unless the pre-use leakage test passes.
2. If the IS Can is not properly attached to the anaesthetic machine, there is no CO₂ absorption occurring.
3. Terminate use immediately and replace the IS Can if at any time during use, the anaesthetic machine indicates expiratory reverse flow message.

When to replace the IS Can absorber (Figure 5).

(Figure 5)

- A new filter or HMEF should be used on the patient Y piece of the breathing system for every patient. When this is done, the IS Can is not restricted to single patient use and can be used until exhaustion is identified. It is advised that a filter or HMEF is also placed at the expiratory port of the machine.

CAUTIONS

1. If the IS Can has been attached to the machine for one month, it is recommended to replace it regardless of extent of usage.
2. The degree of colour change at exhaustion depends on the conditions of use and is an indication only. The IS Can absorber must be used in conjunction with inspiratory FiCO₂ monitoring. The IS Can absorber needs replacing when the inspiratory FiCO₂ concentration exceeds 5mmHg (approximately 0.5% volume FiCO₂).

CAUTION With some absorbents, the colour indicator can sometimes regenerate slowly in between uses and the absorbent may revert back to its original colour. However, this is not a regeneration of the absorption capacity. Therefore, the Intersurgical IS Can should be discarded in the first instance that exhaustion has been determined.

Always refer to the anesthetic machine instructions for use.

SAFETY INFORMATION

The IS Can absorber may be used with oxygen, nitrous oxide, Xenon, Enflurane, Halothane, Isoflurane, Desflurane and Sevoflurane.

CAUTIONS

1. Before using the IS Can it is essential to become familiar with and follow the Instructions For Use of the anaesthetic machine on which it is to be used.
2. The IS Can absorber must be used in conjunction with monitoring of FiCO₂ inspiratory carbon dioxide.
3. Once exhausted, the IS Can cannot be reused. Do not leave attached to the machine once exhaustion is identified.

WARNINGS

1. Do not allow contact with alcohol based cleaning products.
2. Always conduct a pre-use leakage test, according to anaesthetic machine instructions for use.
3. Do not use the IS Can if the pre-use leakage test fails.
4. Do not use with trichloroethylene or chloroform because toxic products may be produced by chemical reactions.
5. Do not use if the scrims/pads are missing or out of place or if granules or dust have escaped beyond these scrims.
6. The carbon dioxide absorbent inside the IS Can absorber has a moisture content set in the manufacturing process. Do not flush the IS Can absorber with dry gas for long periods. This may result in loss of moisture and reduced CO₂ absorption efficiency. Excessive loss of moisture can lead to failure to absorb CO₂ and to clinical problems.
7. Do not sterilize, disinfect or clean the IS Can absorber before or during use. Cleaning can lead to device malfunction and breakage.
8. Do not use the IS Can after expiry date.
9. Never allow pure carbon dioxide to pass through any Carbon Dioxide Absorbents or disposable absorbers, as this will generate high temperatures.
10. Do not use if the IS Can absorber becomes excessively hot.
11. Do not use if it is believed that dry gas has been left running during long periods of non-use at high flows.
12. Do not use if there is a delay in the increase of inspired agent concentration.
13. It is recommended to use water traps in the patient circuit during low flow anaesthesia (fresh gas flows of 1 L/min or less) and when excessive condensed moisture accumulates in the limbs of the breathing system. Use water traps in both the expiratory and inspiratory limbs to remove water condensation.

For more information please refer to the Intersurgical Catalogue or contact Intersurgical.

RISK AND SAFETY PHRASES RELATING TO THE ABSORBENTS WITHIN THE IS CAN ABSORBER.

H315: Causes skin irritation.

H318 Causes serious eye damage.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P302/P352 IF ON SKIN: Wash with plenty of soap and water.

P305/351/338. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P332/313: If skin irritation occurs: Get medical advice/attention.

STORAGE

The IS Can absorber must be stored, sealed, in a clean and dry environment away from direct sunlight and not under heavy load. Recommended storage at room temperature throughout the indicated product shelf-life.

DISPOSAL

Intersurgical advises that after clinical use, the waste IS Can is treated at source with minimal handling; this being by incineration within a mixture of other clinical waste. However, notwithstanding this, the IS Can must be disposed of according to local regulations. Refer to the MSDS for further details.

DATA SHEET

The Material Safety Data Sheet is available on request.

*GE Healthcare is a registered trademark of General Electric Company

Aisys®, Aespire®, Avance® are registered trademarks of Datex-Ohmeda Inc.

ADU™ is a trademark of Datex-Ohmeda Inc.

Datex-Ohmeda, Inc. is a General Electric Company.

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Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK