
The AbCan™

Instructions For Use

en

Instruction for use: The Intersurgical AbCan™ disposable CO₂ absorber.

These instructions must be carefully read before using the Intersurgical AbCan disposable CO₂ absorber. Also refer to the instructions of use provided by the anaesthetic machine manufacturer.

The instruction for use is not supplied individually. Only 1 copy of the instructions for use is provided per box and should therefore be kept in an accessible location for all users. More copies are available upon request.

NOTE: Distribute instructions to all product locations and users. These instructions contain important information for the safe use of the product. Read these instructions for use in their entirety, including warnings and cautions before using this product. Failure to properly follow warnings, cautions and instructions could result in serious injury or death to the patient. Serious incidents to be reported as per the local regulation.

WARNING: A WARNING statement provides important information about a potentially hazardous situation which, if not avoided could result in death or serious injury.

CAUTION: A CAUTION statement provides important information about a potentially hazardous situation which, if not avoided could result in minor or moderate injury.

Technical data

Product code	Absorbent Volume/Weight			
2199001	1.1 L/ 0.9 kg / 2.0 lb			
2199002 (Not available in the US market)	1.1 L/ 0.9 kg / 2.0 lb			
2199003	1.1 L/ 0.8 kg / 1.8 lb			
Resistance to flow of canister at 60 L/min	Casing	Dust filters	Typical pH	Absorbent details
Maximum 2.5 cm H ₂ O	Acrylonitrile Butadiene Styrene (ABS)	Polyester Polypropylene	Used below 10 Unused 12-14 depending on absorbent	For full details please contact Intersurgical

Intended use:

To remove carbon dioxide from anaesthetic and respiratory gases delivered to a patient.

The Intersurgical AbCan disposable CO₂ absorber, is for the absorption of carbon dioxide within an anaesthetic breathing system.

The AbCan is compatible with the Getinge® (Maquet®) Anaesthesia Delivery Systems, Flow-i™ C20, C30, C40, Flow-c™, Flow-e™.

Choices available:

Product codes:

2199001 AbCan, Spherasorb™ disposable CO₂ absorber, white to violet colour change, 1.1L

2199002 (Not available in the US market) AbCan, Spherasorb™ disposable CO₂ absorber, pink to white colour change, 1.1L

2199003 AbCan, LoFloSorb™ disposable CO₂ absorber, green to violet colour change, 1.1L

For full details of Spherasorb and LoFloSorb please contact Intersurgical or refer to the Intersurgical catalogue.

Principle of operation

The absorbent within the AbCan undergoes a chemical reaction with the CO₂ within the breathing system. The hydroxides within the absorbent are converted to carbonates and water in an exothermic process. The CO₂ is therefore, bound within the absorbent.

WARNING If the moisture content of the absorbent is outside that set during the manufacturing process, this will affect the absorbing performance. For more information on this issue please see the Safety Information section or contact Intersurgical.

Fitting and replacing the Intersurgical AbCan to the Getinge® (Maquet®)

Anaesthesia Delivery Systems, Flow-i™ C20, C30, C40, Flow-c™, Flow-e™.

Before attaching the AbCan, it is beneficial to ensure that the absorbent granules are well settled and level. This can be achieved by tapping gently on the side of the AbCan.

WARNING Do not use the AbCan if the absorber is cracked or damaged, or if the green sealing cap is missing. Also do not use if there is any unusual discolouration of the granules.

CAUTION The AbCan must be completely and pass the machine pre-use leakage and compliance tests before use. This must be done in accordance with the instruction for use of the Getinge® (Maquet®)

Anaesthesia Delivery Systems, Flow-i™ C20, C30, C40, Flow-c™, Flow-e™.

CAUTION The AbCan must be correctly fitted to the Getinge® (Maquet®) Anaesthesia Delivery Systems, Flow-i™ C20, C30, C40, Flow-c™, Flow-e™ before starting these tests.

To fit the AbCan into place:

Remove the green sealing cap from the top of the new AbCan absorber.

Complete the following steps, 1 to 5.

Step 1: Turn the locking switch to the unlocked position, as shown in Figure 1. The AbCan can now be fitted to the machine.

Step 2: Ensure that the silicon rubber seals on the machine ports are not damaged. With the AbCan absorber ports facing towards the machine, present the AbCan to the corresponding machine ports. Lift the AbCan absorber ports up into the machine ports and push the base forwards.

WARNING If you encounter resistance when attempting to fit the AbCan, do not use force as this may damage the equipment. Instead, check the alignment of the ports and the position of the AbCan as shown in Figure 2.

Step 3: Turn the locking switch so it is pointing in the upward locked position as shown in Figure 2. The AbCan is now raised into the locked position and is ready for the pre-use leakage and compliance test.

Step 4: Conduct the machine pre-use leakage and compliance tests, according to the Instructions For Use of the Getinge® (Maquet®) Anaesthesia Delivery Systems, Flow-i™ C20, C30, C40, Flow-c™, Flow-e™.

WARNING DO NOT use the AbCan if either the pre-use leakage or compliance tests fail or if there is any pre-use test failure relating to occlusion.

CAUTION Bypass position is not applicable to all machines.

Step 5: If either of the pre-use leakage or compliance tests fail, repeat steps 1 to 4.

Always also refer to and follow the Instructions For Use of the Getinge® (Maquet®) Anaesthesia Delivery Systems, Flow-i™ C20, C30, C40, Flow-c™, Flow-e™.

WARNING DO NOT use the AbCan until both the pre-use leakage and compliance tests are successful. Please note that failure of the pre-use leakage and compliance test may be due to leakage elsewhere in the machine or breathing system. In the event of repeated failure of the pre-use test, also examine other


Figure 1. Unlocked positions


Figure 2. Locked positions

possible causes for leakage in addition to the AbCan.

WARNING When using the AbCan absorber with Flow-i™ machines: C20, C30, C40, Flow-c™, Flow-e™, to avoid the need for frequent draining of water, it is recommended to use a breathing system that includes a water trap on the inspiratory limb. If using a filter on the inspiratory port, to avoid increased resistance to flow, place an additional water trap 1915000 before the filter (Figure 3). To allow easier positioning of the AbCan canister it is recommended to use a spacer, such as 1961000 straight 22M-22F between the inspiratory port and the water trap. Use the filter according to its intended use, and replace as indicated in the filter IFU.

CAUTION It may be necessary to slightly rotate the 1915000 water trap when removing and reinserting the AbCan. However, do not rotate this water trap by more than 45 degrees.

CAUTION The water trap should be emptied as required and in accordance with local hospital infection control recommendations.

CAUTION Use "push and twist" action (Figure 4) during connecting and disconnecting system components.

CAUTION: Laparoscopic procedures use Carbon Dioxide as insufflation gas, resulting in increased levels of end tidal Carbon Dioxide. The absorbent material will be more active, increasing gas humidity. As a result, higher levels of condensate may be experienced in the breathing system and/or breathing filters. Monitor the filters throughout use and change if increased resistance and/or excessive condensate is observed.

For more information please refer to the Intersurgical Catalogue or contact Intersurgical.

Replacing the AbCan:

CAUTION When exhausted, the AbCan absorber must be removed and replaced to ensure that CO₂ absorption in the breathing system is continued. Remove the seal from a new AbCan and complete steps 1 to 5 as above, including the pre-use leakage and compliance tests according to the Instructions For Use of the Getinge® (Maquet®) Anaesthesia Delivery Systems, Flow-i™ C20, C30, C40, Flow-c™, Flow-e™.

CAUTION: CO₂ is not being absorbed when the AbCan absorber is not attached to the Anaesthesia Delivery System.

CAUTION If the AbCan has been attached to the machine for one month, it is recommended to replace it regardless of extent of usage.

CAUTION The degree of colour change at exhaustion depends on the conditions of use and is an indication only. The AbCan must be used in conjunction with inspiratory FiCO₂ monitoring. The AbCan needs replacing when the inspiratory FiCO₂ concentration exceeds 5mmHg (approximately 0.5% volume FiCO₂).

WARNING: If the Flow-i™ machine has a Bypass position (B) on the locking switch, the AbCan will not absorb patient CO₂ with the locking switch in this position.

A new filter or HMEF should be used on the patient Y piece of the breathing system for every patient. When this is done, the AbCan is not restricted to single patient use and can be used until exhaustion is identified. It is advised that a filter is also placed at the expiratory port of the patient cassette.

NOTE: With some absorbents, the colour indicator can sometimes regenerate slowly inbetween uses and the absorbent may revert back to its original colour. However, this is not a regeneration of the absorption capacity. Therefore, the Intersurgical AbCan should be discarded in the first instance that exhaustion has been determined.

Always refer to and follow the Instructions For Use of the Getinge® (Maquet®) Anaesthesia Delivery Systems, Flow-i™ C20, C30, C40, Flow-c™, Flow-e™.

Safety information

NOTE The AbCan must be used in conjunction with monitoring of FiCO₂, inspiratory carbon dioxide.

The AbCan may be used with oxygen, nitrous oxide, Xenon, Enflurane, Halothane, Isoflurane, Desflurane and Sevoflurane.

CAUTION Before using the AbCan it is essential to become familiar with and follow the Instructions For Use of the Getinge® (Maquet®) Anaesthesia Delivery Systems, Flow-i™ C20, C30, C40, Flow-c™, Flow-e™.

CAUTION Once exhausted, the AbCan cannot be reused. Do not leave attached to the machine once exhaustion is identified.

WARNING Do not attempt to remove or attach the AbCan, unless the locking switch is the unlocked position.

WARNING Do not allow contact with alcohol based cleaning products.

WARNING Always conduct a pre-use leak and compliance test, according to the Instructions For Use of the Getinge® (Maquet®) Anaesthesia Delivery Systems, Flow-i™ C20, C30, C40, Flow-c™, Flow-e™. Do not use the AbCan if these pre-use tests fail.

WARNING Do not use with trichloroethylene or chloroform because toxic products may be produced by chemical reactions.

WARNING Do not use if the foam pads or scrims are missing or out of place or if granules or dust have escaped beyond these foam pads.

WARNING The carbon dioxide absorbent inside the AbCan has a moisture content set in the manufacturing process. Do not flush the AbCan with dry gas for long periods. This may result in loss of moisture and reduced CO₂ absorption efficiency. Excessive loss of moisture can lead to failure to absorb CO₂, and to clinical problems.

WARNING Do not sterilize, disinfect or clean the absorber before or during use. Cleaning can lead to device malfunction and breakage.

WARNING Do not use after expiry date.

WARNING Never allow pure carbon dioxide to pass through any Carbon Dioxide Absorbents or disposable absorbers, as this will generate high temperatures.

WARNING Do not use if the absorber becomes excessively hot.

WARNING Do not use if it is believed that dry gas has been left running during long periods of non-use at high flows.

WARNING Do not use if there is a delay in the increase of inspired agent concentration.

Risk and safety phrases relating to the absorbents within the AbCan absorber.

H315: Causes skin irritation.

H318 Causes serious eye damage.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P301/P352 IF ON SKIN: Wash with plenty of soap and water.

P305/351/338. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P332/313: If skin irritation occurs: Get medical advice/ attention.

Storage

The AbCan absorber must be stored, sealed, in a clean and dry environment away from direct sunlight and not under heavy load.

Recommended storage at room temperature during the indicated product shelf-life.

Disposal

Intersurgical advises that after clinical use, the waste AbCan is treated at source with minimal handling by incineration within a mixture of other clinical waste. In any case, the AbCan must be disposed of according to local regulations. Refer to the MSDS for further details.

Standards

This product complies with relevant requirements of the following standard: BS EN ISO 80601-2-13: Particular requirements for basic safety and essential performance of an anaesthetic workstation.

Data sheet

The Safety Data Sheet is available on request.


Figure 3


Figure 4

*Getinge is a registered trademark of Getinge AB. Maquet is a registered trademark of Maquet GmbH. Flow-i, Flow-c, Flow-e are trademarks of Maquet Critical Care AB.

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