O₂/N₂O Breathing Systems

en	O ₂ /N ₂ O BREATHING SYSTEMS
fr	O ₂ /N ₂ O CIRCUIT DE RESPIRATION
de	O ₂ /N ₂ O BEATMUNGSSYSTEM
es	O ₂ /N ₂ O CIRCUITO RESPIRATORIO
pt	O ₂ /N ₂ O SISTEMA RESPRATÓRIO
it	O ₂ /N ₂ O CIRCUITO RESPIRATORIO
nl	O ₂ /N ₂ O BEADEMINGSSYSTEEM
no	O ₂ /N ₂ O RESPIRATOR SLANGESYSTEM
fi	O ₂ /N ₂ O HENGITYSLETKUSTO
sv	O ₂ /N ₂ O ANDNINGSSYSTEM
da	O ₂ /N ₂ O RESPIRATIONSSYSTEM
el	O ₂ /N ₂ O ΑΝΑΠΝΕΥΣΤΙΚΟ ΣΥΣΤΗΜΑ
lt	O ₂ /N ₂ O KŲĖPAVIMO SISTEMA
pl	O ₂ /N ₂ O SYSTEM ODDECHOWY
ru	O ₂ /N ₂ O ДЫХАТЕЛЬНЫЙ КОНТУП
cs	O ₂ /N ₂ O DÝCHACÍ SYSTÉM
hu	O ₂ /N ₂ O LÉLEGEZTETŐ RENDSZER
sl	O ₂ /N ₂ O DIHALNI SISTEM
lv	O ₂ /N ₂ O ELPOŠANAS SISTĒMA
et	O ₂ /N ₂ O HINGAMISKONTUUR
bg	O ₂ /N ₂ O РЕСПИРАТОРНА СИСТЕМА
sr	O ₂ /N ₂ O SISTEM ZA DISANJE
ro	O ₂ /N ₂ O SISTEM DE VENTILAȚIE
sk	O ₂ /N ₂ O DÝCHACÍ SYSTÉM
hr	O ₂ /N ₂ O SUSTAV ZA DISANJE
tr	O ₂ /N ₂ O SOLUNUM SİSTEMİ


Figure 1 / Figure 1. / Abbildung 1. / Figura 1. / Figur 1. / Kuva 1. / Εικόνα 1. / 1 pav. / Rysunek 1. / Рис. 1. / Obrázek 1. / 1. ábra / Slika 1 / 1. attēls. / Joonis 1. / Фигура 1. / Slika 1. / Obrázok 1. / Şekil 1


Figure 2. / Figure 2. / Abbildung 2. / Figura 2. / Figur 2. / Kuva 2. / Εικόνα 2. / 2 pav. / Rysunek 2. / Рис. 2. / Obrázek 2. / 2. ábra / Slika 2 / 2. attēls. / Joonis 2. / Фигура 2. / Slika 2. / Obrázok 2. / Şekil 2

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en O₂/N₂O BREATHING SYSTEMS

The instruction for use (IFU) is not supplied individually. Only 1 copy of the instructions for use is provided per box and should therefore be kept in an accessible location for all users. More copies are available upon request.

NOTE: Distribute instructions to all product locations and users. These instructions contain important information for the safe use of the product. Read these instructions for use in their entirety, including warnings and cautions before using this product. Failure to properly follow warnings, cautions and instructions could result in serious injury or death to the patient. Serious incidents to be reported as per the local regulation.

WARNING: A WARNING statement provides important information about a potentially hazardous situation, which, if not avoided could result in death or serious injury.

CAUTION: A CAUTION statement provides important information about a potentially hazardous situation, which if not avoided could result in minor or moderate injury.

INTENDED USE

To deliver and remove anaesthetic and respiratory gases to and from a patient via a breathing system comprised of tubing and connectors.

INTENDED USER

For use by qualified medical personnel.

INDICATIONS

The breathing system should be connected between the anaesthetic gas source and the patient connection. The product can be used on a single patient up to a maximum of 7 days. The breathing system may be provided in Flextube™ or Compact™ tubing.

CONTRAINDICATIONS

Do not use with equipment that does not comply to international taper connection standards.

PRODUCT USE

The product is intended for a maximum of 7 days use with a variety of anaesthetic equipment meeting international standards. The start of use is defined from the first connection to the equipment. Use in the acute and non-acute environment under the supervision of suitably or appropriately qualified medical personnel. The breathing system is safe for use with atmospheric gases: N₂O and O₂.

PATIENT CATEGORY

The range of O₂/N₂O breathing systems are suitable for use on all patients with intended delivered tidal volumes of >50ml.

CONNECT THE BREATHING SYSTEM

Attach the breathing system to the analgesic gas source using a push and twist action (fig.1).

Remove any red dust cap from the patient connection of the breathing system.

Perform pre-use checks as recommended by the equipment manufacturer.

Attach the breathing system to the patient interface.

ACCESSORIES

Where the breathing system is supplied with additional components the clinician should assess any potential impact upon the patient's respiratory function, and be aware of the subsequent increased resistance to flow, compressible volume and weight of the product.

Product configuration can be identified on the individual product labelling.

Breathing filters (if included in the configuration)

Follow the instructions in individual IFUs supplied.

Reservoir bags (if included in the configuration)

1. Ensure that the reservoir bag is securely attached to the breathing system via a push and twist action.

2. Consider the impact of increased compliance of the reservoir bag upon prescribed ventilation.

3. The reservoir bag allows manual ventilation, and gives a visual indication of rate and depth of breathing.

PRE-USE CHECK

1. Prior to installation, and each use, check that all breathing system components are free of any obstruction or foreign bodies and that there is a patent airway and gas flow to the patient. If the system contains one way valves, check and confirm that they are functioning correctly.

2. Remove any red dust cap before connecting the breathing system between the equipment and the patient.

3. Connect/disconnect the connection of the breathing system to the connection of the analgesic gas source using a push and twist action.

4. Following installation and each use the breathing system and all accessories must be checked for leaks and occlusions and that all connections are secure.

5. Consider any impact that the internal volume of the device may have upon the effectiveness of the patient's ventilation. The clinician should assess this on an individual patient basis.

6. Check that any port cap or accessory is securely attached to prevent leakage.

7. The breathing system should be changed in accordance with the hospital policy (max 7 days), or earlier if it is no longer fit for its intended purpose.

IN-USE CHECK

Monitor the clinical parameters of the patient during the treatment.

WARNING

1. This system is not suitable for use with flammable anaesthetic agents.

2. Monitor the breathing system throughout use and replace if increased resistance or contamination occurs.

3. Refer to equipment manufacturers guidelines for pre-use check requirements.

4. Prior to use, consult the respective documentation and instructions for use of all connected devices, or combination of devices.

5. Do not stretch, cut or reconfigure the tubing.

6. Cleaning can lead to device malfunction and breakage.

7. Adding attachments or other devices to the breathing system can change the pressure gradient across the breathing system and adversely affect the performance.

8. Use of this breathing system with equipment not compliant with appropriate International Standards can result in degraded performance. Prior to use, consult the respective documentation and instructions for use of all connected devices, or combination of devices to ensure compatibility.

9. Consult individual product labelling for presence or absence of phthalates.

10. Consult individual product labelling for presence or absence of latex.

CAUTION

1. Long term exposure to elevated levels of nitrous oxide can be harmful. Ensure adequate room ventilation, gas scavenging and/or environmental monitoring as necessary.

2. The breathing system is only to be used under medical supervision by suitably or appropriately qualified medical personnel.

3. To avoid contamination, the product should remain packaged until ready for use. Do not use the product if the packaging is damaged.

4. When connecting or disconnecting the system, hold connector and use a push and twist action.

5. Do not use product past the expiry date.

FOR USA: Rx Only: Federal law restricts this device to sale by or on the order of a physician.

Working pressure: <60 cmH₂O

STORAGE CONDITIONS

Recommended storage at room temperature over the indicated product shelf-life.

DISPOSAL

Following use, the product must be disposed of in accordance with local hospital, infection control and waste disposal regulations.

NOTE

The length referenced in the product description on the product label denotes the minimum working length of the breathing system. The full breathing system length is shown in the technical data section of the product label.

TECHNICAL AND PERFORMANCE DATA

See product label schematic for full breathing system length (l), compliance (C) and resistance to flow (R) data.

WARNING Ensure that the breathing system technical and performance data are appropriate to the intended patient and treatment.

Patient ventilation and vital signs should be monitored. Please refer to product label for full information.

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