CO₂ Cuvette for Dräger® ventilators


Fig. 1

en WARNING: A WARNING statement provides important information about a potentially hazardous situation which, if not avoided could result in death or serious injury.

CAUTION: A CAUTION statement provides important information about a potentially hazardous situation which, if not avoided could result in minor or moderate injury.

Description

CO₂ Cuvette for Dräger® ventilators.

Intended use

An inline connector to facilitate the measurement of CO₂ levels, which is used in conjunction with an infrared CO₂ sensor.

WARNING: Single use product. Do not reuse. Reuse poses a threat to patient safety in terms of cross infection and device malfunction and breakage. Cleaning can lead to device malfunction and breakage.

Intended user

For use by qualified medical personnel only.

Indications

The cuvette is validated for use with the following Dräger® ventilators:

- Dräger® Oxylog® 3000 plus
- Dräger® Evita® Infinity® V500

CAUTION: The cuvette may only be used with the Dräger 6871950 CO₂ mainstream sensor. The cuvette is intended to be used in:

- Hospital environment
- Ambulances
- Pre-hospital environment

Contraindications

The cuvette is not intended for use in the home-care environment.

Patient category

Intended for use with adult patients.

Instructions for use

Installation and operation

1. Connect the cuvette (A) to the patient connector of the breathing system (B) using a push and twist action. (fig.1).

2. Connect the CO₂ sensor (C) on to the cuvette.

CAUTION: The sensor should fit on top of the cuvette with the cuvette windows positioned side-

ways so condensation runs down and away from the windows. (fig.1). Build up of condensation or secretions in the cuvette can obscure the windows and may result in incorrect CO₂ readings.

CAUTION: The cable of the sensor should point in the direction of the ventilator and away from the patient (fig.1).

WARNING: The cuvette may not function correctly when used with nebulised drugs and could result in incorrect CO₂ readings.

Pre-use Checks

WARNING: Prior to use, consult the respective documentation and instructions for use of all connected devices or combination of devices.

WARNING: Perform a self-test of the ventilator after full assembly of the breathing system and associated accessories (including the cuvette) before use on each patient.

WARNING: Every component of the breathing system must be visually inspected and checked for function, leakage, and occlusions, immediately before use on each patient. Ensure that connections are secure.

Zero calibration of the sensor

Observe the instructions for use of the ventilator.

WARNING: Ensure 'disposable' cuvette type is selected on the ventilator. Selecting the 'reusable' cuvette type will result in incorrect CO₂ readings.

In-use checks

WARNING: Monitor the clinical parameters of the patient during the treatment.

Duration of use

Single use product. Maximum duration of use 7 days. Do not reuse. Reuse poses a threat to patient safety in terms of cross infection and device malfunction and breakage. Cleaning can lead to device malfunction and breakage.

Technical and performance data

Material:	ABS
Dimensions:	56mm x 24mm x 23mm
Internal Volume:	6.6 ml
Resistance to flow at 30L/min:	0.3 mbar
Accuracy:	0.43% by volume +8% of the CO ₂ concentration*

* @ concentrations of 2 - 9% (as per ISO 80601-2-55).

Disposal

Following use, the product must be disposed of in line with the local healthcare, hygiene and waste disposal protocols and regulations.

CAUTION This medical device is not supplied individually. Only one copy of the instructions for use is provided per unit of issue and should therefore be kept in an accessible location for all users. Distribute instructions to all product locations and users.

Standards

This product complies with the following relevant standards:

- BS EN ISO 5356-1: Anaesthetic and respiratory equipment. Conical connectors. Cones and sockets.
- BS EN ISO 18190: Anaesthetic and respiratory equipment. General requirements for airways and related devices.
- ISO 80601-2-55: Particular requirements for the basic safety and essential performance of respiratory gas monitors.

Recommended storage at room temperature over the indicated product shelf-life.

*Dräger, Oxylog, Evita and Infinity are registered trademarks of Drägerwerk AG & Co. KGaA.

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