Reservoir bags

en	Reservoir bags
fr	Ballons réservoirs
de	Reservoirbeutel
es	Balones reservorios
pt	Sacos reservatórios
it	Pallone reservoir
nl	Beademingsballonnen
no	Reservoir bag
fi	Säiliöpussit
sv	Reservoarblåsor
da	Reservoir poser
el	Αποθεματικοί ασκοί
lt	Talpyklos maišeliai
pl	Worki rezerwuaru
ru	Дыхательные мешки
cs	Zásobníkové vaky
hu	Lélegeztető ballon
sl	Zbiralne vrečke
lv	Rezervuāra maisi
et	Reservuaarkotid
bg	Резервоари
sr	Kese rezervoara
ro	Balon rezervor
sk	Dýchacie vaky
hr	Vrećice spremnika
tr	Rezervuar balonları
ja	麻酔バッグ
zh-s	储气囊
zh-t	儲氣袋
ar	أكياس التخزين

en Reservoir bags

The instruction for use (IFU) is not supplied individually. Only 1 copy of the instructions for use is provided per box and should therefore be kept in an accessible location for all users. More copies are available upon request.

NOTE: Distribute instructions to all product locations and users. These instructions contain important information for the safe use of the product. Read these instructions for use in their entirety, including warnings and cautions before using this product. Failure to properly follow warnings, cautions and instructions could result in serious injury or death to the patient. Serious incidents to be reported as per the local regulation.

WARNING: A WARNING statement provides important information about a potentially hazardous situation, which, if not avoided could result in death or serious injury.

CAUTION: A CAUTION statement provides important information about a potentially hazardous situation, which if not avoided could result in minor or moderate injury.

Intended Use

To form a reservoir of respiratory gases as part of a breathing system.

Intended user

For use by appropriately qualified medical personnel.

Indications

The reservoir bag should be connected within the breathing system or the anaesthetic machine. The reservoir bag allows manual ventilation and also provides a visual indication of rate and depth of spontaneous breathing.

Contraindications

Do not use with equipment that does not comply to international taper connection standards.

Product Use:

The product is intended for a maximum of 7 days use with a variety of anaesthetic equipment meeting international standards. The start of use is defined from the first connection to the equipment. This product can be used repeatedly on several patients, up to a maximum of 7 days, if a suitably validated single-use respiratory filter is used between the product and each patient. Use in the acute and non-acute environment under the supervision of suitably or appropriately qualified medical personnel. The reservoir bag is safe for use with atmospheric gases: N₂O, O₂ and CO₂; anaesthetic agents: isoflurane, sevoflurane, desflurane and other comparable volatile agents.

Patient Category

The range of reservoir bags (0.5, 1, 2 & 3L) are suitable for use on all adult, paediatric and neonatal patients with intended delivered tidal volumes of ≤50ml (neonatal), 50<300ml (paediatric) and ≥300ml (adult). Check the size of reservoir bag is appropriate for the size and weight of the patient.

Connect the reservoir bag

Attach the reservoir bag to the anaesthetic machine and/or breathing system using a push and twist action (fig.1). Perform pre-use checks as recommended by the equipment manufacturer.

In-Use Check

Monitor the clinical parameters of the patient during the treatment.


Warning

1. This reservoir bag is not suitable for use with flammable anaesthetic agents.
2. Refer to equipment manufacturers guidelines for machine settings and pre-use check requirements.
3. Monitor the reservoir bag throughout use and replace if increased resistance, damage or contamination occurs.
4. Perform a self-test of the equipment after full assembly of the breathing system, reservoir bag and associated accessories before use on each patient.
5. Prior to use, consult the respective documentation and instructions for use of all connected devices, or combination of devices.
6. Cleaning can lead to device malfunction and breakage.

7. Reuse of a reservoir bag in anaesthesia, that hasn't used a suitably validated single-use respiratory filter between the product and each patient, poses a threat to patient safety in terms of cross infection, device malfunction and breakage.
8. Breathing systems, their parts and accessories are validated for use with specific ventilators and anaesthesia machines. Incompatible parts can result in degraded performance. The responsible organisation is responsible for ensuring the compatibility of the ventilator or anaesthesia machine and all of the parts used to connect to the patient before use.
9. Use of this reservoir bag with equipment not compliant with appropriate international standards can result in degraded performance.

Caution

1. The reservoir bag is only to be used under medical supervision by suitably or appropriately qualified medical personnel.
2. To avoid contamination, the product should remain packaged until ready for use. Do not use the product if the packaging is damaged.
3. The reservoir bag should not be re-used without an Intersurgical filter or other appropriately validated respiratory filter at the patient connection.
4. Consult individual product labelling for presence or absence of latex.
5. When connecting or disconnecting the reservoir bag, use a push and twist action.

FOR USA: Rx ONLY: Federal law restricts this device to sale by or on the order of a physician.

Working pressure: <60 cmH₂O

Storage Conditions

Recommended storage at room temperature over the indicated product shelf-life.

Disposal

Following use, the product must be disposed of in accordance with local hospital, infection control and waste disposal regulations.

Patient category	Intended delivered volume ml	Compliance ml/cmH ₂ O	At Pressure cmH ₂ O
Adult	≥ 300 ml	5	60 ± 3
Paediatric	50 ml < 300 ml	4	60 ± 3
Neonatal	≤ 50 ml	1.5	60 ± 3

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