
30mm Uniflow™ breathing systems

en	Instructions for use	pl	Instrukcja użycia
fr	Mode d'emploi	ru	Инструкции по использованию
de	Gebrauchsanweisung	cs	Návod k použití
es	Instrucciones de uso	hu	Használati utasítás
pt	Instruções de utilização	sl	Navodila za uporabo
it	Istruzioni per l'uso	lv	Lietošanas instrukcija
nl	Gebruiksaanwijzing	et	Kasutusjuhised
no	Brukerveiledning	bg	Инструкции за употреба
fi	Käyttäjän opas	sr	Uputstvo za upotrebu
sv	Bruksanvisning	ro	Instrucțiuni de utilizare
da	Brugsanvisning	sk	Návod na použitie
el	Οδηγίες χρήσης	hr	Upute za uporabu
lt	Naudojimo instrukcija	tr	Kullanma kılavuzu

30 mm Uniflow™ breathing systems
The instruction for use (IFU) is not supplied individually. Only 1 copy of the instructions for use is provided per box and should therefore be kept in an accessible location for all users. More copies are available upon request.

NOTE: Distribute instructions to all product locations and users. These instructions contain important information for the safe use of the product. Read these instructions for use in their entirety, including warnings and cautions before using this product. Failure to properly follow warnings, cautions and instructions could result in serious injury or death to the patient. Serious incidents to be reported as per the local regulation.

WARNING: A WARNING statement provides important information about a potentially hazardous situation, which, if not avoided could result in death or serious injury.

CAUTION: A CAUTION statement provides important information about a potentially hazardous situation, which if not avoided could result in minor or moderate injury.

Intended Use
When used in Critical Care: To deliver and remove respiratory gases to and from a patient via a breathing system comprised of tubing and connectors.

When used in Anaesthesia: To deliver and remove anaesthetic and respiratory gases to and from a patient via a breathing system comprised of tubing and connectors.

Intended user
For use by qualified medical personnel.

Indications
The breathing system should be connected between the ventilator or the anaesthetic machine, and the patient airway connection. In Critical Care the product should only be used on a single patient up to a maximum of 7 days. In Anaesthesia the product can be used repeatedly on several patients up to 7 days, if a suitably validated single-use respiratory filter is used between the product and each patient. The breathing system is provided in Flextube™ tubing. Silver Knight™ breathing systems contain anti-microbial additive in the Flextube.

Contraindications

Do not use with equipment that does not comply to international taper connection standards.

Do not use Silver Knight breathing systems with patients with an intended delivered tidal volume of < 300ml.

Product Use:

The product is intended for a maximum of 7 days use with a variety of anaesthetic and ventilator equipment meeting international taper standards. The start of use is defined from the first connection to the equipment. Use in the acute and non-acute environment under the supervision of suitably or appropriately qualified medical personnel. The breathing system is safe for use with atmospheric gases: N₂O, O₂ and CO₂; anaesthetic agents: isoflurane, sevoflurane, desflurane and other comparable volatile agents.

Patient Category

The range of 30mm Uniflow breathing systems are suitable for use on all adult and paediatric patients with intended delivered tidal volumes of ≥50ml.

All systems which contain Silver Knight anti microbial are only suitable for use with adult patients with intended tidal volume of ≥300ml.

Connect the breathing system

Attach the breathing system to the inspiratory and expiratory ports of the ventilator or anaesthetic machine using a push and twist action (figure 2). The inspiratory port is indicated by (so it should read indicated by a directional label) directional label on the breathing system.

Accessories:

Where the 30mm Uniflow breathing system is supplied with additional components the clinician should assess any potential impact upon the patient's respiratory function, and be aware of the subsequent increased resistance to flow, compressible volume and weight of the product.

When the ventilator and anaesthetic machine parameters are set, any impact that the internal volume and resistance to flow of the device may have upon the effectiveness of the prescribed ventilation, should be assessed on an individual patient basis by the clinician.

Water traps and breathing filters, HMEs, HMEFs if included in the configuration

Follow the instructions in individual IFUs supplied.

Reservoir bags (if included in the system)

1. Ensure that the reservoir bag is securely attached to the breathing system via a push and twist action.
2. Consider the impact of increased compliance of the reservoir bag upon prescribed ventilation.
3. The reservoir bag allows manual ventilation, and gives a visual indication of rate and depth of breathing.

Monitoring lines (if included in the system)

Monitoring lines provide a leak free connection between the breathing system and a respiratory/anaesthetic gas monitor via a luer lock fitting. The connection should be checked prior to use. The line should be checked to confirm it is kink free and a patent pathway exists.

Pre-Use Check - Critical Care Ventilation use (Single patient use)

1. Prior to attachment, check that all breathing system components are free of any obstruction or foreign bodies and that there is a patent airway.
2. Remove any red dust cap before connecting the breathing system between the equipment and the patient.
3. Follow the diagrams in figure 1 to test and confirm the inner lumen of the breathing system is leak free. Further information on how to correctly perform this leak test can be found on the Intersurgical YouTube channel by searching 'Using the Uniflow/Bain Leak Tester'.
4. Connect/disconnect the tapered connections of the breathing system to the tapered connection of the ventilator and/or respiratory equipment using a push and twist action. Perform pre-use checks as recommended by the equipment manufacturer.
5. Following attachment the breathing system and all accessories must be checked for leaks and occlusions prior to use and that all connections are secure.
6. When setting ventilator parameters consider any impact that the internal volume of the breathing system and accessories may have upon the effectiveness of the prescribed ventilation. The clinician should assess this on an individual patient basis.
7. Check that any port cap or accessory is securely attached to prevent leakage.
8. The breathing system should be changed in accordance with the hospital policy (max 7 days), or earlier if it is no longer fit for its intended purpose.

Pre-Use Check - Anaesthesia use (Reusable)

1. Prior to attachment, and each use, check that all breathing system components are free of any obstruction or foreign bodies and that there is a patent airway.
2. Remove any red dust cap before connecting the breathing system between the equipment and the patient.
3. Follow the diagrams in figure 1 to test and confirm the inner lumen of the breathing system is leak free. Further information on how to correctly perform this leak test can be found on the Intersurgical YouTube channel by searching 'Using the Uniflow/Bain Leak Tester'.
4. It is recommended this leak test is carried out prior to use on each individual patient.
4. Connect/disconnect the tapered connection of the breathing system to the tapered connection of the anaesthesia machine using a push and twist action. Perform pre-use checks as recommended by the equipment manufacturer.
5. The breathing system should be protected from cross contamination by an Intersurgical breathing filter, or other appropriately validated respiratory filter, at the patient connection. The filter and patient connection should be changed for each patient to minimise the risk of cross contamination.
6. Following attachment and each use the breathing system and all accessories must be checked for leaks and occlusions and that all connections are secure.
7. When setting ventilator parameters consider any impact that the internal volume of the device may have upon the effectiveness of the prescribed ventilation. The clinician should assess this on an individual patient basis.
8. Check that any port cap or accessory is securely attached to prevent leakage.
9. The breathing system should be changed in accordance with the hospital policy (max 7 days), or earlier if it is no longer fit for its intended purpose.

In-Use Check

Monitor the clinical parameters of the patient during the treatment.

Warning

1. This system is not suitable for use with flammable anaesthetic agents.
2. Monitor the breathing system throughout use and replace if increased resistance or contamination occurs.
3. Refer to equipment manufacturers guidelines for machine settings and pre-use check requirements.
4. Perform a self-test of the equipment after full assembly of the breathing system and associated accessories before use on each patient.
5. Prior to use, consult the respective documentation and instructions for use of all connected devices, or combination of devices.
6. Do not stretch, cut or reconfigure the tubing.
7. Cleaning can lead to device malfunction and breakage.
8. Ventilator accuracy can be affected by use of a gas driven nebuliser.
9. Adding attachments or other devices to the breathing system can change the pressure gradient across the breathing system and adversely affect the ventilator performance. Perform a self-test of the ventilator (or anaesthesia machine) after full assembly of the breathing system and all connected devices before use on each patient.
10. Reuse of a breathing system in anaesthesia, where a suitably validated single-use respiratory filter has not been used between the product and each patient, poses a threat to patient safety in terms of cross infection, device malfunction and breakage.
11. Breathing systems, their parts and accessories are validated for use with

specific ventilators and anaesthesia machines. Incompatible parts can result in degraded performance. The responsible organisation is responsible for ensuring the compatibility of the ventilator or anaesthesia machine and all of the parts used to connect to the patient before use.

12. Use of this breathing system with equipment not compliant with appropriate International Standards can result in degraded performance. Prior to use, consult the respective documentation and instructions for use of all connected devices, or combination of devices to ensure compatibility.

13. The Luer Lock port connector is for the exclusive use of monitoring respiratory and anaesthetic gases. If applicable, check that the cap is securely attached to prevent inadvertent leakage when the port is not in use. Do not introduce gases or liquids via the misconnection of the luer lock port as this can result in serious injury to the patient.

14. Consult individual product labelling to determine if device is MR safe or unsafe.

15. Consult individual product labelling for presence or absence of phthalates.

Caution

1. The breathing system is only to be used under medical supervision by suitably or appropriately qualified medical personnel.
2. To avoid contamination, the product should remain packaged until ready for use. Do not use the product if the packaging is damaged.
3. The breathing system should not be re-used without an Intersurgical filter or other appropriately validated respiratory filter at the patient connection.
4. When connecting or disconnecting the system, use a push and twist action.
5. When the breathing system is not in use, the patient connection which is changed for each patient should be protected by the red dust cap supplied to prevent contamination.
- FOR USA: Rx only: Federal law restricts this device to sale by or on the order of a physician.
- Working pressure: <60 cm H₂O
- Storage Conditions**
Store out of direct sunlight.
Recommended storage at room temperature over the indicated product shelf-life.

Environmental Operating Conditions:
Performance verified at room/ambient temperature at sea level.

Disposal

Following use, the product must be disposed of in accordance with local hospital, infection control and waste disposal regulations.

Note

The length referenced in the product description on the product label denotes the minimum working length of the breathing system. The full breathing system length is shown in the technical data section of the product label.

Standards

This product complies with the following relevant standards:

- BS EN ISO 5356-1 Anaesthetic and respiratory equipment. Conical connectors. Cones and sockets
- BS EN ISO 5367 Anaesthetic and respiratory equipment. Breathing tube and sets
- BS EN ISO 80601-2-12 MEE Critical Care Ventilators
- BS EN ISO 80601-2-13 Particular requirements for basic safety and essential performance of an anaesthetic workstation
- BS EN ISO 80369-7 Small-bore connectors for liquids and gases in healthcare applications. Connectors for intravascular or hypodermic applications


Fig. 2

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UAB Intersurgical
Arnionij g. 60, Pabrad  ,
LT-18170, Lithuania

INTERSURGICAL
COMPLETE RESPIRATORY SYSTEMS

Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK

T: +44 (0)118 9656 300 info@intersurgical.com www.intersurgical.com

Ireland T: +353 (0)1 697 8540 info@intersurgical.ie	Espan�� T: +34 91 665 73 15 info@intersurgical.es	Nederland T: +31 (0)411 243 860 info@intersurgical.nl	Sverige T: +46 08 514 30 600 info@intersurgical.se	Россия T: +7 495 771 6809 info@intersurgical.co.za	South Africa T: +27 (0)11 444 7968 info@intersurgical.co.za	��tland T: +886 4 2380 5430 info@intersurgical.com.tw	Philippines T: +63 2820 4124 info@intersurgical.ph	USA T: +1 800 828 9633 support@intersurgicalinc.com	Canada T: +1 905-319-6500 info@intersurgical.ca	Colombia T: +57 1 742 5161 info@intersurgical.com	Australia T: +61 (0)2 8048 3300 info@intersurgical.com.au
France T: +33 (0)1 48 76 72 30 info@intersurgical.fr	Portugal T: +351 219 108 550 info@intersurgical.pt	Belgi��/Belgique T: +32 5 2456991 info@intersurgical.be	Danmark T: +45 28 69 80 00 info@intersurgical.dk	��esk�� Republika T: +420 272 911 814 info@intersurgical.cz	��tina T: +86 519 6981700 info@intersurgical.cn.com	��tland T: +886 4 2380 5430 info@intersurgical.com.tw	Philippines T: +63 2820 4124 info@intersurgical.ph	USA T: +1 800 828 9633 support@intersurgicalinc.com	Canada T: +1 905-319-6500 info@intersurgical.ca	Colombia T: +57 1 742 5161 info@intersurgical.com	Australia T: +61 (0)2 8048 3300 info@intersurgical.com.au
Deutschland T: +49 (0)2241 25690 info@intersurgical.de	Italia T: +39 0535 20836 info@intersurgical.it	Luxembourg T: +32 5 2456991 info@intersurgical.be	Lietuva T: +370 387 66611 info@intersurgical.lt	T��rkiye T: +90 216 468 88 28 info@intersurgical.com	��tland T: +86 519 6981700 info@intersurgical.cn.com	��tland T: +886 4 2380 5430 info@intersurgical.com.tw	Philippines T: +63 2820 4124 info@intersurgical.ph	USA T: +1 800 828 9633 support@intersurgicalinc.com	Canada T: +1 905-319-6500 info@intersurgical.ca	Colombia T: +57 1 742 5161 info@intersurgical.com	Australia T: +61 (0)2 8048 3300 info@intersurgical.com.au

Distributed in the USA by Intersurgical Incorporated, 6757 Kinne St., East Syracuse, NY 13057
T: (800) 828-9633 F: (315) 451-3696 support@intersurgicalinc.com www.intersurgical.com