

30mm Uniflow™ breathing systems


en	Instructions for use	pl	Instrukcja użycia
fr	Mode d'emploi	ru	Инструкции по использованию
de	Gebrauchsanweisung	cs	Návod k použití
es	Instrucciones de uso	hu	Használati utasítás
pt	Instruções de utilização	sl	Navodila za uporabo
it	Istruzioni per l'uso	lv	Lietošanas instrukcija
nl	Gebruiksaanwijzing	et	Kasutusjuhised
no	Brukerveiledning	bg	Инструкции за употреба
fi	Käyttäjän opas	sr	Uputstvo za upotrebu
sv	Bruksanvisning	ro	Instrucţiuni de utilizare
da	Brugsanvisning	sk	Návod na použitie
el	Οδηγίες χρήσης	hr	Upute za uporabu
lt	Naudojimo instrukcija	tr	Kullanma kılavuzu

1

2

Patient category ¹	Intended delivered volume ² ml	Flow resistance limit ³ cmH ₂ O/l/min	At flow ⁴ l/min
Adult ⁷	≥ 300 ml	0.06	30
Paediatric ⁸	50 ml < 300 ml	0.12	15

Patient category	Intended delivered volume ml	Compliance ⁵ ml/cmH ₂ O	At Pressure ⁶ cmH ₂ O
Adult	≥ 300 ml	≤ 5	60 ± 3
Paediatric	50 ml < 300 ml	≤ 4	60 ± 3

1 Patient category / Catégorie de patients / Patientenkategorie / Categoría de paciente / Categoria de doentes / Categoria di pazienti / Patiëntencategorie / Pasientkategori / Potilastyypit / Patientkategori / Patientkategori / Κατηγορία ασθενών / Pacientų kategorija / Kategorija pacjentów / Категория пациентов / Kategorie pacientů / Betegek köre / Kategorija bolnika / Pacientu kategorija / Patsiendi kategooria / Категория на пациента / Kategorija pacijenata / Categorie de pacienți / Kategoría pacientov / Kategorija pacijenata / Hasta kategorisi

2 Intended delivered volume ml / Volume administré prévu ml / Zuführendes Volumen ml / Volumen de suministro previsto ml / Volume administrado previsto ml / Volume erogato previsto ml / Beoogd volume ml / Tiltenkt leveret volum ml / Aiottu antotilavuus ml / Avsedd levererad volym ml / Tilsigter leveret volumen ml / Προβλεπόμενος χορηγούμενος όγκος ml / Numatomas pristatomas tūris ml / Planowana objętość dostarczana ml / Планируемый дыхательный объем ml / Predpokládaný prijatý objem ml / Leadandó légzési térfogat ml / Predvidena dovedena prostornina ml / Paredzētais piegādes tilpums ml / Kavandatud manustatav maht ml / Планиран осигуряван капацитет ml / Planirana isporučena zapremina ml / Volumul estimat livrat ml / Určený podávaný objem ml / Predvideni isporučeni volumen ml / Amaçlanan uygulama hacmi ml

3 Flow resistance limit cmH₂O/l/min / Limite de résistance à l'écoulement cmH₂O/l/min / Grenzwert des Strömungswiderstands cmH₂O/l/min / Limite de resistencia de flujo cmH₂O/l/min / Limite de resistência ao fluxo cmH₂O/l/min / Limite della resistenza al flusso cmH₂O/l/min / Limiet stromingsweerstand cmH₂O/l/min / Grense for flow-motstand cmH₂O/l/min / Virtausvastuksen raja cmH₂O/l/min / Gräns för flödesmotstånd cmH₂O/l/min / Grænse for flow-motstand cmH₂O/l/min / Όριο αντίστασης στη ροή cm H₂O/λεπτό / Pasipriešinimo srautui riba cmH₂O/l/min / Granica oporu przepływu cmH₂O/l/min / Предел сопротивления потоку cm H₂O/л/мин / Limit odporu proudění cmH₂O/l/min / Áramlási ellenállás határértéke cmH₂O/l/perc / Omejeitev odpornosti za pretok cmH₂O/l/min / Plūsmas pretestības robeža cmH₂O/l/min / Voolutakistuse piirväärtus cmH₂O/l/min / Граница на съпротивление на потока cmH₂O/l/min / Granica vrednost otpornosti protoku cmH₂O/l/min / Limită de rezistență la curgere cmH₂O/l/min / Limit odporu voči prúdeniu cmH₂O/l/min / Granica otpora na protok cmH₂O/l/min / Akış direnci sınırı cmH₂O/l/dakika

4 At flow l/min/ Débit l/min/ Bei Durchfluss l/min/ Flujo l/min/ Com o fluxo l/min/ Al flusso l/min/ Bij het stromen l/min/ Ved flow l/min/ Virtauksella l/min/ Vid flödet l/min/ Ved flow l/min./ Σε ροή l/λεπτό/ Esant srautui l/min/ Przy przepływie l/min/ В потоке л/мин/ Při proudění l/min/ Áramlási sebesség l/perc/ Pri pretoku l/min/ Pie aktīvas plūsmas l/min/ Voolutiirusel l/min/ При поток l/min/ Pri protoku od l/min/ La curgere l/min/ Pri prúdení l/min/ Pri protoku l/min/ Akışta l/dakika

5 Compliance ml/cmH₂O / Tolérance ml/cmH₂O / Compliance ml/cmH₂O / Compliance ml/cmH₂O / Complacência ml/cmH₂O / Compliance ml/cmH₂O / Naleving ml/cmH₂O / Compliance ml/cmH₂O / Kompianssi ml/cmH₂O / Compliance ml/cmH₂O / Overensstemmelse ml/cmH₂O / Συμμόρφωση ml/cm H₂O / Atitiktis ml/cmH₂O / Podatność ml/cmH₂O / Растяжимость мл/см H₂O / Shoda ml/cmH₂O / Megfelelőség ml/cmH₂O / Kompliance ml/cmH₂O / Ättilstība ml/cmH₂O / Vastavus ml/cmH₂O / Еластичност ml/cmH₂O / Komplijansa ml/cmH₂O / Complianță ml/cmH₂O / Poddajnosť ml/cmH₂O / Rastezljivost ml/cmH₂O / Uygunluk ml/cmH₂O

6 At Pressure cmH₂O / Pression cmH₂O / Bei Druck cmH₂O / Presión cmH₂O / Com a pressão cmH₂O / Alla pressione cmH₂O / Onder druk cmH₂O / Ved tryk cmH₂O / Paineeella cmH₂O / Vid tryck cmH₂O / Ved tryk cmH₂O / Σε πίεση cm H₂O / Esant slėgiui cmH₂O / Przy ciśnieniu cmH₂O / Под давлением cm H₂O / Při tlaku cmH₂O / Nyomás cmH₂O / Pri tlaku cmH₂O / Pie spiediena cmH₂O / Rõhul cmH₂O / При налягане cmH₂O / Pri pŕifisku od cmH₂O / Ľa presiune cmH₂O / Pri tlaku cmH₂O / Pri tlaku cmH₂O / Basınçta cmH₂O

7 Adult / Adulte / Für Erwachsene / Adulto / Adulto / Adulto / Volwassene / Voksne / Aikuinen / Vuxen / Voksen / Ενήλικες / Suaugusieji / Dorosly / Для взрослых / Dospělý / Feinótt / Odrasli / Pieaugušajiem / Täiskasvanu / Взрастни / Odrasli / Adulti / Dospelí / Odrasla osoba / Yetişkin

8 Paediatric / Pédiatriques / Pédiatrischen / Pediátrico / Pediátricos / Pediatrico / Pediatricisch / Barn / Lapsille / Barn / Pædiatrisk / Ταϊδιατρικούς / Vaikai / Dzieci / Для педиатрических / Pediatrické / Gyermek / Pediatrični / Pediatrijas / Pediaatriline / Педиатрична / Pedijatrijske / Copii / Pediatrických / Pedijatriju / Pediatrik

30mm Uniflow™ breathing systems

The instruction for use (IFU) is not supplied individually. Only 1 copy of the instructions for use is provided per box and should therefore be kept in an accessible location for all users. More copies are available upon request.

NOTE Distribute instructions to all product locations and users. These instructions contain important information for the safe use of the product. Read these instructions for use in their entirety, including warnings and cautions before using this product. Failure to properly follow warnings, cautions and instructions could result in serious injury or death to the patient. Serious incidents to be reported as per the local regulation.

WARNING A WARNING statement provides important information about a potentially hazardous situation, which, if not avoided could result in death or serious injury.

CAUTION A CAUTION statement provides important information about a potentially hazardous situation, which if not avoided could result in minor or moderate injury.

Intended Use

When used in Critical Care: To deliver and remove respiratory gases to and from a patient via a breathing system composed of tubing and connectors.

When used in Anaesthesia: To deliver and remove anaesthetic and respiratory gases to and from a patient via a breathing system comprised of tubing and connectors.

Intended user

For use by qualified medical personnel.

Indications

The breathing system should be connected between the ventilator or the anaesthetic machine, and the patient airway connection. In Critical Care the product should only be used on a single patient up to a maximum of 7 days. In Anaesthesia the product can be used repeatedly on several patients up to 7 days, if a suitably validated single-use respiratory filter is used between the product and each patient. The breathing system is provided in Flextube™ tubing. Silver Knight™ breathing systems contain anti-microbial additive in the Flextube.

Contraindications

Do not use with equipment that does not comply to international taper connection standards.

Do not use Silver Knight breathing systems with patients with an intended delivered tidal volume of < 300ml.

Product Use

The product is intended for a maximum of 7 days use with a variety of anaesthetic and ventilator equipment meeting international taper standards. The start of use is defined from the first connection to the equipment. Use in the acute and non-acute environment under the supervision of suitably or appropriately qualified medical personnel. The breathing system is safe for use with atmospheric gases: N₂O, O₂ and CO₂; anaesthetic agents: isoflurane, sevoflurane, desflurane and other comparable volatile agents.

Patient Category

The range of 30mm Uniflow breathing systems are suitable for use on all adult and paediatric patients with intended delivered tidal volumes of ≥250ml.

All systems which contain Silver Knight anti microbial are only suitable for use with adult patients with intended tidal volume of ≥300ml.

Connect the breathing system

Attach the breathing system to the inspiratory and expiratory ports of the ventilator or anaesthetic machine using a push and twist action (figure 2).

The inspiratory port is indicated by directional label on the breathing system.

Accessories

Where the 30mm Uniflow breathing system is supplied with additional components the clinician should assess any potential impact upon the patient's respiratory function, and be aware of the subsequent increased resistance to flow, compressible volume and weight of the product.

When the ventilator and anaesthetic machine parameters are set, any impact that the internal volume and resistance to flow of the device may have upon the effectiveness of the prescribed ventilation, should be assessed on an individual patient basis by the clinician.

Water traps and breathing filters, HMEs, HMEFs if included in the configuration

Follow the instructions in individual IFUs supplied.

Reservoir bags (if included in the system)

1. Ensure that the reservoir bag is securely attached to the breathing system via a push and twist action.
2. Consider the impact of increased compliance of the reservoir bag on patient airway.
3. The reservoir bag allows manual ventilation, and gives a visual indication of rate and depth of breathing.

Monitoring lines (if included in the system)

Monitoring lines provide a leak free connection between the breathing system and a respiratory/anaesthetic gas monitor via a luer lock fitting. The connection should be checked prior to use. The line should be checked to confirm it is kink free and a patent pathway exists.

Pre-Use Check - Critical Care Ventilation use (Single patient use)

1. Prior to attachment, check that all breathing system components are free of any obstruction or foreign bodies and that there is a patent airway.
2. Remove any red dust cap before connecting the breathing system between the equipment and the patient.

3. Follow the diagrams in figure 1 to test and confirm the inner lumen of the breathing system is leak free. Further information on how to correctly perform this leak test can be found on the Intersurgical YouTube channel by searching 'Using the Uniflow/Bain Leak Tester'.

4. Connect/disconnect the tapered connections of the breathing system to the tapered connection of the ventilator and/or respiratory equipment using a push and twist action. Perform pre-use checks as recommended by the equipment manufacturer.

5. Following attachment the breathing system and all accessories must be checked for leaks and occlusions prior to use and that all connections are secure.

6. When setting ventilator parameters consider any impact that the internal volume of the breathing system and accessories may have upon the effectiveness of the prescribed ventilation. The clinician should assess this on an individual patient basis.

7. Check that any port cap or accessory is securely attached to prevent leakage.

8. The breathing system should be changed in accordance with the hospital policy (max 7 days), or earlier if it is no longer fit for its intended purpose.

Pre-Use Check - Anaesthesia use (Reusable)

1. Prior to attachment, and each use, check that all breathing system components are free of any obstruction or foreign bodies and that there is a patent airway.

2. Remove any red dust cap before connecting the breathing system between the equipment and the patient.

3. Follow the diagrams in figure 1 to test and confirm the inner lumen of the breathing system is leak free. Further information on how to correctly perform this leak test can be found on the Intersurgical YouTube channel by searching 'Using the Uniflow/Bain Leak Tester'. It is recommended this leak test is carried out prior to use on each individual patient.

4. Connect/disconnect the tapered connection of the breathing system to the tapered connection of the anaesthesia machine using a push and twist action. Perform pre-use checks as recommended by the equipment manufacturer.

5. The breathing system should be protected from cross contamination by an Intersurgical breathing filter, or other appropriately validated respiratory filter, at the patient connection. The filter and patient connection should be changed for each patient to minimise the risk of cross contamination.

6. Following attachment and each use the breathing system and all accessories must be checked for leaks and occlusions and that all connections are secure.

7. When setting ventilator parameters consider any impact that the internal volume of the device may have upon the effectiveness of the prescribed ventilation. The clinician should assess this on an individual patient basis.

8. Check that any port cap or accessory is securely attached to prevent leakage.

9. The breathing system should be changed in accordance with the hospital policy (max 7 days), or earlier if it is no longer fit for its intended purpose.

In-Use Check

Monitor the clinical parameters of the patient during the treatment.

Warnings

1. This system is not suitable for use with flammable anaesthetic agents.
2. Monitor the breathing system throughout use and replace if increased resistance or contamination occurs.
3. Refer to equipment manufacturers guidelines for machine settings and pre-use check requirements.
4. Perform a self-test of the equipment after full assembly of the breathing system and associated accessories before use on each patient.
5. Prior to use, consult the respective documentation and instructions for use of all connected devices, or combination of devices.
6. Do not stretch, cut or reconfigure the tubing.
7. Do not clean. Cleaning can lead to device malfunction and breakage.
8. Ventilator accuracy can be affected by use of a gas driven nebuliser.
9. Adding attachments or other devices to the breathing system can change the pressure gradient across the breathing system and adversely affect the ventilator performance. Perform a self-test of the ventilator (or anaesthesia machine) after full assembly of the breathing system and all connected devices before use on each patient.
10. Reuse of a breathing system in anaesthesia, where suitably validated, may reduce the respiratory filter and is not been used between the product and each patient, poses a threat to patient safety in terms of cross infection, device malfunction and breakage.
11. Breathing systems, their parts and accessories are validated for use with specific ventilators and anaesthesia machines. Incompatible parts can result in degraded performance. The responsible organisation is responsible for ensuring the compatibility of the ventilator or anaesthesia machine and all of the parts used to connect to the patient before use.
12. Use of this breathing system with equipment not compliant with appropriate International Standards can result in degraded performance. Prior to use, consult the respective documentation and instructions for use of all connected devices, or combination of devices to ensure

compatibility.

13. The luer lock port connector is for the exclusive use of monitoring respiratory and anaesthetic gases. If applicable, check that the cap is securely attached to prevent inadvertent leakage when the port is not in use. Do not introduce gases or liquids via the misconnection of the luer lock port as this can result in serious injury to the patient.

14. Consult individual product labelling to determine if device is MR safe or unsafe.

15. Consult individual product labelling for presence or absence of phthalates.

Caution

1. The breathing system is only to be used under medical supervision by suitably or appropriately qualified medical personnel.

2. To avoid contamination, the product should remain packaged until ready for use. Do not use the product if the packaging is damaged.

3. The breathing system should not be re-used without an Intersurgical filter or other appropriately validated respiratory filter at the patient connection.

4. When connecting or disconnecting the system, use a push and twist action.

5. When the breathing system is not in use, the patient connection which is changed for each patient should be protected by the red dust cap supplied to prevent contamination.

FOR USA: Rx ONLY: Federal law restricts this device to sale by or on the order of a physician.

Working pressure: < 60 cm H₂O

Storage Conditions

Store out of direct sunlight.

Recommended storage at room temperature over the indicated product shelf-life.

Environmental Operating Conditions: Performance verified at room/ambient temperature at sea level.

Disposal

Following use, the product must be disposed of in accordance with local hospital, infection control and waste disposal regulations.

Note

The length referenced in the product description on the product label denotes the minimum working length of the breathing system. The full breathing system length is shown in the technical data section of the product label.

Technical and performance data

See product label schematic for full breathing system length (L), compliance (C) and resistance to flow (R) data.

Applicable standards

The product complies with relevant requirements of the following standards:

- BS EN ISO 5356-1 Anaesthetic and respiratory equipment. Conical connectors. Cones and sockets
- BS EN ISO 5367 Anaesthetic and respiratory equipment. Breathing tube and sets
- BS EN ISO 80601-2-12 MEE Critical Care Ventilators
- BS EN ISO 80601-2-13 Particular requirements for basic safety and essential performance of an anaesthetic workstation
- BS EN ISO 80369-7 Small-bore connectors for liquids and gases in healthcare applications. Connectors for intravascular or hypodermic applications

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Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK