

SpO₂ Sensor Instructions for Use

Shanghai International Holding Corp. GmbH (Europe)
Eiffestraße 80, 20537 Hamburg, Germany

mindray

English

Intended Use

Intended Purpose Statement

The SpO₂ Sensor is intended for connecting with Mindray medical devices that support SpO₂ measurements for measuring the arterial oxygen saturation and pulse rate of patients.

Indication for Use

The SpO₂ sensor is intended for noninvasively spotchecking or continuously monitoring the arterial oxygen saturation and pulse rate of adult and pediatric patients.

The SpO₂ sensor must be used in connection with the monitoring device (including accessories) specified by our company, including the MEC, IMEC, PM, VS, and T series patient monitors of Mindray that supports the SpO₂ measurement function. The user shall be responsible for verifying the compatibility of this sensor, the extension cable, and the equipment before use. This sensor must be used by trained clinical professionals. Otherwise, product performance may be compromised or personal injury may occur.

Intended Users

This sensor must be used by trained clinical professionals.

Intended Patient Population

The SpO₂ Sensor is intended for noninvasively spotchecking or continuously monitoring the arterial oxygen saturation and pulse rate of adult and pediatric patients.

Intended Medical Condition

The product is expected to be used in medical institutions.

Contraindications

Not identified yet.

Side effects

According to the conclusion of clinical evaluation and residual risk evaluation, for the intended patients, there is no known side effects that can occur during or after the use of the medical device. And there is no need for the operator to make extra preparations. Thus, no residual risk associated with using the medical device should be disclosed due to the risk management report.

Warnings

- The sensor material that contacts the patients has undertaken the bio-compatibility test and is verified to be in compliance with ISO 10993-1.
- Select the corresponding SpO₂ extension cable based on the type of device connector.
- Check the SpO₂ sensor and the cable before using. Do not use them for SpO₂ monitoring in case of damage, block, deterioration or visible contamination.
- For disposal of accessories, follow local regulations or your hospital's policy regarding disposal of such accessories. Do not dispose randomly.
- Try to keep the patient still (especially the head) and avoid the measured site suffering excessive motion.
- Cover the sensor with opaque material under the condition of strong light. Failure to do so will result in inaccurate measurement.
- Periodically check the circulation of the sensor site, if the sensor site has very low perfusion, the SpO₂ reading may be inaccurate.
- Inspect the application site every two hours and change the application site when taking a long time continuous monitoring. Some patients may require more frequent inspections, such as patients with perfusion disorder or sensitive skin. Prolonged continuous monitoring may increase the risk of undesirable changes in skin characteristics, such as irritation, reddening, blistering, burns, venous pulsation, obstructed blood circulation, pressure marks, pressure necrosis and so on.
- The SpO₂ sensor must be used in connection with the specified equipment only. The user shall be responsible for ensuring the compatibility between this sensor, the extension cable and other equipment by either referring to the instructions for using this sensor or contacting us before use. Otherwise, product performance may be compromised or personal injury may occur.
- Carefully route cables to reduce the possibility of patient entanglement or strangulation.
- To reduce the impact on the measurement of electrical surgical equipment, do not place the sensor cable near to the power lines of surgical equipment.
- Do not bend the cable excessively for a prolonged time.
- Do not use SpO₂ sensors during magnetic resonance imaging (MRI). Induced current could potentially cause burns. The sensor may affect the MRI image, and the MRI unit may affect the accuracy of the oximetry measurements.
- During radiation therapy, avoid exposing the SpO₂ sensor to the radiation field. Otherwise the SpO₂ reading may be inaccurate, or even turns to be zero.
- Do not place the SpO₂ sensor in an environment that exceeds the stated range.
- Do not modify the SpO₂ sensor in any way. Modification may affect performance and accuracy.

Notes

- Remove the SpO₂ sensor from use if the sensor is too loose and cannot be properly applied to the ear lobe.
- To avoid patient discomfort and sensor damage, do not remove the sensor from the patient's ear by pulling the cable.
- Improper sensor placement may affect the measurement accuracy.
- Ear lobes have lower perfusion ability. Before use, rub or warm the ear. This will improve arterial blood flow and measurement accuracy.
- External pressures on the measurement site may affect measurement accuracy.
- The sensor material that contacts the patients has undertaken the bio-compatibility test and is verified to be in compliance with ISO 10993-1.
- A functional tester or SpO₂ simulator can be used to determine the pulse rate accuracy.
- Using host medical equipment and cables other than those specified may result in increased electromagnetic emission or decreased electromagnetic immunity of the equipment.
- The device or its components should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the device or its components should be observed to verify normal operation in the configuration in which it will be used.
- The device needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided below.
- Other devices may interfere with this device even though they meet the requirements of CISPR.
- When the inputted signal is below the minimum amplitude provided in technical specifications, erroneous measurements could result.
- Portable and mobile communication equipment may affect the performance of this device.
- Other devices that have RF transmitter or source may affect this device (e.g. cell phones, PDAs, and PCs with wireless function).

Applying the SpO₂ Sensor

Be sure to understand all warnings listed in the previous section before applying any SpO₂ sensor to the patient. Before use, verify whether the SpO₂ sensor is intact. In case of any damage, do not apply it to the patient.

- Clean the skin of the application site. For patients with ear disease or pierced ears, use the SpO₂ sensor with caution.
- Clip the sensor onto the ear lobe. Ensure that the ear lobe completely cover the sensor's photodetector window, so that light can completely pass through the patient's tissue.
- Secure the cable clip to the proper position, such as clothes or sheets. Cable placement should avoid pulling and twisting the ear lobe.
- Connect the extension cable and the sensor. Then connect the extension cable to the pulse oximetry equipment by referring to the operator's manual of the equipment.

Checking the Functionality of the Sensor

Replace the sensor if the following situations occur when the sensor is connected to the instrument:

- The "Sensor Fault" alarm occurs.
- There is no "Sensor Off" alarm when the sensor is closed but not applied to the ear lobe.
- The "Sensor Off" alarm occurs when the sensor is applied to the ear lobe.

Start SpO₂ monitoring only when you are sure that the sensor functions normally.

Factors affecting measurement accuracy

The following factors may influence measurement accuracy:

- Dysfunctional haemoglobin, such as carboxyhemoglobin (COHb) and methemoglobin (MetHb)
- Drop of arterial blood flow to immeasurable level caused by shock, anemia, low temperature or vasoconstrictor.
- Presence of certain dyes, frostbite and burns on the ears, dirty skin, and significantly pierced ears.
- Physical movement (patient and imposed motion)
- Ambient light
- Diagnostic testing
- Low perfusion
- Electromagnetic interference, such as MRI environment
- Electrosurgical units

Cleaning/Disinfection

Warnings

- Do not immerse the connector of the SpO₂ sensor in any of the cleaning solutions, disinfectants, or other liquids.
- Do not soak the ear clip of SpO₂ sensor in cleaning solutions or disinfectants for a period of time longer than that specified by the cleaning solutions/disinfectants manufacturer.
- Exercise caution during cleaning/disinfection to avoid wetting the pins.
- Do not clean or disinfect the metal parts of this product. Otherwise, metal part corrosion or product reliability deterioration could result.
- We recommend that the sensor be disinfected only when necessary as determined by your hospital's policy, to avoid long term damage to the sensor.
- Never use cleaning agents/disinfectants other than the recommended.

Cleaning

- Clean the sensor with cotton or soft cloth moistened with water. You can also rinse and wipe the ear clip of the SpO₂ sensor in water.
- After cleaning, wipe off the water with a soft cloth.
- Allow the sensor to air dry.

Disinfection

For recommended disinfectants, see the agents with  in the 512E/512G/518B/513A/551B column of the *Applicable Disinfecting Agents List* (P/N: 046-011532-00).

- Clean the sensor as instructed above.
- Disinfect the sensor with cotton or soft cloth moistened with one of the recommended disinfectants. You can also soak the ear clip of the SpO₂ sensor in recommended disinfectants.
- After disinfection, be sure to wipe off the disinfectant left on the sensor with a soft cloth moistened with water.
- Allow the sensor to air dry.

Specifications

Model	513A
Type	Reusable, Ear clip
Patient Category	Adult/Pediatric (weight: >30kg)
Application Site	Ear lobe
Wavelength Range	600nm - 1000nm
Maximum photic output	18 mW
Cable Length	1.6m
Number of usage	1500 (restricted to cleaning and disinfection cycles)
Note: The information about the wavelength range and maximum photic output consumption can be especially useful to clinicians (for example, when photodynamic therapy is performed).	

	SpO ₂	Pulse rate
Measurement range	0% - 100%	20 - 300 bpm
Accuracy	70 to 100%: ±3% 70 to 69%: Not specified.	±3 bpm
Environmental specifications	Temperature	Relative humidity
Operating conditions	0°C - 45°C	15% - 95% (noncondensing)
Storage conditions	-30°C - 70°C	10% - 95% (noncondensing)
		57.0kPa - 107.4kPa 16.0kPa - 107.4kPa

Symbols

	Refer to instruction manual/ booklet		Not made with natural rubber latex
	Date of manufacture		Manufacturer
	Batch code		Protected against the effects of temporary immersion in water
	Caution		Authorized representative in the European Community
Rx Only	Federal law (USA) restricts this device to sale by or on the order of a physician		Comply with the requirements of Directive 2012/19/EU Waste Electrical & Electronic Equipment
UDI	Unique Device Identifier		Medical Device
	This product is provided with a CE marking in accordance with the regulations stated in Regulation (EU) 2017/745 of the European parliament concerning Medical Devices. The number adjacent to the CE marking (0123) is the number of the EU-notified body certified for meeting the requirements of the Regulation. Note: The product complies with the Council Directive 2011/65/EU.		

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