MR 400 Series

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Statement

SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD. (hereinafter called Mindray) owns the intellectual property rights to this product and this manual. Disclosure of the information in this manual in any manner whatsoever without the written permission of Mindray is strictly forbidden. This manual provides the instructions necessary to operate the product in accordance with its function and intended use. Observance of this manual is a prerequisite for proper performance and correct operation, and ensures patient and operator safety. Mindray is responsible for the effects on safety, reliability and performance of this product, only if:

- the product is used in accordance with the instructions for use.
 - this product is not damaged by human factors. Human factors refer to unintentional falling, intentional damaging, etc.
- In the event that it becomes necessary to return a unit to Mindray, please contact the Mindray Service Department and obtain a Mindray Customer Service Authorization Number. The Mindray Customer Service Authorization Number will appear on the outside of the shipping container. Return shipments will not be accepted if the Mindray Customer Service Authorization Number is not clearly visible. Please provide the model number, serial number, and a brief description of the reason for return.The customer is responsible for freight charges when this product is shipped to Mindray for service (including any relevant customs fees or other freight related charges).

mindray  **MINDRAY** and  are the registered trademarks or trademarks owned by Mindray in China.

© 2006-2022 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights reserved.

Notification of Adverse Events

As a health care provider, you may report the occurrence of certain events to SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD., and possibly to the competent authority of the Member state in which the user and / or patient is established. These events, include device-related dead and serious injury or illness. In addition, as part of our Quality Assurance Program, SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD. requests to be notified of device failures or malfunctions. This information is required to ensure that SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD. provides only the highest quality products.

Intended Use

Intended Purpose Statement

The temperature probe is intended for continuous patient temperature measurement and control applications.

Indication for Use

The temperature probe is designed for continuous patient temperature measurement and control applications. The temperature probe must be used in connection with the monitoring device (including accessories) specified by our company, including Mindray BenVision N series patient monitors. The user shall be responsible for ensuring the compatibility between the probe and any monitoring device (including accessories) by either referring to the instructions for using this probe or contacting us before use. The probes must be used by trained clinical professionals. Otherwise, product performance may be compromised or personal injury may occur.

Intended Users

The probe must be used by trained clinical professionals.

Intended Patient Population

The probe can be used in adult, pediatric and neonate populations.

Intended Medical Condition

The probe is expected to be used in medical institutions.

Contraindications

Not identified yet.

Side effects

N/A

Note:

According to the conclusion of clinical evaluation and residual risk evaluation, for the intended patients, there is no known side effects that can occur during or after the use of the medical device. And there is no need for the operator to make extra preparations. Thus, no residual risk associated with using the medical device should be disclosed due to the risk management report.

MR400 Series Probe Model

Temperature Probe			
Type	Model	Applicable to	Application site
Reusable	MR401B	Adult	Ösophageal/Rectal
	MR403B		Skin
	MR401		Ösophageal/Rectal
	MR403		Skin
	MR402B	Pediatric/Infant	Ösophageal/Rectal
	MR404B		Skin
	MR402		Ösophageal/Rectal
	MR404		Skin
Disposable	MR411	Adult/Pediatric/Infant/Neonate	Ösophageal/Rectal
	MR412	Adult/Neonate	Skin
Temperature Probe Extension Cable			
Type	Model	Applied probe	
Reusable	MR420B	MR411, MR412	
	MR420	MR411, MR412	

Warnings

- Federal law (USA) restricts this device to sale by or on the order of a physician.
- Please select temperature probes according to patient category.
- Probes are constructed with the sensors electrically isolated from the other probe surfaces, this isolation could be lost if the probe is damaged by abuse or mishandling.
- Interference of high frequency electricity may affect the accuracy of the measurement.
- When the skin probe is used, please use caution to avoid burns patients. If the probes must be used simultaneously with high-intensity RF equipment, to reduce the possibilities of patient's burning: avoid the cables connecting with the patient from kinking; do not allow the cables to contact the patient directly; the electrical grounding isolation must be provided between the instrument with which the probe is used and the RF equipment.
- Do not place the product in an environment that exceeds the stated range.
- To prevent contamination, the disposable probes shall be recovered for proper disposal.
- When the skin probe is used, please use caution to avoid burns patients. If the probes must be used simultaneously with high-intensity RF equipment, to reduce the possibilities of patient's burning: avoid the cables connecting with the patient from kinking; do not allow the cables to contact the patient directly; the electrical grounding isolation must be provided between the instrument with which the probe is used and the RF equipment.
- Do not use the temperature probes when the high-intensity sources of RF energy (e.g., electro-surgical devices and MRI equipment) are present; otherwise it may cause reading error, probe damage or local heating which burns patients.
- Mishandling of the probes could result in damage to the internal wire and the loss of electrical isolation or improper temperature readings.
- Do not use the temperature probes when the high-intensity sources of RF energy (e.g., electro-surgical devices and MRI equipment) are present; otherwise it may cause reading error, probe damage or local heating which burns patients.
- When the skin probe is used, please use caution to avoid burns patients. If the probes must be used simultaneously with high-intensity RF equipment, to reduce the possibilities of patient's burning: avoid the cables connecting with the patient from kinking; do not allow the cables to contact the patient directly; the electrical grounding isolation must be provided between the instrument with which the probe is used and the RF equipment.
- Do not place the product in an environment that exceeds the stated range.

Applying The Temperature Probe

Be sure to understand all warnings listed in the previous section before applying any probe to a patient.

Applying The Temperature Probe

- Visually inspect the probe for cracks, holes, crazing etc. prior to each use. If any such degradation in the cable jacket is discovered, discard probe according to your hospital's procedure for medical waste. In medical applications, the user must determine that a probe style is suitable and sufficiently flexible for esophageal or rectal use.
- For the Esophageal/Rectal probe: according to the specific measurements, place the probe inside the recta or esophagus.
- For the Skin probe: the auxiliary measurement is recommended; otherwise place the probe on other regions of the body surface, such as arms. Place the probe in the axilla or other regions of the body surface and fix it with adhesive plaster. Note that the metal side shall contact the patient.
- Connect the probe to the instrument with which the probe is used in terms of the operation manual of the instrument.
- The disposable probes cannot be re-disinfected or reused; the reusable probes shall be cleaned and disinfected, but note that the probes cannot be disinfected with steam.

Inspection and Recalibration

Probes cannot be "recalibrated", but should be inspected monthly by the hospital BIOMEDICAL Equipment group to ensure they are working properly. Probes can be tested by plugging into a compatible monitor and looking for an electrical open or short-circuit, intermittent readings or extremely inaccurate readings which would indicate probe wire damage. The sensor stability is well-documented; probe accuracy should not drift out of tolerance over the normal life of the probe. Dispose of any damaged probes as indicated above.

Cleaning/Disinfection

Cautions

- Never immerse or soak the probe.
- Exercise caution during cleaning/disinfection to avoid wetting the metal.
- Do not clean or disinfect the metal parts of this product. Otherwise, metal part corrosion or product reliability deterioration could result.
- We recommend that the probe be disinfected only when necessary as determined by your hospital's policy, to avoid long term damage to the probe.
- Never use cleaning agents/disinfectants other than the recommended.

Cleaning

- Clean the probe with cotton or soft cloth moistened with water.
- After cleaning, wipe off the water with a moistened cloth.
- Allow the probe to air dry.

Disinfection

For recommended disinfectants, see the agents with ☒ in the *Temp Probe* column of the *Applicable Disinfecting Agents List* (P/N: 046-011532-00).

- Clean the probe as instructed above.
- Disinfect the probe with cotton or soft cloth moistened with one of the recommended disinfectants.
- After disinfection, be sure to wipe off the disinfectant left on the probe with a soft cloth moistened with water.

- Allow the probe to air dry.

SPEC/Symbol

The length of the reusable probe is 3.0 m. The length of the disposable probe is 0.6 m. For specifications such as storage temperature, measurement accuracy, etc., refer to the Instructions for Use that accompanies the equipment used in connection with the accessory.

	Not made with natural rubber latex		Refer to instruction manual/booklet
	Manufacturer		Catalogue number
	Authorized representative in the European Community		Batch code
	Federal law (USA) restricts this device to sale by or on the order of a physician		Comply with the requirements of Directive 2012/19/EU Waste Electrical & Electronic Equipment
	Caution		Date of manufacture
	Unique Device Identifier		Medical Device

This product is provided with a CE marking in accordance with the regulations stated in Regulation (EU) 2017/745 of the European Parliament concerning Medical Devices. The number adjacent to the CE marking (0123) is the number of the EU-notified body certified for meeting the requirements of the Regulation.

Note: The product complies with the Council Directive 2011/65/EU.

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© 2006-2021 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. Todos os direitos reservados.