Sidestream CO₂/AG Accessories

Instructions for Use

Sidestream CO₂/AG Accessories

- Sampling line
- Watertrap
- Airway Adapter

Shanghai International Holding Corp. GmbH(Europe)

Eiffestraße 80, 20537 Hamburg Germany

English

SIDESTREAM CO₂/AG ACCESSORIES

INSTRUCTIONS FOR USE

Statement

SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD. (hereinafter called Mindray) owns the intellectual property rights to this product and this manual. Disclosure of the information in this manual in any manner whatsoever without the written permission of Mindray is strictly forbidden. This manual provides the instructions necessary to operate the product in accordance with its function and intended use. Observation of this manual is a prerequisite for proper performance and correct operation, and ensures patient and operator safety. Mindray is responsible for the effects on safety, reliability and performance of this product, only if:

- the product is used in accordance with the instructions for use.
- this product is not damaged by human factors. Human factors refer to unintentional falling, intentional damaging, etc.

In the event that it becomes necessary to return a unit to Mindray, please contact the Mindray Service Department and obtain a Mindray Customer Service Authorization Number. The Mindray Customer Service Authorization Number must appear on the outside of the shipping container. Return shipments will not be accepted if the Mindray Customer Service Authorization Number is not clearly visible. Please provide the model number, serial number, and a brief description of the reason for return. The customer is responsible for freight charges when this product is shipped to Mindray for service (including any relevant customs fees or other freight related charges).

mindray  **MINDRAY** and  are the registered trademarks or trademarks owned by Mindray in China.

The issue date of this manual is 2020-10.

©2013-2020 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights reserved.

Notification of Adverse Events

As a health care provider, you may report the occurrence of certain events to SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD., and possibly to the competent authority of the Member state in which the user and / or patient is established.

These events, include device-related death and serious injury or illness. In addition, as part of our Quality Assurance Program, SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD. requests to be notified of device failures or malfunctions. This information is required to ensure that SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD. provides only the highest quality products.

Product Introduction

Intended Use

Sidestream CO₂/AG sampling lines, water traps and airway adapters are intended for monitoring CO₂/anesthetic gases concentration in patient's respiratory gas.

The CO₂/AG accessories must be used in connection with the monitoring device (including accessories) specified by our company. The user shall be responsible for ensuring the compatibility between the CO₂/AG accessories and the monitoring device (including accessories) by either referring to this instructions for use or contacting us before use. The CO₂/AG accessories must be used by trained clinical professionals. Otherwise, product performance may be compromised or personal injury may occur. The accessory material that contacts the patients has undertaken the bio-compatibility test and is verified to be in compliance with ISO 10993-1.

According to the conclusion of clinical evaluation and residual risk evaluation, for the intended patients, there is no known side effects that can occur during or after the use of the medical device. And there is no need for the operator to make extra preparations. Thus, no residual risk associated with using the medical device should be disclosed due to the risk management report.

Contraindications

None.

Models

Description	Applicable to	Note	Model
Water trap	Adult/pediatric	Reusable	100-000080-00
	Neonate	Reusable	100-000081-00
Sampling line	Adult	Disposable	60-15200-00
	Neonate	Disposable	60-15300-00
Nasal sampling line (only for sidestream CO ₂)	Adult	Disposable	4000
	Pediatric	Disposable	4100
	Infant	Disposable	4200
Airway adapter	Adult/pediatric	Disposable, straight	60-14100-00
	Adult/pediatric	Disposable, elbow	60-14200-00
	Neonate	Disposable, straight	040-001187-00

Warnings


- Federal law (USA) restricts this device to sale by or on the order of a physician.
- Nasal sampling lines shall only be used for sampling from the nose.
- Check the package for damages before use.
- For disposal of accessories, follow local regulations or your hospital's policy regarding disposal of such accessories.

Applying the CO₂/AG Accessories

Be sure to understand all the precautions listed above before applying any CO₂/AG accessory to a patient.

1. Turn on the monitoring device.
2. Install the water trap by aligning it with the receptacle and push gently into place. Make sure that the locking mechanism is fully engaged by pulling the water trap, which should be firmly seated.
3. Choose proper sampling line according to patient category and applicable situation.
4. Attach one end of the sampling line to gas inlet of the water trap, while the other end to the airway adapter.

- (1) Airway connection is shown in figure 1.
- (2) For either an intubated or a non-intubated patient, the connection of the sampling line with the monitoring device is the same, as shown in figure 2.
- (3) The connection of nasal sampling line with the patient is shown in figure 3.

5. Perform CO₂/AG monitoring according to the monitoring device's operation manual.

Notes:


- ☐ Check that the sampling line is not kinked during CO₂/AG monitoring.
- ☐ Ensure that water trap is placed properly. Improper water trap installation may cause airway leakage, resulting in low CO₂/AG readings.

Cleaning/Disinfection

Disposable accessories are not designed to be reused. Never clean and disinfect the disposable accessories for reuse. No parts of the water trap are intended to be cleaned and disinfected.

Precautions for Water Trap Usage

The water trap is used to prevent the monitoring device from being contaminated by moisture, patient secretion, bacteria and dust. The container at lower part of the water trap is removable. Please empty it regularly, as the figure shows.

It is recommended to replace the water trap every month or more often if the monitoring device indicates this. Water trap may be blocked by the filter aid after long term use, resulting in inaccurate readings.

SPEC/Symbol

Maximum emptying interval during normal use (sample gas of 37 °C, room temperature of 23 °C, 100% RH) for adult and pediatric use:

- 45 h/70 ml/min
- 21 h/150 ml/min

Maximum emptying interval during normal use (sample gas of 37 °C, room temperature of 23 °C, 100% RH) for neonate use:

- 45 h/70ml/min
- 31 h/100ml/min

For specifications such as storage temperature, measurement accuracy, etc., refer to the Instructions for Use that accompanies the equipment used in connection with the accessory.

	Output; exit		Input; entrance
	Not for neonatal use.		Baby
	The product bears CE mark indicating its conformity with the provisions of the REGULATION (EU) 2017/745 on medical devices. Note: The product complies with the Council Directive 2011/65/EU.		Catalogue number
Rx ONLY	Federal law (USA) restricts this device to sale by or on the order of a physician.		Batch code
	Manufacturer		Date of manufacture
	Consult instructions for use		Authorized representative in the European Community
	Temperature limit		Humidity limitation
	Do not re-use		Contains no latex.
	Caution		Medical Device

www.mindray.com
Mindray Building, Koji 12th Road South, High-tech Industrial Park,
Nanshan, Shenzhen 518057, P.R. China
Tel: +86 755 81888998
Fax: +86 755 26582680
E-mail: service@mindray.com

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ACESSÓRIOS DE FLUXO LATERAL CO2/AG
©2013-2020 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. Todos os direitos reservados.