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
en Breathing System
Instructions for use
1. Every component of the breathing system must be visually inspected and checked for function, leakage, and occlusions, immediately before use on each patient. Ensure that connections are secure.
2. Remove the safety cap before connecting the system between the equipment and the patient.
3. When connecting or disconnecting the system, use a push and twist action.
4. If the system contains an Adjustable Pressure Limiting (APL) valve, set and adjust it as clinically required throughout the procedure. System pressure is maximised when the APL valve green cap is turned fully clockwise. The 60 cm H₂O pressure release cannot be overridden. It may be necessary to squeeze the reservoir bag to achieve 60 mbar system pressure.
5. If this system is used without a respiratory filter at the patient connection, then it must only be used on a single patient.
6. If this system is protected by an Intersurgical filter, or other appropriately validated respiratory filter, at the patient connection, and the filter and patient connection are changed for each patient, then the system may be used on more than one patient.
7. The system should be changed in accordance with the hospital department cleaning and disinfection procedures, or when it is no longer fit for its intended purpose.
8. When the system is not in use, the patient connection should be protected by the red cap supplied to prevent contamination or damage.

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INTERSURGICAL
COMPLETE RESPIRATORY SYSTEMS

Ireland
T: +353 (0)1 697 8540
info@intersurgical.ie
France
T: +33 (0)1 48 76 72 30
info@intersurgical.fr
Germany
T: +49 (0)2241 25690
info@intersurgical.de
Spain
T: +34 (0)11 665 73 15
info@intersurgical.es
Portugal
T: +351 219 108 550
info@intersurgical.pt
Italy
T: +39 0535 20836
info@intersurgical.it
Nederland
T: +31 (0)413 243 860
info@intersurgical.nl
België/Belgique
T: +32 5 2456991
info@intersurgical.be
Luxembourg
T: +32 5 2456991
info@intersurgical.lu
Sweden
T: +46 08 514 30 600
info@intersurgical.se
Denmark
T: +45 28 69 80 00
info@intersurgical.dk
Lietuva
T: +370 387 66611
info@intersurgical.lt
Russia
T: +7 495 771 6809
info@intersurgical.ru
Czech Republic
T: +420 272 911 814
info@intersurgical.cz
Türkiye
T: +90 216 468 88 28
info@intersurgical.tr
South Africa
T: +27 (0)11 444 7968
info@intersurgical.co.za
China
T: +86 519 69817000
info@intersurgical-cn.com
Japan
T: +81 (0)3 6863 4359
support@intersurgicalinc.com
Canada
T: +1 905 319 6500
info@intersurgical.ca
Colombia
T: +57 1 742 5161
info@intersurgical-ph
Australia
T: +61 (0)2 8048 3300
info@intersurgical.com.au

Distributed in the USA by Intersurgical Incorporated, 6751 Kinne Street, East Syracuse, NY 13057
T: (800) 828-9633 F: (315) 451-3696 support@intersurgicalinc.com www.intersurgical.com



cs	DÝCHACÍ SYSTÉM
hu	LÉLEGEZTETŐ RENDSZER
sl	DIHALNI SISTEM
lv	ELPOŠANAS SISTĒMA
et	HINGAMISAPARAAT
bg	СИСТЕМА ЗА АСИСТИРАНА ВЕНТИЛАЦИЯ
sr	SISTEM ZA DISANJE
ro	SISTEMUL DE RESPIRAT
sk	DÝCHACÍ SYSTÉM
hr	SUSTAV ZA DISANJE
tr	SOLUNUM DEVRESİ
ja	呼吸装置
zh-s	呼吸系統
zh-t	呼吸系統
ar	نظام التنفس

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Rx ONLY 1639

WATER TRAP, SELF SEALING, MEDIUM VOLUME

CE
1639 Rx ONLY

INTERSURGICAL
COMPLETE RESPIRATORY SYSTEMS

Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK

T: +44 (0)118 9656 300 info@intersurgical.com www.intersurgical.com

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Česká Republika
T: +420 272 911 814
info@intersurgical.cz

Türkiye
T: +90 216 468 88 28
info@intersurgical.com.tr

South Africa
T: +27 (0)11 444 7968
info@intersurgical.co.za

中国
T: +86 519 6981 7000
info@intersurgical-cn.com

日本
T: 81 (0)3 6863 4359
info@intersurgical.co.jp

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Filipinas
T: +63 2820 4124
info@intersurgical.ph

USA
T: +1 800 828 9633
support@intersurgicalinc.com

Canada
T: +1 905-319-6500
info@intersurgical.ca

Colombia
T: +57 1 742 5161
info@intersurgical.com

Australia
T: +61 (0)2 8048 3300
info@intersurgical.com.au

Distributed in the USA by Intersurgical Incorporated, 6757 Kinne St., East Syracuse, NY 13057
T: (800) 828-9633 F: (315) 451-3696 support@intersurgicalinc.com www.intersurgical.com

7705 Issue 9 RO 08.21