15mm Single Limb Breathing Systems

en	15MM SINGLE LIMB BREATHING SYSTEMS
fr	CIRCUITS RESPIRATOIRES MONOBRANCHE 15MM
de	15-MM EINSCHLAUCH-BEATMUNGSSYSTEME
es	CIRCUITOS RESPIRATORIOS DE UN RAMAL DE 15 MM
pt	SISTEMAS DE RESPIRAÇÃO DE MEMBRO ÚNICO DE 15MM
it	SISTEMA DI RESPIRAZIONE DA 15MM A LINEA SINGOLA
nl	15 MM ENKELSLANGS BEADEMINGSSYSTEMEN
no	15 MM PUSTESYSTEMER MED ETT LEM LØP
fi	15MM YHDEN LETKUN HENGITYSLETKUSTOT
sv	15 MM ENSLANGS ANDNINGSSYSTEM
da	15 MM ÅNDEDRÆTSSYSTEMER TIL ET ENKELT LEM
el	ΣΥΣΤΗΜΑΤΑ ΑΝΑΠΝΟΗΣ ΕΝΟΣ ΆΚΡΟΥ 15MM
lt	15 MM VIENOS ATŠAKOS KVĖPAVIMO SISTEMOS
pl	15 MM UKŁADY ODDECHOWE JEDNORUROWE
ru	ОДНОШЛАНГОВЫЕ ДЫХАТЕЛЬНЫЕ СИСТЕМЫ 15 MM
cs	15MM DÝCHACÍ SYSTÉMY S JEDNOU VĚTVÍ
hu	15 MM-ES EGY ÁGGAL RENDELKEZŐ LÉGZŐRENDSZEREK
sl	15 MM DIHALNI SISTEM Z ENIM UDOM
lv	15 MM VIENA TRUBAS ELPOŠANAS SISTĒMAS
et	15MM ÜHE HARULINE HINGAMISSÜSTEEM
bg	ДИХАТЕЛНИ СИСТЕМИ С ЕДИН КРАЙНИК 15 MM
sr	15MM SISTEMI ZA DISANJE SA JEDNOM GRANOM
ro	CIRCUIT VENTILAȚIE 15MM CU O TUBULATURĂ
sk	15 MM DÝCHACIE SYSTÉMY S JEDNOU VETVOU
hr	15MM SUSTAVI ZA DISANJE S JEDNIM ODVOJKOM
tr	15MM TEK HATLI SOLUNUM SİSTEMLERİ


fig. 1 / fig. 1 / abb. 1 / fig. 1 / fig. 1 / fig. 1 / aff. 1 / fig. 1 / kuva 1 / bild 1 / fig. 1 / сѣк. 1 / 1 pav. / rys. 1 / рис. 1 / obr. 1 / 1. ábra / sl. 1 / 1. att. / joonis 1 / фвр. 1 / sl. 1 / fig. 1 / obr. 1 / sl. 1 / шек. 1

Patient category ¹	Intended delivered volume ² ml	Flow resistance limit ³ cmH ₂ O/l/min	At flow ⁴ l/min	Compliance ⁵ ml/cmH ₂ O	At Pressure ⁶ cmH ₂ O
Adult ⁷	≥ 300 ml	0.06	30	5	60 ± 3
Paediatric ⁸	50 ml < 300 ml	0.12	15	4	60 ± 3
Neonatal ⁹	≤ 50 ml	0.74	2.5	1.5	60 ± 3

1 Patient category / Catégorie de patients / Patientenkategorie / Categoría de paciente / Categoria de doentes / Categoria de pazienti / Patientenkategorie / Pasientkategorii / Potilastyypii / Patientkategorii / Patientkategorii / Κατηγορία ασθενών / Pacientų kategorija / Kategoria pacjentów / Категория пациентов / Kategorie pacientů / Betegek köre / Kategorija bolnika / Pacientu kategorija / Patsiendi kategooria / Категория на пациента / Kategorija pacijenata / Categorie de pacienți / Categoría pacientov / Kategorija pacijenata / Hasta kategorisi

2 Intended delivered volume ml / Volume administré prévu ml / Zuzuführendes Volumen ml / Volumen de suministro previsto ml / Volume administrado previsto ml / Volume erogato previsto ml / Beoogd volume ml / Tiltenkt levert volum ml / Aiottu antotilavuus ml / Åvsedd levererad volym ml / Tilsigtet leveret volumen ml / Προβλεπόμενος χορηγούμενος όγκος ml / Numatomas pristatomas tūris ml / Planowana objętość dostarczana ml / Планируемый дыхательный объем ml / Předpokládaný přijatý objem ml / Leadadod légzési térfogat ml / Predvidena dovedena prostornina ml / Paredzētais piegādes tilpums ml / Kavandatud manustatav maht ml / Планиран осигуряван капацитет ml / Planirana isporučena zapremina ml / Volumul estimat livrat ml / Určený podávaný objem ml / Predviđeni isporučeni volumen ml / Amaçlanan uygulama hacmi ml

3 Flow resistance limit cmH₂O/l/min / Limite de résistance à l'écoulement cmH₂O/l/min / Grenzwert des Strömungswiderstands cmH₂O/l/min / Limite de resistencia de flujo cmH₂O/l/min / Limite de resistência ao fluxo cmH₂O/l/min / Limite della resistenza al flusso cmH₂O/l/min / Limiet stromingsweerstand cmH₂O/l/min / Grense for flow-motstand cmH₂O/l/min / Virtausvastuksen raja cmH₂O/l/min / Gräns för flödesmotstånd cmH₂O/l/min / Grænse for flow-motstand cmH₂O/l/min. / Όριο αντίστασης στη ροή cm H₂O/l/λεπτό / Pasipriešinio srautui riba cmH₂O/l/min / Granica oporu przepływu cmH₂O/l/min / Предел сопротивления потоку cm H₂O/л/мин / Limit odporu proudění cmH₂O/l/min / Áramlási ellenállás határértéke cmH₂O/l/perc / Omejevit odpornosti za pretok cmH₂O/l/min / Plūsmas pretestības robeža cmH₂O/l/min / Voolutakistuse piirväärtus cmH₂O/l/min / Граница на съпротивление на потока cmH₂O/l/min / Granična vrednost otpornosti protoku cmH₂O/l/min / Limită de rezistență la curgere cmH₂O/l/min / Limit odporu voči prúdeniu cmH₂O/l/min / Granica otpora na protok cmH₂O/l/min / Akış direnci sınırı cmH₂O/l/dakika

4 At flow l/min / Débit l/min / Bei Durchfluss l/min / Flujo l/min / Com o fluxo l/min / Al flusso l/min / Bij het stromen l/min / Ved flow l/min / Virtauksella l/min / Vid flödet l/min / Ved flow l/min. / Σε ροή l/λεπτό / Ésant srautui l/min / Przy przepływie l/min / В потоке л/мин / Při proudění l/min / Áramlási sebesség l/perc / Pri pretoku l/min / Pie aktīvas plūsmas l/min / Voolutkiirusel l/min / При поток l/min / Pri protoku od l/min / La curgere l/min / Pri prūdeni l/min / Pri protoku l/min / Akışta l/dakika

5 Compliance ml/cmH₂O / Compliance ml/cmH₂O / Compliance ml/cmH₂O / Compliance ml/cmH₂O / Complacência ml/cmH₂O / Compliancia ml/cmH₂O / Naleving ml/cmH₂O / Compliance ml/cmH₂O / Komplianssi ml/cmH₂O / Compliance ml/cmH₂O / Overensstemmelse ml/cmH₂O / Συμμόρφωση ml/cm H₂O / Atitiktis ml/cmH₂O / Podatność ml/cmH₂O / Растяжимость мл/см H₂O / Shoda ml/cmH₂O / Megfelelőség ml/cmH₂O / Komplianca ml/cmH₂O / Atbilstība ml/cmH₂O / Vastavus ml/cmH₂O / Эластичность ml/cmH₂O / Komplijansa ml/cmH₂O / Complianță ml/cmH₂O / Poddajność ml/cmH₂O / Rastezljivost ml/cmH₂O / Uygunluk ml/cmH₂O

6 At Pressure cmH₂O / Pression cmH₂O / Bei Druck cmH₂O / Presión cmH₂O / Com a pressão cmH₂O / Alla pressione cmH₂O / Önder druk cmH₂O / Ved trykk cmH₂O / Paineella cmH₂O / Vid tryck cmH₂O / Ved tryk cmH₂O / Σε πίεση cm H₂O / Ésant sléigui cmH₂O / Przy ciśnieniu cmH₂O / Под давлением cm H₂O / Pri tlaku cmH₂O / Nyomatás cmH₂O / Pri tlaku cmH₂O / Pie spiediena cmH₂O / Rõhul cmH₂O / При налягане cmH₂O / Pri pritisku od cmH₂O / La presiune cmH₂O / Pri tlaku cmH₂O / Pri tlaku cmH₂O / Basınçta cmH₂O

7 Adult / Adulte / Für Erwachsene / Adulto / Adulto / Adulto / Volwassene / Voksne / Aikuinen / Vuxen / Voksen / Ενήλικες / Suaugusieji / Dorosly / Dla взрослых / Dospělý / Felnőtt / Odrasli / Pieaugušajiem / Täiskasvanu / Възрастни / Odrasli / Adult / Dospeli / Odrasla osoba / Yetişkin

8 Paediatric / Pédiatriques / Pädiatrischen / Pediátrico / Pediátricos / Pediatrico / Pediatricisch / Barn / Lapsille / Barn / Рæдиатри́ск / Ταϊδριακοῦς / Vajkai / Dzieci / Для пæдиатрических / Pediatrické / Gyermek / Pediatrični / Pediatrijas / Pediatriiline / Пæдиатрична / Pediatrijske / Copii / Pediatrických / Pediatriju / Pedyatrik

9 Neonatal / Nouveau-nés / Neonatalen / Neonatal / Recém-nascidos / Neonatale / Neonataal / Nyfødte / Vastasyntyneille / Nyfödda / Neonatal / Νεογνά / Naujagimiai / Noworodki / Для новорожденных / Novorozenci / Újszülött / Neonatalni / Jaundzimusiem / Vastsūndinute / Неонатална / Neonatalne / Nou-născuți / Novorodenci / Novorođenčad / Yenidoğan

en 15MM SINGLE LIMB BREATHING SYSTEMS

The instruction for use is not supplied individually. Only 1 copy of the instructions for use is provided per box and should therefore be kept in an accessible location for all users. More copies are available upon request.

NOTE Distribute instructions to all product locations and users. These instructions contain important information for the safe use of the product. Read these instructions for use in their entirety, including warnings and cautions before using this product. Failure to properly follow warnings, cautions and instructions could result in serious injury or death to the patient. Serious incidents to be reported as per the local regulation.

WARNING A WARNING statement provides important information about a potentially hazardous situation, which, if not avoided could result in death or serious injury.

CAUTION A CAUTION statement provides important information about a potentially hazardous situation, which if not avoided could result in minor or moderate injury.

INTENDED USE To deliver and remove respiratory gases to and from a patient via a breathing system comprised of tubing and connectors.

INTENDED USER For use by qualified medical personnel.

INDICATIONS The breathing system should be connected between the ventilator or the respiratory equipment, and the patient airway connection. The product should only be used on a single patient up to a maximum of 7 days. The breathing system may be provided in Flextube™ and Smoothbore tubing.

CONTRAINDICATIONS No known contraindications.

PRODUCT USE The product is intended for a maximum of 7 days use with a variety of ventilator and respiratory equipment meeting international standards. The start of use is defined from the first connection to the equipment. Use in the acute and non-acute environment under the supervision of suitably or appropriately qualified medical personnel. The breathing system is safe for use with atmospheric gases.

PATIENT CATEGORY The range of 15mm breathing systems are suitable for use on all paediatric patients with intended delivered tidal volumes of 50ml < 300ml.

CONNECT THE BREATHING SYSTEM Attach the breathing system to the inspiratory port of the ventilator or respiratory equipment using a push and twist action (fig. 1). Remove any red dust cap from the patient end of the breathing system. Perform pre-use checks as recommended by the equipment manufacturer. Attach the patient end to the patient interface.

ACCESSORIES Where the 15mm breathing system is supplied with additional components the clinician should assess any potential impact upon the patient's respiratory function, and be aware of the subsequent increased resistance to flow, compressible volume and weight of the product. When the ventilator and respiratory equipment parameters are set, any impact that the internal volume and resistance to flow of the device may have upon the effectiveness of the prescribed ventilation, should be assessed on an individual patient basis by the clinician.

WATER TRAPS AND BREATHING FILTERS, HMES, HMEFS (IF INCLUDED IN THE CONFIGURATION) Follow the instructions in individual IFUs supplied.

MONITORING LINES (IF INCLUDED IN THE CONFIGURATION)

Monitoring lines provide a leak free connection between the breathing system and a respiratory gas monitor via a luer lock fitting. The connection should be checked prior to use. The line should be checked to confirm it is kink free and a patent pathway exists.

PRE-USE CHECK

1. Prior to attachment, check that all breathing system components are free of any obstruction or foreign bodies and that there is a patent airway.
2. Remove any red dust cap before connecting the breathing system between the equipment and the patient.
3. Connect/disconnect the tapered connections of the breathing system to the tapered connection of the ventilator and/or respiratory equipment using a push and twist action (fig. 1).
4. Following attachment the breathing system must be checked for leaks and occlusions prior to use and that all connections are secure.
5. When setting ventilator parameters consider any impact that the internal volume of the breathing system and accessories may have upon the effectiveness of the prescribed ventilation. The clinician should assess this on an individual patient basis.
6. Check that any port cap or accessory is securely attached to prevent leakage.

IN-USE CHECK The system should be changed if visibly contaminated or if resistance to flow increases. Monitor the clinical parameters of the patient during the treatment.

WARNING

1. This system is not suitable for use with flammable anaesthetic agents.
2. Monitor the breathing system throughout use and replace if increased resistance or contamination occurs.
3. Refer to equipment manufacturers guidelines for machine settings

and pre-use check requirements.

4. Following use, the product must be disposed of in accordance with local hospital, infection control in line with the local healthcare, hygiene and waste disposal protocols and regulations.

5. Perform a self-test of the equipment after full assembly of the breathing system and associated accessories before use on each patient.

6. Prior to use, consult the respective documentation and instructions for use of all connected devices, or combination of devices.

7. Do not stretch, cut or reconfigure the tubing.

8. Do not reuse. Reuse poses a threat to patient safety in terms of cross infection and device malfunction and breakage. Do not clean. Cleaning can lead to device malfunction and breakage.

9. Ventilator accuracy can be affected by use of a gas driven nebuliser.

10. Adding attachments or other devices to the breathing system can change the pressure gradient across the breathing system and adversely affect the ventilator performance. Perform a self-test of the ventilator (or respiratory equipment) after full assembly of the breathing system and all connected devices before use on each patient.

11. Breathing systems, their parts and accessories are validated for use with specific ventilators and respiratory equipment.

Incompatible parts can result in degraded performance. The responsible organisation is responsible for ensuring the compatibility of the ventilator or respiratory equipment and all of the parts used to connect to the patient before use.

12. Use of this breathing system with equipment not compliant with appropriate International Standards can result in degraded performance. Prior to use, consult the respective documentation and instructions for use of all connected devices, or combination of devices to ensure compatibility.

13. Consult individual product labelling for presence or absence of phthalates.

14. If a luer port connector is included in the system, it is for the exclusive use of monitoring respiratory and anaesthetic gases. If applicable, check that the cap is securely attached to prevent inadvertent leakage when the port is not in use. Do not introduce gases or liquids via the misconnection of the luer lock port as this can result in serious injury to the patient.

15. Do not use sterile products if the packaging is damaged, refer to the product label for information.

16. Do not use product past expiry date.

17. Ensure that the breathing system technical and performance data are appropriate to the intended patient and treatment.

18. Patient ventilation and vital signs should be monitored. Please refer to product label for full information.

CAUTION

1. The breathing system is only to be used under medical supervision by suitably or appropriately qualified medical personnel.

2. To avoid contamination, the product should remain packaged until ready for use. Do not use the product if the packaging is damaged.

3. When connecting or disconnecting the system, hold connector and use a push and twist action.

4. When the breathing system is not in use, the patient connection should be protected by the red dust cap supplied to prevent contamination.

FOR USA: Rx ONLY: Federal law restricts this device to sale by or on the order of a physician.

Working pressure: <60 cmH₂O

STORAGE CONDITIONS Store out of direct sunlight. Recommended storage at room temperature over the indicated product shelf-life.

Environmental Operating Conditions: Performance verified at room/ambient temperature at sea level.

DISPOSAL Following use, the product must be disposed of in accordance with local hospital, infection control and waste disposal regulations.

APPROPRIATE STANDARDS This product complies with relevant requirements of the following standards:

- BS EN ISO 5367: Anaesthetic and respiratory equipment. Breathing tube and sets
- BS EN ISO 5356-1: Anaesthetic and respiratory equipment. Conical connectors. Cones and sockets

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