- | | |
|---|---|
| <p>1 CapStrap Headgear
Harnais CapStrap
Dispositivo di fissaggio sul capo
CapStrap
CapStrap Maskenhalterung
Arnés CapStrap</p> <p>2 Four Position Forehead Arm
Bras frontal à quatre positions
Montante del sostegno per la fronte a quattro posizioni
Stirnpolsterarm mit vier Einstellmöglichkeiten
Brazo de soporte con cuatro posiciones</p> <p>3 Blue Standard Elbow (SE)
Raccord coudé bleu standard (SE)
Gomito standard blu (SE)
Blaues gekrümmtes</p> | <p>Standardanschlussstück (SE)
Codo estándar azul (SE)</p> <p>4 Pressure Pick-off Port
Orifice de prise de pression
Porta di pressione
Druckventilanschluss
Orificio de salida de la presión</p> <p>5 Mask Faceplate
Plaque du masque
Flangia della maschera
Gesichtsplatte
Soporte facial de la mascarilla</p> <p>6 Top Headstrap
Harnais supérieur
Cinghia di fissaggio superiore
Oberes Halterungsband
Arnés superior</p> |
|---|---|


Figure 1
Figure 1
Figure 1
Abbildung 1
Figura 1

- | | |
|--|--|
| <p>1 Four Point Headgear
Harnais à quatre points
Dispositivo di fissaggio sul capo a quattro punti
Vierpunkt-Maskenhalterung
Arnés con cuatro posiciones</p> <p>2 Four Position Forehead Arm
Bras frontal à quatre positions
Montante del sostegno per la fronte a quattro posizioni
Stirnpolsterarm mit vier Einstellmöglichkeiten
Brazo de soporte con cuatro posiciones</p> <p>3 Blue Standard Elbow (SE)
Raccord coudé bleu standard (SE)
Gomito standard blu (SE)
Blaues gekrümmtes
Standardanschlussstück (SE)
Codo estándar azul (SE)</p> <p>4 Pressure Pick-off Port
Orifice de prise de pression
Porta di pressione
Druckventilanschluss
Orificio de salida de la presión</p> | <p>5 Mask Faceplate
Plaque du masque
Flangia della maschera
Gesichtsplatte
Soporte facial de la mascarilla</p> <p>6 Strap Connection to Forehead Bracket
Connexion de la sangle à l'appui frontal
Collegamento della cinghia al sostegno per la fronte
Bandbefestigung an Stirnpolsterhalterung
Conexión de la correa a la abrazadera para la frente</p> <p>7 Crown Strap - starts at top of head and around back of head.
Courroie de tête – commence en haut de la tête et entoure l'arrière de la tête.
Cinghia a corona: parte dalla sommità del capo e scende dietro il capo.
Scheitelband – beginnt oben am Kopf und verläuft am Hinterkopf.
Correa de la coronilla: comienza en la parte superior de la cabeza y rodea la parte posterior de la misma.</p> |
|--|--|



Figure 2
Figure 2
Figure 2
Abbildung 2
Figura 2

Specifications / Caractéristiques techniques / Specifiche / Technische Daten / Especificaciones

Deadspace volume: < 350 cc (S), < 525 cc (M), < 650 cc (L)
Volume d'espace mort : < 350 cm³ (S), < 525 cm³ (M), < 650 cm³ (L)
Volume spazio morto: < 350 cc (S), < 525 cc (M), < 650 cc (L)
Totraumvolumen: < 350 cc (S), < 525 cc (M), < 650 cc (L)
Volumen de espacio muerto: < 350 cc (S), < 525 cc (M), < 650 cc (L)

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Intended Use

Medium and Large Size: The AF531 Oro-Nasal SE Face Mask is intended to provide a patient interface for application of noninvasive ventilation. The mask is to be used as an accessory to ventilators which have adequate alarms and safety systems for ventilation failure, and which are intended to administer CPAP or positive pressure ventilation for treatment of respiratory failure, respiratory insufficiency or obstructive sleep apnea. The mask is for single use in the hospital/institutional environment only. The mask is to be used on patients (>30 kg) who are appropriate candidates for noninvasive ventilation.

Small Size: The AF531 Oro-Nasal SE Face Mask is intended to provide a patient interface for application of noninvasive ventilation. The mask is to be used as an accessory to ventilators which have adequate alarms and safety systems for ventilation failure, and which are intended to administer CPAP or positive pressure ventilation for treatment of respiratory failure, respiratory insufficiency or obstructive sleep apnea. The mask is for single use in the hospital/institutional environment only. The mask is to be used on patients 7 years or older (>20 kg) who are appropriate candidates for noninvasive ventilation.

Note: An exhalation port is not built into the AF531 Oro-Nasal SE Face Mask. A separate exhalation device must be used with this mask.

Note: This mask does not contain natural rubber latex or DEHP.

Caution: U.S. Federal law restricts this device to sale by or on the order of a physician.

Warnings:

- This mask is not suitable for providing life support ventilation.
- This mask requires a separate exhalation device.
- This mask is designed for use with CPAP or bi-level systems recommended by your health care professional or respiratory therapist. Do not wear this mask unless the CPAP or bi-level system is turned on and operating properly. **Do not block or try to seal the exhalation port. Explanation of the Warning:** CPAP systems are intended to be used with special masks with connectors which have vent holes to allow continuous flow of air out of the mask. When the CPAP machine is turned on and functioning properly, new air from the CPAP machine flushes the exhaled air out through the attached mask exhalation port. However, when the CPAP machine is not operating, enough fresh air will not be provided through the mask, and exhaled air may be rebreathed. This warning applies to most models of CPAP systems. Rebreathing of exhaled air for longer than several minutes can, in some circumstances, lead to suffocation.
- If oxygen is used with the device, the oxygen flow must be turned off when the device is not operating. **Explanation of the Warning:** When the device is not in operation, and the oxygen flow is left on, oxygen delivered into the ventilator tubing may accumulate within the device enclosure. Oxygen accumulated in the device enclosure will create a risk of fire.
- Oxygen supports combustion. Oxygen should not be used while smoking or in the presence of an open flame.
- At a fixed flow rate of supplemental oxygen flow, the inhaled oxygen concentration will vary, depending on the pressure settings, patient breathing pattern, mask selection, and the leak rate. This warning applies to most types of CPAP and bi-level machines.
- Some users may experience skin redness, irritation, or discomfort. If this happens, discontinue use and contact your healthcare professional.
- The patient's physician should be contacted if the patient experiences the following symptoms while using the masks or after removing it: Unusual chest discomfort, shortness of breath, stomach distension, belching, or severe headache; drying of the eyes, eye pain, or eye infections; blurred vision. (Consult an ophthalmologist if symptoms persist.)
- At low CPAP or EPAP pressures, the flow through the exhalation port may be inadequate to clear all exhaled gas from the tubing. Some rebreathing may occur.
- This mask should not be used for patients who are uncooperative, obtunded, unresponsive, or unable to remove the mask.
- This mask is not recommended for patients who are taking a prescription drug that may cause vomiting.
- Appropriate patient monitoring should be used as medically necessary.
- If an additional exhalation device is added to the patient circuit, you may need to adjust the pressure level to compensate for the additional leak of the exhalation device.
- This mask does not include an anti-asphyxia valve that would allow a patient to breathe if the ventilator failed.
- Hand wash prior to first use. Inspect the mask for damage or wear (cracking, crazing, tears, etc). Discard and replace any components as necessary.
- This mask is for single use only and is not intended to be disinfected or cleaned after patient use. Dispose of the mask after patient use. Multi-patient use without disinfection between patients increases the risk of infection.
- This mask is not designed for reuse. Respiration cannot guarantee the performance specifications once this mask has been disinfected or cleaned after patient use.
- Use of a nasal or full face mask may cause tooth, gum, or jaw soreness or aggravate an existing dental condition. Consult your physician or dentist if symptoms occur.

Contraindications

This mask may not be suitable for patients with the following conditions: hiatal hernia, impaired cardiac sphincter function, excessive reflux, impaired cough reflex; or patients who are unable to remove the mask by themselves.

Before Use

- Read and understand the instructions completely.
- Hand wash the mask.
- Clean the patient's face.
- Verify that the mask and headgear are the correct size.
- If a nasogastric tube or similar device is in place, use the optional nasogastric tube sealing pad. Position the pad so that its flat surface is against the patient's face and the C-shaped opening surrounds the tube.
- Verify that the therapy device, i.e., ventilator, including the alarms and safety systems, has been validated prior to use.
- Check for the blue elbow on the mask, and note that the mask does not include an anti-asphyxia valve that would allow a patient to breathe if the ventilator failed.
- Inspect the mask and replace it if the cushion has hardened or is torn, or if any parts are broken.
- Verify therapy device pressure(s).

Achieving the Right Fit

- Hold the mask lightly against the patient's face and slide the headgear over their head.
 - If using the CapStrap Headgear, connect the bottom clips to the mask with a gentle push (Figure 1).
 - If using the Four Point Headgear strap, connect the bottom clips to the mask with a gentle push (Figure 2).
 - Disengage the headgear tabs from the bottom clips and tighten the straps to achieve a comfortable fit with minimal leak.
- Adjust the Four Position Forehead Arm to ensure minimal mask contact pressure on the nasal bridge while achieving minimal mask cushion leak.
 - The lowest position provides the highest contact pressure at the bridge of the nose and provides the best mask cushion seal.
 - The highest position provides the lowest contact pressure at the bridge of the nose and provides lower mask cushion seal.
- Connect the exhalation device to the blue standard elbow (SE).
- Connect the tubing that is included with the therapy device to the exhalation device.
- Unfasten the headgear tabs on the top strap and gradually adjust to ensure the mask is tight enough to achieve minimal mask cushion leak and a comfortable fit. DO NOT OVERTIGHTEN. Adjust the crown strap as necessary to achieve proper mask support.
- Turn on the noninvasive ventilator. Instruct the patient to breathe normally.
- Make final adjustments to all adjustment points to lessen leaks and to provide a comfortable fit.

Note: Re-adjust the strap tension if leaks occur as the patient's position changes.

Do not overtighten the straps. Overtightening can cause or worsen leaks.

Connecting and Disconnecting the Mask and Headgear

The AF531 Oro-Nasal SE Face Mask can be used with either the CapStrap Headgear or the Four Point Headstrap. Use the following instructions to disconnect and reconnect the headgear to the mask.

Removing the CapStrap Headgear

First, disconnect one or both of the bottom clips. Next, gently twist the CapStrap hook so that it disconnects from the forehead support bracket (Figure 3).

Replacing the CapStrap Headgear

Press the CapStrap hook into the back of the slot on the forehead support bracket and press until they snap together. On the lower straps, press the bottom clips into place on the faceplate.

Removing the Four Point Headstrap

First, undo the bottom clips. Next, open the headgear tabs and pull them out of the slots on the forehead support arm.

Replacing the Four Point Headstrap

On the top strap, thread the headgear tabs through the slots in the forehead support arm (Figure 4). Fold the tabs back so that they catch onto the headgear straps.

Removing the Mask and Headgear from the Patient

Disconnect the bottom clips and swing the mask upward (Figure 5). This allows easy access to the patient. To remove completely, pull the entire assembly over the patient's head.

Replacing the Mask Elbow

To remove: Hold the faceplate hub and gently pull and twist the elbow from the faceplate hub.

To replace: Press the mask elbow into the faceplate hub and slightly twist while holding the faceplate hub. Do not use excessive force or any tools.

Disposal

Dispose of in accordance with local regulations.


Figure 3

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©2011 Koninklijke Philips Electronics N.V.