
RESPIRONICS

Respironics Inc.
1001 Murry Ridge Lane
Murrysville, PA 15668 USA

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PerforMax Pediatric EE Leak 2 Total Face Mask

Intended Use

The PerforMax Pediatric Total Face Mask is intended to provide an interface for application of non-invasive positive airway pressure delivered to patients by a device such as a CPAP or bi-level system. The mask is for multi-patient use in the hospital/institutional environment only. The mask is to be used on patients 1 year or older (> 7 kg) for whom positive airway pressure therapy has been prescribed.

 **Note:** An exhalation port is built into the mask so a separate exhalation port is not required.

 **Note:** This mask is not made with natural rubber latex or DEHP.

 **Caution:** U.S. Federal law restricts this device to sale by or on the order of a physician.

Symbols

Warning or Caution

Note

Tip

Not made with
natural rubber latex

Leak Symbol
and Value

Warnings

- Patients using this mask must remain under constant supervision by trained medical personnel.
- This mask is not suitable for providing life support ventilation.
- Hand wash prior to use. Inspect the mask for damage or wear (cracking, crazing, tears, etc). Discard and replace any components as necessary.
- This mask is designed for use with CPAP or bi-level systems recommended by your health care professional or respiratory therapist. Do not wear this mask unless the CPAP or bi-level system is turned on and operating properly. Do not block or try to seal the exhalation port. Explanation of the Warning: CPAP systems are intended to be used with special masks with connectors which have vent holes to allow continuous flow of air out of the mask. When the CPAP machine is turned on and functioning properly, new air from the CPAP machine flushes the exhaled air out through the attached mask exhalation port. However, when the CPAP machine is not operating, enough fresh air will not be provided through the mask, and exhaled air may be rebreathed. Rebreathing of exhaled air for longer than several minutes can in some circumstances lead to suffocation. This warning applies to most models of CPAP systems.
- If oxygen is used with the device, the oxygen flow must be turned off when the device is not operating. Explanation of the Warning: When the device is not in operation, and the oxygen flow is left on, oxygen delivered into the ventilator tubing may accumulate within the device enclosure. Oxygen accumulated in the device enclosure will create a risk of fire.
- At a fixed flow rate of supplemental oxygen flow, the inhaled oxygen concentration will vary, depending on the pressure settings, patient breathing pattern, mask selection, and the leak rate. This warning applies to most types of CPAP and bi-level machines.

- At low CPAP or EPAP pressures the flow through the exhalation port may be inadequate to clear all exhaled gas from the tubing. Some rebreathing may occur.
- Some patients may experience skin redness, irritation, or discomfort. If this happens, monitor and intervene.
- Monitor and intervene if any of the following symptoms occur: Unusual chest discomfort, shortness of breath, stomach distension, belching, severe headache, blurred vision, drying of the eyes, eye pain or eye infections.
- This mask is not recommended if the patient is taking a prescription drug that may cause vomiting.
- Monitor and intervene if the patient encounters tooth, gum, or jaw soreness. Use of a mask may aggravate a patient's existing dental conditions.
- A minimum of 3cm H₂O pressure must be maintained when using this mask.
- The mask contains small parts which could result in a choking hazard.
- Attaching an exhalation device requires therapy pressure level adjustment to compensate for increased leak.
- Do not overtighten the headgear straps. Watch for signs of overtightening, such as excessive redness, sores or bulging skin around the edges of the mask. Loosen the headgear straps to alleviate symptoms.
- Do not block or seal off the anti-asphyxia valve.


Figure 1
 1 Mask Cushion
 2 Faceplate Hub
 3 Entrainment Elbow with Anti-Asphyxia Valve and Exhalation (Do not block)
 4 Pressure Pick-Off Port
 5 Bonnet
 6 Top Headgear Straps
 7 Bottom Headgear Clips



Contraindications

This mask may not be suitable for use on patients with the following conditions: glaucoma, recent eye surgery or dry eyes, impaired nictitation (blinking), hiatus hernia, impaired cardiac sphincter function, esophageal reflux, impaired cough reflex; or on patients who are not under constant, immediate observation by personnel experienced in noninvasive ventilation management.

Verify the Anti-Asphyxia Valve

The anti-asphyxia valve consists of an air inlet and a flapper. With the airflow turned off, verify that the elbow with the flapper is lying flat (Figure 2A) so that room air can flow in and out through the air inlet. Next, with the airflow on, the flapper should now cover the air inlet and air from the CPAP or bi-level device should flow into the mask (Figure 2B). If the flapper does not close or does not function properly, replace.

⚠ Warning: Do not block or seal the anti-asphyxia valve.

Before Use Read and Understand the Instructions Completely

- Hand wash the mask.
- Wash the patient's face. Do not use moisturizer/lotion on your hands or the patient's face.
- Inspect the mask and replace it if the cushion has hardened or is torn, or if any parts are broken.
- Verify that the therapy device, i.e., ventilator, including the alarms and safety systems, has been validated prior to use.
- Verify therapy device pressure(s).

Pre-cleaning Instructions

Hand wash the mask before first use.

1. Hand wash mask in warm water with liquid dishwashing detergent or wipe with a 70% v/v isopropyl alcohol swab.

2. Rinse thoroughly. Air dry completely before use. Make sure the mask is dry before use.

⚠ Caution: Do not use bleach, cleaning solutions containing bleach, or cleaning solutions containing conditioners or moisturizers.

⚠ Caution: Any deviation from these instructions may impact the performance of the product.

⚠ Caution: Inspect the mask for damage or wear (cracking, crazing, tears, etc). Discard and replace any components as necessary.

Institutional Disinfection

For multi-patient use in the hospital/institutional environment, use the Disinfection Guide to reprocess the PerformMax Pediatric mask between patients. These instructions can be obtained by visiting us online at www.philips.com/NIVmasks or by contacting customer service at 1-800-345-6443 or at 1-724-387-4000.

Leak Symbol and Port Settings

Some ventilators may incorporate the use of a leak symbol and value in the mask selection setup procedures. The leak characteristics of this mask is leak symbol **2**. The leak symbol and value represents the intentional leak characteristics of the interface. On ventilators equipped with a Mask Selection control, enter the leak symbol value **2** that corresponds with the leak symbol value on the mask.

* Comfort Tips

- The most common mistake is overtightening the headgear. The headgear should fit loose and comfortable. If the patient's skin bulges around the mask or if you see red marks on the patient's face, loosen the headgear.
- Adjust the top headgear straps to reduce leaks at the forehead and temples.
- Adjust the bottom headgear straps to reduce leaks at the sides of the patient's mouth.

PerformMax Interchangeable Leak 1 and Leak 2 Elbows

This mask may use the EE Leak 1 and EE Leak 2 elbows interchangeably. The clear EE Leak 1 elbow with an anti-asphyxia valve does not have built in exhalation. A separate exhalation device must be used with the EE Leak 1 elbow. The amber EE Leak 2 elbow with an anti-asphyxia valve includes exhalation. Refer to **F** for Disassembly and **H** for Assembly of the Mask to change the elbow on the mask.

Achieving the Right Fit - Refer to figures A - H on the back cover.

Applying the Mask - A and B

Using the Mask - C

Removing the Mask - D


Disassembly of the Mask - E and F

Assembly of the Mask - G and H

Specifications

⚠ Warning: The technical specifications of the mask are provided to determine if it is compatible with your CPAP or bi-level therapy device. If used outside these specifications, or if used with incompatible devices, the mask may be uncomfortable, the seal of the mask may not be effective, optimum therapy may not be achieved, and leak, or variation in the rate of leak, may affect device function.

Pressure Flow Curve

Resistance with Anti-Asphyxia Valve Closed to Atmosphere

Drop in Pressure at

	50 SLPM	100 SLPM
All sizes	0.7 cm H ₂ O	2.1 cm H ₂ O

Inspiratory and Expiratory Resistance with Anti-Asphyxia Valve Open to Atmosphere

Inspiratory Resistance	1.2 cm H ₂ O
Expiratory Resistance	0.7 cm H ₂ O

Deadspace

XXS	255 mL
XS	284 mL

Sound Levels

A-weighted Sound Power Level	32.7 dBA
A-weighted Sound Pressure Level at 1 m	24.1 dBA

Disposal

Dispose of in accordance with local regulations.

Storage Conditions

Temperature: -4° to 140° F (-20° to +60° C)
Relative Humidity: 15% to 95% non-condensing

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©2012 Koninklijke Philips Electronics N.V. All rights reserved.