Instructions for Use

PerforMax SE Total Face Mask Single Use (SU)

Intended Use

The PerforMax SE Total Face Mask SU is intended to provide a patient interface for application of noninvasive ventilation. The mask is to be used as an accessory to ventilators which have adequate alarms and safety systems for ventilator failure, and which are intended to administer CPAP or positive pressure ventilation for treatment of respiratory failure, respiratory insufficiency, or obstructive sleep apnea. The mask is for single use in the hospital/institutional environment only. The small size mask is to be used on patients 7 years or older (>20 kg) who are appropriate candidates for noninvasive ventilation. The large and extra large size masks are to be used on patients (>30 kg) who are appropriate candidates for noninvasive ventilation.

 **Important Notes:** An exhalation port is not built into the PerforMax SE SU Total Face Mask. A separate exhalation device must be used with this mask.

 **Note:** This mask does not contain natural rubber latex or DEHP.

 **Caution:** U.S. Federal law restricts this device to sale by or on the order of a physician.

Symbols

 *Warning or Caution*  *Note*  *Tip*  *Does Not Contain Natural Rubber Latex*  *Consult Instructions for Use*

Warnings:

- This mask is not suitable for providing life support ventilation.
- Hand wash prior to use. Inspect the mask for damage or wear (cracking, crazing, tears, etc). Discard and replace any components as necessary.
- This mask requires a separate exhalation device.
- This mask is designed for use with CPAP or bi-level systems recommended by your health care professional or respiratory therapist. Do not wear this mask unless the CPAP or bi-level system is turned on and operating properly. **Do not block or try to seal the exhalation port. Explanation of the Warning:** CPAP systems are intended to be used with special masks with connectors which have vent holes to allow continuous flow of air out of the mask. When the CPAP machine is turned on and functioning properly, new air from the CPAP machine flushes the exhaled air out through the attached mask exhalation port. However, when the CPAP machine is not operating, enough fresh air will not be provided through the mask, and exhaled air may be rebreathed. This warning applies to most models of CPAP systems.
- This mask does not include an anti-asphyxia valve that would allow a patient to breathe if the ventilator failed.
- If oxygen is used with the device, the oxygen flow must be turned off when the device is not operating. **Explanation of the Warning:** When the device is not in operation, and the oxygen flow is left on, oxygen delivered into the ventilator tubing may accumulate within the device enclosure. Oxygen accumulated in the device enclosure will create a risk of fire.
- Oxygen supports combustion. Oxygen should not be used while smoking or in the presence of an open flame.
- At a fixed flow rate of supplemental oxygen flow, the inhaled oxygen concentration will vary, depending on the pressure settings, patient breathing pattern, mask selection, and the leak rate. This warning applies to most types of CPAP and bi-level machines.

- System pressure and oxygen therapy may be different for a Respironics PerforMax mask than for a nasal or nasal/oral interface. When changing mask types, determine the appropriate pressure level and oxygen therapy.
- At low CPAP or EPAP pressures the flow through the exhalation port may be inadequate to clear all exhaled gas from the tubing. Some rebreathing may occur.
- If an additional exhalation device is added to the patient circuit, you may need to adjust the pressure level to compensate for the additional leak of the exhalation device.
- Some users may experience skin redness, irritation, or discomfort. If this happens, discontinue use and contact your healthcare professional.
- The patient's physician should be contacted if the patient experiences the following symptoms while using the masks or after removing it: Unusual chest discomfort, shortness of breath, stomach distension, belching, or severe headache; drying of the eyes, eye pain, or eye infections; blurred vision. (Consult an ophthalmologist if symptoms persist.)
- This mask should not be used on patients who are uncooperative, obtunded, unresponsive, or unable to remove the mask.
- This mask is not recommended if the patient is taking a prescription drug that may cause vomiting.
- Appropriate patient monitoring should be used as medically necessary.
- The mask contains small parts which could result in a choking hazard.
- Use of a nasal or full face mask may cause tooth, gum, or jaw soreness or aggravate an existing dental condition. Consult your physician or dentist if symptoms occur.
- This mask is for single use only and is not intended to be disinfected or cleaned after patient use. Dispose of the mask after patient use. Multi-patient use without disinfection between patients increases the risk of infection.
- This mask is not designed for reuse. Respironics cannot guarantee the performance specifications once this mask has been disinfected or cleaned after patient use.

Contraindications

This mask may not be suitable for use on patients with the following conditions: glaucoma, recent eye surgery or dry eyes, hiatal hernia, excessive reflux, impaired cough reflex, impaired cardiac sphincter function; or on patients unable to remove the mask by themselves.

Before Use

- Read and understand the instructions completely.
- Hand wash the mask.
- Clean the patient's face.
- Verify that the mask and headgear are the correct size.
- If a nasogastric tube or similar device is in place, use the optional nasogastric tube sealing pad. Position the pad so that its flat surface is against the patient's face and the C-shaped opening surrounds the tube.
- Verify that the therapy device, i.e., ventilator, including the alarms and safety systems, has been validated prior to use.
- Check for the blue elbow on the mask, and note that the mask does not include an anti-asphyxia valve that would allow a patient to breathe if the ventilator failed.
- Inspect the mask and replace it if the cushion has hardened or is torn, or if any parts are broken.
- Verify therapy device pressure(s).

Cleaning Instructions

1. Hand wash mask in warm water with a mild liquid dishwashing detergent.

 **Caution:** Do not use bleach, alcohol, cleaning solutions containing bleach or alcohol, or cleaning solutions containing conditioners or moisturizers.

2. Rinse thoroughly. Air dry completely before use.

⚠️ Caution: Do not attempt to clean the headstrap. Replace if the headstrap becomes soiled.

Sizing the Mask

To select the appropriate size mask, the mask cushion should surround the patient's face without obstructing the patient's sight.

Figure 1-A: The top of the cushion should rest comfortably against the middle of the forehead above the eye brows.

Figure 1-B: The bottom of the cushion should rest comfortably under the mouth and above the chin.

Connecting the Headgear and Mask

Disconnecting: Undo the top headgear tabs (Figure 2-B), and pull the straps through the slots on the top of the faceplate. Disengage the bottom clips from the mask (Figure 2-D) with a gentle tug. Do not remove the straps from the clips. This will help with the orientation of the headgear when connecting.

Connecting: Position the four-point headgear with the crown strap near the top. Thread the headgear tabs through the slots on the top of the faceplate and fold them back to catch onto the headgear straps. Connect the bottom clips to the mask with a gentle push.

Achieving the Right Fit

⚠️ Warning: Check for the blue elbow on the mask, and note that the mask does not include an anti-asphyxia valve that would allow the patient to breathe if the ventilator failed.

1. Adjust the headgear to a large size setting and disconnect one (or both) of the bottom clips. Slide the headgear and mask on over the patient's head. Hold the mask over the patient's face by the faceplate protrusion. Do not hold the mask by the elbow.
2. Slide the headgear into position. The crown strap will lay at the crown of the patient's head (Figure 2-B). The cross-strap will sit low on the back of the patient's head. Connect the bottom clips and tighten the straps using the headgear tabs. Do not overtighten. Final adjustments will be made when the patient is lying down.
3. The top of the cushion will contact the face just above the eyebrows (Figure 2-A). Adjust the top strap but do not overtighten. The bottom of the cushion will contact the face just above the chin. The entire cushion should fit snugly, but comfortably against the patient's face. Re-evaluate the mask size if you cannot achieve facial contact points described in the *Sizing the Mask* section.
4. The patient should be lying down or the bed should be adjusted to a reclining position (Figure 3).
5. Connect the flexible dual-limb tubing from the noninvasive ventilator or therapy device to the 22 mm blue standard elbow (SE).
6. With the prescribed CPAP or IPAP/EPAP pressure applied, turn on the airflow. Verify there is only minor lifting of the mask from the patient's face (6 mm to 12 mm) during a full respiratory cycle. Have the patient assume various positions (supine, lateral, prone) with the mask in place.
7. Make any final adjustments using the headgear tabs.

Removing the Headgear and Mask

If the patient is lying flat or against the bed, first disconnect one (or both) bottom clips from the mask and slide the mask upward. Gently pull the headgear free of the patient. If the patient can sit up, disconnect one or both of the bottom clips and remove the entire mask and headgear.


Figure 1

Figure 2

A = Mask Cushion
B = Crown Strap
C = Top Headgear Strap
D = Bottom Clip
E = Blue Standard Elbow (SE)


Figure 3

Removing or Replacing the Mask Elbow (SE)

To remove: Hold the faceplate and gently pull and twist the elbow from the faceplate hub (Figure 4-E).

To replace: Press the mask elbow into the faceplate hub and gently twist. Do not use excessive force or any tools. Make sure the elbow is firmly in place and will not separate from the hub during routine use.

Disposal

Dispose of in accordance with local regulations.

Specifications

Deadspace Volume

Small = 375 mL

Large = 550 mL

Extra Large = 717 mL


Storage Conditions

Temperature: -20° to +60° C)

Relative Humidity: 15% to 95% non-condensing

Reorder Numbers

Visit www.philips.com/performax for a complete list of products and reorder information, or contact Philips Respironics Customer Service at 1-800-345-6443 (USA or Canada) or at 1-724-387-4000 or Respironics Deutschland at +49 8152 93060 for additional information.

 **Figure 4**

A = Slot for Headgear Strap

B = Mask Cushion

C = Headgear Clip Post

D = Faceplate

E = Faceplate Hub

F = Blue Standard Elbow (SE)

(No built-in exhalation vent or antiasphyxia valve)

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0123

AUSTRALIAN SPONSOR DETAILS:

Philips Electronics Australia Ltd.
65 Epping Road, North Ryde, NSW 2113
Australia

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© 2014 Koninklijke Philips N.V.